

Pierce[®] NeutrAvidin[®] High Binding Capacity Coated 384-Well Plates

15511 15512 15513

0923.2

Number	Description
15511	Pierce NeutrAvidin High Binding Capacity Coated Plates (clear, 384-well), 5 each
15512	Pierce NeutrAvidin High Binding Capacity Coated Plates (white, 384-well), 5 each
15513	Pierce NeutrAvidin High Binding Capacity Coated Plates (black, 384-well), 5 each

Blocking Buffer: Supplied blocked with SuperBlock[®] Blocking Buffer
Binding Capacity: ~35pmol D-biotin/well
Coating Volume: 50µL
Blocking Volume: 100µL

Storage: Upon receipt store plates at 4°C in unopened pouches. Once opened, place unused plates in a resealable bag with desiccant and store at 4°C.

Introduction

The Thermo Scientific Pierce NeutrAvidin Coated Plates are ideal for capturing biotin-labeled molecules without interference from nonspecific binding. NeutrAvidin Protein is deglycosylated avidin, which reduces lectin binding to undetectable levels while retaining stability and biotin-binding affinity. NeutrAvidin Protein offers the advantages of a near-neutral pI (6.3), to minimize nonspecific adsorption, and the lack of the RYD sequence, which eliminates nonspecific binding to the RGD binding domain of adhesion receptors present in a variety of cells. NeutrAvidin Protein yields the lowest nonspecific binding among the known biotin-binding proteins. The clear, white and black plates can be used with colorimetric, chemiluminescent and fluorescent detection methods, respectively.

Our patented coating technology is used to create high binding capacity (HBC) Pierce NeutrAvidin Coated Plates. These HBC plates provide a wider detection range and better curve linearity for small biotinylated ligands, such as peptides and oligonucleotides, than the standard coated plates.

Example ELISA Protocol using Pierce NeutrAvidin Coated Plates

A. Materials Required

- Wash Buffer: Tris-buffered saline (25mM Tris, 150mM NaCl; pH 7.2; Product No. 28376), 0.1% BSA, 0.05% Tween[®]-20 Detergent; alternatively, use Thermo Scientific Blocker BSA (Product No. 37520) supplemented with 0.05% Tween-20 Detergent
- Biotinylated capture antibody adjusted to 10µg/mL, or other appropriate concentration, with Wash Buffer
- Antigen adjusted to appropriate concentration with Wash Buffer
- Primary antibody adjusted to appropriate concentration with Wash Buffer
- Enzyme-labeled secondary antibody adjusted to appropriate concentration with Wash Buffer
- Appropriate enzyme substrate: example substrates are the Thermo Scientific TMB Substrate Kit (Product No. 34021) for horseradish peroxidase and the Phosphatase Substrate Kit (Product No. 37620) for alkaline phosphatase

B. Method

1. Wash each well three times with 100µL of Wash Buffer. Add 50µL of the biotinylated capture antibody to each well and incubate for 2 hours at room temperature.

2. Wash each well three times with 100µL of Wash Buffer. Make a serial dilution of the antigen and add 50µL to each well. Incubate plate for 30 minutes at room temperature.
3. Wash each well three times with 100µL of Wash Buffer. Add 50µL of the primary antibody to each well and incubate plate for 30 minutes at room temperature.
4. Wash each well three times with 100µL of Wash Buffer. Add 50µL of the enzyme-labeled secondary antibody to each well. Incubate plate for 30 minutes at room temperature.
5. Wash each well three times with 100µL of Wash Buffer.
6. Follow the manufacturer's instructions for the specific detection system.

Related Thermo Scientific Products

37070	SuperSignal® ELISA Pico Chemiluminescent Substrate , 100mL, peroxidase substrate
15169	QuantaBlu™ Fluorogenic Peroxidase Substrate Kit
34028	1-Step™ Ultra TMB-ELISA , 250mL, colorimetric peroxidase substrate
37621	1-Step PNPP , 100mL, colorimetric phosphatase substrate
29339	Biotinylated Alkaline Phosphatase , 1mg
29139	Biotinylated Horseradish Peroxidase , 5mg
15075	Reagent Reservoirs , 200/pkg
15036	Sealing Tape for 96-Well Plates , 100/pkg
15507	Pierce NeutrAvidin Coated High Binding Capacity Plates (clear 96-well) , 5 each
15508	Pierce NeutrAvidin Coated High Binding Capacity Plates (clear 8-well strips) , 5 each
15509	Pierce NeutrAvidin Coated High Binding Capacity Plates (white 96-well) , 5 each
15510	Pierce NeutrAvidin Coated High Binding Capacity Plates (black 96-well) , 5 each

General References

- Denlinger, L.C., *et al.* (2001). Cutting Edge: The nucleotide receptor P2X7 contains multiple protein- and lipid-interaction motifs including a potential binding site for bacterial lipopolysaccharide. *J Immunol* **167**:1871-6.
- Ferre-Aubineau, V., *et al.* (1995). Colorimetric microtiter plate hybridization assay using monoclonal antibody for detection of an amplified human immunodeficiency virus target. *J Virol Meth* **55**:145-51.
- Hiller, Y. *et al.* (1987). Biotin binding to avidin. Oligosaccharide side chain not required for ligand association. *Biochem J* **248**:167-71.
- Holmstrom, K., *et al.* (1993). A highly sensitive and fast non-radioactive method for detection of polymerase chain reaction products. *Anal Biochem* **209**:278-283.
- Simon, M.D., *et al.* (2004). A phage display selection of engrailed homeodomain mutants and the importance of residue Q50. *Nucl Acid Res* **32**(12):3623-31.

Tween is a trademark of Croda International Plc.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2011 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.