

## Thermo Scientific μDrop Plate

The Thermo Scientific μDrop Plate is designed for quick and easy measurement of low sample volumes down to 2 μL. It is an ideal tool for photometric DNA or RNA quantitation and purity analysis.



*The μDrop Plate enables measuring low-sample volumes down to 2 μl and cuvettes*

### Features

- dsDNA detection range from a few micrograms to a few milligrams per milliliter
- 16 sample positions for quick and easy measurement of sample volumes down to 2 μL
- Compatible with an 8-channel pipette for easy sample addition
- Quick and easy to wipe off the samples in serial measurements
- Ready-made Thermo Scientific SkanIt Software sessions for nucleic acid analysis
- Measurement slot for a rectangular cuvette


The μDrop™ Plate low-volume measurement area consists of two quartz slides: the top clear quartz and the bottom Teflon-coated quartz slide. Sample positions are arranged in a 2 x 8 matrix, providing a straightforward way of analyzing up to 16 samples simultaneously.

The fixed light path of the μDrop Plate allows direct calculation of the nucleic acid concentrations of the samples.

Samples are easy to pipette onto the μDrop Plate with a single or an 8-channel pipette. The plate is easily wiped clean, making it convenient to be used in serial measurements.

The μDrop Plate additionally has a slot for a cuvette measurement. It enables endpoint, spectral and kinetic measurements in a standard rectangular cuvette with a stopper.

The μDrop Plate has been designed for use with the Thermo Scientific Multiskan GO UV/Vis microplate and cuvette spectrophotometer, as well as with the Thermo Scientific Varioskan Flash spectral scanning multimode reader. This extends their capabilities to include micro-scale absorbance measurements.



Ready-made SkanIt® Software sessions for DNA analysis



## Thermo Scientific $\mu$ Drop Plate

Specification	
Wavelength range	200–1000 nm
Overall dimensions	127.8 mm (W) x 85.5 mm (D) x 14.5 mm (H) [5.0" (W) x 3.4" (D) x 0.57" (H)]
Weight	0.161 kg / 0.4 lbs
Compatible readers	Multiskan <sup>®</sup> GO, Varioskan <sup>®</sup> Flash
Low-volume function	
Measurement area	16 positions in a 2 x 8 matrix
Sample volume	2 $\mu$ l to 10 $\mu$ l
Cuvette function (Cuvette with stopper)	
Dimensions	12.5 (W) x 12.5 (D) x 47–54 (H) mm
Beam center height	8.5 mm
Beam window	$\geq$ 4 mm

Ordering information	
Product code	Description
N12391	$\mu$ Drop Plate

Contact Information
<a href="http://www.thermoscientific.com/platereaders">www.thermoscientific.com/platereaders</a>

### North America:

USA/Canada  
+1 603 595 0505  
USA toll free  
800 345 0206

### Europe:

Austria +43 1 801 40 0  
Belgium +32 53 73 42 41  
Finland +358 9 3291 0200  
France +33 2 2803 2000  
Germany national toll free 08001-536 376  
Germany international +49 6184 90 6940  
Italy +39 02 95059 552  
Netherlands +31 76 571 4440  
Nordic/Baltic countries +358 9 329 100  
Russia/CIS +7 (495) 739 76 41  
Spain/Portugal +34 93 223 09 18  
Switzerland +41 44 454 12 12  
UK/Ireland +44 870 609 9203

### Asia:

Australia +613 9757 4474  
China +86 21 6865 4588 or +86 10 8419 3588  
China toll free 800-810-5118, 400-650-5118  
India +91 22 6716 2200  
Japan +81-3-5826-1616  
Korea +82 11 796 7771  
Other Asian countries +65 6872 9717

### Countries not listed:

+49 6184 90 6940 or +33 2 2803 2000