

Pierce™ Peroxidase Detection Kit

Pub. No. MAN0011508

Rev. B.0

Pub. Part No. 2161493

36000

Number	Description
36000	<p>Pierce™ Peroxidase Detection Kit, contains sufficient reagents to stain 250-500 tissue sections</p> <p>Kit Contents:</p> <p>DAB/Metal Concentrate (10X), 25 mL, store at -20°C (this solution will not freeze at -20°C)</p> <p><u>Peroxidase Detection Reagent Pack:</u></p> <p>Universal Blocker™ Blocking Buffer in TBS, 250mL, store at 4°C</p> <p>Peroxidase Suppressor, 2 × 100 mL, store at 4°C</p> <p>Stable Peroxide Buffer, 250 mL, store at 4°C</p> <p>BupH™ Tris Buffered Saline, 4 packs (each pack makes 500 mL), store at room temperature</p> <p>Surfact-Amps® 20 (10% Tween®-20 Detergent), 10 mL, store at 4°C</p> <p>Harris Modified Hematoxylin (without acetic acid, mercury-free), 100mL, store at room temperature</p> <p>Mounting Medium, 120 mL, store at room temperature</p>

Storage: Upon receipt store individual components as indicated above. The DAB/Metal Concentrate is shipped with dry ice. The Peroxidase Detection Reagent Pack is shipped at ambient temperature.

Table of Contents

Introduction	1
Important Product Information	2
Procedure for Staining Paraffin-Embedded Tissue Section Slides	2
Troubleshooting	4
Additional Information	4
Related Thermo Scientific Products	5
References	5

Introduction

The Thermo Scientific™ Pierce™ Peroxidase Detection Kit is a complete system for staining, counterstaining and preserving histological and cytological samples using the Metal Enhanced DAB Substrate for colorimetric detection. Immunohistochemical staining makes it possible to visualize the distribution and localization of specific cellular components within a cell or tissue. Detection is achieved using a primary antibody that is specific for the antigen and a horseradish peroxidase-conjugated secondary antibody. Alternatively, a biotinylated secondary antibody and HRP-labeled avidin or streptavidin may be used. The enzyme label is reacted with the Metal Enhanced DAB substrate to yield an intensely colored product that can be analyzed by light microscopy. The Metal Enhanced DAB is a special formulation of cobalt chloride and nickel chloride that produces a dark brown/black precipitate in the presence of HRP. Hematoxylin, a purple nuclear stain, is also included in this kit to counterstain tissue sections, which may improve contrast between background and the colored reaction products.

Important Product Information

- Store the DAB/Metal Concentrate at or below -20°C. To use the solution, remove the quantity required for use and immediately return the bottle to -20°C. Do not allow the solution to equilibrate to room temperature. This solution is packaged under nitrogen for long-term stability; **after every use**, replace the nitrogen by gently bubbling a slow stream of nitrogen into the solution.
- The Stable Peroxide Buffer contains the optimal concentration of hydrogen peroxide. Do not add hydrogen peroxide to the Stable Peroxide Buffer as this will increase background. Storage at -20°C or repeated freeze/thawing will not affect activity of the Stable Peroxide Buffer. For convenience, store the product at 4°C to eliminate the need to thaw.
- To minimize potential microbial contamination, carefully handle the kit reagents and use ultrapure water in all solutions.
- Do not use sodium azide as a preservative for buffers. Sodium azide inhibits HRP activity and could interfere with detection.
- Discard diluted and used solutions along with excess buffer after use.
- Use a humidity chamber set at 20-25°C for all incubations to prevent evaporation. Additionally, completely cover the tissue section with solution during incubations to prevent drying.
- Avoid touching slides and do not allow dust or other debris to contaminate samples, tissues or other material.
- Use acid-cleaned slides. Do not use egg albumin-dipped slides because trace avidin present in the albumin may cause background.
- Adjust the standard protocol when antigen concentrations are too high or low. High antigen concentrations will require less incubation time to obtain optimal staining. To shorten incubation times, increase incubation temperature to 37°C.
- An ABC complex system such as the Thermo Scientific™ ABC Standard Peroxidase Staining Kit (Product No. 32020) or Ultra-Sensitive ABC Standard Peroxidase Staining Kit (Product No. 32050) can increase sensitivity if necessary.

Procedure for Staining Paraffin-Embedded Tissue Section Slides

NOTE: Please read all instructions before beginning the procedure.

A. Materials Required

- Wash Buffer (Tris-buffered saline with 0.05% Tween-20™ Detergent): Dissolve contents of one BupH™ TBS Pack in 500 mL of ultrapure water. Add 2.5mL of 10% Tween-20™ Detergent and mix well. Store reconstituted buffer at 4°C.
- Primary Antibody: Dilute antibody with Universal Blocker Blocking Buffer. For best results, empirically determine the optimal dilution for each specific tissue/antigen type being tested. Typical primary antibody dilutions range from 1:5 to 1:100. Amplification methods, such as the ABC Complex, will typically require less primary antibody.
- Biotinylated Secondary Antibody (an HRP-conjugated secondary antibody also may be used): Dilute antibody with Universal Blocker Blocking Buffer. For best results, empirically determine the optimal dilution for each specific system being tested. Typical secondary antibody dilutions range from 1:100 to 1:1000.
- Avidin-HRP or Streptavidin-HRP: Dilute conjugate with Universal Blocker Blocking Buffer. For best results, empirically determine the optimal dilution for each specific system being tested. Typical conjugate dilutions range from 1:500 to 1:1000.

Note: If using an HRP-conjugated secondary antibody, avidin-HRP or streptavidin-HRP is not needed.

- Coverslips

B. Method

1. Prepare paraffin-embedded tissue section slides according to standard protocols. (See the Additional Information Section for protocol to de-paraffinize embedded tissues.)
2. Quench endogenous peroxidase activity by incubating tissue for 30 minutes in Peroxidase Suppressor.
Note: Omit this step if endogenous activity is not a problem or if the antigen will not survive exposure to H₂O₂.
3. Wash slide two times for 3 minutes with Wash Buffer.
4. Add blocking buffer to the slide and incubate for 30 minutes.
5. Blot excess blocker from the tissue sections. Apply the Primary Antibody and incubate tissue for 30 minutes.
6. Wash slide two times for 3 minutes with Wash Buffer.
7. Apply the Biotinylated Secondary Antibody and incubate tissue for 30 minutes. If using an HRP-conjugated secondary antibody, skip to step 10.
8. Wash slides three times for 3 minutes each with Wash Buffer.
9. Incubate the tissue section with the Avidin/Streptavidin-HRP for 30 minutes.
10. Wash slides three times for 3 minutes each with Wash Buffer.
11. Remove the DAB/Metal Concentrate (10X) from -20°C storage and mix well by inverting the bottle. Remove quantity required for use and immediately return bottle to -20°C.
12. Prepare a 1X working solution of the DAB/Metal Concentrate (10X) by adding the Stable Peroxide Buffer and mixing well. For example, to prepare 5 mL of substrate, add 4.5 mL of the Stable Peroxide Buffer to 500 µL of the DAB/Metal Concentrate.
Note: The 1X substrate solution is stable for several hours at 4°C.
13. Add the Metal Enhanced DAB Substrate Working Solution to the tissue and incubate until the desired staining is achieved. Typical incubations are from 5 to 15 minutes.
Note: Thick sections may require longer staining times.
14. Wash slide two times for 3 minutes each with Wash Buffer. Counterstain if desired (See Section C).
15. Mount slide with Mounting Media (See Section D).
Note: The time required to perform this protocol may be reduced with a slight loss in sensitivity by performing the following alterations:
 - Shorten the incubation steps with the primary antibody, the biotinylated secondary antibody and the streptavidin-HRP to 10 minutes each.
 - Instead of incubating washes, rinse tissue with a gentle stream of Wash Buffer several times at each wash step.
 - Because the primary antibody is prepared in Blocking Buffer, combine the incubation steps with the Blocking Buffer and the Primary Antibody.
 - Increase the concentration of the primary antibody and decrease the incubation times to 2 minutes. Optimization of the concentrations may be required.

C. Counterstaining

1. After immunostaining, rinse the slides with distilled water and drain.
2. Add an adequate amount of the hematoxylin stain to the slide to cover the entire tissue surface.
3. Incubate for 1-2 minutes at room temperature.
4. Drain off the hematoxylin and wash slide several times with distilled water.
5. Wash slide in Wash Buffer for 1 minute.
6. Wash slide with distilled water.
7. Mount slide with Mounting Media.

D. Mounting

1. Make sure tissue sections are dry. Place all slides to be mounted on a flat surface and apply 1-3 drops of Mounting Medium to the tissue section.
2. Apply a coverslip on top of the Mounting Medium.
3. Place the slides on a flat surface to set for 1-2 hours at room temperature.

Troubleshooting

Problem	Possible Cause	Solution
Precipitate is brown instead of black/brown	Cobalt and nickel are heavy metals and will separate during storage	Mix by inverting the bottle before use to obtain a homogeneous solution of DAB and metals
Background is dark, reducing the signal-to-noise ratio	DAB/Metal Concentrate was not maintained at -20°C	Store the DAB/Metal Concentrate at or below -20°C to prevent background problems
The DAB/Metal Concentrate contains a precipitate in the bottle	DAB/Metal Concentrate was not maintained at -20°C	If the precipitate does not go into solution upon mixing, do not use the substrate
High background	Too much HRP in the system	Use less antibody – this substrate is 50 times more sensitive than DAB without metals and requires much less antibody for detection

Additional Information

Procedure for De-paraffinizing Embedded Tissues

1. Incubate slides for 5 minutes in clean xylene;* repeat once.
2. Incubate slides for 5 minutes in 100% ethanol; repeat once.
3. Incubate slides for 5 minutes in 95% ethanol: 5% water; repeat once.
4. Incubate slides for 5 minutes in 85% ethanol: 20% water.
5. Incubate slides for 5 minutes in 70% ethanol: 20% water.
6. Incubate slides for 5 minutes in 50% ethanol: 20% water.
7. Incubate slides for 5 minutes in 30% ethanol: 70% water.
8. Incubate slides for 5 minutes in ultrapure water; repeat once.
9. The slides are now ready for immunohistochemistry. Do not allow slides to dry. Store slides in wash buffer until ready to use.

*For researchers who prefer to avoid organic solvents, use Fisherbrand™ Citrisolv™ Clearing Agent as a substitute.

Related Thermo Scientific Products

35000	Peroxidase Suppressor, 100mL
28376	BupH™ Tris Buffered Saline Packs, 40 packs
28320	Surfact-Amps™-20 (10% Tween-20™ Detergent), 6 × 10mL
34065	Metal Enhanced DAB Substrate Kit, contains 10X Metal Enhanced DAB concentrate and 1X Stable Peroxide Buffer
32020	ABC Standard Peroxidase Staining Kit, contains Avidin and Biotinylated HRP
32050	Ultra-Sensitive ABC Standard Peroxidase Staining Kit, contains Avidin and Biotinylated HRP
32052	Ultra-Sensitive ABC Peroxidase Mouse IgG Staining Kit, contains Avidin, Biotinylated Peroxidase, Blocking Buffer, Biotinylated Affinity-purified Horse Anti-Mouse IgG and mixing bottles
32054	Ultra-Sensitive ABC Peroxidase Rabbit IgG Staining Kit, contains Avidin, Biotinylated Peroxidase, Blocking Buffer, Biotinylated Affinity-purified Goat Anti-Rabbit IgG and mixing bottles

Metal Enhanced DAB References

- Chan, J.R., *et al.* (2000). Progesterone synthesized by Schwann cells during myelin formation regulates neuronal gene expression. *Mol Biol Cell* **11**:2283-95.
- Itoh, H. (2000). Upregulation of HGF activator inhibitor type 1 but not type 2 along with regeneration of intestinal mucosa. *Am J Physiol Gastrointest Liver Physiol* **278**:635-43.
- Kaushic, C., *et al.* (2000). Effects of estradiol and progesterone on susceptibility and early immune responses to *Chlamydia trachomatis* infection in the female reproductive tract. *Infect Immunity* **7**:4207-16.
- Martin, M.N. and Slovin, J.P. (2000). Purified γ -glutamyl transpeptidases from tomato exhibit high affinity for glutathione and glutathione S-conjugates. *Plant Physiol* **122**:1417-26.
- Scidmore, M.A., *et al.* (2003). Restricted fusion of *Chlamydia trachomatis* vesicles with endocytic compartments during the initial stages of infection. *Infect Immunity* **2**:973-84.

General References

- Adams, J.C. (1977). Technical considerations on the use of horseradish peroxidase as a neuronal marker. *Neuroscience* **2**:141-5.
- Adams, J.C. (1981). Heavy metal intensification of DAB-based HRP reaction product. *J Histochem Cytochem* **29**:775.
- De Jong, A.S.H., *et al.* (1985). Sensitivity of various visualization methods for peroxidase and alkaline phosphatase activity in immunoenzyme histochemistry. *J Histochem* **17**:1119-30.
- Graham, R.C. and Karnovsky, M.J. (1966). The early stages of absorption of injected horseradish peroxidase in the proximal tubules of mouse kidney: ultrastructural cytochemistry by a new technique. *J Histochem Cytochem* **14**:291-302.
- Hsu, S-M., *et al.* (1981). A comparative study of the peroxidase-anti-peroxidase method and avidin-biotin complex method for studying polypeptide hormones with radioimmunoassay antibodies. *Am J Clin Pathol* **75**:734-8.
- Hsu, S-M., *et al.* (1981). Use of avidin-biotin-peroxidase complex (ABC) in immunoperoxidase techniques: A comparison between ABC and unlabelled antibody (PAP) procedures. *J Histochem Cytochem* **29**:577-80.
- Hsu, S. and Soban, E. (1982). Color modification of diaminobenzidine (DAB) precipitation by metallic ions and its application for double immunohistochemistry. *J. Histochem. and Cytochem.* **30**:1079-82.

Tween is a trademark of Croda Americas LLC.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at www.thermofisher.com. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2021 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.