MinION Mk1D

Quick start guide

Pre-installation

This Quick start guide contains everything you need to set up your MinION™ Mk1D and to check that the device is ready for use.

Before using the device, familiarise yourself with the following:

IT requirements

community.nanoporetech.com/to/minion-it

MinION Mk1D user manual


community.nanoporetech.com/to/minmk1d

Safety and regulatory information

community.nanoporetech.com/to/safety

For detailed information and troubleshooting, view the user manual.

MinION Mk1D

What's in the box:

MinION Mk1D

Configuration Test Cell

USB-C cable*

*Only USB-C ports are supported. USB-A ports are not supported.

1

Install and log in to MinKNOW

1. Check your computer meets the IT requirements. (see first inside page).
2. Install the latest **MinKNOW™** for the **MinION Mk1D** and **P2 Solo** from community.nanoporetech.com/downloads. Once installed, open MinKNOW.
3. Log in to MinKNOW using your Oxford Nanopore account details.

2

Set up your MinION Mk1D

Connect the MinION to your computer using the [USB-C cable](#)*.


USB Type-C


*Only USB-C ports are supported. USB-A ports are not supported.

3

Insert your CTC

1. Lift the lid of the MinION.

2. Insert the Configuration Test Cell (CTC) **a**. Ensure you press down firmly on the CTC **b**. Then close the MinION lid.

4

Perform a

Hardware check

A [hardware check](#) is needed before your first sequencing run. To run a [hardware check](#), follow the instructions below:

1. Click [Start](#) in the left side bar in MinKNOW.
2. Click [hardware check](#).
3. Click [Start](#) in the bottom right corner of MinKNOW. Wait for the hardware check to complete.

4

Perform a Hardware check

4. Check the flow cell position in MinKNOW shows **Hardware check passed** *.
5. Remove the CTC after you're finished.

*If the hardware check fails, ensure your CTC is fully inserted, then repeat the hardware check. Use the details in Additional information in this guide to contact technical support if your hardware check fails again.

5

Discover the Nanopore Community

Discover the Nanopore Community to find protocols, forums, support and downloads for your experiments.
community.nanoporetech.com

Learn how to analyse your nanopore data at:
nanoporetech.com/analyse

Additional information

Warranty

A license and warranty can be purchased for your device here:

store.nanoporetech.com/device-warranty.html

Flow cell warranty:

community.nanoporetech.com/to/warranty

Recycle used flow cells

Oxford Nanopore is committed to environmental sustainability. You can help by sending your flow cells for recycling.

Find out how:

community.nanoporetech.com/support/returns

Place your next order

Buy more consumables at the Oxford Nanopore Store:

store.nanoporetech.com

Documentation

Documentation for your device is available on the Nanopore Community:

nanoporetech.com/documentation

Support

For all of your customer and technical support needs, visit:

community.nanoporetech.com/support

Technical specification

Model number	MIN-101D (MIN-101D-ELY is used in ElysION MinION)		
Supply voltage (V)	5 DC	Weight (g)	130
Maximum rated current (A)	1.5	Installation ports	1 x USB-C
Maximum rated power (W)	7.5	Software installed	MinION driver
Size (H x W x D) (mm)	13 x 55 x 125	Compute specification	N/A
Environmental conditions	Designed to sequence in environmental temperatures of +10°C to +35°C. Functional range of electronics is within environmental temperatures of +5°C to +40°C. 5 cm clearance around the device is needed (for ventilation and temperature control) when in use. Use within 30%–75% relative non-condensing humidity limits. Intended for indoor use. Can be used up to altitudes of 2,000 m. The device has a Pollution Degree 2.		

Oxford Nanopore Technologies

phone +44 (0)845 034 7900

email support@nanoporetech.com

X @nanopore

Gosling Building, Edmund Halley Road,
Oxford Science Park, Oxfordshire, OX4 4DQ, UK

www.nanoporetech.com

Oxford Nanopore Technologies, the Wheel icon, MinION, and MinKNOW are registered trademarks or the subject of trademark applications of Oxford Nanopore Technologies plc in various countries. Information contained herein may be protected by patents or patents pending of Oxford Nanopore Technologies plc. All other brands and names contained are the property of their respective owners. © 2025 Oxford Nanopore Technologies plc. All rights reserved. Oxford Nanopore Technologies products are not intended for use for health assessment or to diagnose, treat, mitigate, cure, or prevent any disease or condition.

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