

Pierce[®] Antibody Coated 96-Well Plates

0687.3

Number	Description
15134	Pierce Goat Anti-Mouse IgG Coated Plates (clear, 96-well), 5 plates
15234	Pierce Goat Anti-Mouse IgG Coated Plates (white, 96-well), 5 plates
15334	Pierce Goat Anti-Mouse IgG Coated Plates (black, 96-well), 5 plates Binding Capacity: ~7pmol IgG/well
15135	Pierce Goat Anti-Rabbit IgG Coated Plates (clear, 96-well), 5 plates
15136	Pierce Goat Anti-Rabbit IgG Coated Plates (white, 96-well), 5 plates
15137	Pierce Goat Anti-Rabbit IgG Coated Plates (black, 96-well), 5 plates Binding Capacity: ~12pmol IgG/well

Note: Plates are activated to 100µL and supplied pre-blocked with SuperBlock[®] Blocking Buffer.

Storage: Upon receipt store plates at 4°C in unopened pouches. Once opened, place unused plates in a resealable bag and store desiccated at 4°C.

Introduction

The Thermo Scientific Pierce Antibody Coated Plates are ideal for binding assays when available antibodies are in low quantities or denature upon direct adsorption to polystyrene plates.¹⁻⁴ Because these plates are IgG-specific, purified antibodies are not required. Pierce Antibody Coated Plates may be used for direct, indirect, competitive or sandwich assays and are available in clear for colorimetric assays, white for chemiluminescent assays, and black for fluorescent assays.

Note: The given protocols below are general ELISA procedures. Specific applications require optimization.

General Procedure for Direct ELISA

A. Materials Required

- Antigen-specific antibody
- Biotin-labeled antigen
- Dilution Buffer: Thermo Scientific StartingBlock T20 (PBS) Blocking Buffer (Product No. 37539), contains blocking buffer with 0.05% Tween[®]-20 Detergent
- Wash Buffer: Tris-buffered Saline (TBS, Product No. 28376) or Phosphate Buffered Saline (PBS, Product No. 28374) containing 0.05% Tween-20 Detergent (Product No. 28320)
- Detecting Enzyme: Streptavidin or Thermo Scientific NeutrAvidin Biotin-Binding Protein labeled with HRP (Product No. 21124 or 31001) or other enzyme
- Enzyme Substrate (e.g., for HRP use Thermo Scientific TMB Substrate Kit, Product No. 34021)

B. Method

Note: Plates are supplied pre-blocked with Thermo Scientific SuperBlock Blocking Buffer.

1. Rinse each well three times with 200µL of Wash Buffer.
2. Dilute antibody to 0.1-1µg/mL with Dilution Buffer and add 100µL to each well.

3. Incubate plate for 30-60 minutes at room temperature. For all incubation steps use a plate mixer at a speed that creates a vortex in each well.
4. Rinse each well three times with 200µL of Wash Buffer. Add 100µL of the biotinylated antigen to each well (approximately 0.1µg/mL). Incubate plate at 37°C for 1 hour.
5. Rinse each well three times with 200µL of Wash Buffer. Add 100µL of the Detecting Enzyme at 0.1-1µg/mL and incubate at 37°C for 1 hour.
6. Rinse each well three times with 200µL of Wash Buffer.
7. Develop and evaluate plate according to the instructions for the specific enzyme substrate.

General Procedure for Indirect ELISA

A. Materials Required

- Matched antibody pair: antigen-specific capture antibody and antigen-specific enzyme-labeled antibody
- Antigen
- Dilution Buffer: StartingBlock T20 (PBS) Blocking Buffer (Product No. 37539), contains blocking buffer with 0.05% Tween-20 Detergent
- Wash Buffer: Tris Buffered Saline (TBS, Product No. 28376) or Phosphate Buffered Saline (PBS, Product No. 28374) containing 0.05% Tween-20 Detergent (Product No. 28320)
- Enzyme Substrate (e.g., for HRP use TMB Substrate Kit, Product No. 34021)

B. Method

Note: Plates are supplied pre-blocked with SuperBlock Blocking Buffer.

1. Rinse each well three times with 200µL of Wash Buffer.
2. Dilute the capture antibody to 0.1-1µg/mL with Dilution Buffer and add 100µL to each well.
3. Incubate plate for 30-60 minutes at room temperature. For all incubation steps use a plate mixer at a speed that creates a vortex in each well.
4. Rinse each well three times with 200µL of Wash Buffer. Add 100µL (~0.1µg/mL) of antigen to each well. Incubate plate at 37°C for 1 hour.
5. Rinse each well three times with 200µL of Wash Buffer. Add 100µL (0.1-1µg/mL) of the enzyme-labeled detection antibody and incubate at 37°C for 1 hour.
6. Rinse each well three times with 200µL of Wash Buffer.
7. Develop and evaluate plate according to the instructions for the specific enzyme substrate.

Related Thermo Scientific Products

37070	SuperSignal® ELISA Pico Chemiluminescent Substrate , 100mL, peroxidase substrate
15169	QuantaBlu™ Fluorogenic Peroxidase Substrate Kit
34028	1-Step™ Ultra TMB-ELISA , 250mL, colorimetric peroxidase substrate
37621	1-Step PNPP , 100mL, colorimetric phosphatase substrate
15075	Reagent Reservoirs , 200/pkg
15082	Microtube Racked System , 960 tubes
15036	Sealing Tape for 96-Well Plates , 100/pkg
21425	EZ-Link® Sulfo-NHS-Biotinylation Kit
21335	EZ-Link Sulfo-NHS-LC-Biotin , 100mg, biotinylation reagent with 22.4Å spacer arm
21126	Streptavidin, Horseradish Peroxidase Conjugated , 1mg

Cited References

1. Cantarero, L.A., *et al.* (1980). The adsorptive characteristics of protein for polystyrene and their significance in solid-phase immunoassays. *Anal Biochem* **105**:375-82.
2. Dierks, S.E., *et al.* (1986). Altered recognition of surface-absorbed compared to antigen-bound antibodies in the ELISA. *Mol Immunol* **23**:403.
3. Shields, M.J., *et al.* (1991). An appraisal of polystyrene- (ELISA) and nitrocellulose-based (ELIFA) enzyme immunoassay systems using monoclonal antibodies reactive toward antigenically distinct forms of human C-reactive protein. *J Immunol Meth* **141**:253-61.
4. Schwab, C. and Bosshard, H.R. (1992). Caveats for the use of surface-absorbed protein antigen to test the specificity of antibodies. *J Immunol Meth* **147**:125-34.

Tween is a trademark of Croda International Plc.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2011 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.