

# Pierce™ Protein L Coated Plates

15190

0779.2

**Number****Description**

15190

**Pierce Protein L Coated Plates**, 5 plates

**Note:** Plates are activated to 100µL and supplied pre-blocked with SuperBlock™ Blocking Buffer. Use these plates for single-antibody assays only. In multiple antibody assays, such as sandwich ELISAs, the first antibody cannot block all binding sites and subsequent antibodies will bind to the plate, resulting in false positives.

**Storage:** Upon receipt store plates at 4°C in unopened pouches. Place opened, unused plates in a resealable bag and store desiccated at 4°C. Plates are shipped at ambient temperature.

**Introduction**

The Thermo Scientific Pierce Protein L Coated Plates are useful for detecting antibodies using labeled antigens. The pre-coated plates can increase assay sensitivity and are ideal when a low amount of antibody is available. Antibodies bind to Protein L rapidly and dissociate with strong chaotropes and extreme pH.

Protein L is an immunoglobulin-binding protein that was isolated from the bacteria *Peptostreptococcus magnus* and is now produced recombinantly. Protein L binds to immunoglobulin kappa light chains without interfering with the antigen-binding site and binds a wider range of Ig classes and subclasses than other antibody-binding proteins such as Protein A or Protein G. Protein L binds to all classes of Ig (i.e., IgG, IgM, IgA, IgE and IgD). Protein L also binds single chain variable fragments (scFv) and Fab fragments.

**Important Product Information**

- Protein L **only** binds to immunoglobulins containing light chains of type kappa I, III, and IV in human and kappa I in mouse. Protein L also may be specific for certain kappa subgroups in other species. Protein L binds scFv without interfering with antigen binding.
- Protein L binds weakly to rabbit immunoglobulins and does not bind immunoglobulins from bovine, goat or sheep; nor does it bind to lambda light chains.

**Example Procedure for ELISA-based Applications**

The following procedure is an example ELISA protocol. Specific applications and systems will require optimization. Do not use these plates for multiple antibody assays, such as sandwich ELISAs, because the first antibody will not block all binding sites and subsequent antibodies can still bind to the plate, resulting in false positives.

**Note:** For binding information concerning the various proteins, please see the website for Tech Tip #34: Binding characteristics for Immunoglobulin Binding Proteins (Protein A, G, A/G and L).

**Additional Materials Required**

- Capture antibody specific for the antigen of interest
- Wash Buffer: TBS (Product No. 28376) or PBS (Product No. 28374) with added 0.05% Tween®-20 Detergent
- Dilution Buffer: Thermo Scientific SuperBlock Blocking Buffer (Product No. 37535) containing 0.05% Tween®-20 Detergent (Product No. 28320) or other proteinaceous blocking solution
- Antigen labeled with HRP, alkaline phosphatase or biotin

- Appropriate substrate for the enzyme (see Related Thermo Scientific Products Section or our catalog or website for a complete product listing)

**Procedure**

1. Rinse each well three times with 200µL of Wash Buffer.
2. Using the Dilution Buffer dilute antibody to ~1µg/mL and add 100µL to each well.
3. Incubate plate 30-60 minutes at room temperature. For best results, for all microplate incubations use a plate mixer that creates a vortex in each well.
4. Rinse each well three times with 200µL wash buffer.
5. Add the labeled antigen to each well (~0.1µg/mL). Sample may be diluted in Dilution Buffer. Incubate at 37°C for 1 hour.
6. Rinse each well three times with 200µL of Wash Buffer.
7. If using a biotinylated antigen, add enzyme-labeled streptavidin or other biotin-binding protein and incubate at 37°C for 1 hour. Rinse each well three times with 200µL of Wash Buffer.
8. Detect the signal according to the instructions for the specific detection system being used.

**Related Thermo Scientific Products**

<b>34028</b>	<b>1-Step Ultra TMB-ELISA, 250mL</b>
<b>34022</b>	<b>1-Step Turbo TMB, 250mL</b>
<b>37621</b>	<b>1-Step PNPP, 100mL</b>
<b>15075</b>	<b>Reagent Reservoirs, 200/pkg.</b>
<b>15036</b>	<b>Sealing Tape for 96-Well Plates, 100/pkg.</b>
<b>15130</b>	<b>Pierce Protein A Coated Plates (clear, 96-well), 5 each</b>
<b>15131</b>	<b>Pierce Protein G Coated Plates (clear, 96-well), 5 each</b>
<b>15138</b>	<b>Pierce Protein A/G Coated Plates (clear, 8-well strips), 5 each</b>

Tween is a trademark of Croda International PLC.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

**No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).**

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at [www.thermoscientific.com/pierce](http://www.thermoscientific.com/pierce). For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.