

GeneAmp® RNA PCR Core Kit

	Package Contents	Catalog Number N8080143	Size 100 rxns	 Kit Contents
	Storage Conditions	<ul style="list-style-type: none"> Store all contents at -20°C until just prior to use. Reverse Transcriptase and RNase Inhibitor are sensitive to air oxidation. 		
	Required Materials	<ul style="list-style-type: none"> Template: RNA Forward and reverse gene-specific primers DEPC-treated water (Cat. no. AM9916) E-Gel® General Purpose Gels, 1.2% (Cat. no. G5018-01) TrackIt™ 1 kb Plus DNA Ladder (Cat. no. 10488-085) 0.2 or 0.5-mL nuclease-free microcentrifuge tubes <i>Optional:</i> 1 M stock DTT (Cat. no. P2325) 		
	Timing	Varies depending on amplicon length		
	Selection Guides	<p>Go online to view related products. PCR Enzymes and Master Mixes</p>		
	Product Description	<ul style="list-style-type: none"> The GeneAmp® RNA PCR Core Kit is designed for use in detection and analysis of RNA gene expression. Murine Leukemia Virus (MuLV) Reverse Transcriptase is included for reverse transcription of RNA to cDNA and AmpliTaq® DNA Polymerase for subsequent PCR amplification. A recombinant RNase Inhibitor is included for inhibition of certain mammalian RNases. 		
	Important Guidelines	<ul style="list-style-type: none"> Select the correct polymerase, PCR instrument, and cycling conditions for your application. Take precautions to avoid cross-contamination by using aerosol-resistant barrier tips and analyzing PCR products in a separate area from PCR assembly. Dilute the 1 M DTT stock to 100 mM final in water prior to use in your RT reaction. 		
	Online Resources	Visit our product page for additional information and protocols. For support, visit www.lifetechnologies.com/support .		

For Research Use Only. Not for use in diagnostic procedures.



RT Characteristics: MuLV

Polymerase Characteristics

Hot-start:	N/A
Length:	Up to 5 kb
Fidelity vs. Taq:	1X
Format:	Separate components

RT Reaction Setup

Use the measurements below to prepare your RT experiment, or enter your own parameters in the column provided.

Component	20- μL rxn	Custom	Final Conc.
DEPC-treated water	to 20 μL	to μL	—
10X PCR Buffer II	2 μL	μL	1X
25 mM MgCl_2	4 μL	μL	5 mM
10 mM dATP	2 μL	μL	1 μM
10 mM dCTP	2 μL	μL	1 μM
10 mM dGTP	2 μL	μL	1 μM
10 mM dTTP*	2 μL	μL	1 μM
100 mM DTT (optional)**	1 μL	μL	5 mM
50 μM Oligo d(T) ₁₆ 50 μM Random hexamers, or 10 μM Gene-specific reverse primer	1 μL 1 μL 1 μL	μL	2.5 μM 2.5 μM 0.5 μM
Template RNA	varies	μL	< 1 $\mu\text{g}/\text{rxn}$
RNase Inhibitor (20 U/ μL)	1 μL	μL	1 U/ μL
MuLV RT (50 U/ μL)	1 μL	μL	2.5 U/ μL

* dUTP can be substituted for dTTP for UNG (Cat. no. N8080096).

** Dilute the 1 M DTT stock to 100 mM final in water prior to use in your RT reaction.

For qPCR instructions, refer to [AmpErase® Uracil N-Glycosylase \(UNG\) MAN0009788](#).

RT, PCR, and One-Step RT-PCR Protocols

 See pages 2, 3, and 4 to view procedures for preparing and running your RT, PCR, or One-Step RT-PCR experiments.

Optimization Strategies

Refer to the pop-ups below for guidelines to optimize your RNA Samples and RT or PCR reactions.

 [RNA Sample Prep](#)  [RT Guidelines](#)  [PCR Guidelines](#)

 [Limited Warranty, Disclaimer, and Licensing Information](#)


Reverse Transcription (RT) Protocol

The example procedure below shows reagent volumes for a single 20- μ L reaction. For multiple reactions, prepare a master mix of components common to all reactions to minimize pipetting error, and then dispense appropriate volumes into each 0.2–0.5 mL PCR reaction tube prior to adding sample and primers.

		Steps	Procedure Details																								
1		Thaw reagents	Thaw, mix, and briefly centrifuge each component before use. Set up the reactions on ice. Note: Avoid generating bubbles when mixing the enzyme.																								
2		Prepare reverse transcription master mix	Add the following components to each reaction tube. Note: Consider the volumes for all components listed in steps 2 and 3 to determine the correct amount of water required to reach your final reaction volume. <table border="1" data-bbox="892 446 1963 747"> <thead> <tr> <th>Component</th> <th>20-μL rxn</th> <th>Final Concentration</th> </tr> </thead> <tbody> <tr> <td>DEPC-treated water</td> <td>to 20 μL*</td> <td>—</td> </tr> <tr> <td>10X PCR Buffer II</td> <td>2 μL</td> <td>1X</td> </tr> <tr> <td>25 mM MgCl₂</td> <td>4 μL</td> <td>5 mM</td> </tr> <tr> <td>10 mM dNTPs (each)</td> <td>2 μL (each)</td> <td>1 μM (each)</td> </tr> <tr> <td>100 mM DTT (optional)</td> <td>1 μL</td> <td>5 mM</td> </tr> <tr> <td>RNase Inhibitor (20 U/μL)</td> <td>1 μL</td> <td>1 U/μL</td> </tr> <tr> <td>MuLV Reverse Transcriptase (50 U/μL)</td> <td>1 μL</td> <td>2.5 U/μL</td> </tr> </tbody> </table> <p>* Any combination of water and template RNA volumes can be used as long as the final volume is 20 μL. Cap, mix, and briefly centrifuge the components.</p>	Component	20- μ L rxn	Final Concentration	DEPC-treated water	to 20 μ L*	—	10X PCR Buffer II	2 μ L	1X	25 mM MgCl ₂	4 μ L	5 mM	10 mM dNTPs (each)	2 μ L (each)	1 μ M (each)	100 mM DTT (optional)	1 μ L	5 mM	RNase Inhibitor (20 U/ μ L)	1 μ L	1 U/ μ L	MuLV Reverse Transcriptase (50 U/ μ L)	1 μ L	2.5 U/ μ L
Component	20- μ L rxn	Final Concentration																									
DEPC-treated water	to 20 μ L*	—																									
10X PCR Buffer II	2 μ L	1X																									
25 mM MgCl ₂	4 μ L	5 mM																									
10 mM dNTPs (each)	2 μ L (each)	1 μ M (each)																									
100 mM DTT (optional)	1 μ L	5 mM																									
RNase Inhibitor (20 U/ μ L)	1 μ L	1 U/ μ L																									
MuLV Reverse Transcriptase (50 U/ μ L)	1 μ L	2.5 U/ μ L																									
3		Add template RNA and primers	Add your primers and sample to each tube for a final reaction volume of 20 μ L. Note: If using oligo d(T) ₁₆ or reverse gene specific primers, incubate RNA and primer following the first two cycling steps (highlighted below), then add master mix and continue with the next three cycling steps. For random hexamers, combine the master mix, primer, and RNA, and then cycle as noted below. <table border="1" data-bbox="892 982 1963 1161"> <thead> <tr> <th>Component</th> <th>20-μL rxn</th> <th>Final Concentration</th> </tr> </thead> <tbody> <tr> <td>50 μM Oligo d(T)₁₆'</td> <td>1 μL</td> <td>2.5 μM</td> </tr> <tr> <td>50 μM Random hexamers, or</td> <td>1 μL</td> <td>2.5 μM</td> </tr> <tr> <td>10 μM Gene-specific reverse primer</td> <td>1 μL</td> <td>0.5 μM</td> </tr> <tr> <td>Template RNA</td> <td>varies</td> <td>\leq 1 μg total RNA/20-μL rxn</td> </tr> </tbody> </table> <p>Cap each tube, mix, and then briefly centrifuge the contents.</p>	Component	20- μ L rxn	Final Concentration	50 μ M Oligo d(T) ₁₆ '	1 μ L	2.5 μ M	50 μ M Random hexamers, or	1 μ L	2.5 μ M	10 μ M Gene-specific reverse primer	1 μ L	0.5 μ M	Template RNA	varies	\leq 1 μ g total RNA/20- μ L rxn									
Component	20- μ L rxn	Final Concentration																									
50 μ M Oligo d(T) ₁₆ '	1 μ L	2.5 μ M																									
50 μ M Random hexamers, or	1 μ L	2.5 μ M																									
10 μ M Gene-specific reverse primer	1 μ L	0.5 μ M																									
Template RNA	varies	\leq 1 μ g total RNA/20- μ L rxn																									
4		Incubate reactions in a thermal cycler	<table border="1" data-bbox="892 1218 1963 1559"> <thead> <tr> <th>Primer Type</th> <th>Temperature</th> <th>Time</th> </tr> </thead> <tbody> <tr> <td rowspan="5">Oligo d(T)₁₆ or gene-specific primers</td> <td>65°C</td> <td>5 minutes</td> </tr> <tr> <td>4°C</td> <td>2 minutes</td> </tr> <tr> <td>37°C</td> <td>30 minutes</td> </tr> <tr> <td>95°C</td> <td>5 minutes</td> </tr> <tr> <td>4°C</td> <td>indefinitely</td> </tr> <tr> <td rowspan="4">Random hexamers</td> <td>25°C</td> <td>10 minutes</td> </tr> <tr> <td>37°C</td> <td>30 minutes</td> </tr> <tr> <td>95°C</td> <td>5 minutes</td> </tr> <tr> <td>4°C</td> <td>indefinitely</td> </tr> </tbody> </table>	Primer Type	Temperature	Time	Oligo d(T) ₁₆ or gene-specific primers	65°C	5 minutes	4°C	2 minutes	37°C	30 minutes	95°C	5 minutes	4°C	indefinitely	Random hexamers	25°C	10 minutes	37°C	30 minutes	95°C	5 minutes	4°C	indefinitely	
Primer Type	Temperature	Time																									
Oligo d(T) ₁₆ or gene-specific primers	65°C	5 minutes																									
	4°C	2 minutes																									
	37°C	30 minutes																									
	95°C	5 minutes																									
	4°C	indefinitely																									
Random hexamers	25°C	10 minutes																									
	37°C	30 minutes																									
	95°C	5 minutes																									
	4°C	indefinitely																									


PCR Amplification Protocol

The example PCR procedure below shows reagent volumes for a single **100- μ L** reaction, including the 20- μ L reverse transcription (RT) reaction mix described in the RT Protocol (page 2). For multiple reactions, prepare a master mix of components common to all reactions to minimize pipetting error, and then dispense the appropriate volumes into each reverse transcription reaction tube prior to adding the primers.

Steps		Procedure Details																											
1	 Thaw reagents	Thaw, mix, and briefly centrifuge each component before use. Note: Avoid generating bubbles when mixing the enzyme. Note: Consider the volumes for all components listed in steps 2–4 to determine the correct amount of water required to reach your final reaction volume.																											
2	 Prepare PCR master mix	Combine the following components in each reaction tube. <table border="1"> <thead> <tr> <th>Component</th> <th>100-μL rxn</th> <th>Final Concentration</th> </tr> </thead> <tbody> <tr> <td>Autoclaved, distilled water</td> <td>to 76 μL</td> <td>—</td> </tr> <tr> <td>10X PCR Buffer II</td> <td>8 μL</td> <td>1X</td> </tr> <tr> <td>25 mM MgCl₂</td> <td>8 μL</td> <td>2 mM</td> </tr> <tr> <td>AmpliTaq® DNA Polymerase (5 U/μL)</td> <td>0.5 μL</td> <td>0.025 U/μL</td> </tr> </tbody> </table> Mix and briefly centrifuge the components.		Component	100- μ L rxn	Final Concentration	Autoclaved, distilled water	to 76 μ L	—	10X PCR Buffer II	8 μ L	1X	25 mM MgCl ₂	8 μ L	2 mM	AmpliTaq® DNA Polymerase (5 U/ μ L)	0.5 μ L	0.025 U/ μ L											
Component	100- μ L rxn	Final Concentration																											
Autoclaved, distilled water	to 76 μ L	—																											
10X PCR Buffer II	8 μ L	1X																											
25 mM MgCl ₂	8 μ L	2 mM																											
AmpliTaq® DNA Polymerase (5 U/ μ L)	0.5 μ L	0.025 U/ μ L																											
3	 Add PCR master mix to reverse transcription reaction tubes	Dispense 76 μ L of the PCR master mix into each 20- μ L reverse transcription reaction tube. <table border="1"> <thead> <tr> <th>Component</th> <th>100-μL rxn</th> </tr> </thead> <tbody> <tr> <td>PCR master mix (from step 2)</td> <td>76 μL</td> </tr> <tr> <td>RT reaction mix (from RT protocol)</td> <td>20 μL</td> </tr> </tbody> </table>		Component	100- μ L rxn	PCR master mix (from step 2)	76 μ L	RT reaction mix (from RT protocol)	20 μ L																				
Component	100- μ L rxn																												
PCR master mix (from step 2)	76 μ L																												
RT reaction mix (from RT protocol)	20 μ L																												
4	 Add primers	Add your primers to each tube for a final reaction volume of 100 μL . <table border="1"> <thead> <tr> <th>Component</th> <th>100-μL rxn</th> <th>Final Concentration</th> </tr> </thead> <tbody> <tr> <td>10 μM forward primer</td> <td>2 μL</td> <td>0.2 μM</td> </tr> <tr> <td>10 μM reverse primer</td> <td>2 μL</td> <td>0.2 μM</td> </tr> <tr> <td>Total</td> <td>100 μL</td> <td></td> </tr> </tbody> </table> Cap, mix, and then centrifuge the contents.		Component	100- μ L rxn	Final Concentration	10 μ M forward primer	2 μ L	0.2 μ M	10 μ M reverse primer	2 μ L	0.2 μ M	Total	100 μL															
Component	100- μ L rxn	Final Concentration																											
10 μ M forward primer	2 μ L	0.2 μ M																											
10 μ M reverse primer	2 μ L	0.2 μ M																											
Total	100 μL																												
5	 Incubate reactions in a thermal cycler	<table border="1"> <thead> <tr> <th></th> <th>Step</th> <th>Temperature</th> <th>Time</th> </tr> </thead> <tbody> <tr> <td>1 cycle</td> <td>Initial denaturation</td> <td>95°C</td> <td>2 minutes</td> </tr> <tr> <td rowspan="3">35 PCR cycles</td> <td>Denature</td> <td>95°C</td> <td>15 seconds</td> </tr> <tr> <td>Anneal</td> <td>~55°C (depending on primer T_m)</td> <td>30 seconds</td> </tr> <tr> <td>Extend</td> <td>72°C</td> <td>1 minute/kb</td> </tr> <tr> <td>1 cycle</td> <td>Final extension</td> <td>72°C</td> <td>7 minutes</td> </tr> <tr> <td></td> <td>Hold</td> <td>4°C</td> <td>indefinitely</td> </tr> </tbody> </table>			Step	Temperature	Time	1 cycle	Initial denaturation	95°C	2 minutes	35 PCR cycles	Denature	95°C	15 seconds	Anneal	~55°C (depending on primer T _m)	30 seconds	Extend	72°C	1 minute/kb	1 cycle	Final extension	72°C	7 minutes		Hold	4°C	indefinitely
	Step	Temperature	Time																										
1 cycle	Initial denaturation	95°C	2 minutes																										
35 PCR cycles	Denature	95°C	15 seconds																										
	Anneal	~55°C (depending on primer T _m)	30 seconds																										
	Extend	72°C	1 minute/kb																										
1 cycle	Final extension	72°C	7 minutes																										
	Hold	4°C	indefinitely																										
6	 Analyze with gel electrophoresis	Analyze 10 μ L using agarose gel electrophoresis. Use your PCR reaction immediately for down-stream applications, or store it at -20°C.																											

One-Step (endpoint) RT-PCR Protocol

The example procedure below shows reagent volumes for a single 20- μ L reaction. For multiple reactions, prepare a master mix of components common to all reactions to minimize pipetting error, and then dispense appropriate volumes into each 0.2–0.5 mL PCR reaction tube prior to adding sample and primers.

		Steps	Procedure Details																																				
1		Thaw reagents	<p>Thaw, mix, and briefly centrifuge each component before use.</p> <p>Note: Avoid generating bubbles when mixing the enzyme.</p> <p>Note: Consider the volumes for all components listed in steps 2 and 3 to determine the correct amount of water required to reach your final reaction volume.</p>																																				
2		Prepare PCR master mix	<p>Combine the following components in each reaction tube.</p> <table border="1"> <thead> <tr> <th>Component</th> <th>20-μL rxn</th> <th>Final Concentration</th> </tr> </thead> <tbody> <tr> <td>DEPC-treated water</td> <td>14.85 μL</td> <td>—</td> </tr> <tr> <td>10X PCR Buffer II</td> <td>2.0 μL</td> <td>1X</td> </tr> <tr> <td>25 mM MgCl₂</td> <td>1.4 μL</td> <td>1.75 mM</td> </tr> <tr> <td>10 mM dATP</td> <td>0.4 μL</td> <td>0.2 μM</td> </tr> <tr> <td>10 mM dCTP</td> <td>0.4 μL</td> <td>0.2 μM</td> </tr> <tr> <td>10 mM dGTP</td> <td>0.4 μL</td> <td>0.2 μM</td> </tr> <tr> <td>10 mM dTTP</td> <td>0.4 μL</td> <td>0.2 μM</td> </tr> <tr> <td>100 mM DTT (optional)</td> <td>1.0 μL</td> <td>5 mM</td> </tr> <tr> <td>RNase Inhibitor (20 U/μL)</td> <td>0.2 μL</td> <td>0.2 U/μL</td> </tr> <tr> <td>MuLV Reverse Transcriptase (50 U/μL)</td> <td>0.12 μL</td> <td>0.3 U/μL</td> </tr> <tr> <td>AmpliAq® DNA Polymerase (5 U/μL)</td> <td>0.2 μL</td> <td>0.05 U/μL</td> </tr> </tbody> </table> <p>Mix and briefly centrifuge the components.</p>	Component	20- μ L rxn	Final Concentration	DEPC-treated water	14.85 μ L	—	10X PCR Buffer II	2.0 μ L	1X	25 mM MgCl ₂	1.4 μ L	1.75 mM	10 mM dATP	0.4 μ L	0.2 μ M	10 mM dCTP	0.4 μ L	0.2 μ M	10 mM dGTP	0.4 μ L	0.2 μ M	10 mM dTTP	0.4 μ L	0.2 μ M	100 mM DTT (optional)	1.0 μ L	5 mM	RNase Inhibitor (20 U/ μ L)	0.2 μ L	0.2 U/ μ L	MuLV Reverse Transcriptase (50 U/ μ L)	0.12 μ L	0.3 U/ μ L	AmpliAq® DNA Polymerase (5 U/ μ L)	0.2 μ L	0.05 U/ μ L
Component	20- μ L rxn	Final Concentration																																					
DEPC-treated water	14.85 μ L	—																																					
10X PCR Buffer II	2.0 μ L	1X																																					
25 mM MgCl ₂	1.4 μ L	1.75 mM																																					
10 mM dATP	0.4 μ L	0.2 μ M																																					
10 mM dCTP	0.4 μ L	0.2 μ M																																					
10 mM dGTP	0.4 μ L	0.2 μ M																																					
10 mM dTTP	0.4 μ L	0.2 μ M																																					
100 mM DTT (optional)	1.0 μ L	5 mM																																					
RNase Inhibitor (20 U/ μ L)	0.2 μ L	0.2 U/ μ L																																					
MuLV Reverse Transcriptase (50 U/ μ L)	0.12 μ L	0.3 U/ μ L																																					
AmpliAq® DNA Polymerase (5 U/ μ L)	0.2 μ L	0.05 U/ μ L																																					
3		Add primers and template RNA	<p>Add your primers and template RNA to each reaction tube.</p> <table border="1"> <thead> <tr> <th>Component</th> <th>20-μL rxn</th> <th>Final Concentration</th> </tr> </thead> <tbody> <tr> <td>10 μM forward gene-specific primer</td> <td>0.8 μL</td> <td>0.4 μM</td> </tr> <tr> <td>10 μM reverse gene-specific primer</td> <td>0.8 μL</td> <td>0.4 μM</td> </tr> <tr> <td>Template RNA</td> <td>varies</td> <td>< 1 μg</td> </tr> </tbody> </table> <p>Cap, mix, and then centrifuge the contents.</p>	Component	20- μ L rxn	Final Concentration	10 μ M forward gene-specific primer	0.8 μ L	0.4 μ M	10 μ M reverse gene-specific primer	0.8 μ L	0.4 μ M	Template RNA	varies	< 1 μ g																								
Component	20- μ L rxn	Final Concentration																																					
10 μ M forward gene-specific primer	0.8 μ L	0.4 μ M																																					
10 μ M reverse gene-specific primer	0.8 μ L	0.4 μ M																																					
Template RNA	varies	< 1 μ g																																					
4		Incubate reactions in a thermal cycler	<table border="1"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> </tr> </thead> <tbody> <tr> <td>RT step</td> <td>37°C</td> <td>30 minutes</td> </tr> <tr> <td>Initial denaturation</td> <td>95°C</td> <td>2 minutes</td> </tr> <tr> <td rowspan="2">35 PCR cycles</td> <td>95°C</td> <td>20 seconds</td> </tr> <tr> <td>60°C</td> <td>1 minute/kb</td> </tr> <tr> <td>Final extension</td> <td>72°C</td> <td>7 minutes</td> </tr> <tr> <td>Hold</td> <td>4°C</td> <td>indefinitely</td> </tr> </tbody> </table>	Step	Temperature	Time	RT step	37°C	30 minutes	Initial denaturation	95°C	2 minutes	35 PCR cycles	95°C	20 seconds	60°C	1 minute/kb	Final extension	72°C	7 minutes	Hold	4°C	indefinitely																
Step	Temperature	Time																																					
RT step	37°C	30 minutes																																					
Initial denaturation	95°C	2 minutes																																					
35 PCR cycles	95°C	20 seconds																																					
	60°C	1 minute/kb																																					
Final extension	72°C	7 minutes																																					
Hold	4°C	indefinitely																																					
5		Analyze with gel electrophoresis	<p>Analyze 10 μL using agarose gel electrophoresis.</p> <p>Use your PCR reaction immediately for down-stream applications, or store it at –20°C.</p>																																				