

General Recommendations for RNA Electrophoresis

This protocol is for the General Recommendations for RNA Electrophoresis

- RNA ladders, as well as any RNA, are extremely sensitive to degradation by ribonucleases. Use only fresh electrophoresis buffers and freshly poured gels.
- Use clean electrophoresis chambers. For RNA gel analysis, avoid electrophoresis tanks used for DNA miniprep analysis since DNA minipreps may contain RNase A or T1.
- Use the same loading dye for samples and for RNA markers. 2X RNA Loading Dye is available separately and is provided with all RiboRuler™ RNA ladders. It contains ethidium bromide for RNA visualization on denaturing formaldehyde gels. If RNA fragments are separated on native agarose gels or on polyacrylamide/urea gels, additional staining with ethidium bromide is recommended.
- For native gels, add 0.5 µg/ml of ethidium bromide to the agarose gel and to the running buffer.

thermoscientific.com

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

North America
Technical Services:
techservice.genomics@thermofisher.com
Customer Services:
customerservice.genomics@
thermofisher.com
Tel 800 235 9880
Fax 800 292 6088

Europe and Asia
Technical Services:
techservice.emea.genomics@
thermofisher.com
Customer Services:
customerservice.emea.genomics@
thermofisher.com

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific