remelTM

RapIDTM Spot Indole

Reagent

EN

REF R8309002..... 15 ml / Bottle

1. INTENDED USE

RapIDTM Spot Indole is a reagent to be used with the RapID manual identification systems.

The device is for professional use only, is not automated nor is it a companion diagnostic.

2. SUMMARY AND EXPLANATION

In 1963, Vracko *et al.* evaluated the use of *p*-dimethylaminocinnamaldehyde (DMACA) to detect indole production by enteric gram-negative bacilli.¹ Further testing by Lowrance *et al.* demonstrated DMACA to be the most sensitive reagent for rapid indole detection.²

3. PRINCIPLE

Intracellular enzymes collectively called “tryptophanases” mediate the degradation of tryptophan to produce indole and other metabolic end products. Indole can be detected by addition of *p*-dimethylaminocinnamaldehyde, which combines with indole to produce a distinct blue-green colour.

4. REAGENTS

RapID Spot Indole Reagent 15ml/bottle

p-Dimethylaminocinnamaldehyde (CAS 6203-18-5) ...10.0 g

Hydrochloric Acid (CAS 7647-01-0) 100.0 ml

Demineralized Water (CAS 7732-18-5) 900.0 ml

5. PRECAUTIONS

WARNING! Causes eye, skin, and respiratory tract irritation.

This product is for *in vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers and media after use. Directions should be read and followed carefully. Refer to Safety Data Sheet (SDS) on company website for additional information.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

Composition / information on ingredients

Hydrochloric acid 7647-01-0

WARNING

H290	May be corrosive to metals
H314	Causes severe skin burns and eye damage
H335	May cause respiratory irritation
P260	Do not breathe dust/fume/gas/mist/vapors/spray
P264	Wash face, hands and any exposed skin thoroughly after handling
P280	Wear protective gloves/protective clothing/eye protection/face protection
P271	Use only outdoors or in a well-ventilated area
P234	Keep only in original container
P310	Immediately call a POISON CENTER or doctor/physician
P304+P340	IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing
P303+P361 +P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower
P363	Wash contaminated clothing before reuse
P305+P351 +P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
P301+P330 +P331	IF SWALLOWED: Rinse mouth. DO NOT induce vomiting
P390	Absorb spillage to prevent material damage
P405	Store locked up
P403+P233	Store in a well-ventilated place. Keep container tightly closed
P406	Store in corrosive resistant polypropylene container with a resistant inliner
P402	Store in a dry place
P501	Dispose of contents/container to an approved waste disposal plant

6. STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 2-8°C until used. Do not freeze or overheat. Protect from light.

7. PRODUCT DETERIORATION

This product should not be used if (1) the color has changed, (2) the expiration date has passed, or (3) there are other signs of deterioration.

8. SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.⁴

9. MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms.

10. QUALITY CONTROL

All lot numbers of RapID Spot Indole Reagent have been tested and found to be acceptable. Quality control organisms should be tested according to the guidelines in the RapID Identification System in use or according to established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

11. BIBLIOGRAPHY

1. Vracko, R. and J.C. Sherris. 1963. Am. J. Clin. Path. 39:429-432.

2. Lowrance, B.L., P. Reich, and W.H. Traub. 1969. Appl. Microbiol. 17:923-924.

3. MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams and Wilkins, Philadelphia, PA.

4. Carroll, K C, M A Pfaller et al., 2019. Manual of Clinical Microbiology, 12th edition, ASM Press, Washington, US.

12. SYMBOL LEGEND

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Do not use if package is damaged
	Not for near patient testing
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Keep away from sunlight
	Importer
	Authorized representative in the European Community
	UK Conformity Assessed
	European Conformity Assessment
	Manufacturer

CAS (Chemical Abstracts Service Registry No.)

 12076 Santa Fe Drive, Lenexa, KS 66215, USA
www.thermofisher.com/microbiology
Tel: (800) 255-6730 • International: (913) 888-0939
www.oxoid.com/IFU
Europe +800 135 79 135 • US 1 855 2360 190
CA 1 855 805 8539 • ROW +31 20 794 7071

Version	Date of modifications introduced
IFU8309002	March 2023 Addition of Importer symbol to Section 12

Printed in the UK

For technical information contact your local distributor.

©2022 Thermo Fisher Scientific Inc. All rights reserved.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

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12076 Santa Fe Drive, Lenexa, KS 66215, États-Unis

www.thermofisher.com/microbiology

Tél : (800) 255-6730 • International : (913) 888-0939

www.oxid.com/IFU

Europe +800 135 79 135 • États-Unis 1 855 2360 190

CA 1 855 805 8539 • RDM +31 20 794 7071

remelTM RapIDTM Spot Indole Reagent

HR

REF R8309002, RapID Spot Indole reagents.. 15 ml / bočica

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11. BIBLIOGRAFIJA

1. Vracko, R. and J.C. Sherris. 1963. Am. J. Clin. Path. 39:429-432.
2. Lowrance, B. L., P. Reich, and W. H. Traub. 1969. Appl. Microbiol. 17:923-924.
3. MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams and Wilkins, Philadelphia, PA.
4. Carroll, K C, M A Pfaller et al., 2019. Manual of Clinical Microbiology, 12th edition, ASM Press, Washington, US.

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©2022 Thermo Fisher Scientific Inc. Všetky práva vyhradené.