


Power SYBR[®] Green RNA-to-C_T[™] 1-Step Kit

For safety and biohazard guidelines refer to the “Safety” section in the *Power SYBR[®] Green RNA-to-C_T[™] 1-Step Kit Protocol* (PN 4391003). For all chemicals in **bold red** type, read the MSDS and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.


This quick reference card provides abbreviated procedures. For complete procedures, refer to the *Power SYBR[®] Green RNA-to-C_T[™] 1-Step Kit Protocol* (PN 4391003).

1 Prepare the RT-PCR reactions.



Component	Volume for One Reaction		
	10 μ L	20 μ L	50 μ L
Power SYBR[®] Green RT-PCR Mix (2\times)	5.0 μ L	10.0 μ L	25.0 μ L
Forward primer (100 to 200 nM final)	Variable	Variable	Variable
Reverse primer (100 to 200 nM final)	Variable	Variable	Variable
RT Enzyme Mix (125 \times)	0.08 μ L	0.16 μ L	0.4 μ L
RNA template (up to 100 ng)	Variable	Variable	Variable
RNase-free H ₂ O	to 10 μ L	to 20 μ L	to 50 μ L
Total Volume	10.0 μL	20.0 μL	50.0 μL

2 Prepare the reaction plate.



Use a reaction plate appropriate for your real-time PCR system:

- MicroAmp[™] Fast Optical 48-Well Reaction Plate: 20 μ L
- MicroAmp[™] Fast Optical 96-Well Reaction Plate: 20 μ L
- MicroAmp[™] Optical 96-Well Reaction Plate: 50 μ L
- MicroAmp[™] Optical 384-Well Reaction Plate: 10 μ L


3 Run the RT-PCR reaction plate.



- Ramp speed or mode: **Standard**
- Reaction volume (μ L): **10, 20, or 50**
- Thermal cycling conditions:

Stage	Temp	Time
Holding	48 °C	30 min
Holding	95 °C	10 min
Cycling (40 cycles)	95 °C	15 sec
	60 °C	1 min
Melt Curve (optional)	95 °C	15 sec
	60 °C	15 sec
	95 °C	15 sec

4 Analyze the experiment.



- **Standard curve** (standard curve and relative standard curve experiments) – Slope, amplification efficiency, R² values, y-intercept, C_T values, outliers
- **Gene expression plot** (relative standard curve and comparative C_T experiments) – Differences in gene expression, standard deviation in the replicate groups
- **Amplification plots** – Baseline and threshold values, outliers
- **Well table or results table** – C_T values for each well and for each replicate group
- (Optional) **Melt curve** – Number of T_m peaks

Power SYBR® Green RNA-to-C_T™ 1-Step Kit Products

Quantity	Part Number
Reagents sufficient for 40 × 50-µL reactions: <ul style="list-style-type: none"> • Power SYBR® Green RT-PCR Mix (2X), 1 mL • RT Enzyme Mix (125X), 20 µL 	4391178
Reagents sufficient for 200 × 50-µL reactions: <ul style="list-style-type: none"> • Power SYBR® Green RT-PCR Mix (2X), 5 mL • RT Enzyme Mix (125X), 80 µL 	4389986

© Copyright 2007, Applied Biosystems. All rights reserved.

For Research Use Only. Not for use in diagnostic procedures.

NOTICE TO PURCHASER: PLEASE REFER TO THE *Power SYBR® Green RNA-to-C_T™ 1-Step Kit Protocol* (PN 4391003) FOR LIMITED LABEL LICENSE OR DISCLAIMER INFORMATION.

Applied Biosystems, and AB (Design) are registered trademarks and RNA-to-C_T and MicroAmp are trademarks of Applied Biosystems Corporation or its subsidiaries in the U.S. and/or certain other countries.

SYBR is a registered trademark of Molecular Probes, Inc.

12/2007