

JC-1 and JC-9 Mitochondrial Potential Sensors

Table 1. Contents and storage information.

Material	Amount	Storage	Stability
JC-1; CBIC ₂ (3)	5 mg *	<ul style="list-style-type: none"> • 2–6°C • Desiccate • Protect from light 	When stored as directed, product should be stable for 6 months.
JC-9; DiNOC ₁ (3)	5 mg †		

* JC-1 concentrations of 1–5 mg/mL correspond to 1.5–7.7 mM, MW = 652. † JC-9 concentrations of 1–5 mg/mL correspond to 1.9–9.4 mM, MW = 532.

Approximate fluorescence excitation/emission maxima: 514/529 nm, monomer form; 585/590 nm J-aggregate form.

Introduction

JC-1 and JC-9 are cationic dyes (Figure 1) that exhibit potential-dependent accumulation in mitochondria, indicated by a fluorescence emission shift from green (~525 nm) to red (~590 nm). Consequently, mitochondrial depolarization is indicated by a decrease in the red/green fluorescence intensity ratio. The potential-sensitive color shift is due to concentration-dependent formation of red fluorescent J-aggregates.^{1–3} JC-1 can be used as an indicator of mitochondrial potential in a variety of cell types, including myocytes³ and neurons,⁴ as well as in intact tissues⁵ and isolated mitochondria.⁶ JC-1 is more specific for mitochondrial versus plasma membrane potential, and more consistent in its response to depolarization, than other cationic dyes such as DiOC₆(3) and rhodamine 123.⁷

The ratio of green to red fluorescence is dependent only on the membrane potential and not on other factors such as mitochondrial size, shape, and density that may influence single-component fluorescence signals. Use of fluorescence ratio detection therefore allows researchers to make comparative measurements of membrane potential and determine the percentage of mitochondria within a population that respond to an applied stimulus. Subtle heterogeneity in cellular responses can be discerned in this way.^{1,6} For example, four distinct patterns of mitochondrial membrane potential change in response to glutamate receptor activation in neurons have been identified using confocal ratio imaging of JC-1 fluorescence.⁴ The most widely implemented application of JC-1 is for detection of mitochondrial depolarization occurring in the early stages of apoptosis (Figure 2).^{7–10}




Figure 1. Structure of JC-1 (5,5',6,6'-tetrachloro-1,1',3,3'-tetraethylbenzimidazolylcarbocyanine iodide).




Figure 2. Bivariate JC-1 analysis of mitochondrial membrane potential in HL60 cells by flow cytometry. The sensitivity of this technique is demonstrated by the response to K⁺/valinomycin-induced depolarization (panel B). Distinct populations of cells with different extents of mitochondrial depolarization are detectable following apoptosis-inducing treatment with 5 μ M staurosporine for 2 hours (panel C). (Figure kindly supplied by Dr. Andrea Cossarizza, University of Modena and Reggio Emilia.)

Guidelines for Use

Materials Required but Not Provided

- DMSO
- DMF

Preparing the Stock Solutions

Stock solutions can be prepared at 1–5 mg/mL in dimethylsulfoxide (DMSO) or dimethylformamide (DMF). A convenient procedure for storing stock solutions is to divide them into portions, each sufficient for one day of experimental work, and store them in a freezer ($\leq -20^{\circ}\text{C}$) until required for use.³

Fluorescence Microscopy

Staining

Typical staining protocols abstracted from the research literature are summarized in Table 2. Following incubation in dye-containing medium, it is usual to wash the cells before starting experimental observations.

Optical Filters

A number of different optical filter configurations can be used for analysis of JC-1 or JC-9 by fluorescence microscopy (Table 3). For confocal laser scanning microscopy, the monomer and J-aggregate forms can be excited simultaneously by 488 nm argon-ion laser sources. The J-aggregate form can be excited selectively using the 568 nm argon-krypton laser line.

Appearance

Polarized mitochondria are marked by punctate orange-red fluorescent staining. On depolarization, the orange-red punctate staining is replaced by diffuse green monomer fluorescence. Some of the green fluorescence may remain associated with mitochondria, due to potential-independent interactions of the JC-1 monomer with mitochondrial membranes.^{2,3}

Flow Cytometry

Staining

Typical staining protocols abstracted from the re-search literature are summarized in Table 2. Dissociated cells for flow cytometric analysis are diluted to a density of about 1×10^6 cells/mL for staining.

Detector Configuration

When excited simultaneously by 488 nm argon-ion laser sources, the JC-1 monomer and J-aggregate can be detected separately in the conventional flow cytometer FL1 and FL2 channels respectively (Figure 2).

Table 2. JC-1 cell staining conditions.

Cell Type	Adherent or Dissociated	Incubation Conditions			Analysis Method
		Dye Concentration	Temperature	Time	
Neurons (rat) ¹	Adherent	2.0 µg/mL	37°C	20–30 min	Confocal microscope
Neurons (rat) ²	Adherent	1.0 µg/mL	37°C	20 min	Confocal microscope
Human fibroblasts ³	Dissociated	0.3 µg/mL	37°C	1 hour	Flow cytometer
O-2A oligodendrocytes (rat) ⁴	Adherent	10 µg/mL	37°C	10 min	Wide-field microscope
PC12 ⁵	Adherent	10 µg/mL	37°C	10 min	Confocal microscope
Colo-205 ⁶	Dissociated	10 µg/mL	37°C	10 min	Flow cytometer
U937 ⁷	Dissociated	10 µg/mL	22°C	10 min	Flow cytometer
Cardiac myocytes (rat) ⁸	Dissociated	10 µg/mL	37°C	10 min	Wide-field microscope

1. J Neurosci 16, 5688 (1996); **2.** Neuron 15, 961 (1995); **3.** Exp Cell Res 245, 170 (1998); **4.** J Physiol 508, 413 (1998); **5.** Neuronal precursor cell line, J Neurosci 18, 932 (1998); **6.** Human colon adenocarcinoma, J Cell Biol 138, 449 (1997); **7.** Human premonocytic cell line, Proc Natl Acad Sci USA 93, 6458 (1996), Biochem Biophys Res Comm 197, 40 (1993); **8.** J Physiol 486, 1 (1995).

Table 3. Optical filters for fluorescence microscope imaging of JC-1.

Species Detected	Excitation	Dichroic	Emission
Monomer alone	485 ± 11 nm	505 nm	530 ± 15 nm
J-aggregate alone	535 ± 17.5 nm	570 nm	590 ± 17.5 nm
Monomer and J-aggregate, simultaneous	475 ± 20 nm	505 nm	≥510 nm
Monomer and J-aggregate, simultaneous	485 ± 11 nm	505 nm	530 ± 15 AND ≥590 nm

References

1. Proc Natl Acad Sci U S A 88, 3671 (1991); 2. Biochemistry 30, 4480 (1991); 3. J Physiol 486, 1 (1995); 4. J Neurosci 16, 5688 (1996); 5. Methods 18, 104 (1999); 6. Exp Cell Res 222, 84 (1996); 7. FEBS Lett 411, 77 (1997); 8. J Neurosci 18, 932 (1998); 9. J Cell Biol 138, 449 (1997); 10. Exp Cell Res 245, 170 (1998).

Product List

Current prices may be obtained from our website or from our Customer Service Department.

Cat #	Product Name	Unit Size
D22421	3,3'-dimethyl- α -naphthoxycarbocyanine iodide (JC-9; DiNOC ₂ (3))	5 mg
T3168	5,5',6,6'-tetrachloro-1,1',3,3'-tetraethylbenzimidazolylcarbocyanine iodide (JC-1; CBIC ₂ (3))	5 mg

Contact Information

Molecular Probes, Inc.

29851 Willow Creek Road
Eugene, OR 97402
Phone: (541) 465-8300
Fax: (541) 335-0504

Customer Service:

6:00 am to 4:30 pm (Pacific Time)
Phone: (541) 335-0338
Fax: (541) 335-0305
probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209
Order Fax: (800) 438-0228

Technical Service:

8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 335-0353
Toll-Free (800) 438-2209
Fax: (541) 335-0238
probetech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2006, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.