

## Technical Tip: A Better Way to Quantitate the Degree of Biotinylation of Proteins

Why guess? Whether you purchase a commercially available biotinylated reagent or make your own using one of our reactive biotin derivatives, you still need to know the degree of biotinylation before you proceed with your experiment. If you're still relying on the conventional colorimetric HABA-based biotin-binding assay<sup>1</sup>, we've got a better way.

The FluoReporter® biotin quantitation assay is based on the Biotective™ Green reagent (Figure 1). The Biotective™ Green reagent exhibits bright green fluorescence (excitation/emission maxima of 495/519 nm), making this assay compatible with standard fluorescence-based microplate readers. The FluoReporter® Biotin Quantitation Assay Kits offers several important advantages:

### FluoReporter® Biotin Quantitation Assay Kit for proteins (F30751)


Detect from 4 to 80 pmol of biotin, a 50-fold higher sensitivity than the HABA-based assay. In addition, unlike the HABA-based assay, which requires ~1 mg of protein sample, the FluoReporter® Biotin Quantitation Assay Kit for proteins requires <1 µg.

### FluoReporter® Biotin Quantitation Assay Kit for nucleic acid (F30755)

This kit is ideal for determining the degree of biotinylation of cDNA samples used in Affymetrix or RLS microarray protocols and can be applied to as little as 13.2 ng\* of nucleic acid.

In addition to the Biotective™ Green reagent, each FluoReporter® Biotin Quantitation Assay Kit contains a biotin standard for the assay, a biotinylated positive control, protease or nuclease, and the appropriate digestions buffers. For accurate biotinylation data, rely on the FluoReporter® assay technology. Additional information on these products can be found at [probes.invitrogen.com](http://probes.invitrogen.com).

\*Accurate results using 13.2 ng of nucleic acid containing ~1 biotin for every 10 bases or 132 ng of nucleic acid containing ~1 biotin for every 100 bases.



**Figure 1.** Schematic representation of the FluoReporter® biotin quantitation assay. This assay uses the Biotective™ Green reagent, which consists of avidin labeled with a fluorescent dye (D) and with quencher dye ligands (Q) occupying the biotin binding sites. Through fluorescence resonance energy transfer (FRET), the ligand quenches the fluorescence. Biotin (B) attached to a protein or nucleic acid displaces the quencher dye from the Biotective™ Green reagent, yielding fluorescence proportional to the amount of added biotin.

## Reference

---

1. Biochem J 94, 23C (1965).

## Product List **Current prices may be obtained from our website or from our Customer Service Department.**

---

Cat #	Product Name	Unit Size
F30751	FluoReporter <sup>®</sup> Biotin Quantitation Assay Kit *for biotinylated proteins* *5 determinations* .....	1 kit
F30755	FluoReporter <sup>®</sup> Biotin Quantitation Assay Kit *for biotinylated nucleic acids* *10 determinations* .....	1 kit

## Contact Information

---

### Molecular Probes, Inc.

29851 Willow Creek Road  
Eugene, OR 97402  
Phone: (541) 465-8300  
Fax: (541) 335-0504

### Customer Service:

6:00 am to 4:30 pm (Pacific Time)  
Phone: (541) 335-0338  
Fax: (541) 335-0305  
probesorder@invitrogen.com

### Toll-Free Ordering for USA:

Order Phone: (800) 438-2209  
Order Fax: (800) 438-0228

### Technical Service:

8:00 am to 4:00 pm (Pacific Time)  
Phone: (541) 335-0353  
Toll-Free (800) 438-2209  
Fax: (541) 335-0238  
probestech@invitrogen.com

### Invitrogen European Headquarters

Invitrogen, Ltd.  
3 Fountain Drive  
Inchinnan Business Park  
Paisley PA4 9RF, UK  
Phone: +44 (0) 141 814 6100  
Fax: +44 (0) 141 814 6260  
Email: euroinfo@invitrogen.com  
Technical Services: eurotech@invitrogen.com

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Assistance Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

### Limited Use Label License

For research use only. Not intended for any animal or human therapeutic or diagnostic use. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402. Tel: (541) 465-8300. Fax: (541) 335-0504.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation <sup>®</sup> are registered with the U.S. Patent and Trademark Office.

Copyright 2005, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.