

DNA Insert Ligation (sticky-end and blunt-end) into Vector DNA

This protocol is for the DNA Insert Ligation (sticky-end and blunt-end) into Vector DNA.

Sticky-end Ligation

1. Prepare the following reaction mixture:

Linear vector DNA	20-100 ng
Insert DNA	1:1 to 5:1 molar ratio over vector
10x T4 DNA Ligase buffer	2 μ L
Thermo Scientific T4 DNA Ligase (Cat #EL0016)	1 U
Water, nuclease-free	to 20 μ L
Total volume	20 μL

2. Incubate 10 minutes at 22 °C.
3. Use up to 5 μ L of the mixture for transformation of 50 μ L of chemically competent cells or use 1-2 μ L per 50 μ L electrocompetent cells.

Note

- The electrotransformation efficiency may be improved by:
 - heat inactivation of T4 DNA Ligase at 65 °C for 10 minutes or at 70 °C for 5 minutes,
 - purification of DNA, using Thermo Scientific GeneJET™ PCR Purification Kit or by chloroform extraction.
- The overall number of transformants may be increased by extending the reaction time to 1 hour.

- If more than 2 U of T4 DNA Ligase is used in 20 μ L reaction mixture, it is necessary to purify DNA (by spin column or chloroform extraction) before electrotransformation.

Blunt-end Ligation

1. Prepare the following reaction mixture:

Linear vector DNA	20-100 ng
Insert DNA	1:1 to 5:1 molar ratio over vector
Thermo Scientific 10x T4 DNA Ligase buffer	2 μ L
50% PEG 4000 solution	2 μ L
T4 DNA Ligase	5 U
Water, nuclease-free	to 20 μ L
Total volume	20 μL

2. Incubate 1 hour at 22 °C.
3. Use up to 5 μ L of the mixture for transformation of 50 μ L of chemically competent cells. Purify DNA for electrotransformation, using the GeneJET™ PCR Purification Kit or by chloroform extraction. Use 1-2 μ L of DNA solution per 50 μ L of electrocompetent cells.

thermoscientific.com/onebio

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Europe

Customer Service
cs.molbio.eu@thermofisher.com

Technical Support
ts.molbio.eu@thermofisher.com

Tel 00800 222 00 888
Fax 00800 222 00 889

United States

Customer Service
cs.molbio@thermofisher.com

Technical Support
ts.molbio@thermofisher.com

Tel 800 235 9880
Fax 800 292 6088

Canada

Customer Service
cs.molbio@thermofisher.com

Technical Support
ts.molbio@thermofisher.com

Tel 800 340 9026
Fax 800 472 8322

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific