

Hybridoma-Serum Free Media (SFM) and Protein-Free Hybridoma Medium (PFHM-II)

Catalog Numbers 12045084, 12045076, 12300067, 23600042

Pub. No. MAN0007366 Rev. 2.0



WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from [thermofisher.com/support](https://www.thermofisher.com/support).



CAUTION! Human origin materials are non-reactive (donor level) for anti-HIV 1 & 2, anti-HCV, and HBsAg. Handle in accordance with established bio-safety practices.

Product description

The Gibco™ Hybridoma-Serum Free Media (SFM) and Protein-Free Hybridoma Medium (PFHM-II) are serum-free products optimized for growth of hybridomas and monoclonal antibody production. Hybridoma-SFM is a very low (<20 µg/mL) protein medium. The very low protein content facilitates monoclonal antibody purification.

PFHM-II is a protein-free, ready-to-use medium that contains no polypeptide growth or attachment factors, or mediators that may complicate downstream processing and final product purification. PFHM-II also performs well as serum-supplemented media for monoclonal antibody production and also may be used as a growth medium.

Contents and storage

Product	Cat. No.	Amount	Storage	Shelf Life ^[1]
Hybridoma-SFM	12045084	500 mL	2–8°C; Protect from light	12 months
	12045076	1000 mL		
Hybridoma-SFM, powder	12300067	1 × 10 L	2–8°C; Store dark and dry	33 months
PFHM-II, powder	23600042	1 × 10 L	2–8°C; Store dark and dry	24 months

^[1] Shelf Life duration is determined from Date of Manufacture.

Culture conditions

Media: Hybridoma-SFM or PFHM-II

Cell type: Hybridoma

Culture type: Suspension

Culture vessels: Shake flasks, roller bottles or bioreactor

Temperature range: 36°C to 38°C

Incubator atmosphere: Humidified atmosphere of 5–10% CO₂ in air. Ensure proper gas exchange and minimize exposure of cultures to light.

(i.e., Cholesterol Lipid Concentrate (250X)) for growth of cholesterol-dependent cell lines (e.g., NS0 and derivatives).

Details	PFHM-II	Hybridoma-SFM
Protein-free	X	
Contains insulin and transferrin		X
Contains phenol red	X	X
Contains surfactant ^[1]		X
Contains inorganic iron carrier ^[2]	X	

^[1] PFHM-II does not contain a surfactant. If used for agitated suspension culture, supplement with 0.1% Pluronic™ F-68.

^[2] Medium should be pre-screened to determine potential interference of inorganic iron carrier(s) with antibody detection and/or purification method.

Procedural guidelines

- Hybridoma-SFM and PFHM-II require supplementation with a cholesterol supplement or some other source of cholesterol
- In most instances, antibiotics are neither necessary nor advised. However, where antibiotics are required, most general antibiotics are compatible with PFHM-II, including Penicillin-Streptomycin, Gentamicin, anti-PPLO, linocin, and

Reconstitute powdered media

Add powdered media to 9 L room temperature deionized or distilled water. Rinse package to remove all traces of powder.

Hybridoma-SFM

1. Add 2.45 g NaHCO_3 per L of medium. Stir until completely dissolved, do not heat.
2. Adjust to pH 8.0 with 1N NaOH while stirring. Slowly bring the pH to 7.0–7.1 using 1N HCl (Filtration will raise the pH 0.1–0.3 units; final post-filtration pH should be 7.2–7.4).
3. Add distilled water to final volume of 10 L. Check pH and osmolality (Osmolality should be 320–345 mOsm/kg).
4. Stir for 15–20 minutes at room temperature. Filter sterilize using a surfactant-free 0.2 μm filter and dispense into sterile, clean containers. Protect from light.

PFHM-II

1. Adjust pH to 7.0 using either 1N HCl or 1N NaOH while stirring.
2. Add 2.0 g NaHCO_3 per L of medium. Stir until completely dissolved, do not heat.
3. Readjust pH to 7.0, (see step 1). Add distilled water to final volume of 10 L.
4. Stir for 15–20 minutes at room temperature. Filter sterilize using a surfactant-free 0.2 μm filter and dispense into sterile, clean containers. Protect from light.

Prepare media

Hybridoma-SFM and PFHM-II do not require supplementation except as noted for cholesterol dependent cell lines, agitated suspension cultures using PFHM-II, and if antibiotics are desired.

1. Aseptically add 2 mL of Cholesterol Lipid Concentrate (250X) to 500 mL of Hybridoma-SFM or PFHM-II.
2. Aseptically add antibiotics, if required.

Recovery

1. Rapidly thaw (<1 minute) frozen cells in a 37°C water bath.
2. Transfer the entire contents of the cryovial into a tissue culture flask containing 30 mL prewarmed Hybridoma-SFM or PFHM-II without antibiotics.
3. Incubate at 37°C in a humidified atmosphere of 5% CO_2 in air.
4. Subculture cells 3–5 days post thaw.

Adapt hybridoma cells to Hybridoma-SFM or PFHM-II

Successful adaptation will depend upon the particular hybridoma cell line and the culture conditions employed. We recommended that backup cultures in the original medium be maintained until success with the new medium has been achieved.

Note: It is critical that cell viability be at least 90% and cells be in the mid-logarithmic phase of growth prior to adaptation.

Direct adaption

1. Subculture hybridoma cells grown in conventional medium with 5–10% serum or other serum-free medium into prewarmed Hybridoma-SFM or PFHM-II. During the adaptation procedure seeding density should be double the normal seeding density for the cell line.
2. Monitor cell growth using Countess™ 3 Automated Cell Counter (or alternate suitable method) until the viable cell density reaches 1×10^6 viable cells/mL. Subculture the cells to a viable cell density of $1\text{--}2 \times 10^5$ viable cells/mL in fresh prewarmed Hybridoma-SFM or PFHM-II.
3. Continue to monitor and passage cells for 3–5 passages until consistent growth is achieved.

Note: If suboptimal performance is observed over 3–5 passages using the direct adaptation method, use the sequential adaptation method.

Sequential adaption

1. Subculture hybridoma cells grown in conventional medium with 5–10% serum or other serum-free medium into a 25:75 ratio of fresh Hybridoma-SFM or PFHM-II to the original media. During the adaptation procedure seed at double the normal seeding density.
2. Monitor cell growth until the viable cell density reaches 1×10^6 viable cells/mL. Subculture cells (dilute to $1\text{--}2 \times 10^5$ viable cells/mL) into stepwise increasing ratios of fresh PFHM-II or Hybridoma-SFM to original medium with each subsequent passage (50:50, 75:25, 90:10 followed by 100% Hybridoma-SFM or PFHM-II). Multiple passages at each step may be required.
3. Continue to monitor and passage cells until consistent growth is achieved. After several passages of consistent growth and viability in 100% complete Hybridoma-SFM or PFHM-II the culture is considered to be adapted.

Cryopreservation

1. Prepare the desired quantity of cells in a tissue culture flask, harvesting in mid-log phase of growth with viability >90%. Reserve the conditioned medium to prepare cryopreservation medium.
2. Determine the viable cell density and calculate the required volume of cryopreservation medium to give a final cell density of $0.5-1 \times 10^7$ cells/mL.
3. Prepare the required volume of cryopreservation medium of 92.5% medium (50:50 ratio of fresh to conditioned media) +7.5% DMSO on the day of intended use. Filter sterilize and store at 4°C until use.

IMPORTANT! Conditioned medium should be obtained from a high viability, mid-log culture of cells.

4. Harvest cells by centrifugation at $100 \times g$ for 5–10 minutes. Resuspend the pellet in the pre-determined volume of 4°C cryopreservation medium.
5. Dispense aliquots of this suspension into cryovials according to the manufacturer's specifications (i.e., 1 mL in a 2-mL cryovial).
6. Achieve cryopreservation in an automated or manual controlled rate freezing apparatus following standard procedures (1°C decrease per minute).
7. Transfer frozen cells to liquid nitrogen (vapor phase). Storage at –200°C to –125°C is recommended.

Related products

Unless otherwise indicated, all materials are available through [thermofisher.com](https://www.thermofisher.com).

Item	Source
Cholesterol Lipid Concentrate (250X)	12531018
L-Glutamine (200 mM)	25030024
GlutaMAX™ I Supplement	35050061
Water, Distilled	15230170
Water For Injection (WFI) for Cell Culture	A1287301
Pluronic™ F-68	24040032
Penicillin-Streptomycin	15140122
Gentamicin (50 mg/mL)	15750060
Antibiotic-Antimycotic (100X)	15240096
Amphotericin B	15290018
Trypan Blue Solution, 0.4%	15250061
Countess™ 3 Automated Cell Counter Starter Package + REX Extended Warranty	A250298

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Life Technologies Corporation | 3175 Staley Road | Grand Island, New York 14072 USA

For descriptions of symbols on product labels or product documents, go to [thermofisher.com/symbols-definition](https://www.thermofisher.com/symbols-definition).

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. MAN0007366

Revision	Date	Description
2.0	21 March 2022	Updated formatting and legal boilerplates. Removed Cat. No. 12040077.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2022 Thermo Fisher Scientific Inc. All rights reserved. Pluronic is a trademarks of BASF Corporation. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.