



Contents

Catalog Number 12183020

Components	Amount	Storage
Lysis Buffer	20 mL	Room temperature
Wash Buffer I	10 mL	
Wash Buffer II	4 mL	
RNase-free Water	3 mL	
Spin Cartridges (with collection tubes)	10 each	
Collection Tubes	10 each	
Recovery Buffer	10 each	



Product description

- The PureLink™ RNA Mini Kit provides a simple, reliable, and rapid method for isolating high-quality total RNA from a wide variety of samples, including animal and plant cells and tissue, blood, bacteria, yeast, and liquid samples. The purified total RNA is suitable for use in a variety of downstream applications.
- Isolate up to 1 mg of nucleic acid.



Required materials

- 2–mercaptoethanol
- 100% ethanol
- (Optional) PureLink™ DNase Set (Cat. No. 12185010)
- Microcentrifuge capable of 12,000 × g
- 1.5–15 mL round-bottomed RNase-free tubes
- RNase-free pipet tips



Online resources

- Visit our [product pages](#) for protocols, safety, and additional product information.
- Go online to view related [PureLink™ products](#).
- For support, visit [thermofisher.com/support](https://www.thermofisher.com/support).

Before first use of the kit

Prepare Wash Buffer II

Add 16 mL 96–100% ethanol to Wash Buffer II. Check the box on the Wash Buffer II label to indicate that ethanol was added. Store Wash Buffer II with ethanol at room temperature.

(Optional) Prepare PureLink™ DNase Mixture

Add the following components to a clean, RNase-free microcentrifuge tube. Prepare 80 µL for each sample to be processed. Store PureLink™ DNase Mixture at –20°C.

Component	Volume
10X DNase I Reaction Buffer	8 µL
Resuspended DNase (~3U/µL)	10 µL
RNase free water	62 µL

Before each use of the kit

Prepare fresh Lysis Buffer

Add 10 µL of 2–mercaptoethanol for each 1 mL of Lysis Buffer needed for the purification procedure in a RNase-free tube.

Use 0.2 mL of Lysis Buffer for every 0.2 mL of whole blood (or fraction thereof).

Troubleshooting

For detailed troubleshooting instructions see the [PureLink™ RNA Mini Kit User Guide](#) at [thermofisher.com](https://www.thermofisher.com) or contact Technical Support.

Limited product warranty and licensing information

Disclaimer: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Manufacturer: Life Technologies Corporation | 2130 Woodward Street | Austin, TX 78744

©2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.



Protocol for purification of RNA from whole blood

This protocol describes how to purify one sample of total RNA from ≤ 0.2 mL of fresh or frozen whole blood on one column using the PureLink™ RNA Mini Kit. For detailed instructions see the PureLink™ RNA Mini Kit User Guide at [thermofisher.com](https://www.thermofisher.com) or contact Technical Support.




Important guidelines




- Use proper RNA handling techniques when working with RNA.
- Collect whole blood in the presence of anticoagulants such as EDTA or sodium citrate and store at 4°C until ready for use.
- Freshly drawn blood can be used without anticoagulants.
- Processing frozen blood may result in lower yield and quality of RNA.

Step		Action
1	 Lyse cells	<ol style="list-style-type: none">Add ≤ 0.2 mL of whole blood sample to a 1.5 mL RNase-free microcentrifuge tube.Add 0.2 mL of Lysis Buffer with 2-mercaptoethanol to the sample.Vortex the lysate to completely disrupt and lyse the blood cells.Centrifuge at $12,000 \times g$ for 2 min.Transfer the supernatant to a clean 1.5 mL RNase-free microcentrifuge tube.
2	 Bind RNA	<ol style="list-style-type: none">Add 1.5 volumes of 100% ethanol to the sample.Disperse any precipitate by vortexing or pipetting up and down several times.Transfer the sample (including any remaining precipitate) to a Spin Cartridge (with Collection Tube).Centrifuge at $12,000 \times g$ for 15 sec at room temperature. Discard the flow-through, and reinsert the Spin Cartridge into the same Collection Tube.

Protocol for purification of RNA from whole blood

Step		Action	
		No DNase treatment	On-column DNase treatment
3	 <p>Wash RNA</p>	<ul style="list-style-type: none"> a. Add 700 μL Wash Buffer I to the Spin Cartridge. b. Centrifuge at 12,000 \times g for 15 sec at room temperature. c. Discard the flow-through and the Collection Tube. Place the Spin Cartridge into a new Collection Tube. d. Add 500 μL Wash Buffer II with ethanol to the Spin Cartridge. e. Centrifuge at 12,000 \times g for 15 sec at room temperature. f. Discard the flow-through and reinsert the Spin Cartridge in the same Collection Tube. g. Repeat steps d–f one more time. 	<ul style="list-style-type: none"> a. Add 350 μL Wash Buffer I to the Spin Cartridge. b. Centrifuge at 12,000 \times g for 15 sec at room temperature. c. Discard the flow-through and the Collection Tube. Place the Spin Cartridge into a new Collection Tube. d. Add 80 μL PureLink™ DNase Mixture onto the surface of the Spin Cartridge membrane. e. Incubate at room temperature for 15 min. f. Add 350 μL Wash Buffer I to the Spin Cartridge. g. Centrifuge at ~2,600 \times g for 5 min at room temperature. h. Discard the flow-through and the Collection Tube. Place the Spin Cartridge into a new Collection Tube. i. Add 500 μL Wash Buffer II with ethanol to the Spin Cartridge. j. Centrifuge at 12,000 \times g for 15 sec at room temperature. k. Discard the flow-through and reinsert the Spin Cartridge in the same Collection Tube. l. Repeat steps i–k one more time.

Protocol for purification of RNA from whole blood

Step	Action
<p>4</p>  <p>Elute RNA</p>	<p>a. Centrifuge the Spin Cartridge with Collection Tube at $12,000 \times g$ for 1 min at room temperature.</p> <p>b. Discard the Collection Tube and insert the Spin Cartridge into a Recovery Tube.</p> <p>c. Add 30 μL to 3 \times 100 μL RNase-Free Water to the center of the Spin Cartridge.</p> <p>d. Incubate at room temperature for 1 min.</p> <p>e. Centrifuge at $12,000 \times g$ for 2 min at room temperature.</p>
<p>5</p>  <p>Analyze RNA yield and quality</p>	<p>Determine the quantity and quality of the purified total RNA using any of the following techniques (See the PureLink™ RNA Mini Kit User Guide for details).</p> <ul style="list-style-type: none"> ▪ UV absorbance at 260 nm ▪ Fluorescence microplate reader with Quant-iT™ RiboGreen™ RNA Assay Kit
<p>6</p>  <p>Store RNA</p>	<ul style="list-style-type: none"> ▪ Keep purified RNA on ice if using the RNA within a few hours of isolation. ▪ Store purified RNA at -80°C or long-term storage.