PRODUCT INFORMATION SHEET

Pub. No. MAN0000844


Rev. A.0



Contents

 Catalog No.
 Size

 10496016
 200 μL






Storage

- Product is shipped at ambient temperature.
- Store at –20°C.




Product description

- The Invitrogen™ High DNA Mass Ladder is designed for sizing and quantification of double stranded DNA on 0.8% to 1% agarose gels.
- The High DNA Mass Ladder consists of 6 individual chromatography-purified DNA fragments ranging in size from 1,000 bp to 10,000 bp.
- The ladder is supplied with 10X BlueJuice[™] Gel Loading Buffer for sample DNA.





Required materials





Online resources

- Visit our product pages for additional information and protocols.
- Go online to view related DNA ladders and markers.
- For support, visit thermofisher.com/support.



Important guidelines

- Do not heat the High DNA Mass Ladder before loading.
- Load the same volume of DNA sample and DNA ladder.
- For quantification, adjust the concentration of the sample to equalize it approximately with the amount of DNA in the nearest band of the ladder.
- For DNA bands visualization with GelRed[™] use gel staining after electrophoresis to avoid aberrant DNA migration.



Guidelines for agarose gel preparation

 Determine the required agarose concentration for your gel based on the size of DNA fragments to be separated.

Fragment size	Recommended agarose gel %		
	1X TAE	1X TBE	
800-10,000	0.8	0.7	
400-8,000	1.0	0.85	
300-7,000	1.2	1.0	

- Prepare agarose in a flask with 2-4 times the volume of the agarose solution.
- Exercise caution when handling microwaved agarose. The solution may become superheated and foam over when agitated.
- Refer to the product insert for UltraPure[™] Agarose for detailed instructions on agarose preparation.
- Guidelines for staining gels
- Troubleshooting
- Limited product warranty and disclaimer details



Prepare DNA ladders and samples for electrophoresis

Step			Action		
1	THE THE PARTY OF T	Cast agarose gel	a. Prepare agarose solution (w/v) for the gel percentage appropriate for separating your DNA fragments.b. Microwave agarose solution.c. Cast agarose gel.		
2		Prepare DNA ladder	a. Thaw, mix and briefly centrifuge each component before use. b. Add the following components to prepare enough ladder for a single 5 mm well. Component DNA ladder [1] 4 µL (520 ng) 10X BlueJuice™ Gel Loading Buffer Water, nuclease free [1] Scale components up or down depending upon width of wells. Modify volume by 0.2 µL (0.1 µg of DNA) for each 1 mm of width. c. Mix gently. d. Load DNA ladder on gel.		
3		Prepare samples	 a. Dilute your sample with 10X BlueJuice™ Gel Loading Buffer (Cat. no. 10816015): mix 1 volume of loading dye with 9 volumes of the DNA sample. b. Mix gently. c. Load DNA ladder on gel. 		
4		Perform electrophoresis	a. Add appropriate amount of UltraPure TAE or U b. Set appropriate voltage and perform electrophe DNA size		
			1-5 kb 4–10 V/cm >5 kb 1–3 V/cm	TAE or TBE	
5		Stain agarose gel	 a. Incubate gel in staining buffer for 30 minutes. b. Visualize DNA ladder and samples. • Use UV transilluminator to detect DNA bands stained with ethidium bromide. • Use blue light transilluminator to detect DNA bands stained with SYBR™ stains. 		

For support, visit thermofisher.com/support.