

# ElectroMAX™ DH5α-E™ Cells

Cat. No. 11319-019

Size: 0.5 ml

Store at -80°C

(Do not store in liquid nitrogen)

## Description

ElectroMAX™ DH5α-E™ Cells are *E. coli* cells which can be transformed by electroporation (1,2). These cells can only be transformed by electroporation and are **not** transformed by "heat shock". DH5α-E™ Cells are suitable for the generation of cDNA libraries using pUC or pBR322-derived vectors. The  $\phi$ 80*lacZ*ΔM15 marker provides α-complementation of the β-galactosidase gene allowing blue/white screening on agar plates containing X-gal or Bluο-gal. DH5α-E™ allows efficient transformation of large plasmids and can also serve as a host for M13mp cloning vectors if a lawn of DH5αF™, DH5αF'IQ™, JM101, or JM107 is provided to allow plaque formation. We do not recommend using this strain for direct cloning of methylated genomic DNA.

Component	Amount
DH5α-E™ Cells	5 x 100 μl
pUC19 DNA (10 pg/μl)	50 μl
S.O.C. Medium	2 x 6 ml

## Genotype

F<sup>-</sup>  $\phi$ 80*lacZ*ΔM15 Δ(*lacZYA-argF*) U169 *recA1 endA1 hsdR17*

(r<sub>k</sub><sup>-</sup>, m<sub>k</sub><sup>+</sup>) *gal<sup>-</sup> phoA supE44 λ<sup>-</sup> thi<sup>-1</sup> gyrA96 relA1*

## Quality Control

ElectroMAX™ DH5α-E™ Cells are tested for transformation efficiency using the protocol on the next page and the following electroporator conditions: 2.0 kV, 200 Ω, 25 μF. Transformation efficiency should be >1.0 × 10<sup>10</sup> transformants/μg of pUC19 DNA.

Part No. 11319019.pps

Rev. Date: 25 October 2006

## Transformation Procedure

pUC19 control DNA (10 pg/ $\mu$ l) is provided to check transformation efficiency. Use experimental DNA that is free of phenol, ethanol, salts, protein, and detergents to obtain maximum transformation efficiency.

1. Add DNA to microcentrifuge tubes.
  - A. To determine transformation efficiency, add 1  $\mu$ l of the pUC19 control DNA to a microcentrifuge tube.
  - B. For ligation reactions, precipitate the reactions with ethanol and resuspend in TE Buffer (10 mM Tris HCl, pH 7.5; 1 mM EDTA). The concentration of resuspended DNA should not exceed 100 ng/ $\mu$ l. Add 1  $\mu$ l of the DNA to a microcentrifuge tube (see Note 1).
2. Thaw ElectroMAX™ DH5 $\alpha$ -E™ cells on wet ice.
3. When cells are thawed, mix cells by tapping gently. Add 20  $\mu$ l of cells to each chilled microcentrifuge tube.
4. Refreeze any unused cells in a dry ice/ethanol bath for 5 minutes before returning them to the -80°C freezer. Do not use liquid nitrogen. Although the cells are refreezable, subsequent freeze-thaw cycles will decrease transformation efficiency.
5. Pipette the cell/DNA mixture into a chilled 0.1 cm cuvette and electroporate. If you are using the BTX® ECM® 630 or BioRad GenePulser® II electroporator, we recommend using the following electroporator conditions: 2.0 kV, 200  $\Omega$ , 25  $\mu$ F (see Note 2).
6. To the cells in the cuvette, add 1 ml of S.O.C. medium and transfer the solution to a 15 ml snap-cap tube (*e.g.* Falcon™ tube).
7. Shake at 225 rpm (37°C) for 1 hour.

8. Dilute cells transformed with pUC19 control DNA 1:100 with S.O.C. medium. Spread 50  $\mu$ l of the dilution on prewarmed LB plates containing 100  $\mu$ g/ml ampicillin.
9. Dilute experimental reactions as necessary and spread 100-200  $\mu$ l on selective plates.
10. Incubate plates overnight at 37°C.

### **Growth of Transformants for Plasmid Preparations**

Grow ElectroMAX™ DH5 $\alpha$ -E™ Cells which have been transformed with a pUC-based plasmid overnight at 37°C in TB (3). A 100 ml culture in a 500 ml baffled shake flask will yield approximately 1 mg of pUC19 DNA.

### **Notes**

1. Transformation efficiencies will be 10- to 100-fold lower for ligation mixtures and cDNA than for an intact control plasmid such as pUC19. **Salts and buffers severely inhibit electroporation.** Ligation reactions can be diluted 5-fold, and 1  $\mu$ l added to 20  $\mu$ l of cells. For optimal results, precipitate ligation mixtures with ethanol prior to transformation. Use only 1 to 2  $\mu$ l of the resuspended DNA per 20  $\mu$ l reaction. Adding undiluted ligation mixtures or too high a volume of DNA decreases transformation efficiency and increases the risk of arcing.
2. If you are using an electroporator other than a BTX® ECM® 630 or BioRad GenePulser® II electroporator, you may need to vary the settings to achieve optimal transformation efficiency.

3. Transformation efficiency (CFU/ $\mu\text{g}$ ):

$$\frac{\text{CFU on control plate}}{\text{pg pUC19 DNA}} \times \frac{1 \times 10^6 \text{ pg}}{\mu\text{g}} \times \frac{\text{volume of transformants}}{\text{volume plated}} \times \frac{\text{dilution}}{\text{factor}}$$

For example, if 10 pg of pUC19 yields 50 colonies when 50  $\mu\text{l}$  of a 1:100 dilution is plated, then:

$$\text{CFU}/\mu\text{g} = \frac{50 \text{ CFU}}{10 \text{ pg}} \times \frac{1 \times 10^6 \text{ pg}}{\mu\text{g}} \times \frac{1 \text{ ml}}{0.05 \text{ ml plated}} \times 10^2 = 1.0 \times 10^{10}$$

**References**

1. Calvin, N. M., and Hanawalt, P. C. (1988) *J. Bacteriol.* 170, 2796.
2. Dower, William J., *et al.* (1988) *Nucl. Acids Research* 16, 6127.
3. Tartof, K. D. and Hobbs, C. A., (1987) *Focus*<sup>®</sup> 9:2, 12.

©2003 - 2006 Invitrogen Corporation. All rights reserved. For research use only.  
Not intended for any animal or human therapeutic or diagnostic use.