Lovibond® Water Testing

Tintometer® Group

Instruction Manual Addendum

BD 600 GLP

(EN) Instruction manual.....	3	(DE) Gebrauchsanleitung.....	45
(ES) Manual de instrucciones.....	88	(FR) Mode d'emploi.....	132
(PT) Manual de instruções.....	177	(IT) Manuale d'istruzione.....	220
(NL) Gebruiksaanwijzing.....	263		

Table of Contents

BD 600 GLP

1 Introduction	5
1.1 General Information	5
1.1.1 Read instructions before use	5
1.1.2 Purpose of the Product	5
1.1.3 Authorized use	5
1.1.4 Requirements for safe use	5
1.1.5 Keep instructions	5
1.1.6 User qualification	5
1.1.7 Handling of hazardous chemicals	6
1.1.8 Disposal notes	6
1.2 List of all used signs in the document	6
2 Product overview	8
2.1 Guide to symbols	8
2.2 Certification	8
2.3 Features	8
2.4 Product description	9
2.4.1 Description of the keys and switches	9
2.4.2 Interface description	10
3 Commissioning	12
3.1 Operating environment	12
3.2 Contents of delivery	12
3.3 Inserting and replacing the batteries	13
4 Operation	15
4.1 First start up	15
4.2 General operation principles	15
4.3 List of control elements and their function	17
4.4 Main menu	18
4.5.1 Start test series	18
4.5.2 Handling of measurement series	19
4.5.3 Stop measurement series	20
4.5.4 Show current values	21
4.5.5 Display test series	21
4.5.6 Export test series	23
4.5.7 Exporting data files from the rack memory to a PC	25
4.5.8 Data file structure	26
4.5.9 Options	28
4.6 Advanced operations	30
4.7.1 Handling	32
4.7.2 Notice	32
4.7.3 Initial Steps	32

4.7.4 Selection of the sample volume	32
4.7.5 Sample Preparation	32
4.8 Description of setup and handling of important accessories	35
5 Maintenance	37
5.1 Calibration	37
5.2 Decommissioning	37
5.3 Maintenance and cleaning	37
6 Troubleshooting	38
6.1 Error Messages and Notices	38
6.2 Interpretation of BOD curve and Problem Solving Strategies	40
7 Accessories & Replacement parts	41
7.1 List of Accessories	41
8.1 Specifications - Mains adapter	42
8.3 Specifications - BD 600 GLP	43
9 Appendix	44
9.1 List of Trademarks	44

1 Introduction

1.1 General Information

1.1.1 Read instructions before use

This manual provides important information on the safe operation of the product. Please read this manual carefully and familiarize yourself with the product before use.

1.1.2 Purpose of the Product

The Lovibond® sensor system BD 600 is a 6 sample system that allows precise measurements of biological oxygen demand (BOD) based on the manometric principle.

1.1.3 Authorized use

The manufacturer's liability and warranty for damage is voided with improper use, failure to follow this manual, use by unqualified personnel, or unauthorized changes to the product.

The manufacturer is not liable for costs or damages that arise from the user or third parties due to the use of this product, especially in cases of improper use of the product or misuse or faults in the connection of the product.

The manufacturer assumes no liability for print errors.

1.1.4 Requirements for safe use

Note the following points for safe use:

- If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.
- The product may only be used according to the authorized use specified above.
- The product may only be supplied with power by the energy sources mentioned in this operating manual.
- The product may only be used under the environmental conditions mentioned in this operating manual.
- The product must not be opened or modified.

The product must not be used if:

- it is visibly damaged (e.g. after being transported)
- it was stored under adverse conditions for a lengthy period of time (storage conditions, see chapter "Specifications")

1.1.5 Keep instructions

The manual must be kept in the vicinity of the product so you can always find the information you need.

1.1.6 User qualification

The operating personnel must be able to understand and correctly implement the safety labels and safety instructions on the packages and inserts of the products.

The user must be capable and able to read and understand this manual in order to familiarize themselves with the handling and to ensure safe use.

1.1.7 Handling of hazardous chemicals

Chemical and/or biological hazards may exist where this product is used. Abide by all governing laws, regulations and protocols when using this product.

For the development of products, Lovibond® pays close attention to safety. Some hazards from dangerous substances cannot be avoided. If self-produced tests or solutions are used, the responsibility concerning any risks caused by those tests or solutions lies with the user (personal responsibility).

EN

1.1.8 Disposal notes

Dispose of the batteries and electrical devices at a suitable facility according to local legal requirements.

It is illegal to dispose of the batteries with household waste.

Within the European Union, the batteries are removed at a specialized treatment center at the instrument's end of life.

Instruments marked with this symbol must not be disposed of in normal domestic waste.

1.2 List of all used signs in the document

The following symbols are used in this manual to identify sections that require special attention:

Danger!

A hazard exists that will result in death or severe injury if not avoided.

Warning!

Improper handling of certain reagents can cause damage to your health. In any case follow the safety labels on the packing, the safety instructions of the package insert and available SDS. Protective measures specified there have to be followed exactly.

Caution!

A hazard exists that may result in minor or moderate injury.

Notice!

Important information or specific instructions need to be strictly followed.

2 Product overview

2.1 Guide to symbols

Labels attached to the product should be strictly observed to avoid personal injury or damage to the product. Refer to this chapter for information regarding the nature of the danger or risk before taking any action where such label is present.

For professional users in the European Union:

If you wish to discard electrical and electronic equipment (EEE), please contact your dealer or supplier for further information.

For disposal in countries outside of the European Union:

This symbol is only valid in the European Union (EU). If you wish to discard this product please contact your local authorities or dealer and ask for the correct method of disposal.

2.2 Certification

Device/Accessories	Test approvals, EMC
BOD Measurement System	CE, EMC in accordance DIN EN 61326, basic requirements
Stirrer unit	CE, DIN EN 61326:2013
Mains adapter	CE, EN 55024

2.3 Features

The BD 600 is the ideal system for process control or as a secondary test to the dilution method. Using a respirometric system solves many of the issues associated with the dilution method for BOD analysis.

- Significantly reduced sample preparation - it is easy to begin collecting data - simply add sample to the bottle, add nitrification inhibitor and place the sensor head on the bottle. There is no need for sample dilutions or seeding the sample. Because this is a respirometric test, concerns regarding air in the sample are eliminated.
- Operators do not have to be present at completion of tests. Measurements are automatically taken and stored at regular time intervals and can be recalled at your convenience - no more weekend trips to the lab or plant!
- Easy and quick interpretation of measurement values - In addition to displaying results of a sample directly in mg/l BOD, the large, integrated, display is able to show the BOD curve which makes it is easy to see and understand measurement data points and trends. In addition, each sensor head has an integrated LED. This allows operators to visually confirm which sample the data is being displayed for and reduces potential confusion or errors over which sample data is presented.
- Consistent measurement conditions - The automatic start function ensures all measurements are taken at an equal sample temperature. When the measurement conditions are consistent across readings, operators do not have to work about temperature or pressure fluctuations skewing results.

- Sample identification and tracking - The BD 600 interface, makes it easy to assign an identification code or sample name when setting up a new sample. This helps track a sample from the time it is collected and ensures traceability during sample measurement.
- Use the remote control - If using an incubator with a glass door, operators can use the supplied remote control to display measurement data for a sample - without even opening the door! This further ensures that samples stay at a consistent temperature which is essential to accurate results.

2.4 Product description

BOD measurement by means of pressure differential in a closed system (respirometric BOD measurement). The measurement system records a measurement in one selected of the four following modes with the following measurement intervals:

5 days I = every 2 hours, 28 days = every 12 hours, 60 days = every 24 hours and 90 days = every 24 hours the first 30 days and every 48 hours the last 60 days. In addition a separate measuring mode will be used during calibration of sensor heads: 5 days II = every six hours.


BOD5 measurement, BSB7 measurement, OECD 301F / biological degradability testing under GLP / further bio-technical applications in aqueous solutions. For use by qualified personnel only. The measuring time is adjustable in the pre-set configurations mentioned above (5 days, 28 days, 60 days and 90 days).

2.4.1 Description of the keys and switches

Key	Function
	ON/OFF key This key is used for switching the device on and off.
	Function keys The function keys have a different meaning in each menu. Info text above the keys indicates their meaning. If no text is shown above the key, it has no function.
	Backspace Characters which have been entered can be deleted with the backspace key.
	Number keys The number keys are used for entering the date, time, as well as file and head names.

Key	Function
	<p>Quick selection keys</p> <p>The quick selection keys take you directly to the corresponding menu each submenu:</p> <p>Quick selection key 1 ,Start' --> ,Start measurement series'</p> <p>Quick selection key 2 ,List' --> ,Display current values'</p> <p>Quick selection key 3 ,Graph' --> ,Display measurement series'</p>
	
	
	<p>Arrow keys</p> <p>They are provided mainly for navigation through the menus.</p>
	<p>Head keys</p> <p>Individual measurement points and/or heads can be selected in the submenus with the head keys.</p>

2.4.2 Interface description

1. SD card holder: The BD 600 GLP does not allow free access to the SD-card as it stores the life time data memory. The slot is closed with a plastic shield. To remove the card the plastic shield need to be damaged.

2. USB host interface

NOTICE!

The USB host interface is only intended for USB sticks. USB hubs, external hard disks, and adapter sticks (e.g. USB SD card adapters) are not supported.

- Battery operation entails the limitation that the host interface can only provide 200 mA.

3. USB device interface

The USB device interface is located under the USB host interface on the side of the front section of the housing. It makes it possible to access data saved on the SD card with a PC.

For this purpose, connect the device to a PC via the appropriate USB cable. The device must be switched on. If necessary, deactivate the auto-off option. If there is an SD card in the card holder, the device appears as mass storage on the PC.

4. Connection jack for the mains adapter.
5. Window for the IR receiver for the remote control.

3 Commissioning

3.1 Operating environment

Temperature equalisation is essential prior to biological testing and reproducible BOD measurements, as temperature has a major effect on biological activity. BOD measurements, for example, are always performed in a thermostatically controlled cabinet at a temperature of 20 °C. For temperature equalisation, we recommend Lovibond® thermostatically controlled cabinets with a user-selectable temperature from 2 °C to 40 °C.

EN

3.2 Contents of delivery

CAUTION!

Inspect items to ensure no damage has occurred during shipment. If there is damage or something is missing, please contact local distributor immediately.


Scope of delivery

BD 600 or BD 600 GLP	BD 606	
1	2	complete unit with 6 sensors and control unit with batteries (BD 600 GLP with certificate)
1	2	power supply unit incl. Y-cable for common power supply of instrument and stirring unit
1	1	remote control (without batteries)
1	2	inductive stirring unit
6	12	sample bottles
6	12	rubber gaskets
6	12	magnetic stirring rods
1	1	overflow flask, 157 ml
1	1	overflow flask, 428 ml
1	1	bottle, 50 ml potassium hydroxide solution
1	1	bottle, 50 ml Nitrification inhibitor
1	1	instruction manual
1*	2	USB cable

* BD 600 only

3.3 Inserting and replacing the batteries

The device can be supplied with power via the accompanying mains adapter or batteries. If the mains adapter is connected and batteries are inserted, the device is supplied via the mains adapter and the batteries are not used. If the mains adapter is removed during operation, a seamless transition to battery operation takes place automatically.

Remove all bottles from the rack in order to insert the batteries. Open the battery compartment (**A**) and insert three size C alkali-manganese cells (LR14 / ,baby cell') in the tubes provided for this purpose. The tube makes it easier to insert the batteries and prevent the batteries from jumping out again when you attempt to close the compartment.

Otherwise, if the battery polarity is reversed it can be damaged and leak and damage the device. Then insert the batteries with tubes into the battery compartment (ensure correct polarity).

WARNING!

Insert the batteries in the tube such that a positive pole always faces a negative pole.

WARNING!

The device is not designed for rechargeable batteries. Therefore, they must not be used. Rechargeable batteries can be damaged, leak, and damage the device.

CAUTION!

Do not handle depleted batteries or leaking material with bare hands! Wear protective gloves! Avoid contact with the eyes and skin!

**DANGER!**

Only handle the mains adapter with dry hands!

Protect the mains adapter from moisture.

Do NOT open the mains adapter!

Do not subject the mains adapter to great force!

Do NOT use the mains adapter if there is damage to the housing or the connection contacts!

Only use the mains adapter with an appropriate socket intended for this purpose!

Ensure that this socket is in fault-free condition!

The mains voltage and the mains frequency must always lie within the specified range of the mains adapter.

4 Operation


4.1 First start up

1. Place the instrument rack on the inductive stirring plate.
2. Use the Allen key to adjust the height between the rack and the inductive stirring plate. Or adjust the height later if continuous stirring is not achieved.
3. Use the Y-cable to connect the instrument and the inductive stirring plate to the power supply. Also, insert the batteries to prevent the unit from turning off in the event of a power failure.
4. Press the power button to start.
5. Select the desired language using the arrow keys and confirm with F2.
6. Set Date/Time using the arrow keys and confirm with F2.
7. The Lovibond logo appears briefly and will turn to the main menu. The instrument is now ready.

EN

4.2 General operation principles

Measuring Principle

Biochemical Oxygen Demand (BOD)

The biochemical oxygen demand (BOD) in water (e.g. waste water, surface water) is the amount of oxygen that is consumed during the degradation of organic substances through biochemical processes.

Measuring principle

Respirometric methods provide direct measurements of the oxygen consumed by microorganisms in the waste water from an air or an oxygen-enriched environment in a closed vessel under conditions of constant temperature and agitation. Carbon dioxide produced metabolically by the bacteria is chemically bound by the potassium hydroxide solution contained in the sealed cup in the bottle. The result is a pressure drop in the system, which is directly proportional to the BOD value and is measured by the pressure sensor. The BOD level is then displayed directly in mg/l.

General Operation Procedure (Brief Summary)

- Estimate the measurement range of the sample to be tested and use the sample volume as indicated in section "Handling".
- If necessary, pre-treat the sample as indicated in section "Measuring Principle" (e.g. adjust pH value and temperature, filtration, etc.).
- Measure the sample volume precisely with the overflow volumetric flask and pour into the BOD bottle (use a funnel, if necessary).
- If necessary, add the nitrification inhibitor as indicated in section "Handling".
- Place the magnetic stir bar in the BOD bottle.
- Fill the seal cup with 3 - 4 drops of KOH solution and place the seal cup in the test bottle.
- Screw the BOD sensors on the test bottles.
- Place the sample bottle in the bottle rack.
- Start the test (see section "Handling").
- Incubate the sample according to specifications (e.g. BOD₅ at 20 °C).

EN

4.3 List of control elements and their function

Display overview


Description	Function
1 Header	Displays the name of the menu (left) and date/time (right)
2 Display Area	Area that displays application steps of the corresponding menu
3 Footer	Displays the meaning of function keys F1 & F2 (left & right). If nothing is displayed, the corresponding keys have no function. Displays the icon for data storage (SD/USB) and power supply (battery/plug)

Icon overview

Icon	Description	Function
	SD Card*	Displays that SD card is recognized by the system.
	USB drive*	Displays that USB drive is recognized by the system.
	Battery power	Displays that system is powered by batteries and indicates battery life.
	Plug symbol	Displays that system is powered by mains adapter.

*If the icons are faded, the components are not recognized or connected by the system.

4.4 Main menu


Use the up ▲ and down ▼ arrow keys to select a submenu.

Use the right ► arrow key to access a submenu. Alternatively, use the quick selection keys "Start", "List" and "Graphic" to access the first three submenus.

To change settings at the BD 600 GLP that will affect measurements in a measurement series you need to login by using the general password ,2017' and a name (at least two letters) which will appear in the log file afterwards.

4.5.1 Start test series

After sample preparation, enter the submenu "Start test series" to begin a measurement. Use the up ▲ and down ▼ arrow keys to select a setting:

- Bottle position - Name
- Range - Sample volume / ATH
- Measurement duration / Time interval

Function

Description

Bottle Position – Name

When a sensor head is detected by the system, the bottle position and name are displayed on the screen. Use the left ◀ and right ► arrow keys or the head keys to select one of six possible sample bottles. The bottle selection can be tracked by a red LED on the top of the sensor head. The ,Name' of a bottle is just a label to make identification of samples easier. The GLP compliant sample tracing is ensured via the unchangeable head-ID locked in the measurement value data file, which will be saved for each measurement series on the rack. The name of a bottle can only be set prior to starting a new measurement series.

Range – Sample Volume/ATH

Select the expected BOD measuring range of your sample. The required sample volume and the number of drops for the nitrification inhibitor (ATH) are displayed. An overview of all BOD measuring ranges, required sample volumes and number of ATH drops is summarized under section "Handling".

Measurement duration – Time interval

The measurement system records a test series in one of the four selected modes with the following time intervals:

- 5 days l = every 2 hours
- 28 days = every 12 hours

Function	Description
	<ul style="list-style-type: none"> • 60 days = every 24 hours • 90 days = every 24 hours the first 30 days and every 48 hours the last 60 days <p>In addition a separate measuring mode will be used during calibration of sensor heads:</p> <ul style="list-style-type: none"> • 5 days II = every six hours
	The GLP version runs always with deactivated Autostart mode. So the A-symbol is always greyed out and in the submenu ,Options' is no setting for Autostart available.
Back	Press function key F1 to return to the main menu.
Start	Press the function key F2 to start the test series. It is only possible to start a new measurement series on a head if no other measurement series on that head is running and the head is valid for the chosen duration. In case of trying to start a measurement at a head where already a measurement series is pending the info ,On the device there is already a measurement running with the head ID' appears on the display and need to be confirmed to go ahead.

4.5.2 Handling of measurement series

In the BD 600 a sensor head stores all the data of the currently running measurement series. These data can be stored manually on a memory device like an USB stick or SD card. In difference to that the BD 600 GLP rack ensures that all measurement data of a head are stored automatically to the internal memory (internal SD card, that only allows read access and cannot be extracted or changed by user). To give you the most possible flexibility a head/bottle position can be changed on the rack even when a measurement series is active. This will not affect the consistency of data storage. In addition it is possible to start a measurement series at one rack, do several measurements there, put the head/bottle afterwards to different rack and let the measurement series continue on the new rack. To continue the already started measurement on the new rack it is important to tell the new rack that there is a new bottle which measurement series should be continued. To do so go into ,Start test series' menu and select the position of the bottle. The system is showing now the parameters of the running measurement series of that new bottle. Because this measurement series is still active it is not possible to change these parameter. To continue this measurement series at this rack press ,Start' (F2). After this the measurement series will be continued on this rack until it is over. A new data storage file for this measurement series will be created also on this rack. This will contain all setting for the series as well as all already measured values in this series even if those were performed on a different rack before. In such a case two data storage files of the measurement series exist. One series on the rack where the series has been started and one on the rack where the series will be finished. The first one will contain only measurement values until the time of replacing the bottle plus several entries of ,Sensor head not accessible' until it will be finished or it will be stopped manually. The second rack will contain all measurement data including one entry ,The test series has been continued' at that measurement entry that the first measurement on the new rack represents. So in consequence the complete data record of a measurement series

can always be found at that rack where a measurement series has been finished. In opposite to bottles with active measurement series bottles with an active calibration measurement cannot be continued on a different rack. Calibration measurements need to be performed completely at one rack. If a bottle/head with an active calibration measurement is been placed to a different rack the only option there is to stop that calibration series manually (see 5.5.2.2) to be able to use that head for new test series. If a bottle/head with a finalised measurement series has been placed on the rack and the ‚Start‘ (F2) has been pushed in this menu, the user can choose whether he likes to start a new measurement series by overwriting the data inside the head memory (the test series has been stored on the rack, where the series has been finished) or if he likes to download the finalised test series to the memory of this rack. The latter helps if one rack should be used as a data transfer hub to a PC system. After download the test series data file is on this rack and also on the rack where the measurement has been completed.

4.5.3 Stop measurement series

By pressing the backspace key in the menu ‚Start test series‘ the submenu ‚Stop test series‘ appears.

This sub menu displays all active measurement series on the rack. An active measurement series on the rack does not mean that the head/bottle that belongs to that series has to be still on the rack. A measurement series is active until it reaches its scheduled end or it will be stopped manually.

Each measurement series is listed by a number. This number is independent from the number of the bottle position on the rack!

So for example a measurement series started at bottle position 1 can get the number 6 in this list. This is because the position of a bottle on the rack can be changed freely. Because of the Head-ID the rack can identify the measurement series that belongs to each rack bottle position automatically and so it ensures correct data assignment and data storage.


If a bottle has been replaced from the rack a started measurement series will still be active and listed here until the scheduled end of the series. In such a case the data file of the measurement series gets an entry ‚Sensor head not accessible‘ for each measurement where the bottle is not back on the rack. After a replaced bottle has been returned to the rack, the rack will automatically continue the pending measurement series.

By pressing the number keys the details of the active measurement series of the rack can be displayed. If less than six series are active the list shows also empty entries. Only if less than six (measurement or calibration) series are active a new series can be started. It is only possible to run up to six measurement series at the same time. Calibration measurements of sensor heads are listed like a measurement series as well. They use always the time interval ‚5 days II‘ and can be identified by this. A measurement series that is not a calibration can be continued on a different rack (see 5.5.2.1).

To stop a measurement series manually select it by using the number keys. Afterwards press ‚Stop‘ (F2). The BD 600 GLP is asking for the password and your name (see 5.5 addendum). After confirmation the chosen series will be stopped. The data file of this series will contain the entry ‚Stopped by operator‘ plus the entered name.

A measurement series that is complete contains the entry ‚Test series stopped regularly‘ below the last measurement value and stops automatically.


4.5.4 Show current values

The latest measurements can be tracked in the submenu “Show current value”. All six bottle positions are displayed in the center of the screen. The latest value is displayed next to the sensor head if connected to the system.

Function	Description
	Sensor head with a completed test series is recognized.
	Sensor head with an on-going test series is recognized.
	Sensor head is not connected or recognized.
Back	Press function key F1 to return to the main menu.
Update	Press function key F2 to update the current value. Alternatively, press the head keys to update the measurement value of a specific bottle position.

4.5.5 Display test series


A schematic diagram of all bottle positions is displayed on the right site of the screen. The head keys can be used to select a bottle position. After, the measured data of the selected head are listed in mg/L with time stamp.

Use the up ▲ and down ▼ arrow keys to scroll through the measurements.

Use the left ◀ and right ▶ arrow keys to show additional information such as measurement parameters and head name.


Function	Description
	Sensor head is recognized and can be selected.
	Sensor head is recognized and selected.

Function	Description
	Sensor head is not connected or recognized.
x / y / z	<p>Example: "2 / 3 / 39"</p> <p>x: Selected data point in the test series; example "2": second measurement selected in the table.</p> <p>y: Number of data points already recorded during a test series; example "3": three measurements have already been recorded.</p> <p>z: Number of total measurements within the test series; example "39": thirty-nine measurements to be done</p>
Back	Press function key F1 to return to the main menu.
Graphics	Press the function key F2 to plot the measured data in a diagram (see below).

Use the head keys to select and deselect a sensor head to plot the measured data. A maximum of three plots can be displayed simultaneously for comparison.

Possible Errors

Underrange Underrange appears instead of a measurement value, if the recorded measurement value is below the starting measurement value.

Overrange Overrange appears instead of a measurement value, if he recorded measurement value is outside the measuring range.


Fail Fail appears instead of a measurement value, if no measurement could be conducted for the provided time; e.g. the sample was removed before the end of the measurement.

The eye symbol indicates that no usable measurement value was determined in at least one measurement within the test series.

The warning symbol indicates that a test series has been ended, although not all measurements have been performed.


(A) Underrange is shown with zero.

(B) Overrange is shown with the maximum display value within the measuring range.

4.5.6 Export test series

The BD 600 GLP can export measurement series data been stored in heads that are placed on the rack to an external USB memory device. On contrary to a normal BD 600 the GLP version cannot export these data to SD cards as this port is used for the internal life time memory.


The file exported to USB is very similar to the measurement series data file stored automatically inside the rack. But it does not contain the CPU-ID for identification of the rack and further entries that are stored together with a calibration measurement (batch no, test value and tolerance of the used test tablets). Instead of these it tells you the bottle position on the rack.

A schematic diagram of all bottle positions is displayed on the left side of the screen. Use the head keys to select and deselect the sensor heads for data export. Only selected heads are considered for data export. Use the up ▲ and down ▼ arrow keys to select between SD-Card and USB-drive if both are connected. Use the right ► arrow keys to continue with the data export.

EN

Function	Description
	Sensor head is recognized and can be selected.
	Sensor head is recognized and selected.
	Sensor head is not connected or recognized.
M: (SD)	Appears when SD-Card is connected to the instrument.
U: (USB)	Appears when USB-drive is connected to the instrument.
Back	Press function key F1 to return to the main menu.

After selecting a storage medium, existing files and folders are displayed. Use the up ▲ and down ▼ arrow keys to select a folder or a file.

Function	Description
 folder ►	When a folder is selected, use the left ◀ and right ► arrow keys to enter the folder or return to the previous entry.

Function	Description
 bod 05.01.2021.csv▶	When an existing file is selected, use the right ▶ arrow key to enter the status bar. The existing file can now be renamed using the number keys, the backspace key, the left ◀ and right ▶ arrow keys. Press F2 to save and overwrite the existing file. Use the up ▲ and down ▼ arrow keys to exit the status bar.
m:/bod 05.01.2021(1).csv	The destination folder is always displayed at the bottom of the list together with an automatically generated file name. The file name can be modified. Therefore, use the up ▲ and down ▼ arrow keys at the beginning or the end of the list.
	Folder and files names are displayed in lower-case letters. Folder and file names with special characters are not displayed, but are indicated by the warning symbol at the top right of the screen.
Back	Press function key F1 to return to the previous menu.
Save	Press function key F2 to save the measured data in csv-format onto the storage drive. Only data of selected sensor heads are saved.

4.5.7 Exporting data files from the rack memory to a PC

A BD 600 GLP rack has a life time memory to store all measuring data and calibration data produced from heads that have been placed to that rack and the log files of the rack.

To export these data files to a PC the BD 600 GLP and the PC need to be connected via an USB cable. After this it is necessary to set the rack into the ‚data transfer mode‘. Switch the rack off and hold afterwards the list key while pressing the on/off key. The rack will reboot now and display that the ‚data transfer mode‘ is active. When the rack has been connected to a PC via USB interface, the internal memory of the BD 600 GLP now is visible as an external storage device with read access only. All measurement series data files can be found inside the folder ‚data‘. All calibration series data files can be found inside the folder ‚calib‘. And the log files that contains all changes of parameters and options of the rack can be found inside the folder ‚log‘. You can all desired files copy now from the rack to the PC.

As long as the ‚data transfer mode‘ is been active, no measurement will be performed on the rack. After 90 minutes the data transfer mode ends automatically and the system returns back into the normal operation mode. If you want to go back from ‚data transfer mode‘ to normal operation mode manually just press the list key .

To prevent that the system is missing to execute a measurement during being in the ‚data transfer mode‘ the transfer mode should be started at a time where the next scheduled measurement is longer than 90 minutes away. To display when the next measurement will happen on the rack press the on/off key just once before you activate the ‚data transfer mode‘. The data transfer mode is needed to ensure that all data files that you like to copy are not under access of the normal measuring routine at the same time to prevent corrupted data.

4.5.8 Data file structure

The BD 600 GLP stores all values of a measurement series and all values of a calibration series in dedicated files on the internal memory. The measurement series data are stored inside the folder ,data' and the calibration series data are stored inside the folder ,calib'.

The measurement data files are named as followed:

nnnnn_HEAD-ID_name_date_starting-time.csv

nnnnn is a steadily increasing number that is rack dependent

HEAD-ID is the unchangeable identification number of the used sensor head

name is the name that was given to a sample

date and starting time is linked to the start of the measurement series

Example of a measurement data file. *.csv as seen in Excel (semicolon used as separator):

CPU-ID:	335556620 - 2924158247 - 1417818825 - 4110425221
Head ID:	65535-255-65535
Head Name:	Example
Range:	0-400 mg/L
Sample volume / ATH:	157 mL / 5 drops
Measurement Duration:	5 days
Time Interval:	5 days I
Date (format):	31.12.2021
Time (format):	24h

Measurement	Date	Time	BOD \ mg/L
0	17.05.2021	16:15	0
1	17.05.2021	18:15	5
The test series has been continued.	17.05.2021	19:18	-
2	17.05.2021	20:15	8
3	17.05.2021	22:15	10
4	18.05.2021	00:15	13
5	18.05.2021	02:15	17
6	18.05.2021	04:15	19
7	18.05.2021	06:15	19
8	18.05.2021	08:15	19

Measurement	Date	Time	BOD \ mg\L
9	18.05.2021	10:15	20
10	18.05.2021	12:15	20
Stopped by operator.	18.05.2021	13:50	Operator X

The file contains also parameters and measurement values of a measurement series and in addition the unchangeable identification number of the rack (CPU-ID) and messages related to the measuring series.

The calibration data files are named as followed:

nnnnn_calib_HEAD-ID_name_date_starting-time.csv

nnnnn is a steadily increasing number that is rack dependent

calib is a phrase tagged to a calibration data file for easy identification vs. measurement data

HEAD-ID is the unchangeable identification number of the used sensor head

name is the name that was given to a sample (see 5.5.2 and 5.6.4)

date and starting time is linked to the start of the measurement series


Example of a calibration data file. *.csv as seen in Excel (semicolon used as separator):

CPU-ID:	335556620 - 2924158247 - 1417818825 - 4110425221
Head ID:	65535-255-65535
Head Name:	Calibration Example
Range:	0-400 mg/L
Sample volume / ATH:	157 mL / 5 drops
Measurement Duration:	5 days
Time Interval:	5 days II
Date (format):	31.12.2021
Time (format):	24h
Name:	Operator X
Batch No.:	P07A
Test value:	330 mg\L
Tolerance:	30 mg\L

Measurement	Date	Time	BOD \ mg\L
0	12.05.2021	09:01	0
1	12.05.2021	15:01	250
2	12.05.2021	21:01	311
3	13.05.2021	03:01	325
4	13.05.2021	09:01	325
5	13.05.2021	15:01	326
6	13.05.2021	21:01	326

Measurement	Date	Time	BOD \ mg/L
7	14.05.2021	03:01	326
8	14.05.2021	09:01	326
9	14.05.2021	15:01	326
10	14.05.2021	21:01	326
11	15.05.2021	03:01	326
12	15.05.2021	09:01	326
13	15.05.2021	15:01	326
14	15.05.2021	21:01	326
15	16.05.2021	03:01	326
16	16.05.2021	09:01	326
17	16.05.2021	15:01	326
18	16.05.2021	21:01	326
19	17.05.2021	03:01	326
20	17.05.2021	09:01	326
Test series has stopped regularly.	17.05.2021	09:01	-
Head is valid. Next test before: 11.11.2021	17.05.2021	09:01	-

4.5.9 Options

Instrument settings can be changed in the “Options” submenu. Use the up ▲ and down ▼ arrow keys to select a setting. Use the left ◀ and right ▶ arrow keys to change the selected setting. Use function key F2 to accept the new settings or function key F1 to return to the main menu.

Function	Description
Auto start	<p>The autostart function is generally be disabled in the BD 600 GLP version.</p> <p>This means the system will directly begin to save all measurement values and will not wait until pressure drops slightly inside a bottle.</p> <p>In consequence to this a sample need to have the temperature of the incubator already before the measurement starts. If the temperature of a sample would be higher, a pressure decrease caused by cooling down of the gas atmosphere above the sample would create a slight pressure decrease that would create an offset to the measurement values. In opposite a lower temperature of the sample than in the incubator would create a slight expansion of</p>

	the gas atmosphere above the sample, which would lower the result a bit.
Brightness	Adjust the background lighting of the screen in 10 % increments from 0 % to 100 %.
Contrast	Adjust the contrast of the screen in 10 % increments from 0 % to 100 %.
Name heads	The name of a bottle can only be set prior to starting a new measurement series in the ‚Start test series‘ menu. The menu entry to change the name in the ‚Options‘ menu is not available in the GLP version of the BD 600.
Calibration	<p>The BD 600 GLP contains an additional entry in the ‚Options‘ menu to start a calibration procedure. Every sensor head need to be checked regularly (after half a year) if it is working within its specifications. This is done by using the calibration procedure. After selection of this entry in the ‚Options‘ menu choose the bottle with the head that you like to check.</p> <p>To perform a calibration use a test kit (part no. 2418328) and work according to the test kit instructions. In difference to the instruction of the test kit use this calibration mode and not the normal BOD measurement mode.</p> <p>After a sample bottle with the head to be tested has been prepared select the position of that bottle on the rack and press Start (F2).</p> <p>To get access the password and a name that will be documented in the calibration value data file need to be put in. Afterwards the software is requesting the batch number of the test tablets (see tablet packaging), the reference value that should be reached during the test (see also tablet packaging) and the acceptable tolerance value (see the value after the sign at the tablet packaging). Put in all these values and select the correct range that fits for the used test kit tablets (see instruction of the test kit, typically 0-400 mg/l). After confirmation a calibration series will start on the rack.</p> <p>The calibration series can be identified as it uses the ‚5 days II‘ interval (to display data of an active series).</p> <p>After a calibration procedure has been performed successfully (measured value meets the reference value within given tolerances) the sensor head gets a new validity status for half a year (183 days). After that period of time the calibration procedure needs to be performed again.</p> <p>If the calibration procedure fails, the sensor head is no longer valid and gets a fictive validity date that is in the past (31.12.2021). The calibration data file gets the entry ‚Head is invalid. Next test before: 31.12.2021‘</p>

A measurement series can only be started with a sensor head that has got sufficient time left, before its validity status will end.

In contrary to a normal measurement series a calibration series cannot be continued on a different rack and must be finished at the same rack where it has been started.

Language	Change the display language
Date/Time	Change date and time
Date (format)	Change the representation of the date in the header. The following options are available: <ul style="list-style-type: none"> • DD.MM.YYYY / 31.12.2021 • MM.DD.YYYY / 12.31.2021 • YYYY.MM.DD / 2021.12.31
Time (format)	Change the representation of the time in the header. The following options are available: <ul style="list-style-type: none"> • 12 h • 24 h
Auto-Off	Use the auto-off function to save power. Select between 3, 5 or 10 minutes for the unit to switch off automatically at the selected time. When the Auto-Off is deactivated, the unit runs continuously.
Remote Control	Select "on" or "mains supply" to activate the IR interface and enable the use of the remote control. If "mains supply" is selected, the IR interface is only active when the instrument is powered by DC.
Device ID	Assign a Device ID to operate the unit using the remote control
Update	To ensure data consistency and integrity of the BD 600 GLP user cannot perform a software update. This need to be done by the manufacturer or involved dealer


NOTICE!

To prevent loss of stored test results store or print them out before performing an Update. If the update procedure is interrupted (eg. interruption of connection, LoBat., etc.) the instrument isn't able to work (no display). The instrument will only work again after completing the data transfer.

4.6 Advanced operations

Remote control

Purpose

With the remote control, the submenus "Show current values", "Display test series" and "Export test series" can be accessed during operation or when the unit is switched off. This enables remote monitoring of a test series, e.g. when the BD 600 GLP measuring system is placed in an incubator for temperature control.

Initial Set-Up

1. The remote control is delivered without a battery. Before use, insert a lithium button cell (e.g. CR2025).
2. Activate the remote control in the "Options" menu.
3. Assign a Device ID to the instrument in the "Options" menu. If multiple instruments are used, assign different Device IDs.

EN

Device ID

13.00
05.01.2021

01

SID USB 1

How to use the remote control

1. Point the remote control at the instrument. Press the On key to switch on the instrument. The Device ID appears on the screen. (If other instruments are nearby and appropriately configured, their screen will also display the Device ID)
2. Use the number keys to enter the Device ID of the desired instrument. (If no or an incorrect ID is entered, the instrument will return to its previous status)
3. The submenu "Show test series" will appear on the screen.
4. The remote control can now be used to operate the unit in the same way as the keypad.

Press the Off key to switch off the instrument or press the On key to access the Device ID screen in order to continue with another instrument.

4.7.1 Handling

4.7.2 Notice

Notice!

The safety datasheets of the chemicals comprise all instructions on safe handling, occurring hazards, preventive actions and actions to take in hazardous situations.

4.7.3 Initial Steps

Place the BOD system, consisting of rack and inductive stirring plate, in an incubator for temperature control. Set the temperature of the incubator to 20 ± 1 °C as recommended for the determination of BOD according to EN 1899. Connect the system to the mains supply.

4.7.4 Selection of the sample volume

The expected BOD value of the sample defines the necessary sample volume for the experiment. Select a measuring range that results in a BOD value in the upper half of the measuring range, e.g. for a BOD of 150 mg/L select the measuring range 0 - 200 mg/L. If the expected BOD value is unknown, an estimate can be made on the basis of the COD value (COD = chemical oxygen demand), e.g. for domestic waste water the BOD₅ is approx. 80 % of the COD value.

BOD range in mg/L	Sample volume in mL	Nitrification inhibitor ATH dosage
0 – 40	428	10 drops
0 – 80	360	10 drops
0 – 200	244	5 drops
0 – 400	157	5 drops
0 – 800	94	3 drops
0 – 2000	56	3 drops
0 – 4000	21.7	1 drop

4.7.5 Sample Preparation

1. Use clean equipment for sampling.
2. Ensure the use of representative samples.
3. Test the pH value of the sample. The ideal pH value must be in the physiological range of the microorganisms, which is usually between pH 6.5 and 7.5. Any greater deviations provide a lower BOD value. If the sample pH is too high, use diluted hydrochloric acid (1 M) or diluted sulphuric acid (1 M) to adjust the pH. If the pH value is too low, use a sodium hydroxide solution (1 M) to adjust the pH.

4. Prepare the sample in accordance with the respective normative requirements so that it can be used as settled, stirred, filtered or even homogenized with a blender. It is recommended to test each sample two to three times.
5. The sample must be brought to the desired measuring temperature ± 1 °C.
6. Add a clean magnetic stir bar to the BOD sample bottle for continuous agitation.
7. Measure the required sample volume depending on the expected measuring range (see "4.7.3 Selection of the sample volume") with the corresponding overflow flask and fill the exact sample volume into the BOD sample bottle. Use a funnel if necessary.
8. Depending on the sample volume, add the recommended number of drops of nitrification inhibitor (ATH) to the BOD sample bottle (see "4.7.3 Selection of the sample volume").
9. Fill the dry rubber gasket with 3 - 4 drops of 45 % potassium hydroxide solution. Then place the rubber gasket on the BOD sample bottle.
10. Screw on the sensor head hand-tight. The rubber gasket ensures the necessary sealing. (Do not use grease or other lubricants!).
11. Place the prepared BOD sample bottle into the pre-tempered rack inside the incubator at 20 ± 1 °C (see "4.7.2 Initial Steps"). (Note: The BD 600 has an optional autostart function which allows the use of samples with a temperature of 15 to 21 °C. With the auto start function switched on, the system checks within the first three hours if a pressure drop is recognized by the sensor head. If so, the measurement will start immediately or after the three hours.)
12. Check if the sensor head is recognized by the system and if the sample is continuously stirred. If necessary, adjust the 4 ball screws on the rack.
13. Start the test series (see "4.5.1 Start test series")
14. Incubate the sample in accordance with the respective normative requirements (e.g. BOD₅ at 20 °C for 5 days).

Additional Information:

Sample water

- Domestic wastewater usually contains sufficient nutrient salts and suitable microorganisms and no toxic or retarding substances. Therefore, the samples can be used undiluted without the need for additional nutrients or seeding by microorganisms.
- Industrial wastewater may contain inappropriate microorganisms, toxic and inhibitory substances as well as insufficient nutrient content. Such samples must be diluted with dilution water containing a sufficient amount of nutrients to achieve an adequate BOD:N:P ratio and to eliminate the inhibitory effect of toxic substances. Furthermore, an additional seeding by a microbial source (inoculum) is often needed. An application paper on strongly-loaded organic waste water can be downloaded on www.lovibond.com.

Potassium Hydroxide

When oxygen is consumed by the microorganisms to form carbon dioxide during the degradation of organic compounds, there is no direct change in pressure. Potassium hydroxide in the rubber gasket and carbon dioxide react chemically to form potassium carbonate: $2 \text{KOH} + \text{CO}_2 \rightarrow \text{K}_2\text{CO}_3 + \text{H}_2\text{O}$

This removes the carbon dioxide formed from the gas phase creating a measurable negative pressure which correlates with the oxygen consumption of the microorganism and ultimately with the measured BOD.

Nitrification Inhibitor

Nitrifying bacteria also consume oxygen. This consumption can already occur within the first five days, especially in samples with low BOD values. In wastewater monitoring, often only the BOD from the degradation of carbon compounds is determined (carbon BOD), with the addition of allylthiourea (ATH) nitrification is inhibited to prevent oxygen consumption from the degradation of nitrogen compounds. If the oxygen consumption in the course of the nitrification (N-BOD) should be determined, a comparison of the sample with and without nitrification inhibitor can be done. The difference between both BOD values will correspond to the oxygen requirement of nitrifying bacteria.

Biological degradation Test OECD 301F

Sample handling and preparation of biological degradation testing see OECD 301F standard.

As some substances are easy to degrade and may produce a high amount of carbon dioxide during the period of 28 days testing, it is necessary to use more KOH solution as described for BOD measurement. It is also possible to use sodium hydroxide or potassium hydroxide pellets. Use at least 140 mg sodium hydroxide or 200 mg potassium hydroxide per bottle. This is about one to two pellets. In case of use of the KOH solution dose at least 10 drops into the rubber gasket.

Start 28 days or longer testing only by using the 6.5 cm long rubber gaskets instead of the 4.5 cm rubber gaskets that are also available for the BD 600. For the OECD biodegradability testing choose the sample volume according to the ThOD value (theoretical oxygen demand) of the solution that should be tested.


It is only possible to use sensor heads, which take a validity status for the whole planned measuring period. If the validity status ends before planned end of measurement, the sensor head will not be accepted and need to be validated before further use.

It is only possible to start a new measurement at a rack position where no measurement series is currently active. Such measurements need to be finished at first or be stopped manually before a new measurement series can be started at that rack position.

To get the best results it is important that the samples and bottles are at the temperature inside the incubator already before the measurement is started and the rack is been placed into the incubator.

4.8 Description of setup and handling of important accessories

Inductive stirring system

1. Mains adapter
2. DC connection cable
3. Stirrer drive
4. Stirring position
5. Y-cable

EN

Device description and functional description

The inductive stirring system is designed for the stirring of liquids in special BOD bottles. It comprises a super-flate stirrer drive with 6 stirring positions and the mains adapter. The inductive stirrer drive does not have a motor and is thus wear-free. It is especially well-suited for stirring during continuous operation in thermocabinets. It is completely protected from dripping water due to its hermetically sealed encapsulation. It can be used in harsh environmental conditions.

The wide recess between the stirring positions also ensures greater ventilation of the containers in the thermocabinet. The stirrer drive is supplied with power via the mains adapter and has integrated control electronics. When switched on, a soft-start at reduced speed assures an even start-up and centring of the magnetic stir bars.

The electronically controlled automatic monitoring unit reduces the speed about every 40 seconds. In the process, each magnetic stir bar is re-centred in the bottle for a few seconds. Therefore, you can easily exchange individual bottles while the unit is switched on.

Due to the synchronous operation, mutual interference of the magnetic stir bars can be practically ruled out altogether.

Magnetic stir bars

Use the PTFE covered magnetic stir bars included in the scope of delivery.

**NOTICE!**

The length of the magnetic stir bars should not exceed 40 mm.

**NOTICE!****Stirring**

Do not place hot containers on the stirrer drive; maximum temperature: 56 °C.

Fill the BOD bottles.

Place a magnetic stir bar in each bottle.

Place the BOD bottles in the rack.

**Warning due to magnetism!**

The magnetic fields can influence parts that are sensitive to magnetic fields, magnetic parts, or metal parts (e.g. data carrying media, pacemakers, watches, etc.). Keep these parts away from the stirrer drive and magnetic stir bars.

5 Maintenance

5.1 Calibration

A test kit (Art. no. 2418328) is available for testing the BD 600.

The test kit enables the testing of all components. It comprises special reagent tablets which generate a defined vacuum in the closed BOD bottle.

5.2 Decommissioning

Disconnect the external power supply from the unit in order to decommission the system. Disconnect the mains adapter plug from the mains supply. Remove batteries from the measurement system. Remove seal caps from the test bottles and clean in the correct manner. Empty and clean test bottles correctly. Clean sensor heads. Store the stirrer drive and magnetic stir bars in a manner such that no hazards arise for parts which are sensitive to magnetic fields.

CAUTION!

Observe the SAFETY INSTRUCTIONS at the beginning of the manual.
Observe the local legal regulations for all emptying and cleaning.

5.3 Maintenance and cleaning

- High-quality metal alloys are used for the contacts between the BOD sensors and the bottle rack. Carefully clean the contacts as necessary with a soft cloth. To smooth out any unevenness, for optimisation of the contact between the BOD sensor and bottle rack, and for optimisation of the stirrer position, the BOD base unit has 4 adjusting screws on the bottom. Carefully clean the BOD base unit (including bottle rack) and the BOD sensors as necessary with a dry cloth. The parts coming into contact with the samples (BOD bottle, seal cup, magnetic stir bar) must be carefully cleaned after each test. Empty the bottles after the testing is complete – observe local regulations in the process – and rinse out repeatedly with hot water. Rinse thoroughly after using cleansers! Residue from cleansers can destroy the BOD test.
- The stirrer drive is maintenance-free. The magnetic drive is installed inside the housing to ensure that it is waterproof. Clean the surface of the stirrer drive regularly. The stirrer drive can be cleaned with cleansers or di-sinfecting solutions which are suitable for PVC. Wipe off the surface of the mains adapter with a dry cloth.

CAUTION!

The device may only be opened by an authorised service location in the case of a repair. Disconnect the device from the mains before opening!

6 Troubleshooting


6.1 Error Messages and Notices

Message in the display	Meaning
The tolerance may amount to up to 20% of the reference value.	During setting up a calibration measurement a non-valid tolerance value was chosen. The correct tolerance value can be found at the test kit packaging. A typical value is 30 mg/L.
The reference value must be between 0 and 4000 mg/L.	During setting up a calibration measurement a non-valid reference value was chosen. The typical value is close to 300 mg/L and can be found at the test kit packaging.
Setting the language and date/time formats during running test series is not possible.	You can only change the language, date, time and the format of these at a rack where no measurement series or calibration series is active. Setting of these parameters will be logged in the internal memory in the file logfile.txt inside the folder ,log'.
The test series has been stopped.	Occurs after manual stopping of an active measurement series. This line in a measurement or calibration series indicates that the series has been stopped manually. The name of the operator is printed in the last row of this entry.
No more than 6 test series are possible at the same time!	As a rack has got only capacity for six bottles only six measurement or calibration series can be active at the same time. Stop one series first before starting a new one.
On the device there is already a measurement running with the Head ID.	An active measurement or calibration series is already pending on the rack with the same head ID. This could be, if a bottle was replaced and set back to the rack afterwards.
Differing measurement parameters in the head. The measurement has been automatically terminated.	An active measurement or calibration series is already pending on the rack with the same head ID but different setup values. This could be, if a bottle was replaced to another rack and the measurement was stopped manually there. After starting a new measurement with that sensor head and putting it back to the first rack this message will be displayed. To overcome this, at first the active (old) measurement need to be stopped as it is no longer valid. Afterwards the new measurement series can be activated on this rack.
The head will not be valid by the end of the test series.	The validity status of this sensor head expires before the chosen measurement series will be finished.

Message in the display	Meaning
Sensor head not accessible.	If a scheduled measurement is pending and the bottle was replaced this message will be placed in the measurement series data file. After replacing the bottle the measurement series will be continued automatically.
An error occurred during a running measurement!	If no manual bottle replacement during a measurement could be the case control the contact surfaces of the sensor head and the corresponding rack place. If everything is clean and gets contact please get in contact with Lovibond or your local dealer.
No more storage space left on the SD card!	This indicates that the internal memory is full. As about 50.000 measurement series can be stored, this case is very unlikely.
SD card error!	Indicates a problem with the internal memory. Please contact Lovibond or your local dealer.

6.2 Interpretation of BOD curve and Problem Solving Strategies

Example: Measuring Range: 0 – 40 mg/L

Graph	Description	Action
A	Ideal BOD curve	
B	Additional BOD due to nitrification	Add nitrification inhibitor
C	The BOD value is above the measuring range	Increase measuring range or dilute sample
D	Measuring system is leaking	Check BOD bottle for leakage and replace accessories (e.g. rubber gasket) if necessary
E	<ul style="list-style-type: none"> • Not sufficient microorganisms • Not temperature adjusted 	<ul style="list-style-type: none"> • Inoculate sample • Adjust temperature

7 Accessories & Replacement parts

7.1 List of Accessories

6 BOD Sample flasks, brown glass	418645
BD 600 GLP Sensor Certification	999610-GLP
BOD Sample flask, brown glass, 500 ml	418644
BOD sensor GLP	2444470-GLP
Complete Set Overflow measuring flasks	418654
Inductive Stirring Unit	2444456
Magnetic stirring rod	418633
Magnetic stirring rod, 100 pc.	418633-100
Magnetic stirring rod remover	418638
Nitrification inhibitor	2418642
Overflow measuring flask 21.7 ml	418664
Overflow measuring flask 56 ml	418655
Overflow measuring flask 94 ml	418656
Overflow measuring flask 157 ml	418657
Overflow measuring flask 244 ml	418658
Overflow measuring flask 360 ml	418659
Overflow measuring flask 428 ml	418660
Potassium hydroxide solution 45%	2418634
Power supply unit for inductive stirring unit	444454
Remote control	2444481
Rubber gasket, 4.5 cm	418636
Rubber gasket GLP, 6.5 cm	418676
Test set for checking the system (10 tablets)	2418328
USB cable 3 m	2444482
Y-cable	2444475

EN

8.1 Specifications - Mains adapter

Type	SRB1502300P
Design	Mains adapter plug, switching power supply
Input voltage, frequency	100 – 240 V \pm 10 %, 50 / 60Hz
Input current	1000 mA
Primary adapter	Europe, UK, Australia, USA
Protection class	II
Output voltage, frequency	15 V, DC
Max. output current	2300 mA
Output protection	Short-circuit-proof
DC cable length approx.	1800 mm
Environmental conditions	0 °C to 40 °C, with max. 93 % relative air humidity, non-condensing, maximum altitude 2000 m, pollution degree I
Energy efficiency class	VI
Safety standard	EN 60950, EN 62368-1
EMC	IEC 61204-3
Dimensions [mm]	62 x 84 x 53
Weight, with EU adapter	258 g
Approval, EMC	CE, EN 55024

Caution!

Subject to technical modification!
To ensure maximum accuracy of test results, always use the reagent systems supplied by the instrument manufacturer.

8.3 Specifications - BD 600 GLP

Measuring Principle	Manometric; mercury-free; electronic pressure sensor
Measuring Range	0-40, 0-80, 0-200, 0-400, 0-800, 0-2000, 0-4000 mg/l
Display	Backlit
Interfaces	USB
Auto – OFF	No
External Storage	USB
Internal Storage	1 GB or up to 50000 measurement or calibration series
Measurement Time	Selectable: 5, 28, 60 and 90 days possible
Power Supply	<ul style="list-style-type: none"> • 3 alkali-manganese batteries (Baby cells/ size C) • 100 - 240 V / 50-60 Hz
Clock	Real Time Clock and Date
Portability	Benchtop
Autostart	No
Measuring Stations	6
Storage Interval	- every 2 hours (5 days) - every 12 hours (28 days) - every 24 hours (60 days) - every 24 - 48 hours (90 days)
Protection Class	IP 53
Compliance	<ul style="list-style-type: none"> • CE • EMC according to DIN EN 61326
Dimensions	181 x 230 x 375 mm
Weight	(4100 g, unit with bottles & batteries 5775 g, complete with stirring unit)
Weight with Packaging	(4100 g, unit with bottles & batteries 5775 g, complete with stirring unit)

Caution!

Subject to technical modification!
To ensure maximum accuracy of test results, always use the reagent systems supplied by the instrument manufacturer.

9 Appendix

9.1 List of Trademarks

Lovibond® and Tintometer® are registered trademarks of the Tintometer group of companies. All translations and transliterations of Lovibond® and Tintometer® are asserted as trademarks of The Tintometer® Group.

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- Consistent measurement conditions - The automatic start function ensures all measurements are taken at an equal sample temperature. When the measurement conditions are consistent across readings, operators do not have to work about temperature or pressure fluctuations skewing results.
- Sample identification and tracking - The BD 600 interface, makes it easy to assign an identification code or sample name when setting up a new sample. This helps track a sample from the time it is collected and ensures traceability during sample measurement.
- Use the remote control - If using an incubator with a glass door, operators can use the supplied remote control to display measurement data for a sample - without even opening the door! This further ensures that samples stay at a consistent temperature which is essential to accurate results.

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The inductive stirring system is designed for the stirring of liquids in special BOD bottles. It comprises a super-flat stirrer drive with 6 stirring positions and the mains adapter. The inductive stirrer drive does not have a motor and is thus wear-free. It is especially well-suited for stirring during continuous operation in thermocabinets. It is completely protected from dripping water due to its hermetically sealed encapsulation. It can be used in harsh environmental conditions.

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9.1 Hinweis zu Copyright und Warenzeichen
Trademarks of the Tintometer Group
Tintometer Inc.