

# BD™ Cytometric Bead Array (CBA) Human Th1/Th2/Th17 Cytokine Kit Instruction Manual

## Trademarks

FCAP Array is a registered trademark of Softflow, Inc.

Macintosh and Mac are trademarks of Apple Computer, Inc., registered in the US and other countries.

Falcon® is a registered trademark of Corning Incorporated.

*BD flow cytometers are class I (1) laser products*

*© 2019 Becton, Dickinson and Company. All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, stored in retrieval systems, or translated into any language or computer language, in any form or by any means: electronic, mechanical, magnetic, optical, chemical, manual, or otherwise, without prior written permission from BD Biosciences. For research use only. Not for use in diagnostic or therapeutic procedures. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited.*

*BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2019 BD*

# Contents

---

<b>Chapter 1: About this kit</b>	<b>5</b>
Purpose of the kit	6
Limitations	8
Kit contents	9
Storage and handling	11
<b>Chapter 2: Before you begin</b>	<b>13</b>
Workflow overview	14
Required materials	15
<b>Chapter 3: Assay preparation</b>	<b>17</b>
Preparing Human Th1/Th2/Th17 Cytokine Standards	18
Mixing Human Th1/Th2/Th17 Cytokine Capture Beads	20
Diluting samples	21
<b>Chapter 4: Assay procedure</b>	<b>23</b>
Performing the Human Th1/Th2/Th17 Cytokine Assay	24
Data Analysis	27
<b>Chapter 5: Performance</b>	<b>31</b>
Theoretical limit of detection	32
Recovery	33
Linearity	35
Specificity	37
Precision	39
<b>Chapter 6: Reference</b>	<b>41</b>
Troubleshooting	42
References	43
Notes	46



## About this kit

---

This section covers the following topics:

- [Purpose of the kit \(page 6\)](#)
- [Limitations \(page 8\)](#)
- [Kit contents \(page 9\)](#)
- [Storage and handling \(page 11\)](#)

## Purpose of the kit

---

**Use of the kit** The BD™ CBA Human Th1/Th2/Th17 Cytokine Kit can be used to measure Interleukin-2 (IL-2), Interleukin-4 (IL-4), Interleukin-6 (IL-6), Interleukin-10 (IL-10), Tumor Necrosis Factor (TNF), Interferon- $\gamma$  (IFN- $\gamma$ ), and Interleukin-17A (IL-17A) protein levels in a single sample. The kit performance has been optimized for analysis of physiologically relevant concentrations (pg/mL levels) of specific cytokine proteins in tissue culture supernatants, EDTA plasma, and serum samples. The kit provides sufficient reagents for 80 tests.

---

**Principle of CBA assays** BD CBA assays provide a method of capturing a soluble analyte or set of analytes with beads of known size and fluorescence, making it possible to detect analytes using flow cytometry.

Each capture bead in a BD CBA kit has been conjugated with a specific antibody. The detection reagent provided in the kit is a mixture of phycoerythrin (PE)–conjugated antibodies, which provides a fluorescent signal in proportion to the amount of bound analyte.

When the capture beads and detector reagent are incubated with an unknown sample containing recognized analytes, sandwich complexes (capture bead + analyte + detection reagent) are formed. These complexes can be measured using flow cytometry to identify particles with fluorescence characteristics of both the bead and the detector.

## Principle of this assay

The BD CBA Human Th1/Th2/Th17 Cytokine Kit uses bead array technology to simultaneously detect multiple cytokine proteins in research samples. Seven bead populations with distinct fluorescence intensities have been coated with capture antibodies specific for IL-2, IL-4, IL-6, IL-10, TNF, IFN- $\gamma$ , and IL-17A proteins. The seven bead populations are mixed together to form the bead array, which is resolved in a red channel (ie, FL3 or FL4) of a flow cytometer.




Figure 1

During the assay procedure, you will mix the cytokine capture beads with the recombinant standards or unknown samples and incubate them with the PE-conjugated detection antibodies to form sandwich complexes. The intensity of PE fluorescence of each sandwich complex reveals the concentration of that cytokine. After acquiring samples on a flow cytometer, use FCAP Array™ software to generate results in graphical and tabular format.

---

**Advantages over ELISA**

The broad dynamic range of fluorescent detection via flow cytometry and the efficient capturing of analytes via suspended particles enable the BD CBA assay to measure the concentration of an unknown in substantially less time and using fewer sample dilutions compared to conventional ELISA methodology.

- The required sample volume is approximately one-seventh the quantity necessary for conventional ELISA assays due to the detection of seven analytes in a single sample.
  - A single set of diluted standards is used to generate a standard curve for each analyte.
  - A BD CBA experiment takes less time than a single ELISA and provides results that would normally require seven conventional ELISAs.
- 

## Limitations

---

**Assay limitations** The theoretical limit of detection of the BD CBA Human Th1/Th2/Th17 Cytokine Kit is comparable to conventional ELISA, but due to the complexity and kinetics of this multi-analyte assay, the actual limit of detection on a given experiment may vary. See [Theoretical limit of detection \(page 32\)](#) and [Precision \(page 39\)](#).

The BD CBA Kit is not recommended for use on stream-in-air instruments for which signal intensities may be reduced, adversely affecting assay sensitivity. Stream-in-air instruments include the BD FACStar™ Plus, BD FACSVantage™, and BD Influx™ flow cytometers (BD Biosciences).



Quantitative results or protein levels for the same sample or recombinant protein run in ELISA and BD CBA assays may differ. A spike recovery assay can be performed using an ELISA standard followed by BD CBA analysis to assess possible differences in quantitation.

This kit is designed to be used as an integral unit. Do not mix components from different batches or kits.

## Kit contents

### Contents

The BD CBA Human Th1/Th2/Th17 Cytokine Kit contains the following components sufficient for 80 tests.

Vial label	Reagent	Quantity
A1	Human IL-2 Capture Beads	1 vial, 0.8 mL
A2	Human IL-4 Capture Beads	1 vial, 0.8 mL
A3	Human IL-6 Capture Beads	1 vial, 0.8 mL
A4	Human IL-10 Capture Beads	1 vial, 0.8 mL
A5	Human TNF Capture Beads	1 vial, 0.8 mL
A6	Human IFN- $\gamma$ Capture Beads	1 vial, 0.8 mL
A7	Human IL-17A Capture Beads	1 vial, 0.8 mL
B	Human Th1/Th2/Th17 PE Detection Reagent	1 vial, 4 mL
C	Human Th1/Th2/Th17 Cytokine Standards	2 vials lyophilized
D	Cytometer Setup Beads	1 vial, 1.5 mL
E1	PE Positive Control Detector	1 vial, 0.5 mL
E2	FITC Positive Control Detector	1 vial, 0.5 mL

Vial label	Reagent	Quantity
F	Wash Buffer	1 bottle, 130 mL
G	Assay Diluent	1 bottle, 30 mL
H	Serum Enhancement Buffer	1 bottle, 10 mL

**Bead reagents**

**Human Cytokine Capture Beads (A1–A7):** An 80-test vial of each specific capture bead (A1–A7). The specific capture beads, having discrete fluorescence intensity characteristics, are distributed from brightest (A1) to dimmest (A7).

**Cytometer Setup Beads (D):** A 30-test vial of setup beads for setting the initial instrument PMT voltages and compensation settings is sufficient for 10 instrument setup procedures. The Cytometer Setup Beads are formulated for use at 50 µL/test.

**Antibody and standard reagents**

**Human Th1/Th2/Th17 PE Detection Reagent (B):** An 80-test vial of PE-conjugated anti-human IL-2, IL-4, IL-6, IL-10, TNF, IFN $\gamma$ , and IL-17A antibodies that is formulated for use at 50 µL/test.

**Human Th1/Th2/Th17 Cytokine Standards (C):** Two vials containing lyophilized recombinant human cytokine proteins. Each vial should be reconstituted in 2.0 mL of Assay Diluent to prepare the top standard.

**PE Positive Control Detector (E1):** A 10-test vial of PE-conjugated antibody control that is formulated for use at 50 µL/test. This reagent is used with the Cytometer Setup Beads to set the initial instrument compensation settings.

**FITC Positive Control Detector (E2):** A 10-test vial of FITC-conjugated antibody control that is formulated for use at 50  $\mu$ L/test. This reagent is used with the Cytometer Setup Beads to set the initial instrument compensation settings.

---

## Buffer reagents

**Wash Buffer (F):** A 130-mL bottle of phosphate buffered saline (PBS) solution (1X), containing protein and detergent, used for wash steps and to resuspend the washed beads for analysis.

**Assay Diluent (G):** A 30-mL bottle of a buffered protein solution (1X) used to reconstitute and dilute the Human Th1/Th2/Th17 Cytokine Standards and to dilute unknown samples.

**Serum Enhancement Buffer (H):** A 10-mL bottle of a buffered protein solution (1X) used to dilute mixed Capture Beads when testing serum or plasma samples.

**Note:** Source of all serum proteins is from USDA inspected abattoirs located in the United States.

---

## Storage and handling

---

### Storage

Store all kit components at 2 to 8°C. Do not freeze.

**Appearance:** The visual appearance of Assay Diluent may range in color from red to yellow/orange. The visual appearance of ELISA Dilution reagent may range in color from clear to cloudy white.

---

**Warning** Components A1–A7, B, D, E1, E2, F, G, and H contain sodium azide. Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.

Human Th1/Th2/Th17 Cytokine Standards (component 51-9006249) contains 0.02% (w/w) of a CMIT/MIT mixture (3:1), which is a mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC No 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC No 220-239-6] (3:1).

**Hazard statements** May cause an allergic skin reaction.

**Precautionary statements** Wear protective gloves / eye protection.

Wear protective clothing.

Avoid breathing mist/vapours/spray.

If skin irritation or rash occurs: Get medical advice/attention.

IF ON SKIN: Wash with plenty of water.

Dispose of contents/container in accordance with local/regional/national/international regulation

---

# 2

## Before you begin

---

This section covers the following topics:

- [Workflow overview \(page 14\)](#)
- [Required materials \(page 15\)](#)

## Workflow overview

**Workflow**      The overall workflow consists of the following steps:

Step	Description
1	<a href="#">Preparing Human Th1/Th2/Th17 Cytokine Standards (page 18)</a>
2	<a href="#">Mixing Human Th1/Th2/Th17 Cytokine Capture Beads (page 20)</a>
3	<a href="#">Diluting samples (page 21)</a> , if necessary
4	Performing instrument setup with Cytometer Setup Beads (instructions can be found at <a href="http://bdbiosciences.com/cbasetup">bdbiosciences.com/cbasetup</a> )  <b>Note:</b> Can be performed during the incubation in <a href="#">step 5</a> .
5	<a href="#">Performing the Human Th1/Th2/Th17 Cytokine Assay (page 24)</a>
6	Acquiring samples (instructions can be found at <a href="http://bdbiosciences.com/cbasetup">bdbiosciences.com/cbasetup</a> )
7	<a href="#">Data Analysis (page 27)</a>

**Incubation times**      To help you plan your work, the incubation times are listed in the following table:

Procedure	Incubation time
Preparing standards	15 minutes
Preparing mixed capture beads (when analyzing serum or plasma samples only)	30 minutes
Preparing Cytometer Setup Beads	30 minutes
Performing the assay	3 hours

## Required materials

### Materials required but not provided

In addition to the reagents provided in the BD CBA Human Th1/Th2/Th17 Cytokine Kit, the following items are also required:

- A dual-laser flow cytometer equipped with a 488-nm or 532-nm and a 633-nm or 635-nm laser capable of distinguishing 576-nm, 660-nm, and >680-nm fluorescence. The following table lists examples of compatible instrument platforms.

Flow cytometer	Reporter channel	Bead channels
BD FACSAArray™	Yellow	Red
BD FACSCanto™ platform BD™ LSR platform BD FACSAria™ platform	PE	APC
BD FACSCalibur™ (single laser) BD FACSCalibur (dual laser)	FL2	FL3 FL4
<b>Note:</b> Visit <a href="https://bdbiosciences.com/cbasetup">bdbiosciences.com/cbasetup</a> for setup protocols.		

- Falcon® 12 × 75-mm sample acquisition tubes (Catalog No. 352008), or equivalent
- 15-mL conical, polypropylene tubes (Falcon, Catalog No. 352097), or equivalent
- FCAP Array software (Catalog No. 641488 [PC] or 645447 [Mac])

**Materials  
required for  
plate loader-  
equipped flow  
cytometers**

- 
- Millipore MultiScreen<sub>HTS</sub>-BV 1.2 µm Clear non-sterile filter plates [Catalog No. MSBVN1210 (10 pack) or MSBVN1250 (50 pack)]
  - Millipore MultiScreen<sub>HTS</sub> Vacuum Manifold, (Catalog No. MSVMHTS00)
  - MTS 2/4 Digital Stirrer, IKA Works, VWR (Catalog No. 82006-096)
  - Vacuum source
  - Vacuum gauge and regulator (if not using the recommended manifold)
-



## Assay preparation

---

This section covers the following topics:

- [Preparing Human Th1/Th2/Th17 Cytokine Standards \(page 18\)](#)
- [Mixing Human Th1/Th2/Th17 Cytokine Capture Beads \(page 20\)](#)
- [Diluting samples \(page 21\)](#)

# Preparing Human Th1/Th2/Th17 Cytokine Standards


---

**Purpose of this procedure**      The Human Th1/Th2/Th17 Cytokine Standards are lyophilized and should be reconstituted and serially diluted immediately before mixing with the Capture Beads and the PE Detection Reagent.

You must prepare fresh cytokine standards to run with each experiment. Do not store or reuse reconstituted or diluted standards.

- 
- Procedure**      **To reconstitute and serially dilute the standards:**
1. Open one vial of lyophilized Human Th1/Th2/Th17 Standards. Transfer the standard spheres to a 15-mL conical, polypropylene tube. Label the tube “Top Standard.”
  2. Reconstitute the standards with 2 mL of Assay Diluent.
    - a. Allow the reconstituted standard to equilibrate for at least 15 minutes at room temperature.
    - b. Gently mix the reconstituted protein by pipette only. Do not vortex or mix vigorously.
  3. Label 12 × 75-mm tubes and arrange them in the following order: 1:2, 1:4, 1:8, 1:16, 1:32, 1:64, 1:128, and 1:256.
  4. Pipette 300 µL of Assay Diluent in each of the 12 × 75-mm tubes.
  5. Perform serial dilutions:
    - a. Transfer 300 µL from the Top Standard to the 1:2 dilution tube and mix thoroughly by pipette only.

- b. Continue making serial dilutions by transferring 300  $\mu$ L from the 1:2 tube to the 1:4 tube and so on to the 1:256 tube.
- c. Mix thoroughly by pipette only. Do not vortex.



6. Prepare one 12  $\times$  75-mm tube containing only Assay Diluent to serve as the 0 pg/mL negative control.

### Concentration of standards

See the [Procedure for tubes](#) section of [Performing the Human Th1/Th2/Th17 Cytokine Assay \(page 24\)](#) for a listing of the concentrations (pg/mL) of all seven recombinant proteins in each standard dilution.

### Next step

Proceed to [Mixing Human Th1/Th2/Th17 Cytokine Capture Beads \(page 20\)](#).

# Mixing Human Th1/Th2/Th17 Cytokine Capture Beads

---

**Purpose of this procedure**      The Capture Beads are bottled individually (A1–A7). You must pool all seven bead reagents immediately before using them in the assay.

---

- Mixing the beads**      **To mix the Capture Beads:**
1. Determine the number of assay tubes (including standards and controls) required for the experiment (eg, 8 unknowns, 9 cytokine standard dilutions, and 1 negative control = 18 assay tubes).
  2. Vigorously vortex each Capture Bead suspension for 3 to 5 seconds before mixing.  
  
**Note:** The antibody-conjugated beads will settle out of suspension over time. It is necessary to vortex the vial before taking a bead-suspension aliquot.
  3. Add a 10-μL aliquot of each Capture Bead, for each assay tube to be analyzed, into a single tube labeled “mixed Capture Beads” (eg, 10 μL of IL-2 Capture Beads × 18 assay tubes = 180 μL of IL-2 Capture Beads required).
  4. Vortex the bead mixture thoroughly.
- 

**Resuspending the beads**      If you are using serum or plasma samples, you must perform this procedure. This procedure is optional for all other sample types.

- To resuspend the Capture Beads in Serum Enhancement Buffer:**
1. Centrifuge the mixed Capture Beads at 200g for 5 minutes.
  2. Carefully aspirate and discard the supernatant.

3. Resuspend the mixed Capture Beads pellet in Serum Enhancement Buffer (equal to the volume removed in [step 2](#)) and vortex thoroughly.
4. Incubate the mixed Capture Beads for 30 minutes at room temperature, protected from light.

---

**Next step**

The mixed Capture Beads are now ready to be transferred to the assay tubes. Discard excess mixed Capture Beads. Do not store after mixing.

To begin the assay, proceed to [Performing the Human Th1/Th2/Th17 Cytokine Assay \(page 24\)](#). If you need to dilute samples having a high cytokine concentration, proceed to [Diluting samples \(page 21\)](#).

---

## Diluting samples

---

**Purpose of this procedure**

The standard curve for each cytokine covers a defined set of concentrations from 20 to 5000 pg/mL. It might be necessary to dilute samples to ensure that their mean fluorescence values fall within the range of the generated cytokine standard curve. For best results, dilute samples that are known or assumed to contain high levels of a given cytokine. This procedure is not necessary for all samples.

<b>Procedure</b>	<hr/> <p><b>To dilute samples with a known high cytokine concentration:</b></p> <ol style="list-style-type: none"><li>1. Dilute the sample by the desired dilution factor (ie, 1:2, 1:10, or 1:100) using the appropriate volume of Assay Diluent.</li><li>2. Mix sample dilutions thoroughly.</li></ol> <p><b>Note:</b> Optimal recovery from serum samples typically requires a 1:4 dilution.</p>
<b>Next step</b>	<hr/> <p>Perform instrument setup using the Cytometer Setup Beads. For details, go to <a href="https://bdbiosciences.com/cbasetup">bdbiosciences.com/cbasetup</a> and select the appropriate flow cytometer under CBA Kits: Instrument Setup.</p> <p>Or, if you wish to begin staining your samples for the assay, proceed to <a href="#">Performing the Human Th1/Th2/Th17 Cytokine Assay (page 24)</a>, and you can perform instrument setup during the 3-hour staining incubation.</p> <hr/>

## Assay procedure

---

This section covers the following topics:

- [Performing the Human Th1/Th2/Th17 Cytokine Assay \(page 24\)](#)
- [Data Analysis \(page 27\)](#)

# Performing the Human Th1/Th2/Th17 Cytokine Assay

- Before you begin**
- Prepare the standards as described in [Preparing Human Th1/Th2/Th17 Cytokine Standards \(page 18\)](#).
  - Mix the Capture Beads as described in [Mixing Human Th1/Th2/Th17 Cytokine Capture Beads \(page 20\)](#). Be sure to follow the appropriate procedure (cell culture supernatant vs serum/plasma) for your sample type.
  - If necessary, dilute the unknown samples. See [Diluting samples \(page 21\)](#).

**Procedure for tubes**

**To perform the assay:**

1. Vortex the mixed Capture Beads and add 50 µL to all assay tubes.
2. Add 50 µL of the Human Th1/Th2/Th17 Cytokine Standard dilutions to the control tubes as listed in the following table:

Tube label	Concentration (pg/mL)	Cytokine Standard dilution
1	0 (negative control)	no standard dilution (Assay Diluent only)
2	20	1:256
3	40	1:128
4	80	1:64
5	156	1:32
6	312.5	1:16
7	625	1:8
8	1,250	1:4
9	2,500	1:2
10	5,000	Top standard



3. Add 50  $\mu$ L of each unknown sample to the appropriately labeled sample tubes.
4. Add 50  $\mu$ L of the Human Th1/Th2/Th17 PE Detection Reagent to all assay tubes.
5. Incubate the assay tubes for 3 hours at room temperature, protected from light.  
**Note:** If you have not yet performed cytometer setup, you may wish to do so during this incubation.
6. Add 1 mL of Wash Buffer to each assay tube and centrifuge at 200g for 5 minutes.
7. Carefully aspirate and discard the supernatant from each assay tube.
8. Add 300  $\mu$ L of Wash Buffer to each assay tube to resuspend the bead pellet.

---

**Procedure for  
filter plates****To perform the assay:**

1. Wet the plate by adding 100  $\mu$ L of wash buffer to each well.
2. Place the plate on the vacuum manifold.
3. Aspirate for 2 to 10 seconds until the wells are drained.
4. Remove the plate from the manifold, then blot the bottom of the plate on paper towels.
5. Add 50  $\mu$ L of each of the following to the wells in the filter plate:
  - Capture Beads (vortex before adding)
  - Standard or sample (add standards from the lowest concentration to the highest, followed by samples)
  - Human Th1/Th2/Th17 PE Detection Reagent

6. Cover the plate and shake it for 5 minutes at 1,100 rpm on a plate shaker.
  7. Incubate the plate for 3 hours at room temperature on a non-absorbent, dry surface.  
**Note:** Place the plate on a non-absorbent, dry surface during incubation. Absorbent or wet surfaces can cause the contents of the wells to leak.
  8. Remove the cover from the plate and apply the plate to the vacuum manifold.
  9. Vacuum aspirate for 2 to 10 seconds until the wells are drained.
  10. Remove the plate from the manifold, then blot the bottom of the plate on paper towels after aspiration.
  11. Add 120  $\mu$ L of wash buffer to each well to resuspend the beads.
  12. Cover the plate and shake it for 2 minutes at 1,100 rpm before you begin sample acquisition.
- 

### Next step

Acquire the samples on the flow cytometer. For details, go to [bdbiosciences.com/cbasetup](https://bdbiosciences.com/cbasetup) and select the appropriate flow cytometer under CBA Kits: Instrument Setup.

CBA samples must be acquired on the same day they are prepared. Prolonged storage of samples, once the assay is complete, can lead to increased background and reduced sensitivity.

To facilitate the analysis of samples using FCAP Array software, we recommend the following guidelines:

- Acquire standards from lowest (0 pg/mL) to highest (Top Standard) concentration, followed by the test samples.

- If running sample dilutions, acquire sequentially starting with the most concentrated sample.
- Store all FCS files (standards and samples) in a single folder.

When you are finished acquiring samples, proceed to [Data Analysis \(page 27\)](#).

---


## Data Analysis

---


**How to analyze** Analyze Human Th1/Th2/Th17 Cytokine data using FCAP Array software. For instructions on analysis, go to [bdbiosciences.com/cbasetup](https://bdbiosciences.com/cbasetup) and see the *Guide to Analyzing Data from BD CBA Kits Using FCAP Array Software*.

Typical data


The following data, acquired using BD CellQuest software, shows standards and detectors alone.




Negative control (0 pg/mL)



Standard 80 pg/mL




Standard 625 pg/mL




Standard 5000 pg/mL

**Standard curve examples**


The following graphs represent standard curves from the Human Th1/Th2/Th17 Cytokine Standards.




IL-2




IL-4




IL-6




IL-10



TNF



IFN-γ



IL-17A



## Performance

---

This section covers the following topics:

- [Theoretical limit of detection \(page 32\)](#)
- [Recovery \(page 33\)](#)
- [Linearity \(page 35\)](#)
- [Specificity \(page 37\)](#)
- [Precision \(page 39\)](#)

## Theoretical limit of detection

**Experiment details**

The individual standard curve range for a given cytokine defines the minimum and maximum quantifiable levels using the BD CBA Human Th1/Th2/Th17 Cytokine Kit (ie, 20 pg/mL and 5000 pg/mL). By applying the 4-parameter curve fit option, it is possible to extrapolate values for sample intensities not falling within the limits of the standard curve. It is up to the researcher to decide the best method for calculating values for unknown samples using this assay. The theoretical limit of detection for each cytokine using the BD CBA Human Th1/Th2/Th17 Cytokine Kit is defined as the corresponding concentration at two standard deviations above the median fluorescence of 30 replicates of the negative control (0 pg/mL).

**Limit of detection data**

Cytokine	Limit of detection (pg/mL)
IL-2	2.6
IL-4	4.9
IL-6	2.4
IL-10	4.5
TNF	3.8
IFN- $\gamma$	3.7
IL-17A	18.9



## Recovery

---

### Experiment details

Individual cytokine protein was spiked into various matrices at three different levels within the assay range. The spiked samples were assayed and the results were compared with the expected values. The cell culture medium used in these experiments was not diluted before addition of the cytokine protein. Pooled human serum and pooled human plasma samples were diluted 1:4 in Assay Diluent before addition of cytokine protein. The plasma samples in these experiments were EDTA treated.

Recovery data

Cytokine	Matrix	Average % Recovery	Range (%)
IL-2	Media	83	72–95
	Serum	86	83–91
	Plasma	84	78–92
IL-4	Media	81	75–87
	Serum	91	87–94
	Plasma	87	85–88
IL-6	Media	86	79–92
	Serum	90	88–92
	Plasma	93	91–98
IL-10	Media	86	80–92
	Serum	95	93–96
	Plasma	91	89–94
TNF	Media	88	81–95
	Serum	95	93–97
	Plasma	95	92–98
IFN- $\gamma$	Media	80	72–89
	Serum	84	82–87
	Plasma	76	76–77
IL-17A	Media	84	74–93
	Serum	74	60–92
	Plasma	73	55–93

## Linearity

**Experiment details**

In two experiments, the following matrices were spiked with IL-2, IL-4, IL-6, IL-10, TNF, IFN- $\gamma$ , and IL-17A and then were serially diluted with Assay Diluent.

**Linearity data**

Cytokine	Matrix	Sample dilution	Detected (pg/mL)	Average % of expected
IL-2	Media	1:2	1020.8	100
		1:4	469.4	92
		1:8	208.2	82
	Serum	1:2	1161.5	100
		1:4	514.0	89
		1:8	241.0	83
	Plasma	1:2	1001.4	100
		1:4	480.2	96
		1:8	232.9	93
IL-4	Media	1:2	958.4	100
		1:4	464.2	97
		1:8	222.0	93
	Serum	1:2	1119.9	100
		1:4	523.8	94
		1:8	263.7	94
	Plasma	1:2	937.5	100
		1:4	494.6	106
		1:8	250.3	107
IL-6	Media	1:2	1101.5	100
		1:4	498.1	90
		1:8	236.7	86
	Serum	1:2	1203.6	100
		1:4	567.8	94
		1:8	275.7	92
	Plasma	1:2	1063.0	100
		1:4	577.2	109
		1:8	272.5	103

Cytokine	Matrix	Sample dilution	Detected (pg/mL)	Average % of expected
IL-10	Media	1:2	1049.2	100
		1:4	506.0	96
		1:8	240.9	92
	Serum	1:2	1167.6	100
		1:4	543.3	93
		1:8	270.1	93
	Plasma	1:2	1013.6	100
		1:4	506.8	100
		1:8	250.1	99
TNF	Media	1:2	1077.0	100
		1:4	518.4	96
		1:8	238.5	89
	Serum	1:2	1277.8	100
		1:4	574.2	90
		1:8	273.2	86
	Plasma	1:2	1084.1	100
		1:4	573.3	106
		1:8	273.5	101
IFN- $\gamma$	Media	1:2	946.1	100
		1:4	441.6	93
		1:8	219.7	93
	Serum	1:2	1039.1	100
		1:4	466.9	90
		1:8	219.0	84
	Plasma	1:2	857.7	100
		1:4	451.3	105
		1:8	214.6	100


Cytokine	Matrix	Sample dilution	Detected (pg/mL)	Average % of expected
IL-17A	Media	1:2	1106.7	100
		1:4	521.8	94
		1:8	244.8	88
	Serum	1:2	937.3	100
		1:4	496.4	106
		1:8	256.8	110
	Plasma	1:2	752.9	100
		1:4	466.0	124
		1:8	226.4	120

## Specificity

### Experiment details

The antibodies used in the BD CBA Human Th1/Th2/Th17 Cytokine Kit have been screened for specific reactivity with their specific cytokines. Analysis of samples containing only a single recombinant cytokine protein found no cross-reactivity or background detection of cytokine in other Capture Bead populations using this assay.

**Specificity data** Sample data containing only a single recombinant cytokine protein was acquired using BD CellQuest software.



## Precision

### Intra-assay precision

Ten replicates of each of three different levels of IL-2, IL-4, IL-6, IL-10, TNF, IFN- $\gamma$ , and IL-17A were tested.

Cytokine	Sample	Mean (pg/mL)	Standard deviation	%CV
IL-2	Sample 1	67.8	2.9	4
	Sample 2	279.1	14.8	5
	Sample 3	1238.5	46.0	4
IL-4	Sample 1	78.1	3.8	5
	Sample 2	294.5	10.2	3
	Sample 3	1231.2	25.8	2
IL-6	Sample 1	75.8	4.5	6
	Sample 2	303.1	14.1	5
	Sample 3	1284.1	55.5	4
IL-10	Sample 1	76.2	4.1	5
	Sample 2	296.3	11.3	4
	Sample 3	1272.2	46.8	4
TNF	Sample 1	75.1	6.0	8
	Sample 2	302.6	20.4	7
	Sample 3	1316.9	74.5	6
IFN- $\gamma$	Sample 1	69.6	2.4	3
	Sample 2	280.9	11.4	4
	Sample 3	1233.7	53.1	4
IL-17A	Sample 1	76.9	3.6	5
	Sample 2	305.7	7.8	3
	Sample 3	1308.9	52.8	4

Inter-assay  
precision

Three different levels of IL-2, IL-4, IL-6, IL-10, TNF, IFN- $\gamma$ , and IL-17A were tested in four experiments conducted by four different operators.

Cytokine	Sample	Mean (pg/mL)	Standard deviation	%CV
IL-2	Sample 1	71.0	6.3	9
	Sample 2	290.0	20.3	7
	Sample 3	1230.3	82.0	7
IL-4	Sample 1	76.9	8.1	11
	Sample 2	297.7	18.6	6
	Sample 3	1217.6	57.6	5
IL-6	Sample 1	77.7	10.4	13
	Sample 2	297.8	25.3	8
	Sample 3	1254.4	89.0	7
IL-10	Sample 1	77.8	8.3	11
	Sample 2	296.0	18.8	6
	Sample 3	1245.8	93.7	8
TNF	Sample 1	77.8	9.4	12
	Sample 2	300.0	26.9	9
	Sample 3	1263.7	96.8	8
IFN- $\gamma$	Sample 1	74.4	8.3	11
	Sample 2	291.2	26.6	9
	Sample 3	1239.8	95.1	8
IL-17A	Sample 1	75.1	9.4	12
	Sample 2	303.7	22.7	7
	Sample 3	1272.8	76.8	6



## Reference

---

This section covers the following topics:

- [Troubleshooting \(page 42\)](#)
- [References \(page 43\)](#)

## Troubleshooting

**Recommended actions**      These are the actions we recommend you take if you encounter the following problems.

Problem	Recommended actions
Variation between duplicate samples	Vortex Capture Beads before pipetting. Beads can aggregate.
Low bead number in samples	Avoid aspiration of beads during wash step. Do not wash or resuspend beads in volumes higher than recommended volumes.
High background	Test various sample dilutions. The sample might be too concentrated. Remove excess Human Th1/Th2/TH17 II PE Detection Reagent by increasing the number of wash steps, since background may be due to non-specific binding.
Little or no protein detected in sample	Sample may be too dilute. Try various sample dilutions.
Less than seven bead populations observed during analysis, or distribution is unequal	Ensure that equal volumes of beads were added to each assay tube. Vortex Capture Bead vials before taking aliquots. Once Capture Beads are mixed, vortex to ensure beads are distributed throughout solution.
Debris (FSC/SSC) during sample acquisition	Increase FSC threshold or further dilute samples. Increase number of wash steps, if necessary.
Overlap of bead fluorescence (FL3) during acquisition	Samples might have very high cytokine concentration. Ensure instrument settings have been optimized using Cytometer Setup Beads.

Problem	Recommended actions
Standards show low fluorescence or poor standard curve	<p>Ensure that all components are properly prepared and stored. Use new vial of standard with each experiment and once reconstituted, do not use after 12 hours.</p> <p>Ensure that incubation times were of proper length.</p>
All samples are positive or above high standard mean fluorescence value	Dilute samples further. Samples might be too concentrated.
Biohazardous samples	You may treat samples briefly with 1% paraformaldehyde before acquiring samples on the flow cytometer. However, this may affect assay performance and should be validated by the user.

## References

### Related publications

1. Bishop JE, Davis KA. A flow cytometric immunoassay for  $\beta$ 2-microglobulin in whole blood. *J Immunol Methods*. 1997;210:79–87.
2. Camilla C, Defoort JP, Delaage M, Auer R, Quintana J, Lary T, Hamelik R, Prato S, Casano B, Martin M, Fert V. A new flow cytometry-based multi-assay system. 1. Application to cytokine immunoassays. *Cytometry Suppl*. 1998;8:132.
3. Carson R, Vignali D. Simultaneous quantitation of fifteen cytokines using a multiplexed flow cytometric assay. *J Immunol Methods*. 1999;227:41–52.
4. Chen R, Lowe L, Wilson JD, et al. Simultaneous quantification of six human cytokines in a single sample using microparticle-based flow cytometric technology. *Clin Chem*. 1999;45:1693–1694.

5. Collins DP, Luebering BJ, Shaut DM. T-lymphocyte functionality assessed by analysis of cytokine receptor expression, intracellular cytokine expression, and femtomolar detection of cytokine secretion by quantitative flow cytometry. *Cytometry*. 1998;33:249–255.
6. Fulton RJ, McDade RL, Smith PL, Kienker LJ, Kettman JR Jr. Advanced multiplexed analysis with the FlowMetrix system. *Clin Chem*. 1997;43:1749–1756.
7. Kricka LJ. Simultaneous multianalyte immunoassays. In Immunoassay. Diamandis EP, Christopoulos TK, eds. Academic Press. 1996:389–404.
8. Lund-Johansen F, Davis K, Bishop J, de Waal Malefyt R. Flow cytometric analysis of immunoprecipitates: High-throughput analysis of protein phosphorylation and protein-protein interactions. *Cytometry*. 2000;39:250–259.
9. McHugh TM. Flow microsphere immunoassay for the quantitative and simultaneous detection of multiple soluble analytes. *Methods Cell Biol*. 1994;42:575–595.
10. Oliver KG, Kettman JR, Fulton RJ. Multiplexed analysis of human cytokines by use of the FlowMetrix system. *Clin Chem*. 1998;44:2057–2060.
11. Stall A, Sun Q, Varro R, Lowe L, Crowther E, Abrams B, Bishop J, Davis K. A single tube flow cytometric multibead assay for isotyping mouse monoclonal antibodies. Abstract LB77. Experimental Biology Meeting 1998 (late-breaking abstracts).

12. Cook EB, Stahl JL, Lowe L, et al. Simultaneous measurement of six cytokines in a single sample of human tears using microparticle-based flow cytometry: allergics vs. non-allergics. *J Immunol Methods* 2001;254:109–118.
  13. Dotti G, Salvodo B, Takahashi S, et al. Adenovector-induced expression of human-CD40-ligand (hCD40L) by multiple myeloma cells: A model for immunotherapy. *Exp Hematol.* 2001;29:952–961.
-

## Notes

[illegible]



**United States**

877.232.8995

**Canada**

866.979.9408

**Europe**

32.2.400.98.95

**Japan**

0120.8555.90

**Asia/Pacific**

65.66642770

**Latin America/Caribbean**

55.11.5185.9625

Toll Free 0800.771.71.57

**Becton, Dickinson and Company**

**BD Life Sciences - Biosciences**

2350 Qume Drive, San Jose, CA 95131

Ordering (US) 855.236.2772

Technical Service 877.232.8995

Fax 800.325.9637

**answers@bd.com**

**[bdbiosciences.com](http://bdbiosciences.com)**

Catalog No. 560484

23-12381-01 Rev. 03

