

S1 Nuclease

Catalog Number EN0321

Pub. No. MAN0013722 Rev. B00



WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Sheets (SDSs) are available from [thermofisher.com/support](https://www.thermofisher.com/support).

Contents and storage

Contents	Amount	Storage
S1 Nuclease	10000 U, 100 U/ μ L	-25 °C to -15 °C
5X Reaction Buffer	6 x 1 mL	

Description

S1 Nuclease degrades single-stranded nucleic acids, releasing 5'-phosphoryl mono- or oligonucleotides. It is five times more active on DNA than on RNA (1).

S1 Nuclease also cleaves dsDNA at the single-stranded region caused by a nick, gap, mismatch or loop.

S1 Nuclease exhibits 3'-phosphomonoesterase activity. The enzyme is a glycoprotein with a carbohydrate content of 18 %.

Applications

- Removal of single-stranded overhangs of DNA fragments (2).
- S1 transcript mapping (3, 4).
- Cleavage of hairpin loops.
- Creation of unidirectional deletions in DNA fragments in conjunction with Exo III (5).

Source

Aspergillus oryzae cells.

Definition of Activity Unit

One unit of the enzyme produces 1 μ g of acid soluble deoxyribonucleotides in 1 min at 37 °C.

Enzyme activity is assayed in the following mixture: 30 mM sodium-acetate (pH 4.5), 50 mM NaCl, 0.1 mM ZnCl₂, 5 % (v/v) glycerol, 800 μ g/mL heat denatured calf thymus DNA.

Storage Buffer

The enzyme is supplied in: 20 mM Tris-HCl (pH 7.5), 50 mM NaCl, 0.1 mM ZnCl₂ and 50 % (v/v) glycerol.

5X Reaction Buffer

200 mM sodium acetate (pH 4.5 at 25 °C), 1.5 M NaCl and 10 mM ZnSO₄.

Inhibition and Inactivation

- Inhibitors: metal chelators, PP_i, P_i, 5'-ribonucleotides and deoxyribonucleotides.
- Inactivated by heating at 70 °C for 10 min in the presence of EDTA.

Note

S1 Nuclease can introduce breaks into double-stranded DNA, RNA and DNA/RNA hybrids at high enzyme and low salt concentrations (6).

Protocol for Removal of 3'- and 5'-overhangs with S1 Nuclease

S1 Nuclease removes 3' and 5' single stranded DNA overhangs and hairpin loops. The activity of S1 Nuclease is substrate-dependent and the optimal enzyme and DNA amounts for successful blunting should be determined experimentally.

1. Prepare the following reaction mixture:

Components	Volume
DNA	~1 µg
5X Reaction Buffer for S1 Nuclease	6 µL
S1 Nuclease	0.1 µL (10 U)
Water, nuclease-free (#R0581)	to 30 µL
Total volume	30 µL

2. Incubate the mixture at room temperature for 30 min.

3. Stop the reaction by adding 2 µL of 0.5 M EDTA and heating at 70 °C for 10 min.

Note

The S1 Nuclease can be diluted with 1X reaction buffer immediately prior to use.

References

1. Lehman, R.L., Endonucleases specific for single-stranded polynucleotides, *The Enzymes*, 3rd. Ed. (Boyer, P.D., ed.), 4, 193-201, 1981.
2. Roberts T.M., et al., A general method for maximizing the expression of a cloned gene, *Proc. Natl. Acad. Sci. USA*, 76, 760-764, 1979.
3. Berk, A.J., Sharp, P.A., Spliced early mRNAs of simian virus 40, *Proc. Natl. Acad. Sci. USA*, 75, 1274-1278, 1978.
4. Weidle, U., Weissmann, C., The 5'-flanking region of a human IFN-alpha gene mediates viral induction of transcription, *Nature*, 303, 442-446, 1983.
5. Henikoff, S., Unidirectional digestion with exonuclease III creates targeted breakpoints for DNA sequencing, *Gene*, 28, 351-359, 1984.
6. Vogt, V.M., Purification and further properties of single-strand-specific nuclease from *Aspergillus oryzae*, *Eur. J. Biochem.*, 33, 192-200, 1973.

Revision history: Pub. No. MAN0013722

Revision	Date	Description
B00	2025-04-28	New design template, updated buffer quantity, removed COA part

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



Thermo Fisher Scientific Baltics UAB | V.A. Graiciuno 8, LT-02241 Vilnius, Lithuania
For descriptions of symbols on product labels or product documents, go to thermofisher.com/symbols-definition.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2025 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

thermofisher.com/support | thermofisher.com/askaquestion
thermofisher.com

28 April 2025

ThermoFisher
S C I E N T I F I C