

## Lipofectin<sup>®</sup> Reagent

Part no. 18292.lipo.pps	MAN0001376	Rev. Date 14 July 2011
Cat. nos.:	Size:	Store at 4°C (do not freeze)
18292-011	1 mL	
18292-037	4 x 1 mL	

### Description

Lipofectin<sup>®</sup> Reagent is suitable for transfecting DNA (1,2), RNA (3,4), and oligonucleotides (5) into mammalian cells, and DNA into plant protoplasts (6,7). Lipofectin<sup>®</sup> Reagent is specifically recommended for transfecting endothelial cells (8). Refer to the Cell Lines database at [www.invitrogen.com](http://www.invitrogen.com) for a list of other cell types successfully transfected. Lipofectin<sup>®</sup> Reagent is a 1:1 (w/w) liposome formulation of the cationic lipid N-[1-(2,3-dioleoyloxy)propyl]-n,n,n-trimethylammonium chloride (DOTMA) and dioleoyl phosphotidylethanolamine (DOPE) in membrane filtered water.

### Important Guidelines for Transfection

- Form complexes using the recommended amounts of DNA and Lipofectin<sup>®</sup> Reagent. Optimize as necessary (see pages 2–3).  
**Note:** We recommend diluting DNA and Lipofectin<sup>®</sup> Reagent into Opti-MEM<sup>®</sup> I Reduced Serum Medium (Cat. no. 31985-062) before complexing.
- Transfect cells at the recommended confluence or cell density. Optimize as necessary. Maintain the same seeding conditions between experiments.
- *Do not* add antibiotics to media during transfection; this causes cell death.
- For optimal results, transfect cells in medium without serum. Cells may be transfected in the presence of serum, if desired. However, complexes *must* be formed in serum-free medium.
- Test serum-free media for compatibility with Lipofectin<sup>®</sup> Reagent since some serum-free formulations (e.g. CD 293, 293 SFM II, VP-SFM) may inhibit cationic lipid-mediated transfection.

**Intended Use:** For research use only.

Not intended for any animal or human therapeutic or diagnostic use.

## Transfecting Adherent Mammalian Cells

Use the following procedure to transiently or stably transfect mammalian cells. **Note:** All amounts and volumes are given on a per-well basis.

1. One day before transfection, plate cells in growth medium without antibiotics such that they will be at the recommended confluence at the time of transfection.

Condition	Cell no.	Growth med. vol.	Format	Confluence at time of transfection
Transient	$1-2 \times 10^5$	2 mL	6-well	40–60%
Stable	$1-2 \times 10^5$	4 mL	60-mm	30–50%

2. **For each transfection sample**, prepare complexes as follows:
  - a. Dilute 1–2  $\mu\text{g}$  of DNA in 100  $\mu\text{L}$  Opti-MEM<sup>®</sup> I Reduced Serum Medium (or other medium) without serum.
  - b. Mix Lipofectin<sup>®</sup> Reagent before use, then dilute 2–20  $\mu\text{L}$  of Lipofectin<sup>®</sup> Reagent in 100  $\mu\text{L}$  Opti-MEM<sup>®</sup> I Medium (or other medium) without serum. Let stand at room temperature for 30–45 minutes.
  - c. Combine the diluted DNA with diluted Lipofectin<sup>®</sup> Reagent (total volume = 200  $\mu\text{L}$ ). Mix gently and incubate for 10–15 minutes at room temperature (the solution may appear cloudy).
3. Remove the growth medium from the cells and wash the cells once with 2 mL of growth medium without serum. Remove the wash medium.
4. Add 0.8 mL Opti-MEM<sup>®</sup> I Medium (or other medium) without serum (1.8 mL for stable transfection) to the complexes (from step 2c of this procedure); mix gently and add to cells.
5. Incubate cells at 37°C in a CO<sub>2</sub> incubator for 5–24 hours.
6. Replace the medium with 2 mL of growth medium containing serum (4 mL for stable transfection). Alternatively, add 1 mL of growth medium containing 2X the normal concentration of serum (2 mL for stable transfection) without removing the DNA-containing medium.
7. **Transient:** Test for transgene expression 24–72 hours post-transfection.  
**Stable cell lines:** Passage cells at a 1:5 (or higher dilution) into selective medium 48–72 hours post-transfection.

## Transfecting Suspension Mammalian Cells

Use the following procedure to transfect mammalian cells in suspension in a 6-well format. All amounts and volumes are given on a per well basis.

1. On the day of transfection, prepare a single-cell suspension from stock cells. Wash the cells once with serum-free growth medium without antibiotics, and seed cells at a density of  $2\text{--}3 \times 10^6$  cells per well in 0.8 mL of serum-free growth medium without antibiotics.
2. For each transfection sample, prepare complexes as follows:
  - a. Dilute 1–5  $\mu\text{g}$  of DNA in 100  $\mu\text{L}$  of Opti-MEM<sup>®</sup> I Reduced Serum Medium (or other medium) without serum.
  - b. Mix Lipofectin<sup>®</sup> Reagent before use, then dilute 2–25  $\mu\text{L}$  of Lipofectin<sup>®</sup> Reagent in 100  $\mu\text{L}$  of Opti-MEM<sup>®</sup> I Medium (or other medium) without serum. Let stand at room temperature for 30–45 minutes.
  - c. Combine the diluted DNA with diluted Lipofectin<sup>®</sup> Reagent (total volume = 200  $\mu\text{L}$ ). Mix gently and incubate for 10–15 minutes at room temperature (the solution may appear cloudy).
3. Add the 200  $\mu\text{L}$  of complexes to cells. Mix gently by rocking the plate back and forth.
4. Incubate cells at 37°C in a CO<sub>2</sub> incubator for 5–24 hours.
5. The following day, add 4 mL of complete growth medium to the cells.
6. Incubate cells at 37°C in a CO<sub>2</sub> incubator for 24–48 hours prior to testing for transgene expression.

## Scaling Up or Down Transfections

Procedures are provided to transfect cells in a 6-well format (60-mm format for stable mammalian cell transfection). For other formats, vary the amounts of DNA, Lipofectin<sup>®</sup> Reagent, cells, and medium used in proportion to the relative surface area of the tissue culture vessel.

## Optimizing Transfection

To obtain the highest transfection efficiency and low non-specific effects, optimize transfection conditions by varying cell density, DNA and Lipofectin<sup>®</sup> Reagent concentrations, and transfection incubation time.

## Certificate of Analysis

The Certificate of Analysis provides detailed quality control information for each product. Certificates of Analysis are available on our website. Go to [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support) and search for the Certificate of Analysis by product lot number, which is printed on the box.

## References

1. Felgner, P.L., *et al.* (1987) *Proc. Natl. Acad. Sci. USA* 84, 7413.
2. Invitrogen (1989) *Focus*<sup>®</sup> 11, 13.
3. Malone, R.W., *et al.* (1989) *Proc. Natl. Acad. Sci. USA* 86, 6077.
4. Malone, R.W. (1989) *Focus*<sup>®</sup> 11, 61.
5. Chiang, M.-Y., *et al.* (1991) *J. Biol. Chem.* 266, 18162.
6. Sporlein, B. and Koop, H.-U. (1991) *Theor. Appl. Genet.* 83, 1.
7. Zhu, Z., *et al.*, (1990) *Focus*<sup>®</sup> 12, 41.
8. Tilkins, M.L. (1994) *Focus*<sup>®</sup> 16, 117.

## Limited Use Label License

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

©2011 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

For support visit [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support) or email [techsupport@lifetech.com](mailto:techsupport@lifetech.com)

[www.lifetechnologies.com](http://www.lifetechnologies.com)

