

MitoTracker® Mitochondrion-Selective Probes

[Buy Now](#)
[Learn More](#)

Table 1. Contents and storage information.

Material	Amount	Storage	Stability
MitoTracker® probes	20 vials each containing 50 µg lyophilized solid	<ul style="list-style-type: none"> • ≤-20°C * • Desiccate • Protect from light • Avoid freeze-thaw cycles 	When in solid form stored as directed, product is stable for at least 6 months.
Reduced rosamine MitoTracker® probes (M7511, M7513)	20 vials each containing 50 µg lyophilized solid	<ul style="list-style-type: none"> • ≤-20°C * • Store under argon or nitrogen • Desiccate • Protect from light • Avoid freeze-thaw cycles 	

* Before opening a vial, allow the product to warm to room temperature.

Approximate fluorescence excitation and emission maxima, in nm: See Table 2.

Introduction

MitoTracker® Probes

The cell-permeant MitoTracker® probes contain a mildly thiol-reactive chloromethyl moiety for labeling mitochondria (Figure 1). Additionally, some MitoTracker® probes such as MitoTracker® Orange CMTMRos, MitoTracker® Orange CM-H₂TMRos, MitoTracker® Red CMXRos, MitoTracker® Red CM-H₂XRos, and MitoTracker® Deep Red are retained in the mitochondria after fixation. MitoTracker® Green FM and MitoTracker® Red FM are good live cell options.

To label mitochondria, cells are simply incubated with MitoTracker® probes, which passively diffuse across the plasma membrane and accumulate in active mitochondria. Once their mitochondria are labeled, the cells can be treated with an aldehyde-based fixative for samples that need fixation to allow further processing of the sample. Some of the MitoTracker® probes are also retained after permeabilization with some detergents during subsequent processing steps (e.g., immunocytochemistry or *in situ* hybridization). In addition, MitoTracker® probes eliminate some of the difficulties of working with pathogenic cells because once the mitochondria are stained, the cells can be treated with fixatives before the sample is analyzed.

Although conventional fluorescent stains for mitochondria, such as tetramethylrhodamine and rhodamine 123, are readily sequestered by functioning mitochondria, these stains are easily washed out of cells once the mitochondria experience a loss in membrane potential. This characteristic limits the use of such conventional stains in experiments that require cells to be treated with aldehyde fixatives or with other agents that affect the energetic state of the mitochondria. To overcome this limitation, use the MitoTracker® probes—a series of patented

mitochondrion-selective stains that are concentrated by active mitochondria and retained during cell fixation.¹

The MitoTracker® probes differ in spectral characteristics and fixability (Figure 2, Table 2). MitoTracker® probes are provided in specially packaged sets of 20 vials, each containing 50 µg for reconstitution as required.

Rosamine- and carbocyanine-based MitoTracker® dyes: We offer rosamine-based MitoTracker® dyes which include MitoTracker® Orange CMTMRos (Cat. no. [M7510](#)), a derivative of tetramethylrosamine, and MitoTracker® Red CMXRos (Cat. no. [M7512](#)), a derivative of X-rosamine. Reduced MitoTracker® dyes, MitoTracker® Orange CM-H₂TMROS (Cat. no. [M7511](#)) and MitoTracker® Red CM-H₂XROS (Cat. no. [M7513](#)), which are derivatives of dihydrotetramethylrosamine and dihydro-X-rosamine, respectively, are also available. These reduced probes do not fluoresce until they enter live cells, where they are oxidized to the corresponding fluorescent mitochondrion-selective probe and then sequestered in the mitochondria.

The carbocyanine-based MitoTracker® dyes include MitoTracker® Red FM (Cat. no. [M22425](#)), MitoTracker® Green FM dye (Cat. no. [M7514](#)), and MitoTracker® Deep Red FM (Cat. no. [M22426](#)).

MitoTracker® Red CMXRos, MitoTracker® Red FM, and MitoTracker® Deep Red FM are all well suited for multicolor labeling experiments because their red fluorescence is well resolved from the green fluorescence of other probes.

Some of the MitoTracker® probes are well retained after fixation and permeabilization (Table 2). While some mitochondrial labeling may be achieved in fixed cells using MitoTracker® dyes, the signal-to-noise ratios are not generally favorable and these dyes are recommended for use in live cells only. If mitochondrial labeling of fixed cells is desired, using anti-OxPhos antibodies available from Invitrogen (for example, Cat. no. [A21355](#)) is a good option.

Before You Begin

Materials Required but Not Provided

- DMSO
- Suitable buffer or growth medium for live cell imaging
- Aldehyde based fixatives such as formaldehyde for cell fixation
- Aldehyde based detergents such as Triton® X-100




Figure 1. The intracellular reactions of our fixable mitochondrion-selective dye, MitoTracker® Orange CM-H₂TMROS. When this cell-permeant probe enters an actively respiring cell, it is oxidized to MitoTracker® Orange CMTMRos and sequestered in the mitochondria, where it reacts with thiols on proteins and peptides to form an aldehyde-fixable conjugate.

Table 2. MitoTracker® Mitochondrion-Selective Probes

	Product	Catalog Number	Ex (nm) *	Em (nm) *	Oxidation State	Retained after Fixation
Rosamine-based MitoTracker® Probes						
Buy Now	MitoTracker® Orange CMTMRos	M7510	554	576	Oxidized	Yes
Buy Now	MitoTracker® Orange CM-H ₂ TMRos	M7511 ‡	554	576	Reduced	Yes
Buy Now	MitoTracker® Red CMXRos	M7512	579	599	Oxidized	Yes
Buy Now	MitoTracker® Red CM-H ₂ XRos	M7513 ‡	579	599	Reduced	Yes
Carbocyanine-based MitoTracker® Probes						
Buy Now	MitoTracker® Green FM	M7514 †	490	516	NA	No
Buy Now	MitoTracker® Red FM	M22425	581	644	NA	No
Buy Now	MitoTracker® Deep Red FM	M22426	644	665	NA	Yes

* Fluorescence excitation (Ex) and emission (Em) maxima determined in methanol; values may vary somewhat in cellular environments.
† Nonfluorescent in aqueous solution. ‡ Nonfluorescent until oxidized. NA = Not applicable.

Preparing Stock Solutions

Before opening a vial, allow the product to warm to room temperature.

To prepare a stock solution, dissolve the lyophilized MitoTracker® product in high-quality, anhydrous dimethylsulfoxide (DMSO) to a final concentration of 1 mM; the molecular weight (MW) is indicated on the product label. The reduced rosamine MitoTracker® probes (M7511, M7513) are quite sensitive to oxidation, especially in solution, and must be stored under argon or nitrogen, at ≤−20°C and protected from light. It is preferable to use solutions of the dihydro derivatives immediately after they are prepared. Store all other solutions of the MitoTracker® dyes frozen at ≤−20°C and protected from light.

Experimental Protocols**Cell Preparation and Staining**

The concentration of probe for optimal staining varies by application. The initial conditions suggested here are guidelines that may be modified based on the particular cell type or on other factors, such as the permeability of the cells or tissues to the probe. In general, the reduced rosamine MitoTracker® probes (M7511, M7513) are loaded at three- to five-fold higher concentrations than other MitoTracker® probes.

- 1.1 Preparing staining solutions.** Dilute 1 mM MitoTracker® stock solution (see *Preparing Stock Solutions* for preparation) to the final working concentration in appropriate buffer or growth medium. Because the reduced forms of the MitoTracker® probes are susceptible to potential oxidases in serum, we do not recommend using complete media with these dyes.

General Guidelines: Use working concentrations of 25–500 nM. For staining cells that are to be fixed and permeabilized (see *Fixation and Permeabilization after Staining*), use a working concentration of 100–500 nM. To reduce potential artifacts and mitochondrial toxicity from overloading, keep the concentration of dye as low as possible. For the MitoTracker® Green FM probes, use a slightly lower concentration (20–200 nM). At higher concentrations, these probes tend to stain other cellular structures.

- 1.2 Staining adherent cells.** Grow cells on coverslips inside a Petri dish filled with the appropriate culture medium. When cells have reached the desired confluency, remove the media from the dish and add prewarmed (37°C) staining solution containing MitoTracker® probe (prepared in step 1.1). While incubation times vary depending on the model system and probe used, incubation for 15–45 minutes under growth conditions appropriate for the particular cell type is generally sufficient but may need to be optimized. After staining is complete

replace the staining solution with fresh prewarmed media or buffer and observe cells using a fluorescence microscope or fluorescence microplate reader. If the cells are to be fixed and permeabilized, continue to *Fixation and Permeabilization after Staining*.

- 1.3 Staining suspension cells.** Centrifuge to obtain a cell pellet and aspirate the supernatant. Resuspend the cells gently in prewarmed (37°C) staining solution containing the MitoTracker® probe (prepared in step 1.1). While incubation times vary depending on the model system and probe used, incubation for 15–45 minutes under growth conditions appropriate for the particular cell type is generally sufficient but may need to be optimized. After staining is complete, re-pellet the cells by centrifugation and resuspend cells in fresh prewarmed medium or buffer. Cells may be analyzed by flow cytometry, microplate-based analysis, or fluorescence microscopy. If immobilized cells on coverslips are needed, use poly-D-lysine to coat the slides or coverslips before mounting. If the cells are to be fixed and permeabilized, continue to *Fixation and Permeabilization after Staining*.

Optional: Fixation and Permeabilization after Staining

After staining live cells with a MitoTracker® dye, it is often useful to fix and permeabilize the cells for subsequent manipulations. For example, fixation and permeabilization allows you to probe for other intracellular structures by immunocytochemistry. Most of the MitoTracker® dyes are well-retained following fixation and permeabilization using the protocol described below. However, MitoTracker® Green FM and MitoTracker® Red FM are not retained well after fixation.

- 2.1 Washing the cells.** After staining, wash the cells in fresh, pre-warmed buffer or growth medium.
- 2.2 Fixing the cells.** Carefully remove the medium/buffer covering the cells, and replace it with freshly prepared, pre-warmed buffer or growth medium containing 2–4% formaldehyde. For MitoTracker® Red CMXRos, we have found that fixing with 3.7% formaldehyde in complete growth medium at 37°C for 15 minutes works well for endothelial cells.
- 2.3 Rinsing the cells.** After fixation, rinse the cells several times in buffer.
- 2.4 Permeabilization (optional).** When permeabilization is needed for subsequent steps such as immunocytochemistry, incubate fixed cells in buffer containing detergent such as Triton® X-100. Following permeabilization, rinse the cells in buffer and proceed with immunocytochemistry procedure. We found that incubating the endothelial cells for 10 minutes in PBS containing 0.2% Triton® X-100 works well.

Alternatively, the cells may be permeabilized by incubating in ice-cold acetone for 5 minutes, and then washed in PBS. Even when cells are not going to be labeled with an antibody, this acetone-permeabilization step may be useful in improving the signal-to-background ratio.

Fluorescence Microscopy

Spectral characteristics for the MitoTracker® probes are summarized in Table 2.




Figure 2. Fluorescence excitation and emission spectra for MitoTracker® Green FM (A), MitoTracker® Orange (B), MitoTracker® Red CMXRos (C), and MitoTracker® Deep Red (D). All spectra were determined in methanol.

References

1. J Histochem Cytochem 44, 1363 (1996);
2. Cytometry 39, 203 (2000);
3. Mol Cell Biochem 172, 171 (1997).

Product List Current prices may be obtained from our website or from our Customer Service Department.

Cat. no.	Product Name	Unit Size
M22426	MitoTracker® Deep Red FM *special packaging*	20 × 50 µg
M7514	MitoTracker® Green FM *special packaging*	20 × 50 µg
M7511	MitoTracker® Orange CM-H ₂ TMRos *special packaging*	20 × 50 µg
M7510	MitoTracker® Orange CMTMRos *special packaging*	20 × 50 µg
M22425	MitoTracker® Red FM *special packaging*	20 × 50 µg
M7513	MitoTracker® Red CM-H ₂ XRos *special packaging*	20 × 50 µg
M7512	MitoTracker® Red CMXRos *special packaging*	20 × 50 µg
Related Products		
A21355	anti-OxPhos Complex V inhibitor protein, mouse IgG1, monoclonal antibody	100 µg
M36008	MitoSOX™ Red mitochondrial superoxide indicator *for live-cell imaging*	10 × 50 µg

Contact Information

Molecular Probes, Inc.

29851 Willow Creek Road
Eugene, OR 97402
Phone: (541) 465-8300
Fax: (541) 335-0504

Customer Service:

6:00 am to 4:30 pm (Pacific Time)
Phone: (541) 335-0338
Fax: (541) 335-0305
probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209
Order Fax: (800) 438-0228

Technical Service:

8:00 am to 4:00 pm (Pacific Time)
Phone: (541) 335-0353
Toll-Free (800) 438-2209
Fax: (541) 335-0238
probetech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd.
3 Fountain Drive
Inchinnan Business Park
Paisley PA4 9RF, UK
Phone: +44 (0) 141 814 6100
Fax: +44 (0) 141 814 6260
Email: euroinfo@invitrogen.com
Technical Services: eurotech@invitrogen.com

**For country-specific contact information,
visit www.invitrogen.com.**

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Assistance Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402. Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2008, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.