

AmpErase[®] Uracil N-Glycosylase (UNG)

 Package Contents	Catalog Number N8080096	Size 100 Units	 Kit Contents
 Storage Conditions	<ul style="list-style-type: none"> Store at -20°C until just prior to use. 		
 Required Materials	<ul style="list-style-type: none"> Template: cDNA, gDNA, λDNA Forward and reverse gene-specific primers Autoclaved, distilled water E-Gel[®] General Purpose Gels, 1.2% (Cat. no. G5018-01) 0.2 or 0.5-mL nuclease-free microcentrifuge tubes 100-mM dUTP solution (Cat. no. R0133) <p>For PCR: GeneAmp[®] RNA PCR Core Kit (Cat. no. N8080143)</p> <p>For qPCR:</p> <ul style="list-style-type: none"> TaqMan[®] GAPDH Control Reagents (human) (Cat. no. 402869) TaqMan[®] Universal PCR Master Mix, No AmpErase[®] UNG (Cat. no. 4326614) 20X TaqMan[®] Primer/Probe Mix MicroAmp[®] EnduraPlate[™] Optical 384-well Plate (Cat. no. 4483285) 		
 Timing	Varies depending on amplicon length		
 Selection Guides	<p>Go online to view related products. PCR Enzymes and Master Mixes</p>		
 Product Description	<ul style="list-style-type: none"> AmpErase[®] Uracil N-Glycosylase (UNG) is a 26-kDa ultrapure, recombinant enzyme encoded by the <i>E. coli</i> uracil N-glycosylase gene, designed to degrade PCR products from previous amplifications or mis-primed, non-specific products, without degrading native nucleic acid templates. PCR products are made susceptible to degradation by substituting dUTP for dTTP in the reaction mix and testing subsequent mixes with the UNG prior to amplification. 		
 Important Guidelines	<ul style="list-style-type: none"> Take precautions to avoid cross-contamination by using aerosol-resistant barrier tips and analyzing PCR products in a separate area from PCR assembly. Because AmpErase[®] UNG has activity below 55°C, the annealing temperature should be ≥ 55°C to avoid degradation of newly synthesized dU-containing products by residual UNG activity. Dilute 100 mM dUTP to 20 mM prior to using in PCR reaction. For qPCR applications, dilute UNG 1:5 in TE pH 8.0. 		
 Online Resources	<p>Visit our product page for additional information and protocols. For support, visit www.lifetechnologies.com/support.</p>		



PCR Reaction Setup

Use the measurements below to prepare your PCR experiment, or enter your own parameters in the column provided.

Component	100-μL rxn	Custom	Final Conc.
Autoclaved, distilled water	to 100 μL	to μL	–
10X PCR Buffer II	10 μL	μL	1X
10 mM dATP	2.0 μL	μL	0.2 mM
10 mM dCTP	2.0 μL	μL	0.2 mM
10 mM dGTP	2.0 μL	μL	0.2 mM
20 mM dUTP*	1–5 μL	μL	0.2–1.0 mM
25 mM MgCl ₂ **	8.0 μL	μL	2.0 mM
10 μM Forward primer	2.0 μL	μL	0.2 μM
10 μM Reverse primer	2.0 μL	μL	0.2 μM
Template DNA	Varies		< 2.0 μg/100 μL
AmpliAq [®] DNA Polymerase (5 U/μL)	0.5 μL	μL	2.5 U/100 μL
AmpErase [®] UNG (1 U/μL)	1 μL	μL	1 U/100 μL

* Dilute dUTP to 20 mM prior to using it in PCR reaction. Due to the lower efficiency of dUTP incorporation, using any amount of dUTP in the recommended range may result in lower amplification than using the standard dTTP concentration (0.2 mM/50-μL rxn).

** Increase MgCl₂ on an equimolar basis when an increased dUTP concentration is required. The magnesium ion is not required for UNG activity.

PCR and qPCR Protocols

 See pages 2 and 3 to view procedures for preparing and running your PCR and qPCR experiments.

Optimization Strategies



Refer to the pop-ups below for guidelines to optimize your PCR reactions.

 [DUG Guidelines](#)  [PCR Guidelines](#)  [qPCR Guidelines](#)

 [Limited Warranty, Disclaimer, and Licensing Information](#)


PCR Amplification Protocol

The example PCR procedure below shows appropriate volumes for a single 100- μ L reaction. For multiple reactions, prepare a master mix of components common to all reactions to minimize pipetting error, and then dispense appropriate volumes into each 0.2–0.5 mL PCR reaction tube prior to adding template DNA and primers.

Timeline	Steps	Procedure Details																								
1 	Thaw reagents	Thaw, mix, and briefly centrifuge each component before use. Note: Avoid generating bubbles when mixing the enzyme.																								
2 	Prepare PCR master mix	Add the following components to each PCR reaction tube. Note: Consider the volumes for all components listed in steps 2 and 3 to determine the correct amount of water required to reach your final reaction volume.																								
		<table border="1"> <thead> <tr> <th>Component</th> <th>100-μL rxn</th> <th>Final Concentration</th> </tr> </thead> <tbody> <tr> <td>Autoclaved, distilled water</td> <td>to 100 μL</td> <td>–</td> </tr> <tr> <td>10X PCR Buffer II</td> <td>10 μL</td> <td>1X</td> </tr> <tr> <td>10 mM dNTPs (dATP, dCTP, dGTP)</td> <td>2 μL each</td> <td>0.2 mM</td> </tr> <tr> <td>20 mM dUTP*</td> <td>2 μL</td> <td>0.4 mM</td> </tr> <tr> <td>25 mM MgCl₂**</td> <td>8 μL</td> <td>2.0 mM</td> </tr> <tr> <td>AmpliTaq® DNA Polymerase (5 U/μL)</td> <td>0.5 μL</td> <td>2.5 U/100 μL</td> </tr> <tr> <td>AmpErase® UNG (1 U/μL)</td> <td>1 μL</td> <td>1 U/100 μL</td> </tr> </tbody> </table>	Component	100- μ L rxn	Final Concentration	Autoclaved, distilled water	to 100 μ L	–	10X PCR Buffer II	10 μ L	1X	10 mM dNTPs (dATP, dCTP, dGTP)	2 μ L each	0.2 mM	20 mM dUTP*	2 μ L	0.4 mM	25 mM MgCl ₂ **	8 μ L	2.0 mM	AmpliTaq® DNA Polymerase (5 U/ μ L)	0.5 μ L	2.5 U/100 μ L	AmpErase® UNG (1 U/ μ L)	1 μ L	1 U/100 μ L
Component	100- μ L rxn	Final Concentration																								
Autoclaved, distilled water	to 100 μ L	–																								
10X PCR Buffer II	10 μ L	1X																								
10 mM dNTPs (dATP, dCTP, dGTP)	2 μ L each	0.2 mM																								
20 mM dUTP*	2 μ L	0.4 mM																								
25 mM MgCl ₂ **	8 μ L	2.0 mM																								
AmpliTaq® DNA Polymerase (5 U/ μ L)	0.5 μ L	2.5 U/100 μ L																								
AmpErase® UNG (1 U/ μ L)	1 μ L	1 U/100 μ L																								
3 	Add primers and template DNA	* Dilute 100 mM stock to 20 mM prior to use. ** Increase MgCl ₂ on an equimolar basis, when an increased dUTP concentration is required. The magnesium ion is not required for UNG activity. Mix and briefly centrifuge the components.																								
		Add the primers and template DNA to each tube for a final reaction volume of 100 μ L.																								
		<table border="1"> <thead> <tr> <th>Component</th> <th>100-μL rxn</th> <th>Final Concentration</th> </tr> </thead> <tbody> <tr> <td>10 μM forward primer</td> <td>2 μL</td> <td>0.2 μM</td> </tr> <tr> <td>10 μM reverse primer</td> <td>2 μL</td> <td>0.2 μM</td> </tr> <tr> <td>Template DNA</td> <td>varies</td> <td>< 2.0 μg/rxn</td> </tr> </tbody> </table>	Component	100- μ L rxn	Final Concentration	10 μ M forward primer	2 μ L	0.2 μ M	10 μ M reverse primer	2 μ L	0.2 μ M	Template DNA	varies	< 2.0 μ g/rxn												
Component	100- μ L rxn	Final Concentration																								
10 μ M forward primer	2 μ L	0.2 μ M																								
10 μ M reverse primer	2 μ L	0.2 μ M																								
Template DNA	varies	< 2.0 μ g/rxn																								
4 	Incubate reactions in a thermal cycler	Cap, mix, and then centrifuge the contents.																								
		<table border="1"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> </tr> </thead> <tbody> <tr> <td>Pre-incubation</td> <td>50°C</td> <td>2 minutes</td> </tr> <tr> <td>Initial Denaturation</td> <td>95°C</td> <td>2 minutes</td> </tr> <tr> <td rowspan="3">35 PCR Cycles</td> <td>Denature</td> <td>95°C</td> </tr> <tr> <td>Anneal</td> <td>~55°C (depending on Primer T_m)*</td> <td>30 seconds</td> </tr> <tr> <td>Extend</td> <td>37–65°C</td> <td>1 minute/kb</td> </tr> <tr> <td>Hold</td> <td>72°C</td> <td>7 minutes</td> </tr> </tbody> </table>	Step	Temperature	Time	Pre-incubation	50°C	2 minutes	Initial Denaturation	95°C	2 minutes	35 PCR Cycles	Denature	95°C	Anneal	~55°C (depending on Primer T _m)*	30 seconds	Extend	37–65°C	1 minute/kb	Hold	72°C	7 minutes			
Step	Temperature	Time																								
Pre-incubation	50°C	2 minutes																								
Initial Denaturation	95°C	2 minutes																								
35 PCR Cycles	Denature	95°C																								
	Anneal	~55°C (depending on Primer T _m)*	30 seconds																							
	Extend	37–65°C	1 minute/kb																							
Hold	72°C	7 minutes																								
5 	Analyze with gel electrophoresis	* Adjust the temperature according to the primer melting temperature, but ensure it is at least 50°C. Analyze 10 μ L using agarose gel electrophoresis. Use your PCR reaction immediately for down-stream applications, or store it at –20°C.																								

qPCR Protocol

The example qPCR procedure below shows appropriate volumes for a single **20- μ L** reaction. For multiple reactions, prepare a master mix of components common to all reactions to minimize pipetting error, and then dispense appropriate volumes into each 0.2–0.5 mL PCR reaction tube, or well of a MicroAmp® EnduraPlate™ Optical 384-well plate, prior to adding template DNA and primers.

Timeline	Steps	Procedure Details																		
1 	Thaw reagents	Thaw reagents. Set up reactions on ice. For 384-well plates, we recommend a maximum reaction volume of 10 μ L per well. Note: Always prepare no-template control (NTC) reactions to test for DNA contamination of the enzyme/primer mixes. Use TaqMan® GAPDH control reagents to confirm reaction set-up.																		
2 	Prepare qPCR master mix	<table border="1"> <thead> <tr> <th>Component</th> <th>20-μL rxn</th> <th>Final Concentration</th> </tr> </thead> <tbody> <tr> <td>Autoclaved, distilled water</td> <td>to 20 μL</td> <td>–</td> </tr> <tr> <td>TaqMan® Universal PCR Master Mix, No AmpErase® UNG</td> <td>10 μL</td> <td>1X</td> </tr> <tr> <td>AmpErase® UNG (1 U/μL)*</td> <td>0.2 μL</td> <td>0.04 U/20 μL</td> </tr> <tr> <td>20X Fluorescent Primer/Probe Mix</td> <td>1 μL</td> <td>1X</td> </tr> <tr> <td>cDNA (10 to 100 ng)</td> <td>2 μL</td> <td>varies</td> </tr> </tbody> </table>	Component	20- μ L rxn	Final Concentration	Autoclaved, distilled water	to 20 μ L	–	TaqMan® Universal PCR Master Mix, No AmpErase® UNG	10 μ L	1X	AmpErase® UNG (1 U/ μ L)*	0.2 μ L	0.04 U/20 μ L	20X Fluorescent Primer/Probe Mix	1 μ L	1X	cDNA (10 to 100 ng)	2 μ L	varies
Component	20- μ L rxn	Final Concentration																		
Autoclaved, distilled water	to 20 μ L	–																		
TaqMan® Universal PCR Master Mix, No AmpErase® UNG	10 μ L	1X																		
AmpErase® UNG (1 U/ μ L)*	0.2 μ L	0.04 U/20 μ L																		
20X Fluorescent Primer/Probe Mix	1 μ L	1X																		
cDNA (10 to 100 ng)	2 μ L	varies																		
3 	Incubate reactions in a real-time instrument	* Prepare a fresh 1:5 dilution in TE buffer prior to use. It may be necessary to optimize the dilution of UNG (1:5–1:10) for specific targets. We recommend a fresh 1:5 dilution as a starting point. Note: Do not store or reuse diluted enzyme. Cap or seal each PCR tube/plate, gently mix, and centrifuge contents.																		
4 	Collect and analyze data	<table border="1"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> </tr> </thead> <tbody> <tr> <td>Pre-incubation</td> <td>50°C</td> <td>2 minutes</td> </tr> <tr> <td>Heat activation</td> <td>95°C</td> <td>10 minutes</td> </tr> <tr> <td rowspan="2">40 cycles</td> <td>95°C</td> <td>15 seconds</td> </tr> <tr> <td>60°C</td> <td>1 minute</td> </tr> </tbody> </table>	Step	Temperature	Time	Pre-incubation	50°C	2 minutes	Heat activation	95°C	10 minutes	40 cycles	95°C	15 seconds	60°C	1 minute				
Step	Temperature	Time																		
Pre-incubation	50°C	2 minutes																		
Heat activation	95°C	10 minutes																		
40 cycles	95°C	15 seconds																		
	60°C	1 minute																		
		Analyze results following your Real-Time instrument manufacturer's guidelines. Optional: The specificity of the PCR products can be checked by agarose gel electrophoresis.																		