

Radioactive and Non-radioactive DNA Labeling by Nick-translation

This protocol is for the Radioactive and Non-radioactive DNA Labeling by Nick-translation.

I. Radioactive DNA labeling by nick-translation

- Mix the following components:

10× reaction buffer for DNA Polymerase I	2.5 μ L
Mixture of 3 dNTPs, 1 mM* (without the labeled dNTP)	1.25 μ L
[alpha-32P]-dNTP, ~110 TBq/ mmol (3000 Ci/mmol)	1.85-3.7 MBq (50-100 μ Ci)
DNase I, RNase-free freshly diluted to 0.002 u/ μ L**	1 μ L
DNA Polymerase I	0.5-1.5 μ L (5-15 u)
Template DNA	0.25 μ g
Water, nuclease-free	to 25 μ L
Total volume	to 10 μL

- Immediately incubate at 15°C for 15-60 minutes.
- Terminate the reaction by adding 1 μ L of 0.5 M EDTA, pH 8.0.
- Take an aliquot (1 μ L) to determine the efficiency of the label incorporation. A specific activity of DNA at least 108 cpm/ μ g DNA is expected.
- If needed, the labeled DNA may be separated from the unincorporated radioactive precursors on Sephadex G-50 or Bio-Gel P-60 column or using spin column (e.g. Thermo Scientific GeneJET PCR Purification Kit (Cat #K0701)).

Note

* To prepare a mixture of three non-labeled dNTPs (1 mM of each), mix 1 μ L aliquots of stock solutions of each dNTP (100 mM, from Cat #R0181) with 97 μ L of Water, nuclease-free. These dNTP mixes can be stored at -20°C for further use.

** The DNase I, RNase-free can be diluted with the 1× reaction buffer for DNA Polymerase I.

- The reaction volumes can be scaled up or down providing that the final concentrations of the components (DNA, dNTPs, labeled dNTP) are as indicated in the protocol.
- Radioactive DNA probes with higher specific activities can be prepared using two radioactively labeled dNTPs simultaneously. In this case, the composition of the unlabeled dNTP mix should be adjusted accordingly.

II. Non-radioactive DNA labeling by nick-translation

The protocol above can be used for non-radioactive labeling by nick-translation using biotin-11-dUTP, fluorescein-12-dUTP, DIG-dUTP or aminoallyl-dUTP:

- normal dTTP is substituted for labeled-dUTP at a molar ratio of 1:3-1:5,
- reaction time is prolonged to 1-2 hours.

thermoscientific.com

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

North America
Technical Services:
 ts.molbio@thermofisher.com
Customer Services:
 cs.molbio@thermofisher.com
 Tel 800 235 9880
 Fax 800 292 6088

Europe and Asia
Technical Services:
 ts.molbio.eu@thermofisher.com
Customer Services:
 cs.molbio.eu@thermofisher.com
 Tel 00800 2220888
 Fax 00800 2220889

Thermo
 SCIENTIFIC

Part of Thermo Fisher Scientific