Thermo Fisher

Page 1 of 3

OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION

TETRACYCLINE TE10 CT0053B

Product Details

Product Name	TETRACYCLINE
Disc Code	TE10
Disc Concentration	10µg
CT Code	СТ0053В
Shelf Life (years)	3
CLSI Recommended Concentration	No
EUCAST Recommended Concentration	No
Legacy BSAC	Yes

Microbiological Performance Criteria

Legacy BSAC (ISA in air)			
Organism	ATCC/NCTC Number	Range (mm)	
Escherichia coli	ATCC [®] 25922	22 – 28	
Escherichia coli	ATCC®10536	23 – 29	
Staphylococcus aureus	ATCC [®] 9144	31 - 40	
Staphylococcus aureus	ATCC [®] 25923	26 – 35	
Enterococcus faecalis	ATCC [®] 29212	9 – 13	

Legacy BSAC (ISA enriched with 5% v/v horse blood in air)			
Organism	ATCC/NCTC Number	Range (mm)	
Staphylococcus aureus	ATCC [®] 9144	30 – 38	
Staphylococcus aureus	ATCC [®] 25923	28 – 36	

Legacy BSAC (ISA enriched with 5% v/v horse blood in CO ₂)			
Organism	ATCC/NCTC Number	Range (mm)	
Staphylococcus aureus	ATCC [®] 9144	33 – 40	
Staphylococcus aureus	ATCC [®] 25923	27 – 34	
Streptococcus pneumoniae	ATCC [®] 49619	26 – 36	
Pasteurella multocida	NCTC8489	29 – 34	

Thermo Fisher S C I E N T I F I C

Page 2 of 3

OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION

TETRACYCLINE TE10 CT0053B

Legacy BSAC (ISA enriched with 5% v/v horse blood and 20mg/l NAD in CO_2)			
Organism ATCC/NCTC Number Range (mm)		Range (mm)	
Haemophilus influenzae	ATCC [®] 49247	9 – 14	
Haemophilus influenzae	NCTC11931	27 – 35	
Neisseria gonorrhoeae	ATCC [®] 49226	27 – 35	

References

The zone limits are in accordance with the legacy BSAC methods for antimicrobial disc susceptibility testing, version 15 January 2015.

Other additional challenging strains may be employed.

The information given above is believed to be correct, however both the information and the product are offered without warranty for any specific use, other than that specified in the current Oxoid manual.



Page 3 of 3

OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION

TETRACYCLINE TE10 CT0053B

Revision History

Section / Step	Description of Change	Reason for Change	Reference
Entire document/ Microbiological Performance Criteria/ References	Update to current format/ Clarification of testing media/Change of BSAC recommended strains to Legacy BSAC	Change control	BT-CC-2889