

Linker Ligation

This protocol is for the Linker Ligation

Double stranded oligonucleotide linkers are often used to generate compatible overhangs not found in the insert. Linkers normally contain restriction enzyme recognition sequences and are digested after ligation to generate overhangs compatible with cloning vectors. Alternatively, linkers may have overhangs which are ready for ligation with a cloning vector and do not require further manipulation following ligation.

1. Prepare the following reaction mixture:

Linear DNA	100-500 ng
Phosphorylated linkers	1-2 µg
10X T4 DNA Ligase buffer	2 µl
50% PEG 4000 solution	2 µl
T4 DNA Ligase	2 u
Water, nuclease-free	to 20 µl
Total volume	to 20 µl

2. Mix thoroughly, spin briefly and incubate for 1 hour at 22°C.
3. Heat inactivate at 65°C for 10 min or at 70°C for 5 min.

Note

T4 DNA Ligase is active in PCR and restriction digestion buffers (see table below). Therefore, linker ligation reactions can be performed in the restriction enzyme buffer optimal for the subsequent digestion. In this case, the ligation reaction should be supplemented with ATP to a final concentration of 0.5 mM. After inactivation of the T4 DNA Ligase, add the restriction enzyme directly to the reaction mixture and incubate according to the digestion protocol.

Activity in PCR and restriction digestion buffers		
Buffers		Activity*, %
PCR and RT buffers		75
Restriction buffers for restriction enzymes	FastDigest®, 1X / 2X Tango™, B, G, O, R, KpnI, BamHI, EcoRI	75-100
	Ecl136II, SacI	50

* activity of T4 DNA Ligase in various buffers supplemented with 0.5 mM ATP.

thermoscientific.com

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

North America
Technical Services:
 techservice.genomics@thermofisher.com
Customer Services:
 customerservice.genomics@thermofisher.com
 Tel 800 235 9880
 Fax 800 292 6088

Europe and Asia
Technical Services:
 techservice.emea.genomics@thermofisher.com
Customer Services:
 customerservice.emea.genomics@thermofisher.com

Thermo
 S C I E N T I F I C

Part of Thermo Fisher Scientific