
thermo scientific

# GLD Pro Gas Leak Detector User Manual

Version 4

BR21603 October 2020

**ThermoFisher**  
SCIENTIFIC

# Contents

1.0	Introduction	4
2.0	Battery charging	5
3.0	Powering up	6
4.0	Zeroing the unit	6
5.0	Prior to operation	6
6.0	Detecting leaks	7
7.0	Specifications	7
8.0	Maintenance	8
9.0	Troubleshooting	8
10.0	Technology	9
11.0	Interpretation of results	10
12.0	Back label explanation	11
13.0	Case specifications	12
14.0	Service	12

# Operating instructions

## 1.0 Introduction

The Thermo Scientific™ leak detector is specifically designed for use with gas chromatography (GC) systems. It detects minute leaks of any gas with a thermal conductivity different from air. The reference gas inlet (Figure 1) draws in ambient air for comparison to air drawn into the sample probe. A leak is indicated by both an LED light display and an audible alarm.

If this instrument is used in any manner other than described in the manual, the CE declaration is void.

**Leak Detector LED Indicator Display and Interpretation**

	Unit State - Powered On/Off	LED Color	LED State	Interpretation
Unit Battery State	On	Blue > Red & Yellow	On/Blue flashes once, red and yellow cycle for 15 seconds/ Off	Unit is powering on and running startup sequence
	On	Blue	On/Steady	Ready for use, battery does not require charging
	<b>On</b>	<b>Blue</b>	<b>On/Constant Flashing</b>	<b>IMPORTANT: Unit must be charged</b>
	Power Button Depressed	NONE	NONE	Battery has fully discharged, charge immediately
Unit Charging State	On/Off	Green	On/Flashing	Unit is plugged in to the AC adaptor or USB cable and is charging a very low battery (trickle charging)
	On/Off	Green	On/Steady	Unit is plugged in to AC adaptor or USB cable and is in full charge mode
	On/Off	Green	Off	Unit is plugged into AC adaptor or USB cable and the unit is fully charged
Unit Operation	On	Red or Yellow	On/Steady	Unit is indicating a difference in thermal conductivity between probe and reference
	On	Red & Yellow	On/Cycling for 4 seconds/Off	Zero set button has been depressed and unit is rezeroing

## 2.0 Battery charging

**The leak detector should be fully charged prior to use.** Only use the AC adaptor provided. To charge the battery, first install the correct plug for your country's AC outlets onto the provided AC adaptor. Insert the AC adaptor into an electrical outlet, and then insert the barrel plug on the other end of the AC adaptor into the connector on the bottom of the leak detector unit. The green battery charge indicator LED will illuminate. When the battery is fully charged, the green battery charge indicator LED will go out. When the leak detector's charge is low, the BLUE LED located between the red and yellow LEDs will begin to flash.

If unit is off, the BLUE LED may flash when the power button is depressed. If the battery is fully discharged, no LED will illuminate.

**NOTE:** If the battery is fully discharged, no LED will illuminate.

**CAUTION:** For the health of the battery, it is not recommended to continuously operate the leak detector when it is plugged in and charging. However, if your battery has been allowed to fully discharge and will no longer hold a charge, the unit can be operated while plugged in receiving charge directly from the power supply. To operate the leak detector while it is plugged in, follow the instructions for regular operation, and when the blue LED illuminates steady or flashing the unit is ready to use.

**CAUTION:** The leak detector contains a nickel-metal hydride (NiMH) battery. Like other battery-operated devices, if the battery is left for long periods of time without being charged, it can discharge to a point where the battery cannot accept a charge. **We recommend that you charge the battery at least once every 3 months. If your leak detector is used often for routine maintenance, as recommended, more frequent charging may be necessary.**

**CAUTION:** Replacement of the rechargeable battery in this unit is performed at the factory. There are no serviceable parts inside this unit. Opening the case or tampering with the internal parts will void the factory warranty.

**NOTE:** Recharging a low battery can take approximately 6 hours.

**NOTE:** The leak detector may get warm to the touch during charging, but this is normal and not a safety concern.

### 3.0 Powering up

Depress and hold the power button (Figure 1) until the unit responds with the wake-up mode. The leak detector will run through a self calibration sequence for approximately 15 seconds. During this time DO NOT attempt to zero the unit.

### 4.0 Zeroing the unit

After the LED lights stop flashing, the unit is ready for use. The instrument may need to be zeroed periodically between uses, especially if it is moved from room to room or between areas of differing temperature or humidity. Do not attempt to zero the unit while the probe is stored in the holder. The probe MUST be removed from the probe docking station before zeroing the unit. To re-zero, press the zero set button. The unit will run a self-calibration sequence for approximately 4 seconds. When all LED lights stop flashing and the blue LED light is lit, the unit is ready for use.

**NOTE:** To avoid false readings, do not attempt to use or zero the unit while the self-calibration sequence is in progress.


Figure 1 – Leak detector schematic

### 5.0 Prior to operation

Verify the operation of the leak detector before each use by sampling gas from a GC split vent or other source of hydrogen or helium. Also, visually inspect the probe tip, reference gas inlet, and exhaust port for obstructions (Figure 1).

**IMPORTANT:** Fittings being checked must be clean and dry; liquid leak detecting agents, dust, and other debris may damage the leak detector if drawn into the probe.

The leak detector responds to almost any gas you can smell and many gases that you can't smell. Solvent vapors, split vent exhaust, or even strong air currents around the probe or reference inlet can cause instability or false positive readings. Be careful not to breathe into the reference inlet when checking for leaks or to cover/block the inlet with your hand.

## 6.0 Detecting leaks

Slowly move the probe tip around fittings and other potential leak sources. If the leak detector senses a gas other than air, the LED bar graph will begin to light, with more lights indicating a more significant leak. On the 3rd red or 2nd yellow LED, an audible tone will begin to beep. The more LEDs that illuminate, the faster the beep. When the last red or yellow LED illuminates, the beep becomes a steady tone. The red LED lights indicate a helium or hydrogen leak. The yellow LED lights indicate a nitrogen, argon, or carbon dioxide leak. Remove the probe from the vicinity of the leak and allow the unit to return to zero. If a large amount of gas has entered the probe, it may take a few seconds for the instrument to clear itself. Do not attempt to zero the unit while it is clearing out the gas from the probe. This may cause the unit to malfunction. Place the probe near the leak again to confirm its location. The reference gas inlet (Figure 1) must not be restricted or the unit will not operate correctly. Similarly, the exhaust port allows the gas being tested to exit the leak detector and must remain unobstructed. The exhaust port is located in the probe docking station.

**CAUTION:**

This unit is designed to detect TRACE AMOUNTS of hydrogen arising from a small leak in a nonflammable environment, e.g., laboratory room air, etc. This unit is rated for use in a nonflammable atmosphere where the sample gas may become sufficiently high in concentration to become explosive.

**NOTE:**

To disable the audible beep during leak detection, depress and hold the zero set button for 2–3 seconds. After you hear a steady tone for 1 second, release the button; the beep function is disabled. To turn the beep function on again, depress and hold the zero set button. The beep function is always enabled at power up.

**NOTE:**

The leak detector will power down after 6 minutes of operation. This feature prevents excess battery discharge if the unit is accidentally left on.

## 7.0 Specifications

<b>Power Rating:</b>	5 volts DC, 1 A (AC adaptor & USB charging cable supplied)
<b>Battery Rating:</b>	12 hours normal operation
<b>Operating Temp. Range:</b>	50-98.6 °F (10-37 °C)
<b>Humidity Range:</b>	0-97%
<b>Warranty:</b>	1-year limited warranty
<b>Certifications:</b>	CE (EU, Korea, Japan, Australia) CSA/UL tested, not listed
<b>Compliance:</b>	WEEE, RoHS, China RoHS2 
<b>Altitude:</b>	Up to 2000m above mean sea level
<b>Pollution degree:</b>	2
<b>Indoors Use Only</b>	


## 8.0 Maintenance

Avoid spilling liquids onto the unit, or it may malfunction. If a liquid is spilled onto the unit, turn off the power immediately, remove heavy liquids with a dry towel, and let the unit sit until the liquid dries. Dust and debris can enter the probe tip of the leak detector and, over time, can clog the smallbore tubing inside the unit. To prevent this, clean the probe tip periodically. To clean the probe tip, unscrew the cap to expose the brush (Figures 2 and 3). Gently clean the probe using a small brush or your fingers to remove dust and debris, then replace the cap. Do not use liquids to clean the probe. Liquids can damage the leak detector if drawn in through the probe.

Information on where to have the unit sent for maintenance or service\* is listed at the end of this document.

**Figure 2—Cap unscrewed and partially removed**

**Figure 3—Cap removed, exposing probe tip brush for cleaning.**

## 9.0 Troubleshooting


Problem	Possible Cause	Suggested Solution
Sensitivity decreased	Probe clogged	Clean the probe tip to remove any debris
	Probe line punctured	Visually inspect probe line for holes*
Response decreased	Detector not zeroed	Re-zero detector
LED bar graph stays lit during operation	Detector re-zeroed before unit was purged out	Allow adequate time for detector to purge, then re-zero
	Reference gas inlet covered by hand or other object	Remove obstruction
Does not power up	Batteries need to be charged	Charge unit

\*Contact Thermo Fisher Scientific or your local representative for return instructions for servicing a damaged unit. Additional charges may apply if the warranty has expired or the unit is damaged due to misuse.

# 10.0 Technology

The leak detector measurement is based on thermal conductivity comparisons between the probe air and a reference air. The device employs a dual thermistor technology that measures the ratio of [probe]:[reference] heat exchange values and displays the results on an LED scale (Figure 4). Under ideal operating conditions, a ratio of 1:1 indicates identical air samples for both [probe] and [reference], and therefore, no leak is present.

**LEFT:** Dual analysis is achieved with heater elements positioned in separate flow chambers.


**RIGHT:** Probe and reference air streams are simultaneously monitored for thermal conductivity. Differences in air composition are indicated by differences in the heater element currents.

Because of slight differences in air temperature and/or humidity between the reference inlet (Figure 1) and the probe tip, a small response indicated by a single red or yellow LED light is generally insufficient to positively identify a gas leak. Small to moderate leaks are reliably indicated with 3 red or 2 yellow LED lights and an intermittent beep. Larger leaks are indicated with all red or yellow LED lights lit and a continuous audible alarm.

## 11.0 Interpretation of results

Figure 5 illustrates the leak detector's LED light response range. The greater the number of red or yellow LED lights lit correlates in general to the size of the leak. **NOTE:** The leak detector is not a quantitative device, rather it is designed to detect leaks in gas line connections commonly associated with laboratory equipment.

Gas	Minimum Detectable Leak Rate (atm cc / sec)	Indicating LED Light Color
Helium	$1.0 \times 10^{-5}$	Red
Hydrogen**	$1.0 \times 10^{-5}$	Red
Nitrogen	$1.4 \times 10^{-3}$	Yellow
Argon	$1.0 \times 10^{-4}$	Yellow
Carbon dioxide	$1.0 \times 10^{-4}$	Yellow

**Figure 5**–LED light response chart for the leak detector. A 1:1 ratio of I Probe : I Reference indicates no leak present. Red LED lights indicate the presence of helium and/or hydrogen. Yellow LED lights indicate the presence of nitrogen, argon, and/or carbon dioxide


**\*\*CAUTION:** This unit is designed to detect TRACE AMOUNTS of hydrogen and rising from a small leak in a non flammable environment, e.g., laboratory room air, etc. This unit is rated for use in a non flammable atmosphere where the sample gas may become sufficiently high in concentration to become explosive.

### Tip drift

Tip drift is the phenomenon when a false LED light response is registered as the unit is quickly turned or swept in dramatic arc movements. Tip drift is inherent to all dual thermistor leak detector technology and is based in large part on the asymmetry of the flow cells; shaking or tipping the unit influences the air flow profiles, which impacts the rates of heat exchange. If the device is functioning normally, the LED light signal will return to zero in 3–5 seconds after the unit is held still. In extreme cases, the unit may require another “zero” cycle before using. To avoid tip drift, be sure to hold the unit steady while making measurements.


## 12.0 Back label explanation

1. Product name.
2. Product catalog number.
3. Product serial number.
4. Passed California Energy Commission (CEC) Battery Charging System (BCS) testing.
5. CE mark: see Declaration of Conformity
6. KC mark - Korea Certification
7. Electrical parameters.
8. This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
  - 1) This device may not cause harmful interference.
  - 2) This device must accept any interference received, including interference that may cause undesired operation.
9. Manufacturer company name, address, and contact information.
10. AC adaptor polarity
11. RCM (Regulatory Compliance Mark) - Australia certification
12. China RoHS 2
13. Unit is WEEE compliant.
14. This Class A digital apparatus complies with Canadian ICES-003.
15. Units must be sent back to Thermo Fisher Scientific for service.
16. Unit is VCCI compliant
17. Unit is UKCA compliant

For the most up-to-date information, see our website [www.thermofisher.com](http://www.thermofisher.com)

## 13.0 Case specifications

## 14.0 Service

The Thermo Scientific leak detector carries a 1-year limited warranty from time of purchase. Please have the leak detector serial number available when calling Thermo Fisher Scientific with any concerns you may have. Additional charges may apply if the warranty is expired or the damage is due to misuse.

Expected battery lifetime is 2 years from time of purchase. Customers will need to return the unit to Thermo Fisher Scientific for battery replacement. At that time, preventative maintenance services can also be performed on the unit. A fee will be charged for servicing the unit.

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All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless
**For Research Use Only. Not for use in diagnostic procedures.** © 2017 Thermo Fisher Scientific Inc.