T.C.B.S. Cholera Medium

REF CM0333B, CM0333T

EN

Intended Use

T.C.B.S. Cholera Medium (CM0333B, CM0333T) devices are selective isolation medium for the isolation of *Vibrio* species from faecal and clinical samples where *Vibrio* species are suspected to be the cause of infection.

T.C.B.S. Cholera Medium (CM0333B, CM0333T) devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having *Vibrio* species infections

The device is for professional use only and is not intended for self-testing

Summary and Explanation

Vibrio species, members of the Vibrionaceae family, are a well-known genus capable of causing gastrointestinal diseases in humans, with cholera being the most well-known. Medically important *Vibrio* species include *Vibrio cholerae*, the causative agent of cholera, *Vibrio parahaemolyticus*, the causative agent of acute gastroenteritis and, *Vibrio vulnificus*, that causes septicemia. In general, *Vibrio* species are sea dwelling organisms and halophilic conditions are favourable for many species for their growth and survival. They are most frequently isolated from stool, blood, wound infections and vomit¹.

The cultural isolation and selection of *Vibrio* species is carried out using "Thiosulphate Citrate Bile salt Sucrose" (T.C.B.S.) agar¹. T.C.B.S. medium was developed by Kobayashi, Enomoto, Sakazaki and Kuwahara from the selective isolation agar initially developed by Nakanishi^{2,3}. As a modification, noted by Nakanishi and Kobayashi, the Oxoid T.C.B.S. medium specifically contains ox bile. T.C.B.S. agar, whilst inhibiting non-*Vibrio* species, promotes growth of other pathogenic *Vibrio* species after overnight incubation at 35°C and is suitable for the growth of *Vibrio* species including *Vibrio cholerae* and *Vibrio parahaemolyticus*⁴. Consequently, this agar is recommended by UK Standards for Microbiology Investigations UK SMI ID 19¹.

Principle of Method

T.C.B.S. Cholera Medium is designed to be highly selective, to meet the nutritional requirements of *Vibrio* spp. The medium contains yeast extract and peptone to provide the nitrogen and growth promoting vitamins. Sodium citrate, sodium thiosulfate and ox bile in the medium act as selective agents which provide an alkaline pH to inhibit Gram-positive organisms and suppress coliforms. As *Vibrio cholerae* is sensitive to acid environments, the pH of the medium is increased to enhance growth of this organism. When compared to other *Vibrio* species, most of which are halophilic, *Vibrio cholerae* is halotolerant. Therefore, the elevated sodium chloride concentration in this media favours growth of *Vibrio cholerae*. Sucrose is a fermentable carbohydrate, and sodium chloride stimulates growth. Sodium thiosulfate is a sulfur source and acts with ferric citrate as an indicator to detect hydrogen sulfide production. Bromothymol blue and thymol blue are pH indicators.

Typical Formula

| | <u>grams per litre</u> |
|-------------------------|------------------------|
| Yeast extract powder | 5.0 |
| Bacteriological peptone | 10.0 |
| Sodium thiosulphate | 10.0 |
| Tri-sodium citrate | 10.0 |
| Ox bile | 8.0 |
| Sucrose | 20.0 |
| Sodium chloride | 10.0 |
| Iron (III) citrate | 1.0 |
| Bromothymol blue | 0.04 |
| Thymol blue | 0.04 |
| Agar | 14.0 |

Materials Provided

CM0333B: 500g of T.C.B.S. Cholera Medium

CM0333T: 5kg of T.C.B.S. Cholera Medium

500g of T.C.B.S. Cholera Medium powder yields approximately 5.7L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.

- Store away from light.
 - Allow reconstituted product to equilibrate to room temperature before use.
- Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any applicable federal, state and local regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between uses to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following recommended local guidelines, such as the UK Standards for Microbiology (UK SMI ID) 19 and Q5.

Procedure

Suspend 88g in 1 litre of distilled water. With frequent agitation, bring to the boil to dissolve completely. Cool to 50°C. Mix well and pour into sterile Petri dishes. DO NOT AUTOCLAVE. DO NOT OVERHEAT.

Quality Control

It is the responsibility of the user to perform Quality Control testing, taking into account the intended use of the medium, and in accordance with any applicable local regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18h @ 37° ± 2°C aerobic

| | |
|---|--|
| Positive Controls | |
| Inoculum level 10-100 cfu Colony count is ≥ 50% of the control medium count. | |
| <i>Vibrio cholerae</i> NCTC 11348 | 0.75 – 2 mm Yellow colonies |
| Negative Controls | |
| Inoculum level 10 – 1000 cfu | |
| <i>Pseudomonas aeruginosa</i> ATCC® 27853™ | No growth or pinpoint green colonies |
| Inoculum level: 10 ⁴ – 10 ⁶ cfu | |
| <i>Enterococcus faecalis</i> ATCC® 19433™ | No growth or pinpoint – 1 mm colourless or yellow colonies |

Incubation Conditions: 21-27 h @ 37° ± 2°C aerobic

Tested in accordance with ISO 11133:2014:

| | |
|---|-----------------------------|
| Positive Controls | |
| Inoculum level 50-120 cfu Colony count is ≥ 50% of the control medium count. | |
| <i>Vibrio furnissii</i> NCTC 11218 | 0.75 – 2 mm yellow colonies |

| | |
|---|---------------------------|
| <i>Vibrio parahaemolyticus</i> NCTC 10885 | 0.75 – 2mm green colonies |
| Negative Controls Inoculum level: 10 ⁴ – 10 ⁶ cfu | |
| <i>Escherichia coli</i> ATCC® 8739™ | No growth |
| <i>Escherichia coli</i> ATCC® 25922™ | No growth |
| <i>Escherichia coli</i> ATCC® 11775™ | No growth |

Limitations

Identifications are presumptive and should be confirmed using appropriate biochemical and serological methods. Since the nutritional requirements of organisms vary, some strains may grow poorly on this medium. On T.C.B.S. medium, *V. parahaemolyticus* may resemble *Aeromonas hydrophila*, *Plesiomonas shigelloides* and *Pseudomonas* species. Sucrose-fermenting *Proteus* species produce yellow colonies which may resemble those of *Vibrio*. A few strains of *V. cholerae* appear green or colourless on T.C.B.S. medium due to delayed sucrose fermentation. T.C.B.S. medium, is unsatisfactory for oxidase testing of *Vibrio* spp.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of *Vibrio* species is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device. The precision of T.C.B.S. Cholera Medium (CM0333B, CM0333T) devices was demonstrated by an overall pass rate of 100% over 20 months of testing (12-Jul-2021- 27-Apr-2022). This shows that the performance is reproducible.


T.C.B.S. Cholera Medium (CM0333B, CM0333T) devices are tested in-house as part of the QC process since the products were first developed in 1996 and 2004 for variants. For target organisms, when using 10 – 100 cfu inoculum of *Vibrio cholerae* NCTC 11348 incubating the device at 35-39°C for 21-27 hours, the user can recover organisms with colony size and morphology as listed in this document. For non-target organisms, when using 10-1000 cfu of *Pseudomonas aeruginosa* ATCC® 27853™, 10⁴ to 10⁶ cfu of *Enterococcus faecalis* ATCC® 19433™ and incubating the device at 35-39°C for 21-27 hours, the user can recover organisms with colony size and morphology as listed in this document.

Bibliography

1. Public Health England. 2015. 'Identification of *Vibrio* and *Aeromonas* species. UK Standards for Microbiology Investigations. SMI ID19. Issue no: 3 | Issue date: 14.04.15 | <https://www.gov.uk/government/publications/smi-id-19-identification-of-vibrio-species>
2. Kobayashi T., Enomoto S., Sakazaki R. and Kuwahara S. (1963) *Jap. J. Bacteriol.* 18. 10-11, 387-311.
3. Nakanishi Y. (1963) *Modern Media* 9. 246.
4. Furniss A. L., Lee J. V. and Donovan T. J. (1978) *The Vibrios. PHLs Monograph No. 11.*

Symbol Legend

| Symbol | Definition |
|---|---|
|  | Catalogue number |
|  | <i>In Vitro</i> Diagnostic Medical Device |
|  | Batch code |
|  | Temperature limit |
|  | Use-by date |
|  | Keep away from sunlight |
|  | Do not re-use |
|  | Consult instructions for use or consult electronic instructions for use |

| | |
|---|--|
|  | Contains sufficient for <n> tests |
|  | Do not use if packaging damaged and consult instructions for use |
|  | Manufacturer |
|  | Authorized representative in the European Community/ European Union |
|  | European Conformity Assessment |
|  | UK Conformity Assessment |
|  | Unique device identifier |
|  | Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union |
|  | Made in the United Kingdom |

©2022 Thermo Fisher Scientific Inc. All rights reserved.
 ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
 NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.
 All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

| Version | Date of issue and modifications introduced |
|---------|--|
| 2.0 | 2023-12-12 |

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©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien