Thermo Scientific™ Oxoid™ Single Cartridge ONPG Diagnostic Discs DD13

REF DD0013T

Intended Use

The Thermo Scientific™ Oxoid™ ONPG Diagnostic Discs (DD0013T) are impregnated with ONPG (phosphate buffered O-nitrophenyl-beta-D-galactopyranoside) and used in a diagnostic workflow to aid microbiology professionals by indicating the presence of beta-galactosidase via a yellow colour change. To be used with a pure, agar grown culture, the device provides a qualitative result.

The device is for professional use only, is not automated and nor is it a companion diagnostic.

Principle of the Test

The O-nitrophenyl-β-D-galactopyranoside (ONPG) discs are used for the detection of beta-galactosidase activity in specific organisms, allowing the differentiation between lactose fermenting and non-fermenting bacteria. Lactose fermentation requires the presence of two enzymes: permease, that allows the lactose to enter the bacterial cell wall and β-galactosidase and beta-galactosidase that allows the bacteria to metabolise lactose, via hydrolysis of the molecule in glucose plus galactose. True lactose non-fermenters do not possess either of these enzymes. Late lactose fermenting organisms do not have permease but do possess β-galactosidase. ONPG is an artificial substrate structurally similar to lactose. If β-galactosidase is present ONPG, that is normally colourless, is targeted by the enzyme and split into galactose and O-nitrophenol, a yellow compound; this colour change indicates a positive test¹

Components

ONPG Diagnostic Discs consist of 6mm diameter paper discs impregnated with a specific antimicrobial concentration. The discs are marked on both sides with an alpha-numeric code identifying the antimicrobial agent and concentration. ONPG Diagnostic Discs are supplied in cartridges of 50 discs. The cartridge is individually packed in a foil-sealed blister pack with a desiccant.

Materials Provided

- 1 x cartridge with spring cap and plunger
- Desiccant tablet x1
- Foil individually sealing each cartridge with its desiccant
- Single Cartridge ONPG discs x 50

ONPG diagnostic discs contain the active ingredient (O-nitrophenyl-beta-D-galactopyranoside which is an Antimicrobial agent, typically white to light yellow when in dry form.

Materials Required but Not Supplied

Agar plates with appropriate media and test organism, 0.85% sodium chloride (physiological saline), inoculum suspension medium, pipettes, sterile loops and swabs, sterile forceps, McFarland turbidity standards, incubator, modified atmosphere environments; Quality control (QC) strains, and interpretative criteria for local standard methods.

Warnings and Precautions

- Not intended for susceptibility testing.
- This product is for *in vitro* diagnostic use only.
- Follow directions for use.
- Observe aseptic techniques and established precautions against all microbiological hazards throughout all procedures.
- Cultures, containers and other contaminated materials must be sterilized after use in accordance with guidelines for the handling and disposal of biohazardous waste.
- Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (found on www.thermofisher.com).
- Inspect the product packaging before first use. Do not use product if there is any visible damage to the foil seal.
- Do not use beyond the stated expiry date.
- Once the cartridge is open, ensure it is stored in an opaque desiccated environment to prevent degradation.
- If the discs do not produce the expected reaction with recommended control microorganisms, review the entire procedure.
- In the event of malfunction do not use the device.
- The device does not contain any carcinogenic, mutagenic or toxic to reproduction, endocrine disrupting substances or materials likely to cause sensitisation or an allergic reaction under normal use.

Storage and Handling

Unopened cartridges must be stored at 2 to 8°C until required. Allow cartridges to reach room temperature before removing them from the packaging to minimize condensation. The expiry date is valid only for the unopened blister packs stored under correct conditions. Once opened, cartridges should be stored within a dispenser (sold separately) in the container provided (with an unsaturated (orange) desiccant), or other suitable opaque airtight container with a desiccant to protect the discs from light and moisture.

Method

1. Place one disc into a sterile tube.
2. Add 0.1 ml of sterile 0.85% sodium chloride (physiological saline).
3. Pick the colony under test with a sterile loop and emulsify it in the tube containing the disc and physiological saline.
4. Incubate at 35°C.
5. Examine at hourly intervals for up to 6 hours to detect active lactose fermenters.
6. Organisms that are negative after 6 hours should be incubated for up to 24 hours to detect the late lactose fermenters.

Interpretation

Colourless	ONPG Negative
Yellow	ONPG Positive

The reaction speed depends on the size of the inoculum.

User Quality Control Procedures

The use of positive control organisms such as *Escherichia coli* ATCC® 25922™ and *Citrobacter freundii* ATCC® 8090™ is recommended.

Shigella flexneri ATCC® 12022™ and *Salmonella enteritidis* ATCC® 13076™ may be used as a negative control organism. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

Limitations

For professional use only.

This product is for *in vitro* diagnostic use only. It is recommended that the tests are performed on colonies from pure culture.

Failure to observe instructions for use may lead to erroneous results.

Analytical Performance Characteristics

ONPG Diagnostic Discs (DD0013T) QC performance was 96.61% from 2017 to 2022.

The trueness of measurement for ONPG Diagnostic Discs (DD0013T) was demonstrated using historical routine QC performance data.

The precision of the ONPG Diagnostic Discs (DD0013T) was demonstrated by the historical performance achieved through routine QC testing assessed by the manufacturer. This shows that the performance is reproducible.

In conclusion, the submitted evidence of the Analytical Performance of the ONPG Diagnostic Discs (DD0013T) reinforces evidence that the devices meet the requirements required to provide the customer with safe and effective devices to be used in a diagnostic workflow to aid clinicians in treatment options.

Clinical Performance Characteristics

ONPG Diagnostic Discs (DD0013T) device has been evaluated through scientific peer-reviewed literature and published experience gained by routine diagnostic testing, which compared the number of performance-related quality notifications to the total sales within the same time period and demonstrated the performance of the devices in a broad setting, including customers from a clinical setting.

Although a literature review was attempted, no relevant articles which informed on the clinical performance of ONPG Diagnostic Discs (DD0013T) device were identified. As part of the PMS procedures, if any new risks are identified or there is a change in the state of art of these devices, then this will be captured through future literature searches.

Published experience gained by routine diagnostic testing provides evidence of in-use performance and supports ONPG Diagnostic Discs (DD0013T) as a suitable device for the detection of beta-galactosidase activity from specific clinical samples for a diagnostic purpose. There were no reported performance complaints (0%) per unit sold for the ONPG Diagnostic Discs (DD0013T) from 2018 to 2020. This informs on the adequate performance of the product and demonstrates the safety and performance of the device in a clinical setting.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.


References

1. PHE. 2018. "UK Standards for Microbiology Investigations ONPG (β-Galactosidase) Test ONPG (β-Galactosidase) Test- Standards-for-Microbiology-Investigations-Smi-Quality-and-Consistency-in-Clinical-Laboratories PHE Publications Gateway Number: 2018381 UK Standards for M.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limitation (storage temp.)
	Lot / Batch Number
	Catalog Number
	Do not re-use
	Protect from light
	Protect from moisture
	Consult Instructions for Use or consult Electronic instructions for use
	Contains sufficient for <n> tests
	Expiry Date
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Do not use if package is damaged. Consult instructions for use.
	Keep away from sunlight
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union

The ATCC Licensed Derivative® Emblem, the ATCC Licensed Derivative® word mark, and the ATCC catalogue marks are trademarks of ATCC. Thermo Fisher Scientific Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

© 2025 Thermo Fisher Scientific Inc. All rights reserved.
ATCC® is a trademark of American Type Culture Collection.
All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Please contact your local sales representative for all enquiries.

For technical assistance please contact your local distributor.

Revision Information

Version	Date of issue and modifications introduced
3.0	11 July 2025 Typographical error corrected

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©2025 Thermo Fisher Scientific Inc. Všetky práva vyhradené.