	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Mupirocin Discs 200 µg (MUP200)

REF CT0523B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Mupirocin Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent (Product name). The discs are labelled on both sides with details of the antimicrobial (MUP) and amount present (µg): MUP200 (200 µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Mupirocin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Mupirocin AST Discs are used in the semi-quantitative agar diffusion test method for *in vitro* susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which mupirocin has been shown to be active both clinically and *in vitro*. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling, and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with published breakpoints according to current FDA literature:

Gram-Positive

- N/A

Gram-Negative

- N/A

Species with published breakpoints according to current CLSI M100 literature:

Gram-Positive

- *Staphylococcus aureus*

Gram-Negative

- N/A

Species with published breakpoints according to current EUCAST literature:

Gram-Positive

- *Staphylococcus aureus*

Gram-Negative

- N/A

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Mupirocin AST Discs may be used in the semi-quantitative agar diffusion test method *for in vitro* susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of values assigned to calibrators and control materials intended to establish or verify trueness of a method for Mupirocin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Mupirocin AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Mupirocin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Reagent	Description of Function
Mupirocin	Mupirocin antimicrobial agent. Mupirocin, formerly termed pseudomonic acid, is a novel antibacterial agent with a unique chemical structure and mode of action, produced by fermentation using the organism <i>Pseudomonas fluorescens</i> . It primarily works by inhibiting bacterial protein synthesis by inhibiting the activity of bacterial isoleucyl-tRNA synthase.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Mupirocin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Mupirocin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided with an unsaturated (orange) desiccant or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 1. Raw data of zone sizes, in accordance with EUCAST methodology, taken from the 3 most recent batches of Mupirocin (200 µg) MUP200 Antimicrobial Susceptibility Test Discs (CT0523B).

EUCAST Methodology										
Product Batch	Organism	ATCC® number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3345510	<i>Staphylococcus aureus</i>	29213	MHA	31	37	34	31	31	31	-3 ± 0 (CV=0)
3331954	<i>Staphylococcus aureus</i>	29213	MHA	31	37	34	31	31	31	-3 ± 0 (CV=0)
3380463	<i>Staphylococcus aureus</i>	29213	MHA	31	37	34	32	32	32	-2 ± 0 (CV=0)

Clinical Performance Characteristics

A study by Udo et al. (2006) aimed to analyse the resistance trends and genetic location of their resistance determinants in methicillin-resistant *Staphylococcus aureus* (MRSA) isolates. A total of 5,644 consecutive MRSA isolates were obtained for typing and antibiotic susceptibility testing. Isolates were tested for susceptibility to a range of antimicrobial discs including Oxoid Mupirocin (200 µg) MUP200 AST Discs. The Results showed, in isolates obtained in 1996, mupirocin resistance was 6.0%, this increased from 17.0% in 1997 to 28.0% in 1999. Resistance then fell to 22.0% in 2000 followed by a rapid decrease in resistance in 2001 to 2.7%, it then stayed consistently low in the years 2002, 2003 and 2004 with rates of 3.0%, 1.0% and 4.7%, respectively. Oxoid Mupirocin AST Discs proved effective for the in vitro clinical analysis of MRSA isolates given that they displayed good susceptibility rates in comparison to other agents, detected emerging resistance trends in the study area, and aided in the surveillance of treatments for these infections.

A 2009 study by Marais et al. described antibiotic susceptibility patterns of MRSA isolates collected in South Africa. A total of 248 non-duplicated, laboratory confirmed, *mecA*-positive MRSA isolates were collected for antibiotic susceptibility testing. Isolates were tested for sensitivity to a range of antibiotic discs including Oxoid Mupirocin (200 µg) MUP200 AST Discs. Results showed resistance to mupirocin was observed in nine of the 248 isolates (3.6%). The Oxoid AST Discs proved effective for in vitro clinical analysis of MRSA isolates given that they displayed good susceptibility rates in comparison to other agents.

Lorio et al. (2012) aimed to evaluate the characteristics of 35 methicillin-resistant *Staphylococcus epidermidis* (MRSE) isolates such as staphylococcal chromosomal cassette *mec* (SCC*mec*) types, antimicrobial susceptibility profiles, biofilm formation genes, and multilocus sequence types (MLST). Results showed a total of 35 isolates were confirmed as MRSE, with isolates in the main PFGE type displaying resistance to seven of the 13 non-beta lactam antimicrobial agents, while isolates of the sporadic genotypes were resistant to only four agents. The resistance rates to mupirocin were higher in the predominant genotypes than the sporadic one. Resistance to mupirocin was lower than other antimicrobials tested such as gentamicin, erythromycin, sulfamethoxazole-trimethoprim, and ciprofloxacin. The elevated resistance rates observed in this study were most likely due to the isolates having gene characteristics for high resistance to antimicrobials. This study demonstrated that Oxoid Mupirocin AST Discs were effective for in vitro clinical analysis of Gram-positive isolates, given that they displayed acceptable susceptibility by the organisms studied and detected emerging resistance trends in the study area.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Udo, Edet E., Noura Al-Sweih, Eiman Mokaddas, Molly Johny, Rita Dhar, Huda H. Gomaa, Inaam Al-Obaid, and Vincent O. Rotimi. 2006. "Antibacterial Resistance and Their Genetic Location in MRSA Isolated in Kuwait Hospitals, 1994-2004." *BMC Infectious Diseases* 6: 1–8. <https://doi.org/10.1186/1471-2334-6-168>.
5. Marais, Else, N. Aithma, O. Perovic, W. F. Oosthuysen, E. Musenge, and A. G. Dusé. 2009. "Antimicrobial Susceptibility of Methicillin-Resistant *Staphylococcus Aureus* Isolates from South Africa." *South African Medical Journal* 99 (3): 170–73. <https://doi.org/10.7196/SAMJ.2929>.
6. Iorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Goncalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant *Staphylococcus Epidermidis* ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use

	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2022-10-27. New document

Oxoid Limited, Wade Road, Basingstoke, Hampshire RG24 8PW, UK

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Literatura

- Clinical & Laboratory Standards Institute (CLSI). Proučite najnovije izdanje:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Proučite najnoviju verziju:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Udo, Edet E., Noura Al-Sweih, Eiman Mokaddas, Molly Johny, Rita Dhar, Huda H. Gomaa, Inaam Al-Obaid, and Vincent O. Rotimi. 2006. "Antibacterial Resistance and Their Genetic Location in MRSA Isolated in Kuwait Hospitals, 1994-2004." *BMC Infectious Diseases* 6: 1–8. <https://doi.org/10.1186/1471-2334-6-168>.
- Marais, Else, N. Aithma, O. Perovic, W. F. Oosthuysen, E. Musenge, and A. G. Dusé. 2009. "Antimicrobial Susceptibility of Methicillin-Resistant *Staphylococcus Aureus* Isolates from South Africa." *South African Medical Journal* 99 (3): 170–73. <https://doi.org/10.7196/SAMJ.2929>.
- Lorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Goncalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant *Staphylococcus Epidermidis* ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

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Referencer

- Clinical & Laboratory Standards Institute (CLSI). Se den seneste udgave:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Se den seneste version:
 - EUCAST -agar-diffusionsmetode
 - Læsevejledning, EUCAST -agar-diffusionsmetode til antimikrobiel modtagelighedstestning
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Udo, Edet E., Noura Al-Sweih, Eiman Mokaddas, Molly Johny, Rita Dhar, Huda H. Gomaa, Inaam Al-Obaid, and Vincent O. Rotimi. 2006. "Antibacterial Resistance and Their Genetic Location in MRSA Isolated in Kuwait Hospitals, 1994-2004." *BMC Infectious Diseases* 6: 1–8. <https://doi.org/10.1186/1471-2334-6-168>.
- Marais, Else, N. Aithma, O. Perovic, W. F. Oosthuysen, E. Musenge, and A. G. Dusé. 2009. "Antimicrobial Susceptibility of Methicillin-Resistant *Staphylococcus Aureus* Isolates from South Africa." *South African Medical Journal* 99 (3): 170–73. <https://doi.org/10.7196/SAMJ.2929>.
- Lorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Goncalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant *Staphylococcus Epidermidis* ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

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Bibliografia

1. Clinical & Laboratory Standards Institute (CLSI). Consultare l'ultima edizione:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. Comitato europeo sui test di suscettibilità antimicrobica (European Committee on Antimicrobial Susceptibility Testing, EUCAST). Consultare l'ultima versione:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Udo, Edet E., Noura Al-Sweih, Eiman Mokaddas, Molly Johny, Rita Dhar, Huda H. Gomaa, Inaam Al-Obaid e Vincent O. Rotimi. 2006. "Antibacterial Resistance and Their Genetic Location in MRSA Isolated in Kuwait Hospitals, 1994-2004." *BMC Infectious Diseases* 6: 1–8. <https://doi.org/10.1186/1471-2334-6-168>.
5. Marais, Else, N. Aithma, O. Perovic, W. F. Oosthuysen, E. Musenge e A. G. Dusé. 2009. "Antimicrobial Susceptibility of Methicillin-Resistant Staphylococcus Aureus Isolates from South Africa." *South African Medical Journal* 99 (3): 170–73. <https://doi.org/10.7196/SAMJ.2929>.
6. Lorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Goncalves Neves, Regina Maria Cavalcanti Pilotto Domingues e Kátia Regina Netto dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant Staphylococcus Epidermidis ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Regno Unito

	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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
Referências

1. Clinical & Laboratory Standards Institute (CLSI, Instituto de padrões clínicos e laboratoriais dos EUA). Consultar a edição mais recente:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST, Comité Europeu para o Teste à Suscetibilidade Antimicrobiana). Consultar a versão mais recente:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA, Agência dos alimentos e dos medicamentos dos EUA). CFR Title 21, Volume 5, Part 460 (2005).
4. Udo, Edet E., Noura Al-Sweih, Eiman Mokaddas, Molly Johnny, Rita Dhar, Huda H. Gomaa, Inaam Al-Obaid, and Vincent O. Rotimi. 2006. "Antibacterial Resistance and Their Genetic Location in MRSA Isolated in Kuwait Hospitals, 1994-2004." *BMC Infectious Diseases* 6: 1–8. <https://doi.org/10.1186/1471-2334-6-168>.
5. Marais, Else, N. Aithma, O. Perovic, W. F. Oosthuysen, E. Musenge, and A. G. Dusé. 2009. "Antimicrobial Susceptibility of Methicillin-Resistant *Staphylococcus Aureus* Isolates from South Africa." *South African Medical Journal* 99 (3): 170–73. <https://doi.org/10.7196/SAMJ.2929>.
6. Iorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Goncalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. "Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant *Staphylococcus Epidermidis* ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil." *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

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
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Reino Unido

	https://www.thermofisher.com	
	UE	+800 135 79 135
	SUA	1 855 236 0910
	CA	1 855 805 8539
	Restul lumii	+31 20 794 7071

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
Odkazy

1. Clinical & Laboratory Standards Institute (CLSI). Consult latest edition:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Hlava 21, Zväzok 5, Časť 460 (2005).
4. Udo, Edet E., Noura Al-Sweih, Eiman Mokaddas, Molly Johny, Rita Dhar, Huda H. Gomaa, Inaam Al-Obaid, a Vincent O. Rotimi. 2006. „Antibacterial Resistance and Their Genetic Location in MRSA Isolated in Kuwait Hospitals, 1994-2004.“ *BMC Infectious Diseases* 6: 1–8. <https://doi.org/10.1186/1471-2334-6-168>.
5. Marais, Else, N. Aithma, O. Perovic, W. F. Oosthuysen, E. Musenge, a A. G. Dusé. 2009. „Antimicrobial Susceptibility of Methicillin-Resistant Staphylococcus Aureus Isolates from South Africa.“ *South African Medical Journal* 99 (3): 170–73. <https://doi.org/10.7196/SAMJ.2929>.
6. Iorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Goncalves Neves, Regina Maria Cavalcanti Pilotto Domingues, and Ktia Regina Netto Dos Santos. 2012. „Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant Staphylococcus Epidermidis ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil.“ *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

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Referencias

1. Clinical & Laboratory Standards Institute (CLSI). Consultar la última edición:
 - a. M02, Performance Standards for Antimicrobial Disk Susceptibility Tests.
 - b. M100, Performance Standards for Antimicrobial Susceptibility Testing.
 - c. M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consultar la última versión:
 - a. EUCAST Disk Diffusion Method
 - b. «Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing».
3. Food and Drug Administration (FDA). Título 21 del CFR, volumen 5, parte 460 (2005).
4. Udo, Edet E., Noura Al-Sweih, Eiman Mokaddas, Molly Johny, Rita Dhar, Huda H. Gomaa, Inaam Al-Obaid y Vincent O. Rotimi. 2006. «Antibacterial Resistance and Their Genetic Location in MRSA Isolated in Kuwait Hospitals, 1994-2004». *BMC Infectious Diseases* 6: 1–8. <https://doi.org/10.1186/1471-2334-6-168>.
5. Marais, Else, N. Aithma, O. Perovic, W. F. Oosthuysen, E. Musenge y A. G. Dusé. 2009. «Antimicrobial Susceptibility of Methicillin-Resistant *Staphylococcus Aureus* Isolates from South Africa». *South African Medical Journal* 99 (3): 170–73. <https://doi.org/10.7196/SAMJ.2929>.
6. Lorio, Natalia Lopes Pontes, Roberta Ferreira Caboclo, Milena Borgo Azevedo, Ariane Guimares Barcellos, Felipe Piedade Gonçalves Neves, Regina Maria Cavalcanti Pilotto Domingues y Ktia Regina Netto Dos Santos. 2012. «Characteristics Related to Antimicrobial Resistance and Biofilm Formation of Widespread Methicillin-Resistant *Staphylococcus Epidermidis* ST2 and ST23 Lineages in Rio de Janeiro Hospitals, Brazil». *Diagnostic Microbiology and Infectious Disease* 72 (1): 32–40. <https://doi.org/10.1016/j.diagmicrobio.2011.09.017>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien