Blood Agar Base No. 2 (Dehydrated)

REF CM0271B, CM0271K, CM0271R & CM0271T

EN

Intended Use

Blood Agar Base No. 2 (Dehydrated) (CM0271B/K/R/T) devices are nonselective, general purpose media which may be enriched with blood or supplements for the isolation and cultivation of fastidious pathogens and other microorganisms with clearly visible haemolysis from clinical samples (e.g., wounds, throat, genital, nose, groins etc.). Blood Agar Base No. 2 (Dehydrated) (CM0271B/K/R/T) devices are intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The devices are for professional use only, are not automated and nor are they companion diagnostics.

Summary and Explanation

A fastidious organism is any organism that has complex or particular nutritional requirements. There are many different types of fastidious organisms which exhibit different requirements to ensure successful growth. Some of the main fastidious humans pathogens are *Haemophilus influenzae*, *Staphylococcus aureus*, *Streptococcus pyogenes* and *Streptococcus pneumoniae*.

The addition of blood to culture media provides means of identification, as the ability to haemolyse blood can be used to classify certain microorganisms^{1,2}. The presence and pattern of haemolysis, the size of the zone of haemolysis, and colony morphology can be used together to identify isolates cultured using blood agar^{1,2}.

Historically, blood agar consisted of casein hydrolysate or meat infusion media, with the former producing large colonies and the latter producing clearly visible haemolytic reactions³. In 1966, Ellner et al. combined the two blood agar bases to improve identification of *S. pyogenes*, *Diplococcus pneumoniae* and Viridans streptococci³. The resulting optimal formulation, known as CBA, is now widely used as a general-purpose culture media and provides large colonies with clearly defined haemolytic reactions and colony morphology³. Consequently, blood agar is recommended for the primary culture of a number of human pathogens, including some fastidious microorganisms^{1,2,4}.

Principle of Method

Blood Agar Base No. 2 devices contain proteose peptone, liver digest and yeast extract which supply nitrogen, carbon, amino acids and vitamins. Sodium chloride is added to maintain osmotic equilibrium and agar is the solidifying agent. Addition of horse blood or sheep blood adds additional growth factors and enables the determination of haemolytic patterns.

Typical Formula

	<u>grams per litre</u>
Proteose peptone	15.0
Liver digest	2.5
Yeast extract	5.0
Sodium chloride	5.0
Agar	12.0

Materials Provided

CM0271B: 500g of Blood Agar Base No. 2 (Dehydrated)
CM0271K: 25 Kg of Blood Agar Base No. 2 (Dehydrated)
CM0271R: 2.5 Kg of Blood Agar Base No. 2 (Dehydrated)
CM0271T: 5 Kg of Blood Agar Base No. 2 (Dehydrated)

500g of dehydrated Blood Agar Base No. 2 yields approximately 12.5L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Selective supplements (SR0050B/C/D/E/F/R, SR051B/C/D/E)
- Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com). Serious Incidents
Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 4, ID 7, ID 12 and Q 5.

Procedure

Suspend 40g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilise by autoclaving at 121°C for 15 minutes. Cool to 50°C. Mix well and pour into sterile containers. For blood agar, enrich with 7% v/v sterile defibrinated blood.

Interpretation

Once the media is reconstituted, cream colonies indicates *Staphylococcus aureus*. Irregular straw or straw/green colonies indicates *Pseudomonas aeruginosa*. Straw colonies indicates *Escherichia coli*.
When enriched with 7% v/v horse blood pale straw colonies indicates *Streptococcus pyogenes*. Grey/green colonies indicates *Streptococcus pneumoniae*. Colourless colonies indicates *Haemophilus influenzae*

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 h @ 37°C

Positive Controls	
Inoculum level: 10 – 100 cfu Colony count is \geq 70% of the control medium count.	
<i>Staphylococcus aureus</i> ATCC® 25923™	1 - 2mm cream colonies
<i>Staphylococcus aureus</i> ATCC® 9144™	0.5 - 1.5 mm cream colonies
<i>Pseudomonas aeruginosa</i> ATCC® 27853™	1 – 4 mm irregular straw or straw/ green colonies
<i>Escherichia coli</i> ATCC® 25922™	1 – 2 mm straw colonies
Enriched with 7% v/v Horse Blood	
<i>Streptococcus pyogenes</i> ATCC® 19615™	0.25 – 1 mm pale straw colonies, β haemolysis
<i>Streptococcus pneumoniae</i> ATCC® 6303™	1 – 4 mm grey/green colonies α haemolysis
<i>Streptococcus pneumoniae</i> ATCC® 6305™	Pinpoint – 0.25 mm grey/green colonies α haemolysis
<i>Haemophilus influenzae</i> ATCC® 19418™	Pinpoint – 0.25 mm colourless colonies

Limitations

Due to variation in nutritional requirements some strains may grow poorly or fail to grow on this medium. The haemolytic reaction will be strain dependent.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of fastidious pathogens and other microorganisms with clearly visible haemolysis is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the devices, which must meet the defined acceptance criteria. The precision of Blood Agar Base No. 2 (Dehydrated) (CM0271B/K/R/T) was demonstrated by an overall pass rate of 100% obtained for the product over 1.5 years of testing (December 2020 - May 2022; 10 batches). This shows that the performance is reproducible.

Blood Agar Base No. 2 (Dehydrated) (CM0271B/K/R/T) devices are tested in-house as part of the QC process since the products were launched in 1997 (CM0271B/K/R/T). For target organisms, when using 10-100 cfu inoculum of *Staphylococcus aureus*, *Pseudomonas aeruginosa* and *Escherichia coli* and incubating the device as a plain plate at 37°C ± 2°C for 18 hours, the user can recover organisms with colony size and morphology that meets the defined acceptance criteria. For target organisms, when using 10-


100 cfu inoculum of *Streptococcus pyogenes*, *Streptococcus pneumoniae* and *Haemophilus influenzae* and incubating the device enriched with 7% v/v horse blood at 37°C ± 2°C for 18 hours, the user can recover organisms with colony size and morphology as listed in this document.

Bibliography

1. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
2. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms.' UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
3. Ellner, P. D., C. J. Stoessel, E. Drakeford, and F. Vasi. 1966. 'A New Culture Medium for Medical Bacteriology.' *The American Journal of Clinical Pathology* 45 (4): 502-4. https://doi.org/10.1093/ajcp/45.4_ts.502
4. Public Health England. 2021b. 'Identification of Haemophilus species and the HACEK group of organisms.' UK Standards for Microbiology Investigations. ID 12 Issue 4. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use

	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, United Kingdom

For technical assistance please contact your local distributor.

Revision information

Version	Modifications Introduced
2.0	2023-11-29

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1. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
2. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms.' UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
3. Ellner, P. D., C. J. Stoessel, E. Drakeford, and F. Vasi. 1966. 'A New Culture Medium for Medical Bacteriology.' *The American Journal of Clinical Pathology* 45 (4): 502–4. https://doi.org/10.1093/ajcp/45.4_ts.502
4. Public Health England. 2021b. 'Identification of Haemophilus species and the HACEK group of organisms.' UK Standards for Microbiology Investigations. ID 12 Issue 4. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>.

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Literatur

1. Public Health England. 2020. ‚Identification of Staphylococcus species, Micrococcus species and Rothia species.‘ UK Standards for Microbiology Investigations. ID 7 Issue 4. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
2. Public Health England. 2021a. ‚Identification of Streptococcus species, Enterococcus species and morphologically similar organisms.‘ UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
3. Ellner, P. D., C. J. Stoessel, E. Drakeford, and F. Vasi. 1966. ‚A New Culture Medium for Medical Bacteriology.‘ *The American Journal of Clinical Pathology* 45 (4): 502–4. https://doi.org/10.1093/ajcp/45.4_ts.502
4. Public Health England. 2021b. ‚Identification of Haemophilus species and the HACEK group of organisms.‘ UK Standards for Microbiology Investigations. ID 12 Issue 4. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>.

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Bibliografia

1. Public Health England. 2020. "Identification of Staphylococcus species, Micrococcus species and Rothia species". UK Standards for Microbiology Investigations. ID 7 Issue 4. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
2. Public Health England. 2021a. "Identification of Streptococcus species, Enterococcus species and morphologically similar organisms". UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
3. Ellner, P. D., C. J. Stoessel, E. Drakeford, and F. Vasi. 1966. 'A New Culture Medium for Medical Bacteriology.' *The American Journal of Clinical Pathology* 45 (4): 502-4. https://doi.org/10.1093/ajcp/45.4_ts.502
4. Public Health England. 2021b. "Identification of Haemophilus species and the HACEK group of organisms". UK Standards for Microbiology Investigations. ID 12 Issue 4. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>.

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Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, Regno Unito

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Literatūra

1. Public Health England. 2020. „Identification of Staphylococcus species, Micrococcus species and Rothia species.“ UK Standards for Microbiology Investigations. ID 7 Issue 4. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
2. Public Health England. 2021a. „Identification of Streptococcus species, Enterococcus species and morphologically similar organisms.“ UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
3. Ellner, P. D., C. J. Stoessel, E. Drakeford, and F. Vasi. 1966. 'A New Culture Medium for Medical Bacteriology.' *The American Journal of Clinical Pathology* 45 (4): 502–4. https://doi.org/10.1093/ajcp/45.4_ts.502
4. Public Health England. 2021b. „Identification of Haemophilus species and the HACEK group of organisms.“ UK Standards for Microbiology Investigations. ID 12 Issue 4. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>.

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Bibliografi

1. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Utgave 4. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
2. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms.' UK Standards for Microbiology Investigations. ID 4 Utgave 4. <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
3. Ellner, P. D., C. J. Stoessel, E. Drakeford, og F. Vasi. 1966. 'A New Culture Medium for Medical Bacteriology.' *The American Journal of Clinical Pathology* 45 (4): 502–4. https://doi.org/10.1093/ajcp/45.4_ts.502
4. Public Health England. 2021b. 'Identification of Haemophilus species and the HACEK group of organisms.' UK Standards for Microbiology Investigations. ID 12 Utgave 4. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>.

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BIBLIOGRAFIA

1. Public Health England. 2020. „Identification of Staphylococcus species, Micrococcus species and Rothia species”. UK Standards for Microbiology Investigations. ID 7 Issue 4. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
2. Public Health England. 2021a. „Identification of Streptococcus species, Enterococcus species and morphologically similar organisms”. UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
3. Ellner, P. D., C. J. Stoessel, E. Drakeford i F. Vasi. 1966. „New Culture Medium for Medical Bacteriology”. *The American Journal of Clinical Pathology* 45 (4): 502–4. https://doi.org/10.1093/ajcp/45.4_ts.502
4. Public Health England. 2021b. „Identification of Haemophilus species and the HACEK group of organisms”. UK Standards for Microbiology Investigations. ID 12 Issue 4. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>.

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Bibliografia

1. Public Health England. 2020. "Identification of Staphylococcus species, Micrococcus species and Rothia species." UK Standards for Microbiology Investigations. ID 7 Issue 4. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
2. Public Health England. 2021a. "Identification of Streptococcus species, Enterococcus species and morphologically similar organisms." UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
3. Ellner, P. D., C. J. Stoessel, E. Drakeford e F. Vasi. 1966. "A New Culture Medium for Medical Bacteriology." *The American Journal of Clinical Pathology* 45 (4): 502-4. https://doi.org/10.1093/ajcp/45.4_ts.502
4. Public Health England. 2021b. "Identification of Haemophilus species and the HACEK group of organisms." UK Standards for Microbiology Investigations. ID 12 Issue 4. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>.

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Bibliografia

1. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
2. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms.' UK Standards for Microbiology Investigations. <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
3. Ellner, P. D., C. J. Stoessel, E. Drakeford, and F. Vasi. 1966. 'A New Culture Medium for Medical Bacteriology.' *The American Journal of Clinical Pathology* 45 (4): 502–4. https://doi.org/10.1093/ajcp/45.4_ts.502
4. Public Health England. 2021b. 'Identification of Haemophilus species and the HACEK group of organisms.' UK Standards for Microbiology Investigations. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>.

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Referenser

1. Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
2. Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms.' UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.
3. Elner, P. D., C. J. Stoessel, E. Drakeford, and F. Vasi. 1966. 'A New Culture Medium for Medical Bacteriology.' *The American Journal of Clinical Pathology* 45 (4): 502–4. https://doi.org/10.1093/ajcp/45.4_ts.502.
4. Public Health England. 2021b. 'Identification of Haemophilus species and the HACEK group of organisms.' UK Standards for Microbiology Investigations. ID 12 Issue 4. <https://www.gov.uk/government/publications/smi-id-12-identification-of-haemophilus-species-and-the-hacek-group-of-organisms>.

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Övriga varumärken tillhör Thermo Fisher Scientific Inc. och deras dotterbolag.