Key Code TSMX4147D
www.thermofisher.com

Europe 00800 135 79 135 US 1 855 2360 190
CA 1 855 805 8539 ROW +31 20 794 7071

E. coli O157 Latex Test

REF DR0620M.....  100

1. INTENDED USE

A qualitative latex slide agglutination test for the confirmatory identification of *E. coli* serogroup O157 isolates grown on agar. Used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having bacterial infections. The device is not automated, is for professional use only and is not a companion diagnostic.

2. SUMMARY AND PRINCIPLE OF TEST

A latex agglutination test for the identification of *E. coli* serogroup O157.

Certain strains of *Escherichia coli* have been implicated in cases of haemorrhagic colitis (HC) and haemolytic uremic syndrome (HUS).

It has been shown that these strains produce a vero-cytotoxin (VT). The *E. coli* serotype most frequently isolated from HC and HUS cases is O157:H7. Isolation of this serotype from a diarrhoeal stool particularly with the presence of blood is indicative of a vero-cytotoxin producing strain^{1, 2, 3, 4, 5, 6, 7, 8}.

The *E. coli* O157 latex test will demonstrate by slide agglutination *E. coli* strains possessing the O157 serogroup antigen. The test is best used in conjunction with Sorbitol MacConkey Agar (Oxoid CM0813). *E. coli* O157:H7 strains do not ferment Sorbitol and therefore give colourless colonies on this medium. The majority of *E. coli* isolates do ferment Sorbitol and give characteristic pink colonies.

Sorbitol MacConkey Agar should be used as the primary screen. Non-sorbitol fermenting colonies can then be tested with the latex reagents, to determine whether the isolate belongs to the O157 serogroup and therefore a potential VT- producing strain.

3. COMPONENTS OF THE KIT

DR0621M Test Latex

Consists of blue latex particles sensitised with specific rabbit antibody reactive with the O157 somatic antigen. Each kit contains sufficient reagent for 100 tests.

DR0622M Control Latex

Consists of blue latex particles sensitised with pre-immune rabbit globulin. Each kit contains sufficient reagent for 100 tests.

DR0623M Positive Control Suspension

A suspension of inactivated *E. coli* O157 cells in buffer. Sufficient for 25 tests.

DR0624M Negative Control Suspension

A suspension of *E. coli* O116 cells in buffer. Sufficient for 25 tests.

DR0500G Reaction Cards

There are 35 disposable reaction cards provided in the kit.

Instructions for Use (IFU).

4. MATERIALS REQUIRED

The following materials are required but not provided in this kit
Microbiological loop and bunsen burner.

0.85% saline.

Suitable disinfectant e.g. sodium hypochlorite solution >1.3% w/w.

Timer

5. PRECAUTIONS

IVD The reagents are for *in vitro* diagnostic use only.

For professional use only.

In the event of malfunction do not use device.

Please refer to the Material Safety Data Sheet (MSDS), available on the company website, and product labeling for information on potentially hazardous components.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Do not freeze.

Reagents with different lot numbers should not be interchanged.

Reagents contain 0.1% sodium azide as a preservative.

Sodium azide may react with lead or copper plumbing to produce metal azides which are explosive by contact detonation. To prevent azide accumulation in plumbing, flush with copious amounts of water immediately after waste disposal.

Specimen materials may contain pathogenic organisms handle with appropriate precautions.

6. STORAGE

This kit must be stored at 2 to 8°C. Under these conditions the reagents will retain their reactivity until the expiry date shown on the kit box.

7. CONTROL PROCEDURES

The Control Suspensions provided should be used to check the correct working of the latex reagents each day before routine tests are performed.

The Positive Control Suspension must cause visible agglutination with the Test Latex reagent within one minute.

The Negative Control Suspension should cause no agglutination with the Test and Control Latex within one minute.

Do not use the test if reactions with the Control Suspensions are incorrect.

8. IMPORTANT PROCEDURE NOTES

Do not allow the reagents to become contaminated by letting the dropper tip touch the specimen on the reaction card.

Ensure that the caps are securely fitted after each use to prevent contamination and drying out of the reagent.

After use return the kit to the refrigerator ensuring that the bottles are stored in an upright position.

9. CULTURE MATERIAL

Non-sorbitol fermenting (NSF) colonies may be taken from Sorbitol MacConkey Agar (Oxoid CM0813) or alternatively NSF isolates may be inoculated onto non- selective media such as Nutrient Agar for testing. It is necessary to test up to 10 separate

NSF colonies to ensure a high probability of detecting any O157 strains which may be in a mixed culture with NSF *E. coli* of other serotypes. The use of the Control Latex will ensure that the isolate is not an autoagglutinating strain.

10. TEST METHOD

- Allow the reagents to equilibrate to room temperature before use. Make sure the latex suspensions are mixed by vigorous shaking. Expel any latex from the dropper pipette for complete mixing.
- Dispense 1 drop of the Test Latex onto a circle on the reaction card. Place it close to the edge of the circle.
- Add some loopfuls or a pasteur pipette drop of saline to the circle. Ensure that the latex and saline do not mix at this stage.
- Using a loop pick off a portion of the colony to be tested and carefully emulsify in the saline drop. Ensure that the resulting suspension is smooth.
- Mix the Test Latex and suspension together and spread to cover the reaction area using the loop. Flame the loop.
- Rock the card in a circular motion observing for agglutination. Do not rock the card for more than 1 minute and do not use a magnifying glass.
- If no agglutination occurs then proceed to test other NSF colonies if these are present.
- If agglutination with the Test Latex does occur, then it is necessary to test a further portion of the colony with the Control Latex reagent to ensure that the isolate is not an auto-agglutinating strain.
- When finished dispose of the reaction card into disinfectant.

11. INTERPRETATION OF TEST RESULTS

Agglutination of the Test Latex within one minute is a positive result. This indicates the presence of *E. coli* Serogroup O157.

No agglutination occurring within one minute is a negative result. This indicates the absence of *E. coli* Serogroup O157.

Results cannot be interpreted if there is agglutination of both the Test and Control Latex.

Granular or Stringy Reactions

Occasionally granular or stringy reactions may be seen due to the particulate nature of the test material. When such reactions are seen to occur they should be interpreted using the following criteria:

The result is positive when, using the test reagent, greater clearing of the blue background occurs compared with the reaction of the control reagent. The result is negative, when there are no significant changes in the clearing of the blue background using the test and control reagents.

Some strains of *E. coli* are difficult to emulsify in saline and may give a stringy type reaction with both the test and the control reagents. This does not look like true agglutination and should be ignored. If this stringiness is found to be too severe for a correct judgement to be made then the colony should be suspended in 0.3 ml of saline. Allow the lumps to settle and retest this smooth supernatant.

12. LIMITATIONS OF THE TEST

If a positive result is obtained on a colony of unknown species, then biochemical tests should be performed to confirm that the organism is an *E. coli* strain.

Neither the Sorbitol MacConkey Agar nor the *E. coli* O157 latex test will directly confirm the isolate as a toxin-producing strain.

Other serotypes have been found which produce the vero-cytotoxin.

In addition other bacterial strains have been identified with epitopes that mimic those of *E. coli* O157 lipopolysaccharide.¹¹

Some strains of *Escherichia hermannii* may cross react with *E. coli* O157 sera and the latex test due to a shared antigen (Borczyk et al)¹. *E. hermannii* may be differentiated from *E. coli* as they ferment cellobiose, grow in the presence of potassium cyanide (KCN) and produce a yellow pigment, which may be delayed¹⁰.

13. PERFORMANCE CHARACTERISTICS

In clinical comparisons of the performance of the Oxoid *E. coli* O157 latex test and commercially antisera, results demonstrated 99.3% (151/152) agreement between the two methods⁹. The sensitivity of the Oxoid tests was 100% (49/49), and the specificity was 99.0% (101/103). The discrepant sample was further tested and found to be positive for *E. coli* O157. This was in agreement with the Oxoid test result.

The Oxoid reagents did not show cross reactivity with other *E. coli* strains or non sorbitol fermenting organisms.

The Oxoid test was performed directly from a variety of culture media. None of these media were found to cause interference with the latex test.

The performance of the *E. coli* O157 latex test was evaluated at the Central Public Health Laboratory (CPHL), UK, in the Laboratory of Enteric Pathogens (1998), using fully characterised strains of target and non-target organisms.

All strains were tested after sub-culture onto Sorbitol MacConkey agar (SMAC) and Sorbitol MacConkey agar containing Cefixime Tellurite supplement (CT-SMAC). A lower number of organisms were tested from CT-SMAC, as due its selective nature, some organisms could not be cultured on this medium,

Sensitivity was calculated from results of testing known *E. coli* O157 strains. Both sorbitol negative and sorbitol positive strains were tested and used to calculate this. Specificity was calculated from results of testing known non-sorbitol fermenting (NSF), non-*E. coli* O157 strains (i.e. other organisms, including other *E. coli* strains).

Table 1. Results of testing organisms cultured on SMAC and CT-SMAC.

		SMAC (n=120)		CT-SMAC (n=95)	
		Organism			
		+ (i.e. <i>E. coli</i> O157)	- (i.e. non- <i>E. coli</i> O157)	+ (i.e. <i>E. coli</i> O157)	- (i.e. non- <i>E. coli</i> O157)
<i>E. coli</i> O157 Latex Test	+	59	0	52	0
	-	1	60	0	43

Table 2. Statistical analysis of data.

	Sensitivity %	Specificity %
SMAC		
<i>E. coli</i> O157 Latex Test	98.3	100
CT-SMAC		
<i>E. coli</i> O157 Latex Test	100	100

On SMAC agar one culture, an *E. coli* O157:H9, gave a false negative result. This culture was also negative when tested with a competitor kit and was therefore not re-tested. This organism is not a verocytotoxin-producing *E. coli* VTEC) strain. No VTEC strains were misidentified by the kit.

In addition, a total of 65 sorbitol fermenting, non-*E. coli* O157 strains were tested from SMAC and 52 from CT-SMAC. None gave a positive result with the *E.coli* O157 latex test. One of the sorbitol fermenting strain (*Klebsiella oxytoca*) was found to be auto-agglutinating with the kit. On repeat testing this sorbitol positive strain was found to be negative.

14. REFERENCES

1. Borczyk A. Lior H. Crebin B (1987) Int. J. Food. Microbiol. 4, 347-349.
2. Konowalchuk, J., Speirs, J. and Stavric S. (1977) Infect. Immun. 18, 775-779.
3. Scotland S., Day N. and Rowe B. (1980) FEMS. Microbiol. Lett,7, 15-17.
4. Centres for Disease Control (1982) Morbid Mortal Wkly 31 580-585.
5. Karmali M, Still, B., Petrick M., and Lim C. (1983),Lancet i,619-620.
6. Johnson,W., Liroh, H. and Bezanson. (1983) Lancet i, 76.
7. March S., and Ratnam (1986) J.Clin. Microbiol.23. - 869-872.
8. Krishnan, C., Fitzgerald V., Dakin S., and Behme R. (1987) J. Clin. Microbial. 25, 1043-1047.
9. P.A.Chapman (1989) Evaluation of Commercial latex slide test for identifying Escherichia coli O157.42:1109-1110.
10. Brenner et al (1982) Journal of Clinical Microbiology Vol. 15 No.4 p.703
11. Westerman, R. B., Yongsheg, H. et al. (1997). Production and Characterisation of Monoclonal Antibodies Specific for the Lipopolysaccharide of Escherichia coli O157. J. Clin. Microbiol. 35, 679–684

15. SYMBOL LEGEND

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
	Unique Device Identifier
	Authorized representative in the European Community
	UK Conformity Assessed
	European Conformity Assessment
	Manufacturer

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For Technical assistance see, www.thermofisher.com

Version	Date of modifications introduced
X4147D	April 2024
	Updated to meet IVDR requirements

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Oxoid Limited, Wade Road, Hampshire