Key Code TSMX7382A
www.oxid.com/ifu

Europe + 800 135 79 135
CA 1 855 805 8539

US 1 855 236 0910
ROW +31 20 794 7071

Oxoid Biochemical Identification System (O.B.I.S.) campy

REF ID0800M 60

INTENDED USE

The Oxoid Biochemical Identification System (O.B.I.S.) campy is a rapid, qualitative, colorimetric test for the detection of L-alanyl aminopeptidase. The test is used for the presumptive differentiation of *Campylobacter*, *Helicobacter* and *Arcobacter* species grown on plated media from other Gram-negative organisms. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

PRINCIPLE OF THE TEST

The O.B.I.S. campy test will differentiate species of *Campylobacter*, *Helicobacter* and *Arcobacter* from all other Gram negative organisms¹. Unlike other Gram-negative organisms, *Campylobacteraceae* do not possess the enzyme L-alanyl aminopeptidase. The O.B.I.S. campy test (L-ALA) incorporates a rapid test to detect this enzyme and a Gram-lysis reagent to rapidly determine the Gram status.

First, the Gram-lysis test (or a Gram stain) must be carried out. This test differentiates between Gram-positive and Gram-negative bacteria². The test is carried out on a glass slide. Sodium Hydroxide (0.5M) is used to lyse the cell wall of Gram-negative organisms and release the DNA. The DNA forms a 'string' which can be seen when the loop is raised from the surface of the slide. This reaction does not occur with Gram-positive organisms.

Once the organism has been identified as Gram-negative, the L-ALA test can be carried out.

Each O.B.I.S. reaction card has been impregnated with the L-ALA substrate (L-alanyl-7-amido-4-methylcoumarin) in each of the six reaction zones. An acidic solution of dimethylaminocinnamaldehyde (DMAC) is used as a colour developer. If the substrate has been hydrolysed by the organism, the free 7-amino-4-methylcoumarin will combine with the developer to produce a purple Schiff's base.

COMPONENTS OF THE O.B.I.S. CAMPY KIT (ID0800M)

Each kit contains the following reagents with enough material for 60 tests:

- | | |
|----------------|---|
| ID0803M | O.B.I.S. campy Test Cards – one resealable pouch containing 10 cards, each with six reaction zones.. |
| ID0804M | O.B.I.S. campy Buffer – one white capped dropper bottle containing 7 ml of solution. |
| ID0221M | O.B.I.S. DMAC Developer – one purple capped dropper bottle containing 7 ml of 0.5% w/v dimethylaminocinnamaldehyde in 1M hydrochloric acid. |

ID0802M 0.5M Sodium Hydroxide – one flat white capped bottle containing 6 ml of 0.5M NaOH (sufficient for 600 Gram-lysis screening tests).

ID0898 Paddle pastettes – pack of 60.
Instruction leaflet.

Materials required but not included

Clean glass slides.

Sterile 10 µl plastic disposable inoculating loops.

Positive and negative quality control organisms.

Scissors

PRECAUTIONS

This product is for *in vitro* diagnostic use only.

Do not use the O.B.I.S. campy reagents beyond the stated expiry date.

The DMAC Developer contains a weak acid and will stain the skin. The 0.5M Sodium Hydroxide Solution is classified per applicable European Community (EC) Regulation as corrosive. The following are the appropriate Hazard (H) and Precautionary (P) statements.

DANGER

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P301+P330 +P331	IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
P303+P361 +P353	IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/shower.
P305+P351 +P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Used Test Cards and inoculating loops should be disposed of as biohazardous waste. This should be incinerated or autoclaved at 121°C for at least 15 minutes.

Campylobacter are pathogens. A low infective dose can cause gastroenteritis with potentially serious complications. Take appropriate precautions when handling potentially contaminated samples.

Please refer to the Safety Data Sheet (SDS) on company website and product labelling for information on potentially hazardous components.

Directions should be read and followed carefully.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

STORAGE AND OPENING

The O.B.I.S. campy kit must be stored at 2°C to 8°C. Allow the pouches to equilibrate to room temperature before use to prevent condensation forming on the Test Cards.

Open the pouches by cutting at the notch between the end seal and the clip-lock opening.

Remove the number of Test Cards required and reseal the pouch.

Use within 60 minutes.

If fewer tests are required than the number on the Test Card, cut the card and return the unused portion to the pouch. Do not return used Test Cards to the pouch.

Precipitation will occur after long-term storage of the Sodium Hydroxide – this does not affect performance. A sterile loop should always be used. Discard if there is any sign of contamination.

QUALITY CONTROL PROCEDURE

The following procedure should be performed each time the kit is used:

Gram-lysis Test

- Positive Control** – Use a known Gram-negative (Gram-lysis positive) organism, such as *Campylobacter jejuni* ATCC™33291™ (Thermo Scientific™ Culti-Loops™ R4601400). Follow the method given in the test procedure.
- Negative Control** – Use a known Gram-positive (Gram-lysis negative) organism such as *Bacillus cereus* ATCC™11778™ (Thermo Scientific™ Culti-Loops™ R4601220). Follow the method given in the test procedure.

L-ALA Test

- Positive Control** – use a known L-ALA positive organism, such as *Pseudomonas aeruginosa* ATCC™27853™ (Thermo Scientific™ Culti-Loops™ R4607060). Follow the method given in the test procedure. Ensure that a purple colour forms around the colony material within 20 seconds.
- Negative Control** – Use a known L-ALA negative organism such as *Campylobacter jejuni* ATCC™33291™ (Thermo Scientific™ Culti-Loops™ R4601400). Follow the method given in the test procedure. Ensure that no purple colour forms within 20 seconds.

SPECIMENS

The test is designed for use from purity plates such as Columbia Blood Agar.

Primary isolation media should not be used as the colonies may be too small or too few to carry out an effective test. Pick colonies which have typical *Campylobacter* morphology from selective *Campylobacter* isolation media and streak onto Columbia Blood Agar. Incubate in a microaerobic atmosphere for 48 hours then conduct the O.B.I.S. test..

TEST PROCEDURE AND INTERPRETATION OF RESULTS

Gram-lysis Test

- Take a 10 µl loopful of 0.5M Sodium Hydroxide solution (white flat-capped bottle) and place onto a clean, glass slide.
- Using a sterile loop, take a small amount of material from a purity plate.
- Mix the sample into the NaOH on the slide for up to one minute.
- At intervals, carefully raise the loop from the mixture to check for the presence of a 'string' between the loop and the mixture.
- Record the result:

- A positive result is characterised by a viscous appearance to the mixture and the presence of the DNA string. This indicates that the organism is **Gram-negative**.
- A negative result is characterised by the formation of a cell suspension with no string and indicates that the organism is **Gram-positive**.

This test may also be carried out using a 0.5M potassium hydroxide (KOH) solution. Alternatively, a traditional Gram stain may be used to ascertain the Gram status of the organism.

L-ALA Test

Once the organism has been identified as Gram-negative, the L-ALA portion of the test can be carried out.

- Remove one of the Test Cards from the pouch.
- Using the paddle end of an unused plastic pastette, transfer the equivalent of 5x1 mm diameter colonies from a purity plate to the test area.
- Spread the sample across the reaction zone (inside a circle) of a Test Card.
- Dispense one drop of O.B.I.S. campy buffer (white capped dropper bottle) onto each of the inoculated reaction zones.
- Wait for 30 seconds.
- Dispense one drop of O.B.I.S. DMAC Developer (purple cap) onto each of the inoculated reaction zones.
- The appearance of a purple colour around the original colony material within 20 seconds is a positive L-ALA reaction. A positive reaction indicates the organism is not a *Campylobacter*, *Helicobacter* or *Arcobacter* species. If no colour develops around the original colony material within 20 seconds this is a negative reaction. A negative reaction indicates the organism is a presumptive *Campylobacter*, *Helicobacter* or *Arcobacter* species.


Interpretation Chart

Typical reactions:

Organism	0.5M NaOH	L-ALA
<i>Campylobacteraceae</i>	+ (i.e. Gram negative)	–
Gram-positives	– (i.e. Gram positive)	Not Tested
Other Gram-negatives	+ (i.e. Gram negative)	+

Venn Diagram

The Venn Diagram demonstrates how *Campylobacter* and *Campylobacter*-related organisms can be differentiated from other common bacteria using the Gram-lysis and L-ALA tests of O.B.I.S. campy.

LIMITATIONS OF THE TEST

O.B.I.S. campy is intended for the detection of L-alanyl aminopeptidase in Gram-negative organisms. It can be used for the presumptive identification of *Campylobacter* and related organisms from pure culture.

The reaction with the O.B.I.S. campy L-ALA test is a marker for enzyme activity and atypical strains may occasionally occur. *Bacteroides ureolyticus* can give the same reactions as *Campylobacter* species. However, their colonial morphology and anaerobiosis aid differentiation.

Over time, the DMAC Developer may give a slightly pink reaction with the campy buffer even in negatives. However, this is easily differentiated from the clearly purple reaction seen around the colony material in a positive.

Do not use nichrome wire loops to inoculate cards as this material can interfere with the test.

PERFORMANCE (DATA ON FILE)

In an internal study, 46 *Campylobacter* species and 6 *Arcobacter* species were tested. All yielded a Gram-negative, L-ALA negative reaction. 252 other species (non-*Campylobacter*, *Arcobacter* or *Helicobacter*) were tested. Only one organism gave a result similar to *Campylobacter*, and was due to an atypical Gram-lysis reaction. This resulted in a sensitivity and specificity of 100% and 99.6% respectively³.

REFERENCES

1. Hoosain, N. and Lastovica A. J. (2005). Evaluation of the Oxoid Biochemical Identification System (O.B.I.S.) for the differentiation of *Campylobacter* and *Arcobacter* from other Gram-negative organisms. In: Abstracts of CHRO 2005. 13th International Workshop on *Campylobacter*, *Helicobacter* and related organisms. Sept 4-8, 2005, Gold Coast, Queensland, Australia. Griffith University.
2. Carlone, G. M., Valadez, M. J. and Pickett, M. J. (1982) Methods for distinguishing Gram-positive from Gram-negative bacteria. *Journal of Clinical Microbiology* 16 (6), 1157-1159.
3. Smith, C. M., Colborne, N. R., Stephens, P. J. and Druggan, P. (2006) a simple and rapid biochemical screening test for the differentiation of *Campylobacter* spp. from other contaminating micro-organisms. In: Abstracts of Emerging *Campylobacter* spp. in the food chain, CAMPYCHECK. Feb 8th 2006, Croke Park Conference Centre, Dublin, Ireland.

SYMBOL LEGEND

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Do not re-use
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
	Unique Device Identifier
	Authorized representative in the European Community
	UK Conformity Assessed
	European Conformity Assessment

	Manufacturer
---	--------------

IFU X7382C Revised May 2025

Printed in the UK

Oxoid Ltd
Wade Road, Basingstoke
Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

©2025 Thermo Fisher Scientific Inc. All rights reserved.

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

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REFERENZEN

1. Hoosain, N. and Lastovica A. J. (2005). Evaluation of the Oxoid Biochemical Identification System (O.B.I.S.) for the differentiation of *Campylobacter* and *Arcobacter* from other Gram-negative organisms. In: Abstracts of CHRO 2005. 13th International Workshop on *Campylobacter*, *Helicobacter* and related organisms. 4–8. Sept, 2005, Gold Coast, Queensland, Australien. Griffith University.
2. Carlone, G. M., Valadez, M. J. and Pickett, M. J. (1982) Methods for distinguishing Gram-positive from Gram-negative bacteria. Journal of Clinical Microbiology 16 (6), 1157-1159.
3. Smith, C. M., Colborne, N. R., Stephens, P. J. and Druggan, P. (2006) a simple and rapid biochemical screening test for the differentiation of *Campylobacter* spp. from other contaminating micro-organisms. In: Abstracts of Emerging *Campylobacter* spp. in the food chain, CAMPYCHECK. 8. Feb. 2006, Croke Park Conference Centre, Dublin, Irland.

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REFERENCIAS

1. Hoosain, N. and Lastovica A. J. (2005). Evaluation of the Oxoid Biochemical Identification System (O.B.I.S.) for the differentiation of *Campylobacter* and *Arcobacter* from other Gram-negative organisms. In: Abstracts of CHRO 2005. 13th International Workshop on *Campylobacter*, *Helicobacter* and related organisms. Sept 4-8, 2005, Gold Coast, Queensland, Australia. Griffith University.
2. Carlone, G. M., Valadez, M. J. and Pickett, M. J. (1982) Methods for distinguishing Gram-positive from Gram-negative bacteria. Journal of Clinical Microbiology 16 (6), 1157-1159.
3. Smith, C. M., Colborne, N. R., Stephens, P. J. and Druggan, P. (2006) a simple and rapid biochemical screening test for the differentiation of *Campylobacter* spp. from other contaminating micro-organisms. In: Abstracts of Emerging *Campylobacter* spp. in the food chain, CAMPYCHECK. Feb 8th 2006, Croke Park Conference Centre, Dublin, Ireland.

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IFU X7382C Mayo 2025

 Oxoid Ltd, Wade Road, Basingstoke, Hampshire, RG24 8PW, Reino Unido

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BIBLIOGRAFIA

1. Hoosain, N. and Lastovica A. J. (2005). Evaluation of the Oxoid Biochemical Identification System (O.B.I.S.) for the differentiation of *Campylobacter* and *Arcobacter* from other Gram-negative organisms. In: Abstracts of CHRO 2005. 13th International Workshop on *Campylobacter*, *Helicobacter* and related organisms. Sept 4-8, 2005, Gold Coast, Queensland, Australia. Griffith University.
2. Carlone, G. M., Valadez, M. J. and Pickett, M. J. (1982) Methods for distinguishing Gram-positive from Gram-negative bacteria. Journal of Clinical Microbiology 16 (6), 1157-1159.
3. Smith, C. M., Colborne, N. R., Stephens, P.J. and Druggan, P. (2006) a simple and rapid biochemical screening test for the differentiation of *Campylobacter* spp. from other contaminating micro-organisms. In: Abstracts of Emerging *Campylobacter* spp. in the food chain, CAMPYCHECK. Feb 8th 2006, Croke Park Conference Centre, Dublin, Ireland.

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IFU X7382B

Maggio 2025

Oxoid Ltd, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

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©2025 Thermo Fisher Scientific Inc. Tutti i diritti riservati.
Tutti gli altri marchi sono di proprietà di Thermo Fisher Scientific Inc. e delle sue consociate.