remel

HIPPURATE DISK

REF R21085.....  25

1. INTENDED USE

The Hippurate Disk is a reagent-impregnated disk for use in qualitative procedures to detect the hydrolysis of sodium hippurate by hippuricase. The device is used as an aid to the presumptive identification of Group B streptococci (*Streptococcus agalactiae*), *Gardnerella vaginalis* and *Campylobacter jejuni*. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

2. SUMMARY AND EXPLANATION

The hippurate hydrolysis test determines the ability of Group B streptococci, as well as other bacteria, to enzymatically hydrolyze sodium hippurate.^{1,2} Hwang and Ederer employed the ninhydrin reagent in a rapid method to detect hippurate hydrolysis by the end product glycine.³ Jolly later outlined minimal criteria for identification of *G. vaginalis* including hippurate hydrolysis.⁴ In 1989, Cacho et al. further evaluated the use of the hippurate disk in the identification of *Campylobacter* spp.⁵

3. PRINCIPLE

Hippuric acid is hydrolyzed by the enzyme hippuricase to glycine and benzoic acid. Ninhydrin evokes a five-step reaction beginning with the deamination of glycine to form hydrindantin (a reduced form of ninhydrin), carbon dioxide, and ammonia. A condensation reaction occurs with hydrindantin, ammonia, and residual ninhydrin to produce the final purple colored complex.

4. REAGENTS

Hippurate Disk (25 disks/vial)

Reactive Ingredient: Sodium Hippurate

5. PRECAUTIONS

This product is for in vitro diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after their use. Please refer to the Safety Data Sheet (SDS) on company website and product labelling for information on potentially hazardous components.

Directions should be read and followed carefully.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction do not use device.

Individual test disks are for single use and must be disposed of after use.

6. STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 2-8°C until used. Allow product to equilibrate to room temperature before use. Do not incubate prior to use.

7. PRODUCT DETERIORATION

This product should not be used if (1) the color has changed from white, (2) the expiration date has passed, (3) the desiccant has changed from blue to pink, or (4) there are other signs of deterioration. Protect disks from moisture by removing from the vial only those disks necessary for testing. Promptly replace the cap and return the vial to 2-8°C.

8. SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.⁶

9. MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Plastic test tubes, (7) Forceps, (8) Ninhydrin Reagent (R21238) or BactiDrop™ Ninhydrin (R21534), (9) Demineralized water.

10. PROCEDURE

1. Add 0.1 ml of sterile demineralized water to a small plastic test tube.
2. Take 1-3 colonies of a fresh isolate and emulsify in the water.
3. Using forceps drop a Hippurate Disk into the suspension.
4. Incubate aerobically for 2 hours at 35-37°C.
5. Dispense 2 drops of Ninhydrin Reagent or 3 drops BactiDrop Ninhydrin into each test tube, mix, and reincubate aerobically for 30 minutes at 35-37°C.
6. Observe for a blue to purple color development.

11. INTERPRETATION

Positive Test -	Blue to purple color development within 30 minutes
Negative Test -	Light gray color development or no color change within 30 minutes

12. QUALITY CONTROL

All lot numbers of Hippurate Disk have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory

quality control procedures. If aberrant quality control results are noted, patient results should not be reported. Quality control testing should be run with each shipment and new lot number received.

CONTROL	INCUBATION	RESULTS
<i>Gardnerella vaginalis</i> ATCC® 14018	Aerobic, 2.5 h @ 35-37°C	Positive
<i>Streptococcus agalactiae</i> ATCC® 12386	Aerobic, 2.5 h @ 35-37°C	Positive
<i>Streptococcus pyogenes</i> ATCC® 19615	Aerobic, 2.5 h @ 35-37°C	Negative

13. LIMITATIONS

1. Positive and negative controls should be included with each group of tests to better correlate the color changes produced by both positive and negative tests.
2. Plastic tubes should be used, as glass may yield a false-negative result.

14. BIBLIOGRAPHY

1. Greenwood, J.R. and M.J. Pickett. 1979. J. Clin. Microbiol. 9:200-204.
2. MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Lippincott Williams & Wilkins, Philadelphia, PA.
3. Hwang, M. and G.M. Ederer. 1975. J. Clin. Microbiol. 1:114-115.
4. Jolly, J.L.S. 1983. J. Clin. Pathol. 36:476-478.
5. Cacho, J.B., P.M. Aguirre, A. Hernanz, and A.C. Velasco. 1989. J. Clin. Microbiol. 27:359-360.
6. Jorgensen J.H., M.A. Pfaller et al. 2015. Manual of Clinical Microbiology. 11th edition. Washington ASM Press.

15. PACKAGING

REF R21085, Hippurate Disk 25 Disks/Vial

16. SYMBOL LEGEND

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Batch Code (Lot Number)
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Do not use if package is damaged
	Do not re-use
	Use By (Expiration Date)
	Authorized representative in the European Community
	Importer
	UK Conformity Assessed
	European Conformity Mark
	Manufacturer

For technical information contact your local distributor.
 ©2022 Thermo Fisher Scientific Inc. All rights reserved.
 ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
 All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

 Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, USA

www.thermofisher.com/microbiology
 Tel: (800) 255-6730 • International: (913) 888-0939

www.oxid.com/IFU
 Europe +800 135 79 135 • US 1 855 2360 190
 CA 1 855 805 8539 • ROW +31 20 794 7071

Version	Date of modifications introduced
IFU21085	March 2023 Addition of Importer symbol to Section 16

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