Lactophenol Cotton Blue Droppers

Blue stain for microscopic examination of fungi,
in dropper bottles.

DESCRIPTION

Lactophenol Cotton Blue Droppers is a mounting medium and a staining agent used in the preparation of slides for microscopic examination of fungi from clinical specimens, food and environmental samples.

KIT CONTENTS

- 10 bottles of Lactophenol Cotton Blue Droppers.
- 1 instruction sheet.

METHOD PRINCIPLE

Upon the addition of Lactophenol Cotton Blue, fungi stain blue allowing for easier visualization and examination under microscope.

COMPOSITION

Aqueous solution of the following ingredients (per 100 ml):

Glycerol 50 mg, Phenol 25 mg, Lactic acid 25 mg, Methyl blue (Cotton blue) 50 mg.

TEST PROCEDURE

1. Dispense a drop of Lactophenol Cotton Blue stain on a slide.
2. Using an inoculating needle, carefully spread the fungal culture into a thin preparation.
3. Place a coverslip edge on the drop and slowly lower it. Avoid trapping air bubbles under the coverslip. Wait for about 5 minutes.
4. Observe under microscope.

INTERPRETATION OF RESULTS

Yeast cells, mycelia and fruiting structures are stained a delicate blue colour, while the background appear a faint, pale blue.

QUALITY CONTROL FOR THE USER

Prepare a wet mount with the stain according to laboratory methodology and carry out microscopic analysis.

Control strains

Aspergillus niger

ATCC® 16404

Trichophyton mentagrophytes

ATCC® 9533

Appearance

Delicate blue hyphae and fruiting structures with a pale blue background

Delicate blue hyphae and fruiting structures with a pale blue background

PRECAUTIONS

Lactophenol Cotton Blue Droppers is classifiable as hazardous under current legislation; it is recommended that the Safety Data Sheet be consulted on its use. The product is intended for *in vitro* diagnostic use and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the storage or transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, do not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

ELIMINATING USED MATERIAL

After use, used Lactophenol Cotton Blue Droppers and the material that has come into contact with the sample must be decontaminated and disposed of in accordance with the laboratory procedures for the decontamination and disposal of potentially infected material.


BIBLIOGRAPHY

1. Larone D.H. (1995) Medically important fungi: a guide to identification. Washington DC: ASM Press.
2. Lennette, E.H., A. Balows, W.J. Hausler, Jr., and H.J. Shadomy (ed.). (1985) Manual of clinical microbiology, 4th ed. American Society for Microbiology, Washington, D.C.
3. Finegold, S.M., and W.J. Martin. (1982) Bailey and Scott's diagnostic microbiology, 6th ed. The C.V. Mosby Co., St. Louis.

PRESENTATION

Product	Ref.	Content
Lactophenol Cotton Blue Droppers	87008	10 x 10 ml bottles

TABLE OF SYMBOLS

 Batch code	 <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
 Catalogue number	 Temperature limitation	 Contains Sufficient for <n> tests	 Caution, consult accompanying documents	 Do not reuse

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