Legionella Latex Test

REF	DR0801M Legionella pneumo group 1	50
REF	DR0802M Legionella pneumo groups 2-14	50
REF	DR0803M Legionella species test kit.....	50

1. INTENDED USE

This qualitative latex agglutination test is for confirmatory identification of *Legionella* species from isolates grown on agar:

- DR0801M for *L. pneumophila* serogroup 1;
- DR0802M for *L. pneumophila* serogroup 2-14;
- DR0803M for non-pneumophila *Legionella* species (see below).

Used in a diagnostic workflow to aid clinicians in treatment options for patients suspected of having bacterial infections. The device is not automated, is for professional use only and is not a companion diagnostic.

Note. DR0800M is a combination test which can detect all these *Legionella* groups.

The tests may also be used for the identification of *Legionella* species isolates from environmental sources.

2. INTRODUCTION

Legionnaires' disease, named after the outbreak in 1976 at the American Legion Convention in Philadelphia, is caused by *Legionella pneumophila* and other *Legionella* species. It is characterised as an acute febrile respiratory illness ranging in severity from mild illness to fatal pneumonia. Since that time, it has been recognised that the disease occurs in both epidemic and endemic form and that the sporadic cases are not readily differentiated from other respiratory infections by clinical symptoms.

It is estimated that worldwide about 25,000 cases of *Legionella* infections occur annually. Known risk factors include immunosuppression, cigarette smoking, alcohol consumption and concomitant pulmonary disease. The mortality rate, which can be as high as 25% in untreated immunosuppressed patients, can be lowered if the disease is diagnosed rapidly and appropriate antimicrobial therapy started earlier.

Legionella pneumophila has been shown to be a major cause of both pneumonia and an acute self limiting febrile disease called Pontiac Fever. *L. pneumophila* strains and other *Legionella* species are isolated from patients with pneumonia and from the environment (mainly water).

Rare isolations have also been made in cases other than pneumonia, such as wound abscesses. The major reservoir of *Legionella* species appears to be fresh water sites, air-conditioning units and various water plumbing fixtures.

L. pneumophila is the most common cause of Legionnaires'

disease. At present, 14 different serotypes exist of which *L. pneumophila* serogroup 1 accounts for 90% of cases.

The *Legionella* Latex Test uses antibody sensitised blue polystyrene 'latex' particles which will agglutinate in the presence of specific *Legionella* cell wall antigens to form visible clumps. This provides a fast and simple screening procedure for predominant pathogenic *Legionella* species and serotypes.^{1,2}

3. COMPONENTS OF THE KIT

This kit contains one of the following test reagents:

DR0801 *Legionella pneumophila* serogroup 1 Test Reagent consists of a suspension of blue polystyrene 'latex' particles sensitised with specific rabbit antibody reactive with *Legionella pneumophila* serogroup 1 antigen. Each kit contains sufficient reagent for 50 tests.

OR

DR0802 *Legionella pneumophila* serogroups 2–14 Test Reagent consists of a suspension of blue polystyrene 'latex' particles sensitised with specific rabbit antibody reactive with *Legionella pneumophila* serogroups 2–14 antigen and is also able to detect the *Legionella pneumophila* serogroup 15 (ATCC® 35251). Each kit contains sufficient reagent for 50 tests.

OR

DR0803 *Legionella* species Test Reagent consists of a suspension of blue polystyrene 'latex' particles sensitised with specific rabbit antibody reactive with the following species and serotypes:

- L. longbeachae* 1 and 2
- L. bozemanii* 1 and 2
- L. dumoffii*
- L. gormanii*
- L. jordanis*
- L. micdadei*
- L. anisa*

Each kit contains sufficient reagent for 50 tests.

This kit also contains:

DR0806 *Legionella* Control Latex consists of blue polystyrene 'latex' particles sensitised with non-reactive rabbit globulin. Each kit contains sufficient reagent for 50 tests.

DR0807 *Legionella* Suspension Buffer

A phosphate buffered saline solution. pH 7.3.

Instructions for Use

4. MATERIALS REQUIRED

The following materials are required but not provided in this kit:

Microbiological loop

Bunsen burner

0.85% saline (for optional tube method)

Suitable laboratory disinfectant e.g. Sodium hypochlorite solution >1.3% w/v.

DR0500G Disposable Reaction Card

5. PRECAUTIONS

IVD These reagents are for *in vitro* diagnostic use only.

Do not freeze.

Reagents contain 0.1% sodium azide as a preservative.

In the event of malfunction do not use device.

Sodium azide may react with lead or copper plumbing to produce metal azides which are explosive by contact detonation. To

prevent azide accumulation in plumbing flush with copious amounts of water immediately after waste disposal.

Specimen materials may contain pathogenic organisms, handle with appropriate precautions.

Aerosol formation should be avoided. Particular care should be taken during vortexing.

Please refer to the safety data sheet, available on the ThermoFisher website, and the product labelling for information on potentially hazardous components.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

6. STORAGE

This kit must be stored 2-8°C. Under these conditions the reagents will retain their reactivity until the expiry date shown on the kit box.

7. IMPORTANT PROCEDURE NOTES

Do not allow the reagents to become contaminated by letting the dropper tip touch the specimens on the reaction card. Ensure that the caps are securely fitted after use to prevent contamination and drying out of reagents. After use, return the kit to the refrigerator ensuring that the bottles are stored in an upright position.

8. SPECIMEN COLLECTION AND PREPARATION

Isolates derived from environmental and clinical samples may be cultured on standard non-selective or selective *Legionella* culture media. *Legionella* species on primary isolation have an absolute requirement for L-cysteine. To ensure that an isolate is a *Legionella* it is necessary to show that it cannot grow on any medium which does not contain L-cysteine. This confirmation may be performed prior to or after the latex test.

Cultures may be tested at any stage of growth providing that the colonies are of sufficient size. Older cultures, however, may produce stringy reactions making interpretation more difficult.

9. TEST METHODS

There are two test methods which may be used. Both direct and tube methods give reliable results. If an isolate has a stringy consistency it is recommended that the tube method is used.

(a) Direct Test

1. Bring the latex reagents to room temperature. Make sure the latex suspensions are mixed by vigorous shaking. Expel any latex from the dropper pipette for complete mixing.
2. Dispense 1 drop of *Legionella* latex Test Reagent within and close to the edge of a circle on a reaction card (DR0500G).
3. Dispense 1 drop of *Legionella* Control Latex within and close to the edge of a second circle on the reaction card.
4. Add one drop of diluent buffer to each circle. Ensure that the latex and buffer do not mix at this stage.
5. Using a loop, pick off a colony of at least 1 mm (use 2 or more if the colonies are smaller) and carefully emulsify in the first drop of buffer. For optimal results ensure that the suspension is smooth. Repeat for similar colonies with the second drop of buffer.

6. Mix the latex reagents and suspensions together and spread to cover the reaction areas using the loop. Flame the loop.
7. Gently rock the card in a circular motion and look for agglutination. Do not rock the card for more than 1 minute and do not use a magnifying glass to aid reading the result.
8. When finished, dispose of the reaction card into a suitable disinfectant.
9. Recap the bottles and return to the refrigerator.

(b) Tube Method

1. Bring the latex reagents to room temperature. Make sure the latex suspensions are mixed by vigorous shaking. Expel any latex from the dropper pipette for complete mixing.
2. Label test tubes appropriately and dispense 0.4 ml of 0.85% saline into each tube.
3. Select 4-10 colonies of a similar colonial appearance with a loop and emulsify in the saline.
4. Vortex the cell suspension for 5 seconds.
5. Dispense 1 drop of the latex Test Reagent within and close to the edge of a circle on a reaction card (DR0500G) and repeat for the Control Latex.
6. Using a Pasteur pipette add 1 drop of cell suspension to each of the 2 circles, and mix this into the latex reagents. Spread to cover the reaction areas.
7. Gently rock the card in a circular motion and look for agglutination. Do not rock the card for more than 1 minute and do not use a magnifying glass to aid reading of the result.
8. When finished, dispose of the reaction card into a suitable disinfectant.
9. Recap the bottles and return to the refrigerator.

10. CONTROL PROCEDURES

Quality control testing should be run with each shipment and new kit lot number received. Each laboratory should follow their state and local requirements. In addition, known stock cultures of *Legionella* should be used periodically as controls.

11. READING AND INTERPRETATION OF RESULTS

Positive Results

A result is positive if agglutination of the blue polystyrene 'latex' Test Reagent particles occurs within 1 minute and with no agglutination in the Control Latex circle. A positive reaction indicates that antigens to that serogroup of *Legionella* species have been detected in the sample.

Negative Results

A negative result is obtained if no agglutination occurs and a smooth blue suspension remains after 1 minute in the test Reagent circles.

Uninterpretable Result

The test is uninterpretable if the Control Latex shows agglutination. This indicates that the culture causes autoagglutination.

Granular or Stringy Reactions

Occasional granular or stringy reactions may be seen due to the particulate nature of the test material. When such reactions are seen to occur they should be interpreted using the following criteria:

The result is positive when there is noticeable clearing of the blue background in the Test Reagent.

12. LIMITATIONS

1.

The latex agglutination test is presumptively diagnostic. Confirm positive results using biochemical tests.⁵

2.

A negative latex agglutination test does not mean that the culture is not a *Legionella* species. It only indicates that the culture is not a *Legionella* species detected by the test reagent used.

3.

A cross reaction may occur between *L. pneumophila* serogroup 1 and serogroup 9 due to naturally occurring group antigens. If both the *L. pneumophila* serogroup 1 and 2-14 reagents agglutinate with the isolate then this cross reaction should be suspected.

4.

Cross reactions with the Legionella Species Test Reagent have been reported to occur occasionally with certain serotypes of other *legionellae* (e.g. *L. parisiensis*, *L. saintelensi*, *L. steigerwaltii*, *L. wadsworthii*, *L. santicrucis*, *L. tusconensis*, *L. gratiana*, *L. cincinnatiensis*).⁶

5.

The test is designed to differentiate between different species and serotypes of *Legionella*. Cultures should be confirmed as Gram-negative rods which do not grow on cysteine-deficient media.

For further information please refer to instruction leaflet in the Legionella Latex Test (DR0800M).

13. PERFORMANCE CHARACTERISTICS

The reagents in the Oxoid Legionella Latex Test have been tested for cross-reactivity against a panel of organisms listed below. No cross-reactivity was observed with any of the organisms.

Aeromonas hydrophila
Bacillus subtilis
Citrobacter freundii
Escherichia coli
L. birminghamensis
L. brunensis
L. cherri
L. erythra
L. fairfieldensis
L. feelei
L. hackeliae
L. israeliensis
L. jamestowniensis
L. maceachernii
L. moravica
L. oakridgensis
L. quinlivanii
L. rubrilucens
L. spiritensis
Pseudomonas aeruginosa
Pseudomonas cepacia
Pseudomonas fluorescens
Serratia marcescens
Stenotrophomonas maltophilia

The Legionella Latex Test was evaluated in a clinical and an environmental laboratory. A total of 40 clinical isolates, 318 environmental isolates and 3 additional isolates of unconfirmed origin were tested, covering *Legionella pneumophila* serogroups 1–14 and non-*pneumophila* *Legionella*. Each isolate was confirmed by serology. The performance of the kit was also compared against other commercially available Legionella Latex Reagent Kits. The results of the trial are summarised opposite.⁷

L. pneumophila serogroup 15 has not currently been isolated from clinical or environmental samples in Europe⁹ and has only been isolated once in the USA.⁸

This most recent serogroup to be designated contains only one strain (Lansing-3 (ATCC® 35251). A 16th serogroup was proposed from studies of the Jena-1 isolate^{10,11} until further analysis showed that the strain did not form a unique serogroup but was a member of serogroup 4 *L. pneumophila* (monoclonal group Portland 1).^{12,13}

Additional internal data has been collected which demonstrates that the Legionella Latex Kits (DR0800M and DR0802M) are able to detect the *Legionella pneumophila* serogroup 15 (ATCC® 35251) as well as detecting *L. pneumophila* 1-14 and other pathogenic non-*L. pneumophila*, *Legionella* strains. However, as only one serogroup 15 strain is available for testing, it is not thought appropriate to rename the serogroup 2-14 latex reagent based on the results from a single isolate.

The Legionella kits benefit the user by allowing discrimination of samples into three groups: *L. pneumophila* serogroup 1, *L. pneumophila* serogroups 2-15 (with the 2-14 reagent) and other *Legionella* species in a fast and simple screening procedure.

	Legionella Latex Test	Latex Test/Serology
	Number	%
<i>Legionella pneumophila</i> Serogroup 1	59/59	100
<i>Legionella pneumophila</i> Serogroup 2-14	134/134	100
Other <i>Legionellae</i> included in the kit	63/65	97
Other <i>Legionellae</i> not included in the kit 0/60	0/93	100
Other Organisms	0/10	100

The overall sensitivity of the Legionella Latex Test was 99%.

The overall specificity of the Legionella Latex Test was 100%.

14. REFERENCES

1.

Sedgwick, A. K. and Tilton, R. C. (1983). J. Clin. Microbiol., 17: 365–368.

2.

Ciesielski, C. A., Blaser, M. J. and Wang, W. L. (1986). Infect. Immun., 51: 397–404.

3.

Dennis, P. J. L. (1988). Isolation of Legionellae from Environmental Specimens p. 31–44. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.

4.

Dournon, E. (1988). Isolation of Legionellae from Clinical Specimen p.13–30. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.

5.

Cowan and Steel’s Manual for the identification of Medical Bacteria 3rd Ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993) p. 161–163. University Press, London.

6.

Harrison, T. G. and Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.

7.

Data on file Oxoid Ltd.

8.

Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumonia* subsp. *fraseri* subsp. nov., and *L. pneumonia* subsp. *pascalii* subsp. nov. Journal of Clinical Microbiology; 1988; 26: 1695-1703.

9.

Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires’ Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710-716.

10.

Lück, C et al. Isolation of a Legionella pneumophila Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35-39.
11.

Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. Molecular Identification and Epidemiology; 1998; 47: 667-678.
12.

Lück, P. C et al. DNA Polymorphisms in Strains of Legionella pneumophila Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. Applied and Environmental Microbiology; 1995; 61: 2000-2003.
13.

Helbig, J. H. 2003. Personal Communication.

15. SYMBOL LEGEND

	Catalogue Number
	In Vitro Diagnostic Medical Device
	Consult Instructions for Use (IFU)
	Temperature Limitations (Storage temp.)
	Contains sufficient for <N> tests
	Not for near patient testing
	Batch Code (Lot Number)
	Use By (Expiration Date)
	Importer
	Unique Device Identifier
	Authorized representative in the European Community
	UK Conformity Assessed
	European Conformity Assessment
	Manufacturer

Version	Date of modifications introduced
X5058F	February 2025 Updated to meet IVDR requirements

Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, UK

For technical assistance www.thermofisher.com

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5. Cowan and Steel's Manual for the identification of Medical Bacteria 3rd Ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993) p. 161–163. University Press, London.
6. Harrison, T. G. and Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
7. Data on file Oxoid Ltd.
8. Brenner et al. Legionella pneumophila Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of L. pneumonia subsp. pneumophila subsp. nov., L. pneumophila subsp. fraseri subsp. nov., and L. pneumophila subsp. pascullei subsp. nov. Journal of Clinical Microbiology; 1988; 26: 1695-1703.
9. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires' Disease: Distribution of Legionella pneumophila Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710-716.
10. Lück, C et al. Isolation of a Legionella pneumophila Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35-39.
11. Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. Molecular Identification and Epidemiology; 1998; 47: 667-678.
12. Lück, P C et al. DNA Polymorphisms in Strains of Legionella pneumophila Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. Applied and Environmental Microbiology; 1995; 61: 2000-2003.
13. Helbig, J. H. 2003. Personal Communication.

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8. Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumonia* subsp. *fraseri* subsp. nov., and *L. pneumonia* subsp. *pascullei* subsp. nov. *Journal of Clinical Microbiology*; 1988; 26: 1695-1703.
9. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires’ Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. *European Journal of Clinical Microbiology and Infectious Disease*; 2002; 21: 710-716.
10. Lück, C et al. Isolation of a *Legionella pneumophila* Strain Serologically Distinguishable from all known Serogroups. *Zentralblatt fuer Bakteriologie*. 1995; 282: 35-39.
11. Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of *Legionella* species. *Molecular Identification and Epidemiology*; 1998; 47: 667-678.
12. Lück, P. C et al. DNA Polymorphisms in Strains of *Legionella pneumophila* Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. *Applied and Environmental Microbiology*; 1995; 61: 2000-2003.
13. Helbig, J. H. 2003. Personlig kommunikation.

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14. LITERATUR

- 1. Sedgwick, A. K. und Tilton, R. C. (1983). J. Clin. Mikrobiol., 17: 365–368.
- 2. Ciesielski, C. A., Blaser, M. J. und Wang, W. L. (1986). Infect. Immun., 51: 397–404.
- 3. Dennis, P. J. L. (1988). Isolation of Legionellae from Environmental Specimens p. 31–44. In Harrison, T. G. und Taylor, A. G. (Hrsg.). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
- 4. Dournon, E. (1988). Isolation of Legionellae from Clinical Specimen p.13–30. In Harrison, T. G. und Taylor, A. G. (Hrsg.). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
- 5. Cowan und Steel's Manual for the Identification of Medical Bacteria 3. Ausgabe. Barrow, G. I. und Feltham, R. K. A. (Hrsg.) (1993) p. 161–163. University Press, London.
- 6. Harrison, T. G. und Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. und Taylor, A. G. (Hrsg.). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
- 7. Archivierte Daten Oxoid Ltd.
- 8. Brenner et al. Legionella pneumophila Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of L. pneumonia subsp. pneumophila subsp. nov., L. pneumophila subsp. fraseri subsp. nov., and L. pneumophila subsp. pascullei subsp. nov. Journal of Clinical Microbiology; 1988; 26: 1695-1703.
- 9. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires' Disease: Distribution of Legionella pneumophila Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710-716.
- 10. Lück, C et al. Isolation of a Legionella pneumophila Strain Serologically Distinguishable from all known Serogroups. Zentralblatt für Bakteriologie. 1995; 282: 35-39.
- 11. Fry, N. K und Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. Molecular Identification and Epidemiology; 1998; 47: 667-678.
- 12. Lück, P. C et al. DNA Polymorphisms in Strains of Legionella pneumophila Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. Applied and Environmental Microbiology; 1995; 61: 2000-2003.
- 13. Helbig, J. H. 2003. Personal Communication.

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1. Sedgwick, A. K. και Tilton, R. C. (1983). J. Clin. Microbiol., 17: 365–368.
2. Ciesielski, C. A., Blaser, M. J. και Wang, W. L. (1986). Infect. Immun., 51: 397–404.
3. Dennis, P. J. L. (1988). Isolation of Legionellae from Environmental Specimens p. 31–44. In Harrison, T. G. και Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, HB.
4. Dournon, E. (1988). Isolation of Legionellae from Clinical Specimen σελ.13–30. In Harrison, T. G. και Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, HB.
5. Cowan and Steel's Manual for the identification of Medical Bacteria 3rd Ed. Barrow, G. I. και Feltham, R. K. A. (eds) (1993) σελ. 161–163. University Press, London.
6. Harrison, T. G. και Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. και Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, HB.
7. Data on file Oxoid Ltd.
8. Brenner et al. Legionella pneumophila Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of L. pneumonia subsp. pneumophila subsp. nov., L. pneumonia subsp. fraseri subsp. nov., and L. pneumonia subsp. pascullei subsp. nov. Journal of Clinical Microbiology, 1988, 26: 1695-1703.
9. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires' Disease: Distribution of Legionella pneumophila Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease, 2002, 21: 710-716.
10. Lück, C et al. Isolation of a Legionella pneumophila Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35-39.

11. Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. Molecular Identification and Epidemiology, 1998, 47: 667-678.
12. Lück, P. Cet al. DNA Polymorphisms in Strains of Legionella pneumophila Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. Applied and Environmental Microbiology, 1995, 61: 2000-2003.
13. Helbig, J. H. 2003. Personal Communication.

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14. REFERENCIAS

- 1. Sedgwick, A. K. y Tilton, R. C. (1983). J. Clin. Microbiol., 17: 365–368.
- 2. Ciesielski, C. A., Blaser, M. J. y Wang, W. L. (1986). Infect. Immun., 51: 397–404.
- 3. Dennis, P. J. L. (1988). Isolation of Legionellae from Environmental Specimens págs. 31–44. En Harrison, T. G. y Taylor, A. G. (editores) A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, Reino Unido.
- 4. Dournon, E. (1988). Isolation of Legionellae from Clinical Specimen págs. 13–30. En Harrison, T. G. y Taylor, A. G. (editores) A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, Reino Unido.
- 5. Cowan and Steel’s Manual for the identification of Medical Bacteria 3.^ª ed. Barrow, G. I. y Feltham, R. K. A. (editores) (1993) págs. 161–163. University Press, Londres.
- 6. Harrison, T. G. y Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. En Harrison, T. G. y Taylor, A. G. (editores) A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, Reino Unido.
- 7. Datos en archivo de Oxoid Ltd.
- 8. Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumonia* subsp. *fraseri* subsp. nov., and *L. pneumonia* subsp. *pascullei* subsp. nov. Journal of Clinical Microbiology; 1988; 26: 1695-1703.
- 9. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires’ Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710-716.
- 10. Lück, C et al. Isolation of a Legionella pneumophila Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35-39.
- 11. Fry, N. K y Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. Molecular Identification and Epidemiology; 1998; 47: 667-678.
- 12. Lück, P. C et al. DNA Polymorphisms in Strains of Legionella pneumophila Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. Applied and Environmental Microbiology; 1995; 61: 2000-2003.
- 13. Helbig, J. H. 2003. Comunicación personal.

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 Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Reino Unido

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- 6. Harrison, T. G. and Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
- 7. Data on file Oxoid Ltd.
- 8. Brenner et al. Legionella pneumophila Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of L. pneumonia subsp. pneumophila subsp. nov., L. pneumonia subsp. fraseri subsp. nov., and L. pneumonia subsp. pascuelli subsp. nov. Journal of Clinical Microbiology; 1988; 26: 1695-1703.
- 9. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires’ Disease: Distribution of Legionella pneumophila Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710-716.
- 10. Lück, C et al. Isolation of a Legionella pneumophila Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35-39.
- 11. Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. Molecular Identification and Epidemiology; 1998; 47: 667-678.
- 12. Lück, P. C et al. DNA Polymorphisms in Strains of Legionella pneumophila Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. Applied and Environmental Microbiology; 1995; 61: 2000-2003.
- 13. Helbig, J. H. 2003. Personal Communication.

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14. RÉFÉRENCES

- Sedgwick, A. K. et Tilton, R. C. (1983). J. Clin. Microbiol., 17: 365–368.
- Ciesielski, C. A., Blaser, M. J. et Wang, W. L. (1986). Infect. Immun., 51: 397–404.
- Dennis, P. J. L. (1988). Isolation of *Legionellae* from Environmental Specimens p. 31–44. Dans Harrison, T. G. et Taylor, A. G. (eds.) A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, Royaume-Uni.
- Dournon, E. (1988). Isolation of *Legionellae* from Clinical Specimen p.13–30. Dans Harrison, T. G. et Taylor, A. G. (eds.) A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, Royaume-Uni.
- Cowan and Steel's Manual for the identification of Medical Bacteria 3e éd. Barrow, G. I. et Feltham, R. K. A. (eds) (1993) p. 161–163. University Press, Londres.
- Harrison, T. G. et Taylor, A. G. (1988). Identification of *Legionellae* by Serological Methods. Dans Harrison, T. G. et Taylor, A. G. (eds.) A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, Royaume-Uni.
- Données archivées, Oxoid Ltd.
- Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumonia* subsp. *fraseri* subsp. nov., and *L. pneumonia* subsp. *pascullei* subsp. nov. Journal of Clinical Microbiology; 1988; 26: 1695-1703.
- Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires' Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710-716.
- Lück, C et al. Isolation of a *Legionella pneumophila* Strain Serologically Distinguishable from all known Serogroups. Zentralblatt für Bakteriologie. 1995; 282: 35-39.
- Fry, N. K et Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of *Legionella* species. Molecular Identification and Epidemiology; 1998; 47: 667-678.
- Lück, P. C et al. DNA Polymorphisms in Strains of *Legionella pneumophila* Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. Applied and Environmental Microbiology; 1995; 61: 2000-2003.
- Helbig, J. H. 2003. Communication personnelle.

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Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, Royaume-Uni

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5. Cowan and Steel’s Manual for the identification of Medical Bacteria 3rd Ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993) p. 161–163. University Press, London.

6. Harrison, T. G. and Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.

7. Podaci društva Oxoid Ltd.

8. Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumophila* subsp. *fraseri* subsp. nov., and *L. pneumophila* subsp. *pasculli* subsp. nov. Journal of Clinical Microbiology; 1988; 26: 1695-1703.

9. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires’ Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710-716.

10. Lück, C et al. Isolation of a *Legionella pneumophila* Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35-39.

11. Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of *Legionella* species. Molecular Identification and Epidemiology; 1998; 47: 667-678.

12. Lück, P. C et al. DNA Polymorphisms in Strains of *Legionella pneumophila* Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. Applied and Environmental Microbiology; 1995; 61: 2000-2003.

13. Helbig, J. H. 2003. Privatna komunikacija.

15. OBJAŠNJENJE SIMBOLA

	Kataloški broj
	Medicinski proizvodi za dijagnostiku in vitro
	Pogledajte upute za upotrebu
	Ograničenja temperature (Temp. skladištenja)
	Sadržava dovoljnu količinu za <N> testova
	Brzi test namijenjen isključivo za upotrebu u laboratorijskom okruženju
	Kôd šarže (Serijski broj)
	Upotrijebiti do (Rok trajanja)
	Ovlašteni zastupnik u Europskoj zajednici
	Proizvođač

IFU X5058F Revidirano u Veljača 2025.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

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14. IRODALOMJEGYZÉK

- 1. Sedgwick, A. K. and Tilton, R. C. (1983). J. Clin. Microbiol., 17: 365–368.
- 2. Ciesielski, C. A., Blaser, M. J. and Wang, W. L. (1986). Infect. Immun., 51: 397–404.
- 3. Dennis, P. J. L. (1988). Isolation of *Legionellae* from Environmental Specimens p. 31–44. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.
- 4. Dournon, E. (1988). Isolation of *Legionellae* from Clinical Specimen p.13–30. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.
- 5. Cowan and Steel's Manual for the identification of Medical Bacteria 3rd Ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993) p. 161–163. University Press, London.
- 6. Harrison, T. G. and Taylor, A. G. (1988). Identification of *Legionellae* by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for *Legionella*. John Wiley & Sons Ltd., Chichester, UK.
- 7. Data on file Oxoid Ltd.
- 8. Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumonia* subsp. *fraseri* subsp. nov., and *L. pneumonia* subsp. *pascullei* subsp. nov. Journal of Clinical Microbiology; 1988; 26: 1695-1703.
- 9. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires' Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710-716.
- 10. Lüch, C et al. Isolation of a *Legionella pneumophila* Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35-39.
- 11. Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of *Legionella* species. Molecular Identification and Epidemiology; 1998; 47: 667-678.
- 12. Lüch, P. C et al. DNA Polymorphisms in Strains of *Legionella pneumophila* Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. Applied and Environmental Microbiology; 1995; 61: 2000-2003.
- 13. Helbig, J. H. 2003. Personal Communication.

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14. BIBLIOGRAFIA

- Sedgwick, A. K. and Tilton, R. C. (1983). J. Clin. Microbiol., 17: 365–368.
- Ciesielski, C. A., Blaser, M. J. and Wang, W. L. (1986). Infect. Immun., 51: 397–404.
- Dennis, P. J. L. (1988). Isolation of Legionellae from Environmental Specimens p. 31–44. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
- Dournon, E. (1988). Isolation of Legionellae from Clinical Specimen p.13–30. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
- Cowan and Steel's Manual for the identification of Medical Bacteria 3rd Ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993) p. 161–163. University Press, London.
- Harrison, T. G. and Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, UK.
- Data on file Oxoid Ltd.
- Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of L. pneumonia subsp. pneumophila subsp. nov., L. pneumophila subsp. fraseri subsp. nov., and L. pneumophila subsp. pascullei subsp. nov. Journal of Clinical Microbiology; 1988; 26: 1695-1703.
- Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires' Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710-716.
- Lück, C et al. Isolation of a Legionella pneumophila Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35-39.
- Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. Molecular Identification and Epidemiology; 1998; 47: 667-678.
- Lück, P. C et al. DNA Polymorphisms in Strains of Legionella pneumophila Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. Applied and Environmental Microbiology; 1995; 61: 2000-2003.
- Helbig, J. H. 2003. Personal Communication.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Regno Unito

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14. BIBLIOGRAFIA

- Sedgwick, A. K. i Tilton, R. C. (1983). J. Clin. Microbiol., 17: 365–368.
- Ciesielski, C. A., Blaser, M. J. i Wang, W. L. (1986). Infect. Immun., 51: 397–404.
- Dennis, P. J. L. (1988). Isolation of Legionellae from Environmental Specimens str. 31–44. Harrison, T. G. i Taylor, A. G. (red.). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, Wielka Brytania.
- Dournon, E. (1988). Isolation of Legionellae from Clinical Specimen str. 13–30. Harrison, T. G. i Taylor, A. G. (red.). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, Wielka Brytania.
- Cowan and Steel’s Manual for the identification of Medical Bacteria 3rd Ed. Barrow, G. I. i Feltham, R. K. A. (red.) (1993) str. 161–163. University Press, London.
- Harrison, T. G. i Taylor, A. G. (1988). Identification of Legionellae by Serological Methods. Harrison, T. G. i Taylor, A. G. (red.). A Laboratory Manual for Legionella. John Wiley & Sons Ltd., Chichester, Wielka Brytania.
- Dane dostępne w firmie Oxoid Ltd.
- Brenner et al. Legionella pneumophila Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia i Descriptions of L. pneumonia subsp. pneumophila subsp. nov., L. pneumophila subsp. fraseri subsp. nov. oraz L. pneumophila subsp. pascullei subsp. nov. Journal of Clinical Microbiology; 1988; 26: 1695–1703.
- Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires’ Disease: Distribution of Legionella pneumophila Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710–716.
- Lück, C et al. Isolation of a Legionella pneumophila Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35-39.
- Fry, N. K i Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. Molecular Identification and Epidemiology; 1998; 47: 667–678.
- Lück, P. C et al. DNA Polymorphisms in Strains of Legionella pneumophila Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. Applied and Environmental Microbiology; 1995; 61: 2000–2003.
- Helbig, J. H. 2003. Personal Communication.

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14. REFERÊNCIAS

- Sedgwick, A. K. and Tilton, R. C. (1983). *J. Clin. Microbiol.*, 17: 365–368.
- Ciesielski, C. A., Blaser, M. J. and Wang, W. L. (1986). *Infect. Immun.*, 51: 397–404.
- Dennis, P. J. L. (1988). Isolation of *Legionellae* from Environmental Specimens p. 31–44. In Harrison, T. G. and Taylor, A. G. (eds). *A Laboratory Manual for Legionella*. John Wiley & Sons Ltd., Chichester, UK.
- Dournon, E. (1988). Isolation of *Legionellae* from Clinical Specimen p.13–30. In Harrison, T. G. and Taylor, A. G. (eds). *A Laboratory Manual for Legionella*. John Wiley & Sons Ltd., Chichester, UK.
- Cowan and Steel's Manual for the identification of Medical Bacteria 3.^a Ed. Barrow, G. I. and Feltham, R. K. A. (eds) (1993) p. 161–163. University Press, London.
- Harrison, T. G. and Taylor, A. G. (1988). Identification of *Legionellae* by Serological Methods. In Harrison, T. G. and Taylor, A. G. (eds). *A Laboratory Manual for Legionella*. John Wiley & Sons Ltd., Chichester, UK.
- Dados guardados de Oxoid Ltd.
- Brenner et al. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumonia* subsp. *fraseri* subsp. nov., and *L. pneumonia* subsp. *pascullei* subsp. nov. *Journal of Clinical Microbiology*; 1988; 26: 1695-1703.
- Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires' Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. *European Journal of Clinical Microbiology and Infectious Disease*; 2002; 21: 710-716.
- Lück, C et al. Isolation of a *Legionella pneumophila* Strain Serologically Distinguishable from all known Serogroups. *Zentralblatt fuer Bakteriologie*. 1995; 282: 35-39.
- Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of *Legionella* species. *Molecular Identification and Epidemiology*; 1998; 47: 667-678.
- Lück, P. C et al. DNA Polymorphisms in Strains of *Legionella pneumophila* Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. *Applied and Environmental Microbiology*; 1995; 61: 2000-2003.
- Helbig, J. H. 2003. Personal Communication.

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7. Údaje v súbore Oxoid Ltd.
8. Brenner a kol. *Legionella pneumophila* Serogroup Lansing 3 Isolated from a Patient with Fatal Pneumonia, and Descriptions of *L. pneumonia* subsp. *pneumophila* subsp. nov., *L. pneumonia* subsp. *fraseri* subsp. nov., and *L. pneumonia* subsp. *pascuelli* subsp. nov. Journal of Clinical Microbiology; 1988; 26: 1695 - 1703.
9. Helbig, J. H a kol. Pan-European Study on Culture Proven Legionnaires’ Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710 - 716.
10. Lück, C a kol. Isolation of a *Legionella pneumophila* Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35-39.
11. Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of *Legionella* species. Molecular Identification and Epidemiology; 1998; 47: 667 - 678.
12. Lück, P. Ca kol. DNA Polymorphisms in Strains of *Legionella pneumophila* Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. Applied and Environmental Microbiology; 1995; 61: 2000 – 2003.
13. Helbig, J. H. 2003. Personal Communication.

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- 9. Helbig, J. H et al. Pan-European Study on Culture Proven Legionnaires’ Disease: Distribution of *Legionella pneumophila* Serogroups and Monoclonal Subgroups. European Journal of Clinical Microbiology and Infectious Disease; 2002; 21: 710-716.
- 10. Lück, C et al. Isolation of a Legionella pneumophila Strain Serologically Distinguishable from all known Serogroups. Zentralblatt fuer Bakteriologie. 1995; 282: 35-39.
- 11. Fry, N. K and Harrison, T. G. An evaluation of intergenic rRNA gene sequence length polymorphism analysis for the identification of Legionella species. Molecular Identification and Epidemiology; 1998; 47: 667-678.
- 12. Lück, P. C et al. DNA Polymorphisms in Strains of Legionella pneumophila Serogroups 3 and 4 Detected by Macrorestriction Analysis and Their Use for Epidemiological Investigation of Nosocomial Legionellosis. Applied and Environmental Microbiology; 1995; 61: 2000-2003.
- 13. Helbig, J. H. 2003. Personlig kommunikation.

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- 7. Registrerade data hos Oxoid Ltd.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien