Principle of Method SCIENTIFIC

The CO₂Gen (CD0020C and CD0025A) are gas generating sachets suitable for the generation of optimal atmospheric conditions for microorganisms requiring an enriched CO₂ environment when used with the related plastic pouches, boxes or jars, that are supplied separately.

The paper sachets contain ascorbic acid which reacts with air to produce an atmosphere containing approximately 6% carbon dioxide, while the final concentration of oxygen will diminish accordingly. The expected percentage of CO₂ after the reaction do not vary depending on the sachets size and/or the accessories part of the AnaeroGen™ range they are used with. For both CO₂Gen Compact Sachet (CD0020C) and CO₂Gen™ Sachet (CD0025A) the expected percentage of CO₂ after reaction is around 6%.

CO₂Gen™ Compact Sachet & CO₂Gen™ Sachet

REF	CD0020C.....	20
REF	CD0025A.....	10

Intended Use

CO₂Gen sachets (CD0020C, CD0025A) are accessory devices intended to be used with the *in vitro* diagnostic devices part of the CO₂Gen products range, for the generation of optimal atmospheric conditions for microorganisms requiring an enriched carbon dioxide atmosphere.

These accessories are for professional use only, are not automated and nor are they companion diagnostics.

Summary and Explanation

The different metabolisms and oxygen needs represent one of the main differences between bacterial species and one of the main categorisations used in microbiology. In particular, different oxygen requirements represent different mechanisms that bacteria use to satisfy their energy needs, and basing on different oxygen requirements, microorganisms can be classified into categories that include aerobes and anaerobes (obligate or facultative), as well as microaerophilic and capnophilic. Aerobic and anaerobic bacteria are characterised by the requirement for the presence or absence of oxygen respectively in order to grow, and the prerequisite can be mandatory (= obligate) or optional (= facultative). Additionally, microaerophilic organisms require low concentration of oxygen to be able to grow, while capnophilic microorganisms need high concentration of carbon dioxide¹.

The most important microorganisms classified as capnophilic bacteria include clinically relevant species responsible for pathogenic conditions in humans, such as *Neisseria gonorrhoeae* and *Haemophilus influenzae*. The *N. gonorrhoeae* bacterium is responsible for a sexually transmitted disease called gonorrhea, that starts with an infection affecting genital tissues in both men and women, and it is in most cases asymptomatic; however, it can lead to permanent health problems if left untreated². *H. influenzae* is the causing agent for many mild infections, such as ear infections, but they can easily turn into more severe complications when affecting immunocompromised individuals or if not treated promptly³.

Having gas generating systems that offer effective performance in creating the optimal atmospheric conditions for the growth of microorganisms requiring a carbon-dioxide enriched atmosphere is therefore essential, given the clinical significance of several capnophilic species that can pose great risks to health when affecting immunocompromised patients. Therefore, it is extremely beneficial to have a simple and easy-to-use system to generate the ideal CO₂ enriched conditions necessary to grow capnophilic microorganisms such as CO₂Gen Compact Sachet (CD0020C) and CO₂Gen Sachet (CD0025A), assisting the clinicians' workflow, reducing costs and, ultimately, giving results more efficiently, thus allowing early diagnosis and prompt treatments for patients.

Description of Components

CO₂Gen Compact Sachet (CD0020C) is a pack of 20 sachets to be used with W-Zip seal pouches (with integral seal) or plastic pouches with clips used for the generation of optimal atmospheric conditions for microorganisms requiring an enriched carbon dioxide atmosphere.

CO₂Gen Sachet (CD0025A) is a pack of 10 sachets to be used with the 2.5L jar or box supplied as an accessory for the CO₂Gen range for the generation of optimal atmospheric conditions for microorganisms requiring an enriched carbon dioxide atmosphere.

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms
- AnaeroGen Compact Sealing Clip (AN0005C) A pack of 5 sealing clips intended to be used with the pouches included in AnaeroGen Compact (AN0010C).
- AnaeroGen Compact Plastic Pouch (AG0020C) A pack of 20 plastic pouches (without clips).
- AnaeroGen W-Zip Seal Pouches (AG0060C) A pack of 20 W-Zip seal pouches (with integral seal).
- AnaeroJar™ 2.5L Jar (AG0025A) A 2.5L jar supplied with 2.5L capacity polycarbonate jar base and lid secured to base by four self-venting clips, intended to be used with AnaeroGen 2.5L Sachet (AN0025A) CampyGen 2.5L Sachet (CN0025A), or CO₂Gen Sachet (CD0025A).
- Anaerobox™ 2.5L box with lid and seal (AB0025A). A 2.5L box intended to be used with AnaeroGen 2.5L Sachet (AN0025A), CampyGen 2.5L Sachet (CN0025A), or CO₂Gen Sachet (CD0025A).

Storage

- Store product at 2 - 25°C.
- The product may be used until the expiry date stated on the label.

Warnings and Precautions

- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their

nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

- As soon as the CO₂Gen™ sachet is exposed to the air, the reaction will start. It is therefore essential that the paper sachet is placed in the pouch and sealed within one minute.
- The reaction of the ascorbic acid with oxygen is exothermic. However, the temperature of the CO₂Gen™ sachet will not exceed 65°C.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Procedure

1. For the CO₂Gen Compact Sachet (CD0020C) place 2 inoculated media plates in the plastic pouches (AG0020 or AG0060).

NOTE: If only one plate is to be inoculated, an uninoculated plate should also be placed in the plastic pouch.

Or

For the CO₂Gen Sachet (CD0025A) place up to 12 inoculated media plates onto the plate carrier and placed into the Anaerobar 2.5L (AG0025) or Anaerobic Box (AB0025A).

Disposable plastic Petri dishes should be of the vented variety to aid gas transfer between the interior and exterior of the plates.

2. Tear open an CO₂Gen foil sachet at the tear-nick indicated and remove the CO₂Gen sachet from within.
3. Immediately place the CO₂Gen sachet in the plastic pouch or plate carrier.

NOTE: The CO₂Gen sachet will become warm to touch on exposure to air.

4. Expel excess air from the plastic pouch. Seal the plastic pouch immediately with the Anaerobar Compact Clip (AN0005C).

OR

Close the jar or box lid immediately.

NOTE: The time taken between opening the foil sachet and sealing the plastic pouch should not exceed one minute. Extended exposure will result in loss of reactivity.

5. Incubate appropriately.
6. After the incubation period remove the plates or ID panel and examine for the presence of colonies. If the plates require re-incubation, then a fresh CO₂Gen Compact sachet must be used following steps 2-5 described above.

NOTE: The plates may be initially inspected through the transparent plastic pouch. If the pouch is opened, a fresh CO₂Gen Compact sachet is required for re-incubation.

7. After incubation, the exhausted CO₂Gen Compact sachet and plastic pouch or CO₂Gen and W-zip seal pouch should be sterilised and discarded with the non-hazardous laboratory waste.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the

product, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The user should periodically check their system for its ability to provide adequate conditions for the growth of appropriate bacteria.

Performance Characteristics


Evidence of the repeatability and reproducibility of the CO₂Gen (CD0020C, CD0025A) are demonstrated by the performance data generated from In-Process QC testing performed during the manufacturing of the anaerobic gas generating sachets and by the historical performance achieved through routine QC assessed by the manufacturer.

Bibliography

1. MicrobeOnline. (2022). Oxygen Requirements for Pathogenic Bacteria • Microbe Online. Accessed on 14.09.2022. <https://microbeonline.com/oxygen-requirements-for-pathogenic-bacteria/>
2. CDC. (2022). Detailed STD Facts - Gonorrhea. Accessed on 14.09.2022. <https://www.cdc.gov/std/gonorrhea/stdfact-gonorrhea-detailed.htm>
3. Centre for Disease Control and Prevention. (2022). Haemophilus influenzae: Types of Infection | CDC. Accessed on 14.09.2022. <https://www.cdc.gov/hi-disease/about/types-infection.html>

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment

	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in Japan

©2022 Thermo Fisher Scientific Inc. All rights reserved.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke,
RG24 8PW, United Kingdom

For technical assistance please contact your local distributor.

Revision Information

Version	Date of modifications introduced
2.0	2023-12-08

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© 2022 Thermo Fisher Scientific Inc. Med ensamrätt. Alla
övriga varumärken tillhör Thermo Fisher Scientific Inc. och