www.thermofisher.com

Wilkins-Chalgren Anaerobe Agar (Dehydrated)

REF CM0619B, CM0619T & CM0619W

Intended Use^{IVD}

Wilkins-Chalgren Anaerobe Agar (Dehydrated) (CM0619B, CM0619T and CM0619W) devices are media for the growth of anaerobic microorganisms.

The devices are for professional use only, are not automated, nor are they companion diagnostics.

Summary and Explanation

Anaerobic bacteria are important pathogens that can cause a variety of infections in humans. The site of anaerobic infection is commonly the site of normal colonization. The spectrum of infections ranges from superficial abscesses to life-threatening infections. Anaerobic microorganisms, including pathogens such as *Clostridium perfringens*, *Bacteroides fragilis* and *Fusobacterium nucleatum*, can be highly infectious and can cause life-threatening diseases. *Fusobacterium necrophorum* has been linked to infections such as necrobacillosis and Lemierre's disease.¹ Anaerobes are endogenous in nature, as they are a primary component of the flora of mucus membranes and gastrointestinal tract. However, they are often difficult to isolate and recognise due to their fastidious nature and overgrowth of facultative anaerobes such as coliforms. This can lead to delays in diagnosis, and therefore delays in determining appropriate therapy which may lead to clinical failures.² Selective supplements (e.g. Neomycin Selective Supplement (SR0163H)) can be added to the medium as required.

Principle of Method

Wilkins Chalgren Agar contains peptones derived from the single protein sources casein and gelatin to improve standardisation of the medium. Salt is present to maintain the osmotic equilibrium and agar as a solidifying agent. Yeast extract provides vitamins and other growth factors such as purines and pyrimidines that are necessary for good growth of *Peptostreptococcus anaerobius* and *Prevotella melaninogenica*. Arginine ensures that sufficient amino acid is available for the growth of *Eubacterium lentum*. Pyruvate is an energy source for asaccharolytic cocci such as *Veillonella* and also acts similarly to catalase, degrading traces of hydrogen peroxide. Hydrogen peroxide may be produced by the action of molecular oxygen on medium constituents and, without pyruvate, can interfere with the metabolism of anaerobes. Haemin is essential for the growth of *Bacteroides* species and menadione for *Prevotella melaninogenica*. Horse blood is added to supply extra nutrients.

Typical Formula

	grams per litre
Tryptone	10.0
Gelatin peptone	10.0
Yeast extract	5.0
Glucose	1.0
Sodium chloride	5.0
L-Arginine	1.0
Sodium pyruvate	1.0

Thermo
S C I E N T I F I C
0.0005
0.005
10.0

Menadione
Haemin
Agar

Materials Provided

CM0619B: 500g of Wilkins-Chalgren Anaerobe Agar, or
CM0619T: 5kg of Wilkins-Chalgren Anaerobe Agar, or
CM0619W: 7 x 1.85kg packs of Wilkins-Chalgren Anaerobe Agar

500g of dehydrated Wilkins-Chalgren Anaerobe Agar yields approximately 11.6L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Selective supplements (e.g. SR0163H)
- Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water for at least 15 minutes. Get medical attention immediately if symptoms occur.
- If in eyes, rinse cautiously with water for at least 15 minutes. Get medical attention.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, remove to fresh air, get medical attention if symptoms occur.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Strict anaerobic conditions must be maintained in order to produce the optimum recovery for organisms.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to

minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B11, B 17, ID 14, ID 25 and Q 5.

Procedure

Suspend 43g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilise by autoclaving at 121°C for 15 minutes. Cool to 50°C. Mix well and pour into sterile containers.

Interpretation

Once the medium is reconstituted, *Clostridium* spp. and *Bacteroides fragilis* will show as straw colonies.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 hours @ 37°C, anaerobic

Positive Controls		
Inoculum level:	10-100	cfu
Colony count is ≥ 50% of the control medium count		
<i>Clostridium perfringens</i> ATCC® 13124™	1-4mm irregular, straw colonies	
<i>Clostridium perfringens</i> NCTC 10613	1-4mm irregular, straw colonies	
<i>Clostridium tetani</i> ATCC® 9441™	0.25-1.5 mm straw colonies	
<i>Clostridium sporogenes</i> ATCC® 19404™	0.25-1.5 mm straw colonies	
<i>Clostridium sphenoides</i> ATCC® 19403™	0.25-1.5 mm straw colonies	
<i>Bacteroides fragilis</i> ATCC® 25285™	0.25-1.5 mm straw colonies	

Limitations

Since Wilkins-Chalgren Anaerobe Agar is not a selective medium, it allows the growth of obligate and facultative anaerobes alike. Some anaerobes with specific growth requirements may not grow on this medium. Strict anaerobic conditions must be maintained in order to produce the optimum recovery for organisms. Some anaerobes with specific growth requirements may not grow on this medium, so it should not be used as the sole isolation medium for anaerobic microorganisms.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of anaerobic microorganism strains is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Wilkins-Chalgren Anaerobe Agar (Dehydrated) (CM0619) was demonstrated by an overall pass rate of 100% obtained for the product over 5 years of testing (03.08.2017 – 16.03.2022; 10 batches). This shows that the performance is reproducible.

Wilkins-Chalgren Anaerobe Agar (CM0619B/T/W) device is tested in-house as part of the QC process since the products were launched in 1997. For target organisms, when using 10-100 cfu of *Bacteroides fragilis*, *Clostridium perfringens*, *Clostridium tetani*, *Clostridium sporogenes*, and *Clostridium sphenoides* and incubating the device anaerobically at 37°C for 18 hours, the user can recover organisms with colony size and morphology as listed in this document.

Bibliography

1. Public Health England. 2015b. 'Identification of Anaerobic Gram-Negative Rods'. UK Standards for Microbiology Investigations ID 25 (2). <https://www.gov.uk/government/publications/smi-id-25-identification-of-anaerobic-gram-negative-rods>.
2. Brook, I. and Frazier, E. H., 1998. 'Aerobic and anaerobic microbiology of infection after trauma'. American Journal of Emergency Medicine 16(6):585-591. [https://doi.org/10.1016/S0735-6757\(98\)90225-X](https://doi.org/10.1016/S0735-6757(98)90225-X)

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use

	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.
CLSI is a trademark of the Clinical Laboratory and Standards Institute.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Modifications Introduced
1.0	2022-11-01. Original Document

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Bibliographie


- Public Health England. 2015b. 'Identification of Anaerobic Gram-Negative Rods'. UK Standards for Microbiology Investigations ID 25 (2).
<https://www.gov.uk/government/publications/smi-id-25-identification-of-anaerobic-gram-negative-rods>.
- Brook, I. and Frazier, E. H., 1998. 'Aerobic and anaerobic microbiology of infection after trauma'. American Journal of Emergency Medicine 16(6):585-591.
[https://doi.org/10.1016/S0735-6757\(98\)90225-X](https://doi.org/10.1016/S0735-6757(98)90225-X)

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Royaume-Uni

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	<u>Gramm pro Liter</u>
Trypton	10,0
Gelatine-Pepton	10,0
Hefeextrakt	5,0
Glukose	1,0
Natriumchlorid	5,0

L-Arginin	1,0
Natriumpyruvat	1,0
Menadion	0,0005
Hämin	0,005
Agar	10,0

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Reino Unido
©2022 Thermo Fisher Scientific Inc. Todos los derechos reservados.