Sorbitol MacConkey Agar

REF CB0813G, CM0813B & CM0813R

EN

Intended Use

Sorbitol MacConkey Agar (CB0813G, CM0813B/R) devices are selective and differential medium for the detection of *Escherichia coli* O157 in faecal and urine samples.

These devices are intended for use in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having enteric infections.

The devices are for professional use only, are not automated, and are not a companion diagnostic.

Summary and Explanation

Escherichia coli O157 is recognised as a cause of haemorrhagic colitis, an illness characterised by bloody diarrhoea, severe abdominal pain and haemolytic uraemic syndrome (HUS)^{1, 2} and as such, is a significant human pathogen. The intestinal tract of ruminants is the prime reservoir for *E. coli* O157 and other enterohaemorrhagic *E. coli* (EHEC) strains, therefore meat derived from cattle, sheep, goat and deer can be expected to be contaminated. Foods implicated in human illness related to *E. coli* include meats, dairy products, vegetables, salads, apple juice and water³. Sorbitol MacConkey Agar is recommended for the isolation of pathogenic *E. coli* O157. The formulation, based on that described by Rappaport and Henig⁴, is identical to MacConkey Agar No. 3, except that lactose has been replaced with sorbitol. *E. coli* O157 does not ferment sorbitol and, therefore, produces colourless colonies. In contrast, most *E. coli* strains ferment sorbitol and form pink colonies. The efficiency of Sorbitol MacConkey Agar has been confirmed by March and Ratnam⁵. These workers reported a sensitivity of 100% and a specificity of 85%, and recommended the medium as a simple, inexpensive, rapid and reliable means of screening for *E. coli* O157.

Principle of Method

Peptone is present in the medium as a nitrogen and nutrient source. Sodium chloride is added to maintain the osmotic balance and agar is a solidifying agent. Gram-positive microorganisms are inhibited by the combination of bile salts and crystal violet. Neutral red is added as a pH indicator. Sorbitol is a fermentable carbohydrate; sorbitol non-fermenters produce straw/colourless colonies and sorbitol fermenters produce red colonies.

Typical Formula

	grams per litre
Peptone	20.0
Sorbitol	10.0
Bile salts No.3	1.5
Sodium chloride	5.0
Neutral red	0.03
Crystal violet	0.001
Agar	15.0

Physical Characteristics of Reconstituted Media

Colour	Dark Red
Clarity	Clear
Fill weight	19.0 ± 2.0g
pH	7.1 ± 0.2

Materials Provided

- CB0813G: 10kg dehydrated culture media
- CM0813B: 500g dehydrated culture media
- CM0813R: 2.5Kg dehydrated culture media

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dishes
- Supplements (SR0172E)

Storage

- Store product in its original packaging at 10–30°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.

- Do not use the product if there is any visible damage to the packaging or tubes.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any applicable federal, state and local regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 22.

Procedure

Suspend 51.5g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilize by autoclaving at 121°C for 15 minutes. Cool to 50°C, mix well and pour into sterile containers. Alternatively, cool to 50°C and aseptically add the contents of 2 vials of Cefixime-Tellurite supplement (SR0172E) reconstituted as directed. Mix well and pour into sterile containers.

Quality Control

It is the responsibility of the user to perform Quality Control testing, taking into account the intended use of the medium, and in accordance with any applicable local regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 21 ± 3 hours @ 37 ± 2°C

Positive Controls Inoculum level: 10-100 cfu Colony count is ≥ 50% of the control medium count.	
<i>Proteus mirabilis</i> NCTC 10975	0.5 – 2 mm straw colonies, no swarming
Inoculation with mixed culture using a diminishing sweep technique.	
Differentiation between sorbitol positive and sorbitol negative strains shall be comparable to the standard after incubation at 37°C for 24 hours.	
Escherichia coli ATCC® 25922™ (sorbitol positive) and Escherichia coli NCTC 12900 (sorbitol negative).	
Negative Controls Inoculum level: 104-106 cfu	
<i>Staphylococcus aureus</i> ATCC® 25923™	No growth
<i>Staphylococcus aureus</i> ATCC® 6358™	No growth

Testing performed in accordance with ISO11133: 2014- Tested with the addition of Cefixime-Tellurite Supplement SR0172E.

Incubation Conditions: Incubation Conditions: 21 ± 3 hours @ 37 ± 2°C

Positive Controls Inoculum level: 50-120 cfu Colony count is ≥ 70% of the control medium count.	
<i>Escherichia coli</i> NCTC 12900	1-2 mm straw colonies
Inoculum level: >10 ⁴ cfu	
<i>Escherichia coli</i> ATCC® 25922™	No growth or pinpoint to 0.25mm pink colonies
<i>Escherichia coli</i> ATCC® 8739™	No growth or pin point to 0.25mm pink colonies
Negative Controls Inoculum of 10 ⁴ – 10 ⁶ colony forming units (cfu).	

<i>Staphylococcus aureus</i> ATCC® 25923™	No growth
<i>Staphylococcus aureus</i> ATCC® 6538™	No growth
Tested with the addition of Cefixime-Tellurite Supplement SR0172E	
<i>Proteus mirabilis</i> NCTC 10975	No growth
Medium is challenged with $\geq 10^2$ cfu	
<i>Escherichia coli</i> ATCC®11775	No growth to pinpoint to 0.25 mm pink colonies

Limitations

Delay in reading plates beyond 24 hours should be avoided because the colour intensity of sorbitol-fermenting colonies fades, reducing the contrast with non-fermenting colonies. Strains of other organisms that do not ferment sorbitol (such as *Escherichia hermannii*) may grow on Sorbitol MacConkey Agar. Other Gram-negative organisms including *Pseudomonas*, *Proteus* and *Klebsiella* species are able to grow on Sorbitol MacConkey Agar (CM0813B, CM0813R, CB0813G and PO0232A) devices but may generally be differentiated by the appearance of their colonies. *E. coli* O157 may be overgrown by other organisms and the more selective medium, Cefixime-Tellurite Sorbitol MacConkey Agar (CT-SMAC), may be inoculated in parallel. Although most *Escherichia coli* O157 strains are typical in appearance some strains are atypical. Sorbitol MacConkey Agar (CM0813B, CM0813R, CB0813G and PO0232A) devices cannot be used on its own to detect VTEC strains of *Escherichia coli* as some non-toxicogenic strains will not ferment sorbitol. Identification is presumptive and suspect colonies should be confirmed using appropriate biochemical, molecular or serological methods.

Performance

Accuracy has been demonstrated through review of the QC data. Correct detection of *Escherichia coli* O157 is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device. The precision of Sorbitol MacConkey Agar (CM0813) was demonstrated by an overall pass rate of 100% over two years of testing (15-Apr- 2020 – 08-Jun-2022). This shows that the performance is reproducible.



Sorbitol MacConkey Agar (CM813) is tested in-house as part of the QC process since the formulation was developed in 1995. For target organisms, *Escherichia coli* NCTC12900/ ATCC® 700728™ when incubating the device at $37 \pm 2^\circ\text{C}$ for 21 ± 3 hours, the user can recover organisms with colony size and morphology that meets the defined acceptance criteria. For non-target organisms, when incubating the device at $37 \pm 2^\circ\text{C}$ for 21 ± 3 hours, the user can recover organisms with colony size and morphology that meets the defined acceptance criteria

Bibliography

- (1) Karmali M A, Petric M, Lim C, Fleming P C, Arbus G S, Lior H. (1985) *Infect Dis.* 151(5):775-82.
- (2) Karmali M A, Steele B T, Petric M, Lim C. 1983. *Lancet* ; 19(1): 619-620.
- (3) Desmarchelier P M, Grau, F H (1997) *Escherichia coli* in Foodborne microorganisms of public health significance. 5th ed, pp.231–264. A D Hocking (Ed) AIFST (NSW Branch) Food Microbiology Group Australia.
- (4) Rappaport F and Henig E (1952) *J. Clin. Path.* 5. 361-362.
- (5) March S B and Ratnam S (1986) *J. Clin. Microbiol.* 23. 869-872.
- (6) Karmali M A (1988) *Culture* 9. 2.
- (7) Lior H and Borczyk A A (1987) FALSE POSITIVE IDENTIFICATIONS OF ESCHERICHIA COLI 0157. *Lancet.* 329(8528): 333.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date

	Keep away from sunlight
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
 ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
 NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.
 All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, United Kingdom

For technical assistance please contact your local distributor.

Revision Information

Version	Date of issue and modifications introduced
4.0	2024-01-25

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©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien