www.thermofisher.com

Alkaline Peptone Water

REF **CM1028B**

Intended Use

Alkaline Peptone Water (CM1028B) is an enrichment medium for the isolation of *Vibrio* species from faecal samples. Can also be used for testing food and water samples.

Alkaline Peptone Water (CM1028B) device is used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having *Vibrio* species infections

The devices are for professional use only, are not automated, nor are they companion diagnostics.

Summary and Explanation

Vibrio species, members of the Vibrionaceae family, are a well-known genus capable of causing gastrointestinal diseases in humans, with cholera being the most well-known. Medically important *Vibrio* species include *Vibrio cholerae*, the causative agent of cholera, *Vibrio parahaemolyticus*, the causative agent of acute gastroenteritis and, *Vibrio vulnificus*, that causes septicemia. In general, *Vibrio* species are sea dwelling organisms and halophilic conditions are favourable for many species for their growth and survival. They are most frequently isolated from stool, blood, wound infections and vomit¹.

Principle of Method

The 2% (w/v) sodium chloride incorporated in this medium promotes the growth of *Vibrio cholerae*, while the alkalinity of this medium inhibits most of the unwanted background flora.

Typical Formula

	grams per litre
Peptone	10.0
Sodium chloride	20.0

Materials Provided

CM1028B: 500g of Alkaline Peptone Water

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 25°C.

Warnings and Precautions

- If inhaled: remove to fresh air. Get medical attention immediately if symptoms occur.
- If ingested: clean mouth with water and drink afterwards plenty of water. Get medical attention if symptoms occur.
- If contact with skin: wash off immediately with plenty of water for at least 15 minutes. Get medical attention immediately if symptoms occur.
- If contact with eyes: rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Get medical attention.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 1, ID 7, ID 16 and Q 5.

Procedure

Suspend 30g in 1 litre of distilled water and distribute into final containers. Sterilize by autoclaving at 121°C for 15 minutes.

Interpretation

Turbidity indicates bacterial growth.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 hours @ 37°C

Positive Controls

Inoculum level: 10 – 100 cfu

<i>Vibrio cholerae</i> NCTC 11348	Turbid growth
<i>Vibrio parahaemolyticus</i> NCTC 10885	Turbid growth
<i>Vibrio furnissii</i> NCTC 11218	Turbid growth
<i>Vibrio vulnificus</i> ATCC® 29307™	Turbid growth

Limitations

Vibrio spp. are oxidase positive and ferment glucose with the production of acid only. As oxidase testing may lead to false negative results on media containing carbohydrates (such as TCBS) subculture to nutrient or blood agar before testing.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of *Vibrio* species is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Alkaline Peptone Water (CM1028B) devices was demonstrated by an overall pass rate of 100% obtained for the product over one year of testing (13.APR.2021 – 17.AUG.2022; 10 batches). This shows that the performance is reproducible.

Alkaline Peptone Water (CM1028B) device is tested in-house as part of the QC process since the products were launched in 2001. For target organisms, when using 10-100 cfu inoculum of *Vibrio cholerae* NCTC 11348, *Vibrio parahaemolyticus* NCTC 10885, *Vibrio furnissii* NCTC 11218 or *Vibrio vulnificus* ATCC® 29307™ and incubating the device at 37°C for 18 hours, the user can recover organisms with growth as listed in this document.

Bibliography

- Public Health England. 2015. 'Identification of *Vibrio* and *Aeromonas* species. UK Standards for Microbiology Investigations. SMI ID19. Issue no: 3 | Issue date: 14.04.15 | <https://www.gov.uk/government/publications/smi-id-19-identification-of-vibrio-species>

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight

	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Made in the United Kingdom

ATCC Licensed Derivative

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
1.0	2022-11-22 Original Document

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© 2022 Thermo Fisher Scientific Inc. Todos los derechos reservados.
Todas las demás marcas comerciales son propiedad de Thermo Fisher Scientific Inc. y sus filiales.
Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, Reino Unido