Heart Infusion Broth

EN

REF CM1032B

Brain Heart Infusion Broth

REF CM1135B, CM1135K, CM1135R, CM1135T

Brain Heart Infusion

REF BO1230D, BO1230M

Brain Heart Infusion Broth with Horse Serum

REF BO1340E

Intended Use

Brain Heart Infusion Broth are highly nutritious, general purpose media for the growth of fastidious microorganisms from clinical samples (e.g. wounds, genital swabs). These devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The devices are for professional use only, are not automated, nor are they companion diagnostics.

Summary and Explanation

Fastidious microorganisms require special nutrition and environmental conditions to grow but are able to thrive in their ecological niche which can be replicated by some culture media. Unlike non fastidious microorganisms, they cannot synthesise all organic molecules required for growth.¹ Prominent fastidious human pathogens include *Streptococcus pyogenes*, *Streptococcus pneumoniae*, *Neisseria gonorrhoeae*, *Neisseria meningitidis*, *Bacteroides fragilis* and *Haemophilus influenzae*, among others.² These organisms may be associated with severe infection, especially on immunocompromised individuals. Due to the fastidiousness of such microorganisms, there exists a possibility of false negative results, which highlights the importance of culture media being able to grow them.³ Heart Infusion Broth and Brain Heart Infusion Broth devices are therefore key tools to support the clinician for the growth of fastidious microorganisms from clinical samples.

Principle of Method

Heart Infusion Broth and Brain Heart Infusion Broth contain peptones and glucose which form the nutritional sources within the media. Sodium chloride maintains the osmotic equilibrium of the media whilst disodium phosphate is added to maintain the pH level. The formulation for Heart Infusion devices has been developed to give equivalent performance to Brain Heart Infusion devices, whilst carrying a lower regulatory burden as a result of the exclusion of brain products.

Typical Formula

CM1032B

	grams per litre
Beef heart infusion solids	17.5
Protease peptone	10.0
Glucose	2.0
Sodium chloride	5.0
Disodium phosphate	2.5

*Adjusted as required to meet performance standards

*CM1135B, CM1135K, CM1135R, CM1135T, BO1230D, BO1230M, BO1340**

	grams per litre
Brain infusion solids	12.5
Beef heart infusion solids	5.0
Protease peptone	10.0
Glucose	2.0
Sodium chloride	5.0
Disodium phosphate	2.5

*Adjusted as required to meet performance standards.

** BO1340E – Contains 100 mL Horse Serum

Materials Provided

Product trade name	Product Code	
Heart Infusion Broth	CM1032B	500 g dehydrated
Brain Heart Infusion Broth	CM1135B	500 g dehydrated

Brain Heart Infusion Broth	CM1135R	2.5 g dehydrated
Brain Heart Infusion Broth	CM1135T	5 kg dehydrated
Brain Heart Infusion Broth	CM1135K	25 kg dehydrated
Brain Heart Infusion	BO1230D	24 x 30 mL universal bottles
Brain Heart Infusion	BO1230M	10 x 125 mL sirop bottles
B.H.I. Broth with Horse Serum	BO1340E	24 x 30 mL universal bottles

500 g of dehydrated broth powder yields approximately 13.5 litres after reconstitution.

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

Storage

- Store CM1032B, CM1135B, CM1135K, CM1135R, CM1135T between 10°C and 30°C.
- Store BO1230D, BO1230M between 2°C and 25°C.
- Store BO1340E between 2°C and 10°C.
- Keep out of direct sunlight.
- Keep upright whenever possible.

Warnings and Precautions

- Do not inhale.
- Causes serious eye irritation.
- May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove victim to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Materials of Animal Origin

The devices contain yeast extract manufactured from microbial raw materials, peptone manufactured from bovine and porcine raw materials, brain substance from porcine sources (CM1135 devices only), and serum from equine sources (BO1340 only). Country of Origin details are available for the bovine, porcine, and microbial raw materials.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local the recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 2, B 28 and B 29^{4,5,6}.

Procedure for CM1135 and CM1032

1. Add 37 g to 1 litre of distilled water.
2. Heat to dissolve.
3. Cool, mix well and distribute into final containers.
4. Sterilize by autoclaving at 121° for 15 minutes.

Procedure for BO1230D, BO1230M

Refer to any local protocols and guidelines.

- Inoculate the specimen directly into the broth.
- Incubate in racks, aerobically at 37°C ± 2°C for 22 to 26 hours and 48 hours.
- Examine for the presence of turbid growth.

Procedure for BO1340E

Refer to any local protocols and guidelines.

1. Inoculate the specimen directly into the broth.
2. Incubate in racks, aerobically at 37°C ± 2°C for 18 to 24 hours and 30°C ± 2°C for 48 hours.
3. Examine for the presence of turbid growth.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).
The performance of this medium can be verified by testing the following reference strains.

Heart Infusion Broth CM1032B and Brain Heart Infusion Broth CM1135B, CM1135K, CM1135R, and CM1135T

Incubation Conditions: 37 ± 2°C for 18 hours under aerobic conditions Inoculum level: 10-100 cfu	
<i>Streptococcus pyogenes</i> ATCC® 19615™	Turbid Growth
<i>Streptococcus pneumoniae</i> ATCC® 6303™ and ATCC® 6305™	Turbid Growth
<i>Enterococcus faecalis</i> ATCC® 19433™	Turbid Growth
<i>Pseudomonas aeruginosa</i> ATCC® 27853™	Turbid Growth
Incubation Conditions: 37 ± 2°C for 18 hours under anaerobic conditions Inoculum level: 10-100 cfu	
<i>Streptococcus pneumoniae</i> ATCC® 6305™	Turbid Growth
Incubation Conditions: 37 ± 2°C for 48 hours under aerobic conditions Inoculum level: 10-100 cfu	
<i>Candida albicans</i> ATCC® 10231™	Turbid Growth
Incubation Conditions: Addition of 0.1% w/v bacteriological agar and incubation at 37 ± 2°C for 48-72 hours under aerobic conditions Inoculum level: 10-100 cfu	
<i>Bacteroides fragilis</i> ATCC® 25285™	Turbid Growth
Incubation Conditions: Enrichment of 1 x 10 mL volume with two X factor discs (DD0003) and two V factor discs (DD0004) and aerobic incubation at 37 ± 2°C for 24 hours Inoculum level: 10-100 cfu	
<i>Haemophilus influenzae</i> ATCC® 33391™	Turbid Growth
Incubation Conditions: Enrichment with 10% v/v horse blood, aerobic incubation at 37 ± 2°C for 18 hours (CM1032) or 48 hours (CM1135), and subculture onto chocolate agar plates and aerobic incubation at 37 ± 2°C for 48 hours (CM1032) or 24-48 hours (CM1135) in CO2 atmosphere Inoculum level: 5-10 cfu	
<i>Nisseria meningitidis</i> ATCC® 13077™	1-2 mm grey/brown colonies
<i>Nisseria gonorrhoeae</i> NCTC 11148 and ATCC® 19424™	
Incubation Conditions: Tube coagulase test - aerobic incubation at 37 ± 2°C for four hours Inoculum level: 10 ⁴ -10 ⁶ cfu	
<i>Staphylococcus aureus</i> ATCC® 9144™	Coagulase positive
<i>Staphylococcus epidermidis</i> ATCC® 14990™	Coagulase negative
Incubation Conditions: Testing performed in accordance with current CLSI M22 A Inoculum level: 10-100 cfu	
<i>Staphylococcus aureus</i> ATCC® 25923™	Turbid Growth
<i>Escherichia coli</i> ATCC® 25922™	Turbid Growth
Incubation Conditions: Testing in accordance with ISO 11133:2014 - Incubation at 37 ± 2°C for 24 ± 2 hours (CM1135 only). Inoculum level: 10-100 cfu	
<i>Staphylococcus aureus</i> ATCC® 25923™	Turbid Growth

Brain Heart Infusion BO1230D, BO1230M

Incubation Conditions: Testing in accordance with ISO 11133 - Incubation at 35-39°C for 22-26 hours Inoculum level: ≤ 100 cfu	
<i>Staphylococcus aureus</i> ATCC® 25923™	Turbid Growth
Incubation Conditions: 35-39°C for 22-26 hours under aerobic conditions Inoculum level: ≤ 100 cfu	
<i>Streptococcus pneumoniae</i> ATCC® 6305™	Turbid Growth
Incubation Conditions: 35-39°C for 48 hours under aerobic conditions Inoculum level: ≤ 100 cfu	
<i>Candida albicans</i> ATCC® 10231™	Turbid Growth

B.H.I Broth with Horse Serum (BO1340E)

Incubation Conditions: 35-39°C for 18-24 hours under aerobic conditions Inoculum level: 10 – 100 cfu	
<i>Streptococcus pneumoniae</i> ATCC® 6305™	Turbid Growth

Incubation Conditions: 35-39°C for 36-48 hours under aerobic conditions Inoculum level: 10-100 cfu	
<i>Haemophilus influenzae</i> ATCC® 35056™	Turbid Growth
<i>Neisseria gonorrhoeae</i> ATCC® 49226™	Turbid Growth
Incubation Conditions: 28-32°C for 48 hours under aerobic conditions Inoculum level: 10-100 cfu	
<i>Candida albicans</i> ATCC® 10231™	Turbid Growth

Limitations

Atypical organisms may give anomalous reactions. A small number of atypical strains may give a weak reaction or fail to grow, especially when low numbers are present in the sample.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct growth of microorganisms is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the devices.

The precision of Heart Infusion Broth (CM1032B), Brain Heart Infusion Broth (CM1135B, CM1135R, CM1135T, CM1135K), Brain Heart Infusion (BO1230D and BO1230M) and B.H.I Broth with Horse Serum (BO1340E) devices was demonstrated by an overall pass rate of 100% obtained for the product over three years (CM1032: 19.06.2018 – 31.01.2022; 10 batches), 11 months (CM1135: 01.04.2021 – 02.03.2022; 10 batches), six months (BO1230D: 05.07.2021 – 31.01.2022; 10 batches), two years (BO1230M: 23.03.2020 – 21.02.2022; 10 batches) and a year (BO1340E: 09.11.2020 – 20.12.2021; 10 batches) of testing. This shows that the performance is reproducible.


These devices have been tested in-house as part of the QC process since they were launched. This data shows that when using target organisms and following manufacturer instructions, users can recover organisms with growth, colony size and/or morphology that meets the defined acceptance criteria.

Bibliography

1. Haden, RL. Elective localization in the eye of bacteria from infected teeth. *Arch Intern Med.* American Medical Association; 1923 Dec 1;32(6):828.
2. Hitchens AP. Advantages of Culture Media Containing Small Percentage of Agar. *J Infect Dis.* Oxford University Press; 1921;29(4):390-407.
3. Eley A, Greenwood D, O'Grady F. Comparative growth of *Bacteroides* species in various anaerobic culture media. *J Med Microbiol.* 1985 Apr;19(2):195-201.
4. Public Health England. 2017. 'Investigation of bacterial eye infections'. UK Standards for Microbiology Investigations. B 2 Issue 6.1. <https://www.gov.uk/government/publications/smi-b-2-investigation-of-eye-swabs-and-canalicular-pus>.
5. Public Health England. 2017a. 'Investigation of genital tract and associated specimens'. UK Standards for Microbiology Investigations. B 28 Issue 4.6. <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
6. Public Health England 2019. 'Investigation of infectious causes of dyspepsia'. UK Standards for Microbiology Investigations. B 55 Issue 7. <https://www.gov.uk/government/publications/smi-b-55-investigation-of-gastric-biopsies-for-helicobacter-pylori>.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limitation (storage temp.)
	Use by (expiration date) YYYY-MM
	Keep away from sunlight
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Do not use if packaging damaged

	Authorized representative in the European Community
	Manufacturer
	European Conformity Assessment
	UK Conformity Assessment
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2023 Thermo Fisher Scientific Inc. All rights reserved.

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

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- Haden, RL. Elective localization in the eye of bacteria from infected teeth. *Arch Intern Med.* American Medical Association; 1923 Dec 1;32(6):828.
- Hitchens AP. Advantages of Culture Media Containing Small Percentage of Agar. *J Infect Dis.* Oxford University Press; 1921;29(4):390-407.
- Eley A, Greenwood D, O'Grady F. Comparative growth of Bacteroides species in various anaerobic culture media. *J Med Microbiol.* 1985 Apr;19(2):195-201.
- Public Health England. 2017. 'Investigation of bacterial eye infections'. UK Standards for Microbiology Investigations. B 2 Issue 6.1. <https://www.gov.uk/government/publications/smi-b-2-investigation-of-eye-swabs-and-canalicular-pus>.
- Public Health England. 2017a. 'Investigation of genital tract and associated specimens'. UK Standards for Microbiology Investigations. B 28 Issue 4.6. <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
- Public Health England 2019. 'Investigation of infectious causes of dyspepsia'. UK Standards for Microbiology Investigations. B 55 Issue 7. <https://www.gov.uk/government/publications/smi-b-55-investigation-of-gastric-biopsies-for-helicobacter-pylori>.

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Literatura

1. Haden, RL. Elective localization in the eye of bacteria from infected teeth. *Arch Intern Med.* American Medical Association; 1923 Dec 1;32(6):828.
2. Hitchens AP. Advantages of Culture Media Containing Small Percentage of Agar. *J Infect Dis.* Oxford University Press; 1921;29(4):390-407.
3. Eley A, Greenwood D, O'Grady F. Comparative growth of Bacteroides species in various anaerobic culture media. *J Med Microbiol.* 1985 Apr;19(2):195-201.
4. Public Health England. 2017. 'Investigation of bacterial eye infections'. UK Standards for Microbiology Investigations. B 2 Issue 6.1. <https://www.gov.uk/government/publications/smi-b-2-investigation-of-eye-swabs-and-canalicular-pus>.
5. Public Health England. 2017a. 'Investigation of genital tract and associated specimens'. UK Standards for Microbiology Investigations. B 28 Issue 4.6. <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
6. Public Health England 2019. 'Investigation of infectious causes of dyspepsia'. UK Standards for Microbiology Investigations. B 55 Issue 7. <https://www.gov.uk/government/publications/smi-b-55-investigation-of-gastric-biopsies-for-helicobacter-pylori>.

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
Literaturverzeichnis

1. Haden, RL. Elective localization in the eye of bacteria from infected teeth. *Arch Intern Med.* American Medical Association; 1923 Dec 1;32(6):828.
2. Hitchens AP. Advantages of Culture Mediums Containing Small Percentage of Agar. *J Infect Dis.* Oxford University Press; 1921;29(4):390-407.
3. Eley A, Greenwood D, O'Grady F. Comparative growth of Bacteroides species in various anaerobic culture media. *J Med Microbiol.* 1985 Apr;19(2):195-201.
4. Public Health England. 2017. 'Investigation of bacterial eye infections'. UK Standards for Microbiology Investigations. B 2 Issue 6.1. <https://www.gov.uk/government/publications/smi-b-2-investigation-of-eye-swabs-and-canalicular-pus>.
5. Public Health England. 2017a. 'Investigation of genital tract and associated specimens'. UK Standards for Microbiology Investigations. B 28 Issue 4.6. <https://www.gov.uk/government/publications/smi-b-28-investigation-of-genital-tract-and-associated-specimens>.
6. Public Health England 2019. 'Investigation of infectious causes of dyspepsia'. UK Standards for Microbiology Investigations. B 55 Issue 7. <https://www.gov.uk/government/publications/smi-b-55-investigation-of-gastric-biopsies-for-helicobacter-pylori>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, VK

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©2023 Thermo Fisher Scientific Inc. Med ensamrätt.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien