	https://www.thermofisher.com	
	EU	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

Oxoid™ Bacitracin Discs 10 µg (B10)

REF CT0005B

Antimicrobial Susceptibility Test Discs

NOTE: This IFU should be read in conjunction with the AST generic instructions for use supplied with the product and available online.

Bacitracin Antimicrobial Susceptibility Test (AST) Discs are 6 mm paper discs which contain specific amounts of the antimicrobial agent bacitracin. The discs are labelled on both sides with details of the antimicrobial (B) and amount present (µg): B10 (10µg).

Discs are supplied in cartridges containing 50 discs. There are 5 cartridges per pack. Each cartridge is individually sealed together with a desiccant tablet in a foil covered transparent blister pack. Bacitracin AST Discs can be dispensed using an Oxoid Disc Dispenser (sold separately). Each individual disc should only be used once.

Intended Use

Bacitracin Antimicrobial Susceptibility Test Discs are used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having a microbial infection, these discs are intended to determine susceptibility against microorganisms for which bacitracin has been shown to be active both clinically and in vitro. To be used with a pure, agar grown culture.

The device is not automated, is for professional use only and is not a companion diagnostic.

The test provides information to categorise organisms as either resistant, intermediate or susceptible to the antimicrobial agent.

Further requirements for specimen collection, handling and storage of samples can be found in local procedures and guidelines. There is no specified testing population for these discs.

Published clinical breakpoints in the current version of the FDA,³ CLSI M100^{1ac} or EUCAST² breakpoint tables must be used to interpret the zone size result.

Species with internally derived breakpoints according to historical data:

Gram-Positive

- *Staphylococcus aureus*

Gram-Negative

- N/A

Each disc is single use only. The pack contains sufficient test devices for multiple single-use tests.

Principle of Method

Oxoid Bacitracin AST Discs may be used in the semi-quantitative agar diffusion test method for in vitro susceptibility testing. For full instructions relating to the generation and interpretation of the results according to CLSI^{1bc}/EUCAST^{2ab} methodology please refer to the relevant current standards. Tables showing CLSI^{1bc}/EUCAST² compound/concentrations can be found in their documentation referenced below. Pure cultures of clinical isolates are inoculated onto the test medium and the AST disc placed on the surface. The antibiotic within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific antimicrobial agent(s)/organism(s) under test.

Metrological Traceability of Calibrator and Control Material Values

Metrological traceability of Bacitracin values assigned to calibrators and control materials intended to establish or verify trueness of a method for Oxoid Bacitracin AST Discs is based on internationally recognized procedures and standards.

For recommended concentrations, the zone limits are in accordance with the current performance standards for antimicrobial disc susceptibility tests as detailed by CLSI^{1bc} and/or EUCAST².

Materials Provided

Bacitracin AST Discs consist of 6 mm diameter paper discs impregnated with a specific amount of antimicrobial agent. The discs are marked on both sides to indicate the agent and amount. Bacitracin AST Discs are supplied in cartridges of 50 discs. There are 5 cartridges in each pack. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

See Table 1 below for a description of components associated with the device. For a description of active reagents that influence the result of the device, refer to Table 2.

Table 1. Materials Provided with CT0005B	
Component Description	Material Description
Cartridge with spring, cap, and plunger (x5)	Assembly components and plastic cartridge containing 50x AST discs.
Desiccant tablet (x5)	Light beige to brown, small lozenge shaped tablets. 1 supplied with each cartridge.
Foil	Foil individually sealing each cartridge with its desiccant.
Susceptibility Test Discs (x250)	Absorbent paper individual discs. 6mm. 50 in each cartridge. 5 cartridges per pack.

Table 2. Description of Bacitracin AST Disc Reagents	
Reagent	Description of Function
Bacitracin	White to pale buff powder. Bacitracin is a cyclic polypeptide antimicrobial agent. Bacitracin exerts its mechanism of action through binding to divalent metal ions and inhibiting the dephosphorylation of C ₅₅ -isoprenyl pyrophosphate, a molecule required for the synthesis of the peptidoglycan cell wall of bacteria. This subsequently inhibits the synthesis of the bacterial cell wall, ultimately leading to cell lysis.

The concentration of antibiotic on the AST disc is analysed for every batch and is controlled using internal and external specifications (e.g. FDA³). The actual concentration is detailed on the Certificate of Analysis.

Bacitracin AST Discs can be dispensed using an Oxoid Disc Dispenser, which is not included with the device.

Shelf-life and Storage Conditions

Unopened cartridges of Bacitracin AST discs have a shelf-life of 36 months if stored under the recommended conditions. Unopened cartridges must be stored at -20°C to 8°C until required.

Once opened, cartridges should be stored within the dispenser in the container provided (with an unsaturated (orange) desiccant), or other suitable opaque airtight container with a desiccant to protect the discs from moisture. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU.

Analytical Performance Characteristics

Table 3. Raw data of zone sizes, in accordance with internally derived breakpoints, taken from the 3 most recent batches of Bacitracin (10 units) B10 Antimicrobial Susceptibility Test Discs (CT0005B).

Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference between the experimental values and reference mid-point values (mm) + SD (Coefficient of variance) (%)
							1	2	3	
3331895	<i>Escherichia coli</i>	25922	MHA	0	0	0	0	0	0	0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	0	0	0	0	0	0	0 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	11	17	14	14	14	14	0 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	0	0	0	0	0	0	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3230859	<i>Escherichia coli</i>	25922	MHA	0	0	0	0	0	0	0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	0	0	0	0	0	0	0 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	11	17	14	14	14	14	0 ± 0 (CV=0)

	<i>Pseudomonas aeruginosa</i>	27853	MHA	0	0	0	0	0	0	0 ± 0 (CV=0)
Product Batch	Organism	ATCC® Number	Test Media	Lower Limit	Upper Limit	Calculated reference mid-point value	Reading (mm)			Mean difference (mm) + SD (CoV)(%)
							1	2	3	
3164316	<i>Escherichia coli</i>	25922	MHA	0	0	0	0	0	0	0 ± 0 (CV=0)
	<i>Haemophilus influenzae</i>	49247	HTM	0	0	0	0	0	0	0 ± 0 (CV=0)
	<i>Staphylococcus aureus</i>	25923	MHA	11	17	14	14	14	14	0 ± 0 (CV=0)
	<i>Pseudomonas aeruginosa</i>	27853	MHA	0	0	0	0	0	0	0 ± 0 (CV=0)

Clinical Performance Characteristics

A cross-sectional study was held to determine the prevalence and antimicrobial susceptibility patterns of *Salmonella* isolates causing typhoid fever among febrile patients in Ethiopia.⁴ Blood and stool samples were collected from febrile patients with enteric fever visiting the Ambo Hospital, Ethiopia, from February to December 2016. Isolates were identified using standard microbiological techniques and antimicrobial susceptibility testing was performed through the disc diffusion method. The susceptibility testing was conducted in accordance with the National Committee for Clinical Laboratory Standards (NCCLS), now CLSI, guidelines for a variety of antibiotics, including the Oxoid Bacitracin AST Discs. The following reference strains were utilised as quality control, *S. aureus* (ATCC® 25923) and *E. coli* (ATCC® 25923).

The study screened a total of 372 febrile patients and stool cultures yielded seven *Salmonella paratyphi* isolates and three *Salmonella typhi* isolates, giving a typhoid fever prevalence of 2.7%. All 10 isolates were found to be resistant to bacitracin, amoxicillin, and erythromycin, and 80.0% of isolates were resistant to cefotaxime and streptomycin. Much higher susceptibility was seen among the isolates to chloramphenicol and ciprofloxacin, with all isolates found to be susceptible to amikacin, ceftriaxone and gentamicin.

This study found the *Salmonella* isolates to be frequently resistant to many routinely used antibiotics, including bacitracin. The study hypothesised that high antimicrobial resistance among typhoid patients may be due to the cheap and wide availability of antibiotics, the increased and indiscriminate prescribing of antimicrobials, high frequency use, patient self-medicating and sub-standard antimicrobial preparation. The study did not find any inconsistent or out of specification results for the Oxoid Bacitracin AST Discs. Therefore, this study demonstrated that the Oxoid Bacitracin AST Discs are an effective method for in vitro analysis of *S. typhi* and *S. paratyphi* and for detecting emerging resistance.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

References

- Clinical and Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
- Deksissa, Tolera, and Endrias Zewdu Gebremedhin. 2019. "A Cross-Sectional Study of Enteric Fever among Febrile Patients at Ambo Hospital: Prevalence, Risk Factors, Comparison of Widal Test and Stool Culture and Antimicrobials Susceptibility Pattern of Isolates." *BMC Infectious Diseases* 19 (1): 288. <https://doi.org/10.1186/s12879-019-3917-3>.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device

	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-13

	https://www.thermofisher.com	
	EC	+800 135 79 135
	US	1 855 236 0910
	CA	1 855 805 8539
	ROW	+31 20 794 7071

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- Clinical and Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method
 - Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
- Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).

4. Deksissa, Tolera, and Endrias Zewdu Gebremedhin. 2019. "A Cross-Sectional Study of Enteric Fever among Febrile Patients at Ambo Hospital: Prevalence, Risk Factors, Comparison of Widal Test and Stool Culture and Antimicrobials Susceptibility Pattern of Isolates." *BMC Infectious Diseases* 19 (1): 288. <https://doi.org/10.1186/s12879-019-3917-3>.

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1. Clinical Laboratory Standards Institute (CLSI).
 - a. Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard – M02 - Latest Edition.
 - b. Performance Standards for Antimicrobial Susceptibility Testing – M100 – Latest Edition.
 - c. Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria; Approved Guideline – M45 – Latest Edition.
2. European Committee on Antimicrobial Susceptibility Testing (EUCAST).
 - a. Antimicrobial Susceptibility Testing, EUCAST disk diffusion method – Latest version.
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing – Latest version.
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Ahman J, Matuschek E, Kahlmeter G. The quality of antimicrobial discs from nine manufacturers-EUCAST evaluations in 2014 and 2017. *Clin Microbiol Infect.* 2019 Mar;25(3):346–52.
5. O'Halloran C, Walsh N, O'Grady MC, Barry L, Hooton C, Corcoran GD, et al. Assessment of the comparability of CLSI, EUCAST and Stokes antimicrobial susceptibility profiles for Escherichia coli uropathogenic isolates. *Br J Biomed Sci.* 2018 Jan;75(1):24–9.

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Referințe

- Clinical and Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.

2. European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - a. EUCAST Disk Diffusion Method
 - b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Deksissa, Tolera, and Endrias Zewdu Gebremedhin. 2019. "A Cross-Sectional Study of Enteric Fever among Febrile Patients at Ambo Hospital: Prevalence, Risk Factors, Comparison of Widal Test and Stool Culture and Antimicrobials Susceptibility Pattern of Isolates." *BMC Infectious Diseases* 19 (1): 288. <https://doi.org/10.1186/s12879-019-3917-3>.

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Referencias

- Clinical and Laboratory Standards Institute (CLSI). Consult latest edition:
 - M02, Performance Standards for Antimicrobial Disk Susceptibility Tests
 - M100, Performance Standards for Antimicrobial Susceptibility Testing
 - M45, Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria.
- European Committee on Antimicrobial Susceptibility Testing (EUCAST). Consult latest version:
 - EUCAST Disk Diffusion Method

- b. Reading guide, EUCAST disk diffusion method for antimicrobial susceptibility testing
3. Food and Drug Administration (FDA). CFR Title 21, Volume 5, Part 460 (2005).
4. Deksissa, Tolera, and Endrias Zewdu Gebremedhin. 2019. "A Cross-Sectional Study of Enteric Fever among Febrile Patients at Ambo Hospital: Prevalence, Risk Factors, Comparison of Widal Test and Stool Culture and Antimicrobials Susceptibility Pattern of Isolates." *BMC Infectious Diseases* 19 (1): 288. <https://doi.org/10.1186/s12879-019-3917-3>.

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