remel

BactiDrop™ Desoxycholate (10%) EN

REF R21508

INTENDED USE

BactiDrop Desoxycholate (10%) is a reagent used in qualitative procedures for the detection of bile solubility for the differentiation of *Streptococcus pneumoniae* from other alpha-hemolytic streptococci. The device is used in a diagnostic workflow to aid clinicians in the treatment options for patients suspected of having bacterial infections.

The device is not automated, is for professional use only and is not a companion diagnostic.

SUMMARY AND EXPLANATION

S. pneumoniae possesses an autocatalytic enzyme that lyses the organism's own cell wall during cell division. In 1900, Neufeld used sodium desoxycholate, a bile salt, to detect autolysis, which differentiated bile-soluble *S. pneumoniae* from bile-insoluble alpha-hemolytic streptococci.¹ Evans employed 10% bile to lyse broth cultures of pneumococci.² Hawn and Beebe evaluated the effectiveness of applying bile salts directly to a blood agar plate culture of pneumococci.³

PRINCIPLE

The bile solubility test is based on the observation that *S. pneumoniae* cells lyse when 10% desoxycholate is applied, while other streptococci and gram-positive cocci do not. The actual mechanism for bile solubility is not fully understood. Bile salts lower surface tension at the medium-membrane interface and also cause a cell membrane derangement.^{4,5} *S. pneumoniae* produces an autolytic, intracellular enzyme that causes the organism to undergo rapid autolysis when cultivated on artificial medium.⁶ Goebel and Avery state that the action of bile salts is independent of the autolytic enzyme.⁷ Dubos and Hirsch found that bile salts inactivate the naturally occurring autolytic enzyme either by the alteration or the removal of a normal autolytic inhibitor from the pneumococcal cells.⁸

REAGENTS (CLASSICAL FORMULA)*

Desoxycholic Acid Sodium (CAS 302-95-4).....	100.0 g
Ethanol (CAS 64-17-5).....	200.0 ml
Sodium Merthiolate, 1:1000 (CAS 54-64-8)	0.15 ml
Demineralized Water (CAS 7732-18-5).....	800.0 ml

*Adjusted as required to meet performance standards.

PRECAUTIONS

WARNING! Flammable liquid and vapor. Causes eye, skin, and respiratory tract irritation. This substance has caused adverse reproductive and fetal effects in humans. May cause central nervous system depression. May cause liver, kidney, and heart damage.

This product is for *In Vitro* diagnostic use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards

by properly sterilizing specimens, containers and media after use. Directions should be read and followed carefully. Refer to Safety Data Sheet on company website for detailed information on chemical reagents.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

In the event of malfunction, do not use the device

STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 20-25°C until used. Do not freeze or overheat.

PRODUCT DETERIORATION

This product should not be used if (1) there is evidence of dehydration, (2) the color has changed, (3) the expiration date has passed, or (4) there are other signs of deterioration. The expiration date applies to the product in its intact container when stored as directed.

Note: A flocculant precipitate may form at temperatures below 15°C. This flocculation will not reduce the effective-ness of the reagent.

SPECIMEN COLLECTION, STORAGE, TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{9,10}

MATERIALS PROVIDED

50 BactiDrop Desoxycholate (10%) – glass ampule enclosed in a protective plastic sheath with a dropper dispenser tip (0.75/ampule)
1 Ampule crusher (reusable)

MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Saline, (7) Test tubes, (8) Phenol Red Indicator (R21250), (9) 0.1 N NaOH.

PROCEDURE

Place dropper in the assembled, reusable ampule crusher provided. Hold the dropper/crusher in an upright position and lightly tap the bottom to dislodge any bubbles that may have formed. Grasp the middle of the dropper/ crusher with the thumb and forefinger, and with the tip pointing away, press gently to crush the ampule. Invert dropper and squeeze slightly to dispense in a dropwise fashion. Discard remaining portion of partially used ampule at end of workday.

Spot Test:

1. To a well-isolated colony on sheep blood agar, add 1 drop of reagent and incubate the plate in an upright position at 35-37°C for 30 minutes.
2. Examine for colony disintegration.

Tube Test:

1. Prepare a heavy suspension of a pure culture in 1 ml of normal saline.
2. Add 1 drop of Phenol Red Indicator and adjust pH to 7.0 (pink) by adding 0.1 N NaOH if necessary.

3. Divide this suspension into two tubes, labeling one Test and the other Control.
4. Add 0.5 ml (17 drops) of the reagent to the tube marked Test.
5. Add 0.5 ml of normal saline to the tube marked Control.
6. Gently agitate the tubes to suspend the bacteria and incubate at 35-37°C for up to 3 hours.
7. Examine hourly for clearing. The Control will remain turbid.

INTERPRETATION

Spot Test:

Positive Test* - Colony disintegration within 30 minutes

Negative Test - Colony remains intact within 30 minutes

Tube Test:

Positive Test* - Clearing of Test suspension within 3 hours; Control suspension remains turbid

Negative Test - Test and Control suspensions remain turbid within 3 hours

*Organism is bile-soluble.

QUALITY CONTROL

All lot numbers of BactiDrop Desoxycholate (10%) have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

CONTROL

Streptococcus pneumoniae
ATCC® 6305

RESULTS

Positive

Streptococcus sanguis
ATCC® 10556

Negative

LIMITATIONS

1. When using the spot test, read a positive test as colony disintegration and not a floating-away of the colony on the surface of the reagent.¹⁰
2. Use bile solubility only to differentiate alpha-hemolytic streptococci from *S. pneumoniae*.¹⁰
3. Older cultures of *S. pneumoniae* may give negative results due to loss of enzyme activity.⁹
4. Some strains of pneumococci are insoluble in bile, whereas some strains of alpha-hemolytic streptococci are soluble. Further testing is required for definitive identification.⁹

BIBLIOGRAPHY

1. Neufeld, F. 1900. Z. Hyg. Infektionskr. 34:454-464.
2. Evans, A.C. 1936. J. Bacteriol. 31:423-437.
3. Hawn, C.V.Z. and E. Beebe. 1965. J. Bacteriol. 90:549.
4. Downie, A.W., L. Stent, and S.M. White. 1931. Br. J. Exp. Pathol. 12:1-9.
5. Oginsky, E.L. and W.W. Umbreit. 1959. An Introduction to Bacterial Physiology. 2nd ed. W.H. Freeman and Co., San Francisco, CA.

6. Cowan, S.T. 1974. Cowan & Steel's Manual for the Identification of Medical Bacteria. 2nd ed. Cambridge University Press, Cambridge, MA.
7. Goebel, W.F. and O.T. Avery. 1929. J. Exp. Med. 49:267-286.
8. Dubos, R.J. and J.G. Hirsch. 1965. Bacterial and Mycotic Infections of Man. 4th ed. J.B. Lippincott Co., Philadelphia, PA.
9. Versalovic, J., K.C. Carroll, G. Funke, J.H. Jorgensen, M.L. Landry, and D.W. Warnock. 2011. Manual of Clinical Microbiology. 10th ed. ASM Press, Washington, D.C.
10. Garcia, L.S. 2010. Clinical Microbiology Procedures Handbook. 3rd ed. ASM Press, Washington, D.C.

Symbol Legend

	Catalog Number
	In Vitro Diagnostic Medical Device
	Batch Code (Lot Number)
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
	Use By (Expiration Date)
	Contains sufficient for <n> tests
	Do not use if packaging damaged
	Authorized European Representative
	Manufacturer
	UK Conformity Assessment
	European Conformity Assessment

Remel Inc.
12076 Santa Fe Trail Drive
Lenexa, KS 66215, USA
www.thermofisher.com/microbiology
(800) 255-6730
International: (913) 888-0939

For technical information contact your local distributor.

©2022 Thermo Fisher Scientific Inc. All rights reserved
BactiDrop is a trademark of Thermo Fisher Scientific and its subsidiaries.

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

CAS (Chemical Abstracts Service Registry No.)

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Remel Inc., 12076 Santa Fe Trail Drive, Lenexa, KS 66215, EE. UU.
©2022 Thermo Fisher Scientific Inc. Reservados todos los derechos.
Todas las demás marcas comerciales son propiedad de Thermo Fisher Scientific Inc. y sus filiales.