www.thermofisher.com

Thermo
SCIENTIFIC

XLT-4 Agar Base

[REF] CM1061B, CM1061T

Intended Use^[VD]

XLT-4 Agar Base (CM1061B, CM1061T) when supplemented with XLT-4 Selective Supplement (SR0237C) is a selective medium intended to be used for the isolation and identification of *Salmonella* species from faecal samples.

XLT-4 Agar Base (CM1061B, CM1061T) is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having infections caused by *Salmonella* species. The device is for professional use only, is not automated and is not a companion diagnostic.

Summary and Explanation

Salmonella species are rod shaped motile bacteria with peritrichous flagella with 2 to 5µm length and 0.7 to 1.5µm diameter of *Enterobacteriaceae* family¹. In humans *Salmonella* species can cause serious health problems ranging from fatal typhoid fever, acute gastroenteritis and enterocolitis. The annual health care burden within the US is estimated to be 25 million infections and 200,000 deaths². Human salmonellosis is notifiable as second most common gastrointestinal infection to campylobacteriosis, leading to 5,000 annual cases in Canada³. *Salmonella* species are broadly classified into two major groups that are *Salmonella enterica* and *Salmonella bongori* and the former subdivided into serovars Typhoid and Paratyphoid. Almost 70% of the *Salmonella* infections in the UK are caused by the non-typhoidal serovar *Salmonella* Enteritidis¹. *Salmonella enterica* often causes gastroenteritis in humans and the mode of transmission can be food like undercooked meat, egg, dairy products, contaminated water, infected people and animals⁴. Non-typhoid *Salmonella* (NTS) are well known to cause gastroenteritis, focal infections, bacteraemia and occasionally urinary tract infections⁵.

Principle of Method

The formulation of XLT-4 Agar Base (CM1061B, CM1061T) Agar Base (CM1061B, CM1061T) includes proteose peptone and yeast extract which provide a source of amino nitrogen, along with essential nutrients and vitamins; this ensures optimal growth of salmonellae. The identification *Salmonella* species is determined by the fermentation of xylose, lactose and sucrose along with the decarboxylation of L- lysine and the production of hydrogen sulphide. Phenol red is the pH indicator and agar is the solidifying agent. XLT-4 Agar Base (CM1061B, CM1061T) is supplemented with XLT-4 Selective Supplement (SR0237C). XLT-4 Selective Supplement (SR0237C) is a selective agent that contains Tergitol™ 4 (also known as Niaproof 4 or sodium tetradecylsulfate) which is an anionic surfactant when added to CM1061T/B formulation largely inhibits or reduces the growth of unwanted background flora.

Typical Formula

	grams per litre
Proteose peptone	1.6
Yeast extract	3.0
L-Lysine	5.0
Xylose	3.75
Lactose	7.5
Sucrose	7.5
Ammonium iron (III) citrate	0.8
Sodium thiosulphate	6.8
Sodium chloride	5.0
Phenol red	0.08
Agar	18.0

Materials Provided

CM1061B: 500 g of XLT-4 Agar Base

CM1061T: 5k g of XLT-4 Agar Base

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Selective supplements (SR0237C)
- Sterile containers

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 24 and Q5.

Procedure

Suspend 59 g in 1 litre of distilled water. Add 4.6 mL of XLT-4 Selective supplement (SR0237C). With frequent agitation, bring to the boil to dissolve completely. Cool to 50°C. Mix well and pour into sterile containers. DO NOT AUTOCLAVE. DO NOT OVERHEAT

Interpretation

Once the medium is reconstituted, *Salmonella abony*, *Salmonella* Enteritidis, *Salmonella* Typhimurium, and *Salmonella arizonae* will show as 1-2 mm red colonies with black centres. *Citrobacter freundii* will show as 0.5-2 mm yellow colonies

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 ± 2 hours @ 37°C aerobic

Positive Controls	
Inoculum level: 10-100 cfu Colony count is ≥ 70% of the control medium count	
<i>Salmonella abony</i> NCTC 6017	1-2 mm red colonies, black centre
<i>Salmonella</i> Enteritidis ATCC® 13076™	1-2 mm red colonies, black centre
<i>Salmonella</i> Typhimurium ATCC® 14028™	1-2 mm red colonies, black centre
Colony count is ≥ 50% of the control medium count	
<i>Salmonella arizonae</i> ATCC® 13314™	1-2 mm red colonies, black centre
<i>Citrobacter freundii</i> ATCC® 8090™	0.5–2 mm yellow colonies
Inoculum level: 10-100 cfu Colony count is ≤ 110% of the control medium count	
<i>Pseudomonas aeruginosa</i> ATCC® 9027™	No growth or pinpoint – 1mm red colonies

Negative Controls	
Inoculum level: 10 ³ -10 ⁵ cfu Negative strains are inhibited or shall produce at least a 2 log (10) reduction	
<i>Escherichia coli</i> ATCC® 25922™	No growth or 1-2 mm yellow colonies
Inoculum level: 10 ⁴ -10 ⁶ cfu	
<i>Proteus mirabilis</i> ATCC® 29906™	No growth
<i>Enterococcus faecalis</i> ATCC® 29212™	No Growth

Limitations

Identifications are presumptive and should be confirmed using appropriate biochemical and serological methods. XLT-4 Agar is not suitable for use in the isolation of *Salmonella Typhi* or *Salmonella Paratyphi* type A as these strains will appear as atypical colonies. Lactose positive and Hydrogen Sulphide (H₂S) negative *Salmonella* may produce colonies with an atypical appearance. Non-*Salmonella* strains that are not completely inhibited on XLT-4 Agar may be encountered and must be differentiated from *Salmonella*. The specified plating medium for the detection of *Salmonella* is Xylose Lactose Desoxycholate Agar (XLD). This medium, like XLT-4 Agar relies on lactose, sucrose and xylose fermentation, lysine decarboxylation and the production of Hydrogen Sulphide for the presumptive identification of *Salmonella* colonies. Therefore, XLT-4 Agar is not recommended as the second plating medium when following the ISO protocol. No single plating medium may be totally reliable therefore it is best to combine with a second *Salmonella* selective plating medium which utilises a different indicator and inhibitor system.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of *Salmonellae* species is confirmed by the inclusion of a well-characterised isolate in the QC processes performed as part of the manufacture of each batch of the device. The precision of XLT-4 Agar (CM1061B/T) was demonstrated by an overall pass rate of 100% over 16 months of testing (22.05.2020 - 14.09.2021). This shows that the performance is reproducible.

XLT-4 Agar Base (CM1061B/T) devices are tested in-house as part of the QC process since the products were launched in 2004. For target organisms, when using 100 cfu inoculum of *Salmonella* Enteritidis, *Salmonella* Typhimurium, *Salmonella arizonae*, and *Citrobacter freundii* and incubating the device at 37 °C for 18± 2 hours, the user can recover organisms with colony size and morphology as listed in this document. For non-target organisms, when using 10³ cfu of *Pseudomonas aeruginosa*, 10⁵ cfu of *Escherichia coli*, and 10⁶ cfu of *Proteus mirabilis* and *Enterococcus faecalis* and incubating the device at 37 °C for 18± 2 hours, the user can recover organisms with colony size and morphology as listed in this document.

Bibliography


1. Public Health England. 2021. 'Standards for Microbiology Investigations Identification of *Salmonella* species.' Issue no: 4 | Issue date: 12.03.21. <https://www.gov.uk/government/publications/smi-id-24-identification-of-salmonella-species>
2. Dekker, J. P., & Frank, K. M. (2015). *Salmonella*, *Shigella*, and *Yersinia*. *Clinics in Laboratory Medicine*, 35(2), 225–246.
3. Nesbitt, A., Ravel, A., Murray, R., McCormick, R., Savelli, C., Finley, R., ... & Canadian Public Health

Laboratory Network. (2012). Integrated surveillance and potential sources of Salmonella Enteritidis in human cases in Canada from 2003 to 2009. *Epidemiology & Infection*, 140(10), 1757-1772.

4. Ford, L., Glass, K., Veitch, M., Wardell, R., Polkinghorne, B., Dobbins, T., ... & Kirk, M. D. (2016). Increasing incidence of Salmonella in Australia, 2000-2013. *PLoS one*, 11(10), e0163989.

5. Gulcan, A., Bayram, P., Levent, B., & Gulcan, E. (2012). A case of urinary tract infection due to Salmonella enterica serovar Virchow and review of the related literature. *Acta microbiologica et immunologica Hungarica*, 59(1), 85-89.


Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier

Made in the United Kingdom	Made in the United Kingdom
----------------------------	----------------------------

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection. NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
1.0	2022-09-01 Original Document

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien