	Document Owner Department: QC	BT-SPEC-0001
		Page 1 of 4
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
SIGNAL BLOOD CULTURE SYSTEM BC0100M		

SIGNAL BLOOD CULTURE SYSTEM

BC0100M

The OXOID SIGNAL Blood Culture System is used to culture samples of blood collected from patients where the condition of bacteraemia is suspected.

Each kit (20 tests per pack) to contain the following:-

- BC102 A sealed blood culture bottle containing 84 ± 2 ml of broth medium.
BC101 A growth indicator device which has been given a dose of gamma-irradiation (25-40kGy), and is vented through a 0.2 micron hydrophobic membrane.
Instruction leaflet.

Formula

Tryptone Soya Broth	grams per litre	10.0
Gelatin peptone		10.0
Yeast extract		5.0
Meat extract		5.0
Sodium chloride		8.0
Potassium nitrate		2.0
Glucose		1.0
L-arginine		1.0
Sodium pyruvate		1.0
Gelatin		1.0
Sodium thioglycollate		0.5
Cysteine HCl		0.4
Sodium bicarbonate		0.4
Phosphate buffer		0.3
Sodium polyanethol sulphonate		0.3
Dithiothreitol		0.2
Adenine sulphate		0.01
Sodium succinate		0.01
Ammonium chloride		0.008
Magnesium sulphate		0.008
Menadione		0.005

Physical Characteristics

Straw-coloured medium
Clarity - passes test
Volume check - passes test
Vacuum check - passes test
Sterility check - passes test
Dead cell count - passes test
pH - 7.0 ± 0.2 at 25°C

ThermoFisher SCIENTIFIC	Document Owner Department: QC	BT-SPEC-0001 Page 2 of 4
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
SIGNAL BLOOD CULTURE SYSTEM BC0100M		

Buffering capacity - passes test

Constituents of Blood Culture Media

Constituents of Blood Culture Media are pre-tested with the following quality control organisms and must perform in accordance with Product QC Specification.

Microbiological Tests using Optimum Inoculum Dilution

The total inoculum challenge for each test organism per bottle is 10-50 colony-forming units (cfu).

<i>Clostridium perfringens</i>	ATCC®13124
<i>Bacillus cereus</i>	ATCC®10876
<i>Bacteroides fragilis</i>	ATCC®25285
<i>Clostridium novyi</i>	ATCC®27606
<i>Escherichia coli</i>	ATCC®25922
<i>Fusobacterium nucleatum</i>	ATCC®10953
<i>Haemophilus influenzae</i>	ATCC®19418
<i>Klebsiella pneumoniae</i>	ATCC®29665
<i>Neisseria meningitidis</i>	ATCC®13077
<i>Pseudomonas aeruginosa</i>	ATCC®27853
<i>Staphylococcus aureus</i>	ATCC®25923
<i>Staphylococcus epidermidis</i>	ATCC®14990
<i>Streptococcus pneumoniae</i>	ATCC®6303
<i>Streptococcus mutans</i>	ATCC®25175
<i>Candida albicans</i>	ATCC®10231
<i>Prevotella bivia</i>	ATCC®29303

Performance

Each lot of Signal Blood Culture System is tested using the following quality control organisms and must perform in accordance with Product QC Specification.

Control Organisms

<i>Clostridium perfringens</i>	ATCC®13124
<i>Streptococcus pneumoniae</i>	ATCC®6303
<i>Candida albicans</i>	ATCC®10231

Precautions - refer to product leaflet

All blood samples should be considered infectious and handled with care.

ThermoFisher SCIENTIFIC	Document Owner Department: QC	BT-SPEC-0001 Page 3 of 4
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
SIGNAL BLOOD CULTURE SYSTEM BC0100M		


Although the majority of blood cultures are negative after incubation, they should be sterilized by autoclaving at 121°C for 15 minutes before being discarded. Do not separate the bottle and the growth indicator device before sterilization.

User Quality Assurance - refer to product leaflet

1. Examine the bottles of broth for turbidity and/or change of colour before adding any blood. Discard the bottles showing abnormal characteristics.
2. If further user quality control is required, it is recommended that 3 aerobes and 1 anaerobe from the above list be used.

Inoculation Procedure - refer to product leaflet

1. Inoculate up to 10ml of blood.
2. Place inoculated bottle at $36 \pm 1^\circ\text{C}$ for 1 hour before inserting Signal device.
3. Continuously shake for 24 hours at $36 \pm 1^\circ\text{C}$.
4. Incubate at $36 \pm 1^\circ\text{C}$ for at least 7 days or until positive. Terminal subculture is recommended.
5. Examine for positive and negative bottles. Subculture positives.

	Document Owner Department: QC	BT-SPEC-0001
		Page 4 of 4
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
SIGNAL BLOOD CULTURE SYSTEM BC0100M		

Revision History

Section / Step	Description of Change	Reason for Change	Reference
Kit contents	Clarity of irradiation dose for BC101	Change control	BT-CC-3063