www.thermofisher.com

SS Agar

REF CM0099B & CM0099K

EN

Intended Use

SS Agar (CM0099B, CM0099K) devices are a selective and differential medium for the isolation of *Salmonella* and some *Shigella* species from faecal samples. Can also be used for testing food samples.

The devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having *Salmonella* and *Shigella* infections.

The devices are for professional use only, are not automated and are not companion diagnostics.

Summary and Explanation

Salmonella and *Shigella* are major causes of bacterial enteric illness with transmission occurring via the fecal-oral route, person-to-person contact, or by the ingestion of contaminated food or water. *Salmonella Shigella* (SS) Agar is a selective and differential medium for the isolation, cultivation and differentiation of *Salmonella* spp. and *Shigella* spp.

Principle of Method

The inclusion of bile salts, sodium citrate and brilliant green serve to inhibit gram-positive organisms, coliforms and inhibit swarming *Proteus* spp., while allowing *Salmonella* spp. to grow. Included peptones provide sources of nitrogen, carbon, and vitamins required for organism growth with lactose providing the carbohydrate source. Sodium thiosulfate and ferric citrate permit detection of hydrogen sulfide by the production of colonies with black centers. Neutral red turns red in the presence of an acidic pH, thus showing fermentation has occurred.

Typical Formula

	<u>grams per litre</u>
'Lab Lemco' Powder	5.0
Peptone	5.0
Lactose	10.0
Bile Salts	8.5
Sodium Citrate	10.0
Sodium Thiosulphate	8.5
Iron (III) Citrate	1.0
Brilliant Green	0.00033
Neutral Red	0.025
Agar	15.0

Materials Provided

CM0099B: 500g of SS Agar.
CM0099K: 25kg of SS Agar.

500g of SS Agar (Modified) yields approximately 7.9L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dishes

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 12°C.

Warnings and Precautions

- Do not autoclave.
- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) S 7 and Q 5.

Procedure

Suspend 63g in 1 litre of distilled water. Bring to the boil with frequent agitation and allow to simmer gently to

dissolve the agar. **DO NOT AUTOCLAVE**. Cool to about 50°C mix well and pour into sterile Petri dishes.

Inoculate the medium heavily with the specimen, spreading a portion of the original inoculum in order to obtain well separated colonies on some part of the plate. Incubate for 18-24 hours at 37°C. Visually inspect plates to assess colony growth and colour under good lighting.

Interpretation

Salmonella form transparent colonies with black centres.
Shigella species form transparent colonies.
Proteus species/*Citrobacter* species will form transparent colonies with grey-black centres.

Non-lactose fermenters form colourless colonies.
Occasional resistant coliforms and other lactose fermenters produce pink or red colonies.

In mixed culture, using the diminishing sweep technique, a satisfactory result is represented by diagnostic reactions of salmonellae and shigellae strains and *Escherichia coli*. Clear differentiation must be seen and is based on the colour and morphology of the colonies.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18-24 hours @ 37°C.

Positive Controls Inoculation with mixed cultures Medium is challenged with 10 ³ to 10 ⁵ cfu of <i>Salmonella</i> and <i>Shigella</i> spp. and 10 ³ to 10 ⁵ cfu for <i>Escherichia coli</i> ATCC® 8739.	
<i>Salmonella enteritidis</i> ATCC® 13076™	0.5-1.5mm straw colonies with or without black centre
<i>Salmonella typhimurium</i> ATCC® 14028™	0.5-2mm straw colonies with black centre
<i>Salmonella virchow</i> NCTC 5742	0.5-1.5mm straw colonies with black centre
<i>Shigella sonnei</i> ATCC® 25931™	1-3mm irregular, translucent pink colonies
<i>Shigella flexneri</i> ATCC® 12022™	1-3mm translucent pink colonies
<i>Shigella boydii</i> NCTC 11462	1-3mm irregular, translucent pink colonies
Inoculation with pure cultures Medium is challenged with 10-100 colony-forming units Colony count is ≤ 100% of the control medium count.	
<i>Pseudomonas aeruginosa</i> ATCC® 27853™	0.5-1mm straw colonies
<i>Proteus mirabilis</i> ATCC® 12453™	0.25-2.5mm straw colonies, with or without grey/black centre, no swarming
<i>Proteus mirabilis</i> ATCC® 29906™	0.25-2.5mm straw colonies, with or without grey/black centre, no swarming
Medium is challenged with 10 ² to 10 ³ cfu Colony count is ≥ 10% of the control medium	
<i>Shigella dysenteriae</i> NCTC 9721	1-3mm translucent pink colonies

Negative controls	
<i>Escherichia coli</i> ATCC® 11775™	No growth or 0.5-2mm pink colonies

Incubation Conditions: 35°C for 18-24 hours

Testing performed in accordance with current CLSI M22 A Positive Controls Medium is challenged with 10-100 colony-forming units Colony count is ≥ 70% of the control medium count.	
<i>Salmonella typhimurium</i> ATCC® 14028™	0.5-2mm straw colonies with black centre
<i>Shigella flexneri</i> ATCC® 12022™	1-3mm translucent pink colonies
Negative controls Inoculum level: 10-100 cfu	
<i>Escherichia coli</i> ATCC® 25922™	No growth or 0.5-1mm pink colonies
Inoculum level: 10 ⁴ -10 ⁶ cfu	
<i>Enterococcus faecalis</i> ATCC® 29212™	No growth

Limitations

Identifications are presumptive and colonies should be confirmed using appropriate methods.
Occasional resistant coliforms and other lactose fermenters produce pink or red colonies.
This medium is highly selective and sensitive strains of *Shigellae* will not grow.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of *Salmonella* and *Shigella* species is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the device. The precision of SS Agar (CM0099B/K) was demonstrated by an overall pass rate of 100% over 1 year 10 months of testing (02.07.2020 – 06.05.2022). This shows that the performance is reproducible.

SS Agar (CM0099B/K) device is tested in-house as part of the QC process since the products were launched in 1972. For target organisms, when using 10³ to 10⁵ cfu inoculum of *Salmonella enteritidis* ATCC® 13076™, *Salmonella typhimurium* ATCC® 14028™, *Salmonella virchow* NCTC 5742, *Shigella sonnei* ATCC® 25931™, *Shigella flexneri* ATCC® 12022™, *Shigella boydii* NCTC 11462 or 10² to 10³ cfu of *Shigella dysenteriae* NCTC 9721 and incubating the device at 37°C for 18-24 hours, the user can recover organisms with colony size and morphology as listed in this document. For pure cultures, when using 10-100 cfu inoculum of *Pseudomonas aeruginosa* ATCC® 27853™, *Proteus mirabilis* ATCC® 12453™ or *Proteus mirabilis* ATCC® 29906™ and incubating the device at 37°C for 18-24 hours, the user can recover organisms with colony size and morphology as listed in this document. For non-target organisms, when using 10-100 cfu inoculum of *Escherichia coli* ATCC® 11775™ and incubating the device at 37°C for 18-24 hours, the user can expect colony morphology as listed in this document.

For testing in accordance to CLSI M22 A, when using 10-100 cfu inoculum of *Salmonella typhimurium* ATCC® 14028™ or *Shigella flexneri* ATCC® 12022™ and incubating the device at 35°C for 18-24 hours, the user can recover organisms with colony size and morphology as listed in this document. For non-target organisms, when using 10-100 cfu of *Escherichia coli* ATCC® 25922™ or 10⁴-10⁶ cfu of *Enterococcus faecalis* ATCC® 29212™ and incubating the device at 35°C for 18-24 hours, the user can recover organisms with colony size and morphology as listed in this document.

©2022 Thermo Fisher Scientific Inc. All rights reserved.

ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
3.0	2023-07-19

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© 2022 Thermo Fisher Scientific Inc. Todos los derechos reservados.
Todas las demás marcas comerciales son propiedad de Thermo Fisher Scientific Inc. y sus filiales.
Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, Reino Unido