www.thermofisher.com

Thermo
SCIENTIFIC

Tryptose Blood Agar Base

REF CM0233B, CM0233T

EN

Intended Use

Tryptose Blood Agar Base (CM0233B and CM0233T) devices are highly nutritive media which, when supplemented with blood, support the growth of fastidious organisms from clinical samples (e.g. mouth, nasal, and genital swabs).

The devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The devices are for professional use only, are not intended for self-testing, and are not companion diagnostics.

Summary and Explanation

A fastidious organism is any organism that has complex or particular nutritional requirements. There are many different types of fastidious organisms which exhibit different requirements to ensure successful growth. Some of the main fastidious humans' pathogens are *Haemophilus influenzae*, *Streptococcus pyogenes* and *Streptococcus pneumoniae*.

The addition of blood to culture media provides means of identification, as the ability to haemolyse blood can be used to classify certain microorganisms^{1,2}. The presence and pattern of haemolysis, the size of the zone of haemolysis, and colony morphology can be used together to identify isolates cultured using blood agar^{1,2}.

Principle of Method

Tryptose Blood Agar Base (CM0233B and CM0233T) devices, when supplemented with blood, facilitate the growth of fastidious organisms found in clinical samples. Tryptose and 'Lab-Lemco' powder provide carbon, nitrogen, vitamins, and minerals necessary for bacterial growth. Sodium chloride provides essential electrolytes and maintains osmotic equilibrium within the medium. Agar is the solidifying agent. Blood can be added to the medium to provide additional nutrients and aid in the identification of haemolytic organisms.

Typical Formula

	<u>grams per litre</u>
Tryptose	10.0
"Lab-Lemco" powder	3.0
Sodium chloride	5.0
Agar	12.0

Materials Provided

500g of dehydrated Tryptose Blood Agar powder that yields approximately 16.7L after reconstitution.

5kg of dehydrated Tryptose Blood Agar powder that yields approximately 167L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Defibrinated horse blood (SR0050) (optional)
- Bacitracin disks (DD0002) (optional)
- Petri dish

Storage

- Store product in its original packaging between 10°C and 30°C.

- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 8°C.

Warnings and Precautions

- Do not inhale. May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- Causes serious eye irritation.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove subject to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 29 (Public Health England, 2020).

Procedure

Suspend 30g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilize by autoclaving at 121°C for 15 minutes. Cool to 50°C. Mix well and pour into sterile Petri dishes. For blood agar, enrich with 7% v/v sterile defibrinated horse blood (SR0050C).

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 h @ 37°C aerobic

Positive Controls Inoculum level: 10 – 100 cfu Plain plates	
<i>Staphylococcus aureus</i> ATCC® 25923™	1-2mm straw colonies
<i>Escherichia coli</i> ATCC® 25922™	1-2mm straw colonies
<i>Pseudomonas aeruginosa</i> ATCC® 27853™	1-3mm straw colonies with or without green pigmentation
Positive Controls Inoculum level: 10 – 100 cfu Enriched with 7% v/v horse blood A satisfactory result is represented by recovery of positive strains equal to or greater than 70% of the control medium.	
<i>Streptococcus pyogenes</i> ATCC® 19615™*	0.25-0.75mm pale straw colonies, β haemolysis
<i>Streptococcus pneumoniae</i> ATCC® 6303™	1-2mm grey/green colonies, α haemolysis
<i>Streptococcus pneumoniae</i> ATCC® 6305™	1-2mm grey/green colonies, α haemolysis
<i>Haemophilus influenzae</i> ATCC® 19418™	0.25-0.75mm straw/grey colonies

*Zones of inhibition with Bacitracin discs (DD0002) shall be 10-22mm when 7% v/v horse blood plates are inoculated with *Streptococcus pyogenes* ATCC® 19615™ after incubation at 37°C for 18 hours.

Limitations

As with all media, atypical organisms may give anomalous reactions. Some strains may be encountered that grow poorly or fail to grow.

Performance Characteristics

Accuracy has been demonstrated through a review of the QC data. Correct detection of fastidious and non-fastidious organisms is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the devices. The precision of Tryptose Blood Agar Base (CM0233B and CM0233T) was demonstrated by an overall pass rate of 100% obtained for the product over 2 years of testing (26.06.2020 – 29.04.2022; 10 batches). This shows that the performance is reproducible.

Tryptose Blood Agar Base (CM0233B and CM0233T) devices are tested in-house as part of the QC process since the products were launched in at least 1972, with manufacturing data available since 1996.. For plain plates, when using 10-100 cfu inoculum of *Staphylococcus aureus* (ATCC®25923™), *Escherichia coli* (ATCC®25922™), and *Pseudomonas aeruginosa* (ATCC®27853™), and incubating the device at 37°C for 18 hours under aerobic conditions, the user can recover organisms with colony size and morphology as listed in this document. For plates enriched with 7% v/v horse blood, when using 10-100 cfu inoculum of *Streptococcus pyogenes* (ATCC®19615™), *Streptococcus pneumoniae* (ATCC®6303™), *Streptococcus pneumoniae* (ATCC®6305™), and


Haemophilus influenzae (ATCC®19418™), and incubating the device at 37°C for 18 hours under aerobic conditions, the user can recover organisms with colony size and morphology as listed in this document. For plates enriched with 7% v/v horse blood, inoculated with *Streptococcus pyogenes* (ATCC®19615™) and with bacitracin discs applied to the agar surface, after incubating the device at 37°C for 18 hours under aerobic conditions, the user can obtain inhibition zones as listed in this document.

Bibliography

- Public Health England. 2020. 'Identification of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
- Public Health England. 2021a. 'Identification of Streptococcus species, Enterococcus species and morphologically similar organisms.' UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use

	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, UK

For technical assistance please contact your local distributor.

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
Bibliographie

- Public Health England. 2020. "Identification of *Staphylococcus* species, *Micrococcus* species and *Rothia* species." UK Standards for Microbiology Investigations. ID 7 Issue 4. <https://www.gov.uk/government/publications/smi-id-7-identification-of-staphylococcus-species-micrococcus-species-and-rothia-species>.
- Public Health England. 2021a. "Identification of *Streptococcus* species, *Enterococcus* species and morphologically similar organisms." UK Standards for Microbiology Investigations. ID 4 Issue 4. <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.

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Oxoid Limited, Wade Road, Basingstoke,
RG24 8PW, Royaume-Uni

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©2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.