Pseudomonas Agar Base (Dehydrated)*

EN

REF CM0559B, CM0559R & CM0559T

*This instructions for use (IFU) document is intended to be read in conjunction with the IFU for Oxoid C-N Supplement (SR0102E) or C-F-C Supplement (SR0103E).

Intended Use

Pseudomonas Agar Base (Dehydrated) (CM0559B, CM0559R, CM0559T), when supplemented with C-N Supplement (SR0102E) or C-F-C Supplement (SR0103E), is a selective medium for the isolation of *Pseudomonas* spp. from a wide range of clinical samples including urine, wound swabs, and sputum. The devices can also be used for testing food, water, and environmental samples. When supplemented with Pseudomonas C-N Supplement (SR0102E) it is recommended for the selective isolation of *Pseudomonas aeruginosa* and when supplemented with Pseudomonas C- F-C Supplement (SR0103E) it is recommended for the selective isolation of *Pseudomonas* spp. generally.

Pseudomonas Agar Base with C-N Supplement or C-F-C Supplement added, is used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having Pseudomonas infections.

The device is for professional use only, is not automated, and nor is it a companion diagnostic.

Summary and Explanation

Pseudomonas species are motile, gram negative bacteria that are straight or curved in appearance¹. *Pseudomonas* species can cause a variety of infections in humans, ranging from lungs (pneumonia), wound infections, septicaemia, urinary tract and gastrointestinal infections².

Pseudomonas aeruginosa is a non-fermentive, strictly aerobic, catalase and oxidase positive gram-negative rod. It is widely distributed in soil and water and may cause disease in animals, plants and insects. *P. aeruginosa* is known as the key respiratory pathogen in patients with cystic fibrosis³. In addition to it, people with, chronic obstructive pulmonary disease (COPD), pregnant ladies and people with a weakened immune system are most affected by the *Pseudomonas* infection². *P. aeruginosa* are known to cause a condition known as pyelonephritis particularly in patients with catheters⁴. In general infections caused by *Pseudomonas aeruginosa* are treated with antibiotics. However, with an increase in the number of multi drug resistant strains of *Pseudomonas aeruginosa*, presumptive identification of the strains is highly important. For the identification of best antibiotics to treat the infection clinical specimens are sent to a laboratory and are grown against a set of antibiotics like fluoroquinolones².

Principle of Method

Contains gelatin peptone and casein hydrolysate to supply nitrogen and growth factors, glycerol as a carbon and energy source and agar as the solidifying agent. Magnesium chloride and potassium sulphate are present to amplify pigment production through the enhanced production of pyocyanin by *P. aeruginosa*. The colour produced by the presence of pyoverdine is yellow/green and the presence of pyocyanin produces blue/green colonies. Identification of *P. aeruginosa* is achieved through green pigmentation of colonies, produced through the combination of these pigments.

Typical Formula

	grams per litre
Gelatin peptone	16.0
Casein hydrolysate	10.0
Potassium sulphate	10.0
Magnesium chloride	1.4
Agar	11.0

Materials Provided

CM0559B: 500g dehydrated Pseudomonas Agar Base

CM0559R: 2.5kg dehydrated Pseudomonas Agar Base

CM0559T: 5kg dehydrated Pseudomonas Agar Base

24.2g of dehydrated Pseudomonas Agar Base yields approximately 500ml after reconstitution

Materials Required but Not Supplied

- (1) Inoculating loops, swabs, collection containers
- (2) Incubators
- (3) Quality control organisms
- (4) Selective supplements (SR0102E, SR0103E)
- (5) Petri dish

Storage

Store dehydrated media in its original packaging between 10°C and 30°C. Keep out of direct sunlight. Keep away from moisture.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- Do not inhale.

- Causes serious eye irritation.
- May cause allergy or asthma symptoms or difficulty breathing if inhaled.
- May cause an allergic skin reaction.
- If on skin wash with plenty of soap and water.
- If in eyes, rinse cautiously with water for several minutes.
- Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists, seek medical advice/attention.
- If inhaled, if breathing is difficult, remove person to fresh air and keep in a position comfortable for breathing. If experiencing respiratory symptoms, call a POISON CENTER or doctor/physician.
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local the recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 17, B 11, B 41 and B 57.

Procedure

Suspend 24.2g in 500ml of distilled water. Add 5ml of glycerol. Bring to the boil to dissolve completely. Sterilize by autoclaving at 121°C for 15 minutes. Cool to 50°C and aseptically add the contents of 1 vial of either *Pseudomonas* CN selective supplement (SR0102E) or *Pseudomonas* CFC Selective Supplement (SR0103E) reconstituted as directed. Mix well and pour into sterile Petri dishes.

Interpretation

Once reconstituted the presence of straw colonies may indicate *Pseudomonas* species. Other identification procedures should be performed to confirm findings.

More details at: www.thermoscientific.com/microbiology

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium with C-N Supplement SR0102E can be verified by testing the following reference strains.

Incubation Conditions: 18 ± 2 h @ 36 ± 2°C

Positive Controls Colony count is ≥ 70% of the control medium count Inoculum level 10-100 cfu & 10 ⁴ - 10 ⁶ cfu	
<i>Pseudomonas aeruginosa</i> ATCC® 25668	1-4mm straw colonies, green pigmentation
Negative Controls Inoculum level 10 ⁴ - 10 ⁶ cfu	
<i>Burkholderia cepacia</i> ATCC® 25416	No growth or pinpoint to 0.25mm straw colonies
<i>Pseudomonas putida</i> NCIMB 9034	No growth
<i>Pseudomonas fluorescens</i> NCTC 3756	No growth
<i>Proteus hauseri</i> ATCC® 13315	No growth
<i>Klebsiella pneumoniae</i> ATCC® 29665	No growth
<i>Staphylococcus aureus</i> ATCC® 25923	No growth

<i>Citrobacter freundii</i> ATCC® 8090	No growth
---	-----------

The performance of this medium with C-F-C Supplement SR0103E can be verified by testing the following reference strains.

Incubation Conditions: 44 ± 4 h @ 30 ± 2°C

Positive Controls Colony count is ≥ 70% of the control medium count Inoculum level 10-100 cfu & 10 ⁴ - 10 ⁶ cfu	
<i>Pseudomonas aeruginosa</i> ATCC® 27853	1-4mm straw colonies, green pigmentation
<i>Burkholderia cepacia</i> ATCC® 25416	1-4mm straw colonies
Negative Controls Inoculum level 10 ⁴ – 10 ⁶ cfu	
<i>Klebsiella pneumoniae</i> ATCC® 29665	No growth
<i>Staphylococcus aureus</i> ATCC® 25923	No growth
Inoculum level 10 ⁴ – 10 ⁵ cfu	
<i>Proteus hauseri</i> ATCC® 13315	No growth
Inoculum level 10 ² - 10 ⁴ cfu	
<i>Citrobacter freundii</i> ATCC® 8090	No growth or 1-2mm straw colonies

Incubation Conditions: 44 ± 4 h @ 25 ± 2°C

Positive Controls Colony count is ≥ 70% of the control medium count Inoculum level 10-100 cfu & 10 ⁴ - 10 ⁶ cfu	
<i>Pseudomonas putida</i> NCIMB 9034	1-4mm straw colonies

Limitations

Organisms, other than *P. aeruginosa* that are resistant to the antimicrobials used in the medium may be able to grow; this includes some Enterobacteriaceae, other non- fermenters and *Bacillus* spp. Some non-target strains may also produce pigment. The presence of blue-green or brown pigmentation or fluorescence may be taken as presumptive evidence of *P. aeruginosa* but further tests must be carried out to confirm the identity of the organism. Some strains with particular growth requirements may not grow; pigmentation and fluorescence is also variable. *Burkholderia cepacia*, formerly part of *Pseudomonas*, may grow on the media and further tests must be carried out to confirm identification.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of *Pseudomonas aeruginosa* OR *Pseudomonas* species strains is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the device. The precision of Pseudomonas Agar Base (Dehydrated) (CM0559) was demonstrated by an overall pass rate of 100% obtained for the product over nine months of testing (24.05.2021 – 22.02.2022; 10 batches).

Pseudomonas Agar Base (Dehydrated) (CM0559) with C- N Supplement (SR0102E) device is tested in-house as part of the QC process. For target organisms, when using 10-


100 cfu and 10⁴-10⁶ cfu inoculum of *Pseudomonas aeruginosa* ATCC® 25668 and incubating the device at 36±2°C for 18±2°C hours, the user can recover organisms with colony size and morphology as described in this document. For non-target organisms, when using 10⁴-10⁶ cfu inoculum of *Burkholderia cepacia* ATCC® 25416 and incubating the device at 36±2°C for 18±2°C hours, the test medium must show no growth to pinpoint to 0.25mm straw colonies. For non-target organisms, when using *Pseudomonas putida* NCIMB 9034, *Pseudomonas fluorescens* NCTC 3756, *Klebsiella pneumoniae* ATCC®29665, *Proteus hauseri* ATCC® 13315, *Staphylococcus aureus* ATCC® 25923 and *Citrobacter freundii* ATCC® 8090, and incubating the device at 36±2°C for 18±2°C hours, the test medium must show no growth.

Pseudomonas Agar Base (Dehydrated) (CM0559) with C- F-C Supplement (SR0103E) device is tested in-house as part of the QC process. For target organisms, when using 10-100 cfu and 10⁴-10⁶ cfu inoculum of *Pseudomonas aeruginosa* ATCC® 27853 and *Pseudomonas putida* NCIMB 9034 and *Burkholderia cepacia* ATCC® 25416 and incubating the device at 30±2°C for 44±4 hours and 25±2°C for 44±4 hours (*Pseudomonas putida*), the user can recover organisms with the appropriate morphology as listed in this document. For non-target organisms, when using 10²-10⁴ cfu inoculum of *Citrobacter freundii* ATCC® 8090 and incubating the device at 30±2°C for 44±4 hours, the test medium must show no growth or 1-2mm straw colonies. For non-target organisms, when using 10⁴-10⁶ cfu inoculum of *Klebsiella pneumoniae* ATCC® 29665 and *Staphylococcus aureus* ATCC® 25923, 10⁴-10⁵ cfu inoculum of *Proteus hauseri* ATCC® 13315, and incubating the device at 30±2°C for 44±4 hours, the test medium must show no growth.

Bibliography

1. Public Health England 'Identification of Pseudomonas species and other NonGlucose Fermenters.' iD17i3 24.04.2015. <https://www.gov.uk/government/publications/smi-id-17-identification-of-glucose-non-fermenting-gram-negative-rods>
2. Centers for Disease control and Prevention (CDC). 2019. 'Pseudomonas aeruginosa in Healthcare Settings' <https://www.cdc.gov/hai/organisms/pseudomonas.html>
3. Langan, K. M., Kotsimbos, T., & Peleg, A. Y. (2015). Managing Pseudomonas aeruginosa respiratory infections in cystic fibrosis. *Current Opinion in Infectious Diseases*, 28(6), 547–556.
4. Newman, J. W., Floyd, R. v, & Fothergill, J. L. (2017). The contribution of Pseudomonas aeruginosa virulence factors and host factors in the establishment of urinary tract infections. *FEMS Microbiology Letters*, 364(15).

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.

NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, RG24
8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Date of issue and modifications introduced
2.0	2023-11-10

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
- Centers for Disease control and Prevention (CDC). 2019. 'Pseudomonas aeruginosa in Healthcare Settings' <https://www.cdc.gov/hai/organisms/pseudomonas.html>
- Langan, K. M., Kotsimbos, T., & Peleg, A. Y. (2015). Managing Pseudomonas aeruginosa respiratory infections in cystic fibrosis. *Current Opinion in Infectious Diseases*, 28(6), 547–556.
- Newman, J. W., Floyd, R. v, & Fothergill, J. L. (2017). The contribution of Pseudomonas aeruginosa virulence factors and host factors in the establishment of urinary tract infections. *FEMS Microbiology Letters*, 364(15).

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Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.
Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, UK