Yersinia Selective Agar Base (Dehydrated)*

REF CM0653B, CM0653T

EN

*This instructions for use (IFU) document is intended to be read in conjunction with the IFU for Oxoid Yersinia Selective Supplement (product code: SR0109E).

Intended Use **IVD**

Yersinia Selective Agar Base (Dehydrated) (CM0653B and CM0653T), when supplemented with Yersinia Selective Supplement (SR0109E), is a selective medium for the isolation of *Yersinia enterocolitica* from faecal samples.

Yersinia Selective Agar Base (Dehydrated), with Yersinia Selective Supplement added, is intended to be used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having yersiniosis.

The device is for professional use only, is not automated nor is it a companion diagnostic..

Summary and Explanation

Yersinia enterocolitica are Gram-negative, facultatively anaerobic, coccoid shaped bacteria and are the main causative agent of the foodborne zoonotic disease yersiniosis.^{1, 2} In humans, *Y. enterocolitica* have primarily been isolated from faecal samples, but they have also been found in blood and lymph node tissues.¹

Y. enterocolitica infections are recognized as a main cause of yersiniosis.² In 2019, yersiniosis was the fourth most reported foodborne zoonotic disease in Europe, with more than 7,000 confirmed cases, with 99% of these infections caused by *Y. enterocolitica*.³ Infections present with a range of manifestations, including gastroenteritis, acute diarrhoea, mesenteric adenitis, terminal ileitis, 'pseudo- appendicitis', and in rare cases sepsis.^{1, 2}

Consequently, it is important to isolate *Yersinia enterocolitica* in clinical samples as it is a considerable human pathogen associated with yersiniosis.^{1, 2}

Isolation and presumptive identification of *Y. enterocolitica* is conventionally carried out using a combination of microscopic appearance, cultural appearance on *Y. enterocolitica* selective agar, growth requirements and serotyping.

Principle of Method

Isolation of *Yersinia enterocolitica* is achieved through the inclusion of a combination of selective agents. Special peptone provides nitrogen and amino acids, and yeast extract provides vitamins and other nutrients necessary for bacterial growth. Neutral red acts as a pH indicator. Mannitol is fermented by *Y. enterocolitica*, producing an acidic pH which causes the colonies to develop a red 'bull's eye' surrounded by a transparent border. Most other microorganisms that are able to grow will produce larger colonies with diffuse pinkish centres and opaque outer zones. Sodium deoxycholate, crystal violet and novobiocin inhibit the growth of most gram-positive microorganisms, as well as some strains of *Proteus*. Cefsulodin is a bactericidal antimicrobial agent with activity against *Pseudomonas aeruginosa*, whilst irgasan is a broad-spectrum antimicrobial agent.

Typical Formula

	<u>grams per litre</u>
Special peptone	20.0
Yeast Extract	2.0
Mannitol	20.0
Sodium pyruvate	2.0
Sodium chloride	1.0
Magnesium sulphate	0.01
Sodium deoxycholate	0.5
Neutral red	0.03
Crystal violet	0.001
Agar	12.5

Materials Provided

CM0653B: 500g of dehydrated Yersinia Selective Agar Base
CM0653T: 5kg of dehydrated Yersinia Selective Agar Base

29g of dehydrated Yersinia Selective Agar Base yields approximately 500mL after reconstitution.

Materials Required but Not Supplied

- (1) Inoculating loops, swabs, collection containers
- (2) Incubators
- (3) Quality control organisms
- (4) Selective supplements (SR0109E)
- (5) Petri dish

Storage

- Store product in its original packaging at 10– 30°C until used.
- The product may be used until the expiry date stated on the label.

- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local the recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 21 and S 7.

Procedure

Suspend 29g in 500ml of distilled water. Bring to the boil to dissolve completely. Sterilize by autoclaving at 121°C for 15 minutes. Cool to 50°C and aseptically add the contents of 1 vial of Yersinia Selective Supplement (SR0109E) reconstituted as directed. Mix well and pour into sterile Petri dishes.

Interpretation

Once reconstituted, the presence of transparent colonies with pink centres indicates the presence of *Yersinia enterocolitica*.

Most non-target organisms will be inhibited. Microorganisms that are able to grow will produce larger colonies (>2 mm in diameter) with diffuse pinkish centres and opaque outer zones.

Microorganisms that do not metabolize mannitol form colourless, translucent colonies.

Identification is presumptive and should be confirmed.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium* can be verified by testing the following reference strains.

Incubation Conditions: 21 ± 3 h @ 30° ± 2°C

Positive Controls	
Colony count is ≥ 70% of the control medium count Inoculum level 10 – 100 cfu	
<i>Yersinia enterocolitica</i> ATCC® 27729	0.5-2mm transparent, red, bull's-eye colonies
<i>Yersinia enterocolitica</i> NCTC 10460	0.5-2mm transparent, red, bull's-eye colonies
Colony count is ≥ 50% of the control medium count Inoculum level 10 – 100 cfu	
<i>Citrobacter freundii</i> ATCC® 8090	1-3mm pink colonies, with or without opaque zones
Negative Controls	
Inoculum level 10 ⁴ - 10 ⁵ cfu	
<i>Proteus mirabilis</i> ATCC® 29906	No growth or 1-2mm straw colonies
Tested in accordance with current CLSI M22 A	
Positive Controls	
Inoculum level 10 ³ - 10 ⁴ cfu A satisfactory result is represented by a positive diagnostic reaction	

<i>Yersinia enterocolitica</i> ATCC® 9160	0.5-2mm transparent, red, bull's-eye colonies
Negative Controls Inoculum level 10 ⁴ – 10 ⁶ cfu	
<i>Escherichia coli</i> ATCC® 25922	No growth
<i>Pseudomonas aeruginosa</i> ATCC® 27853	No growth
<i>Enterococcus faecalis</i> ATCC® 29212	No growth

*Tested with the addition of Yersinia Selective Supplement (SR0109E)

Limitations

Some strains of *Y. enterocolitica* may grow poorly or not at all.

Identification is presumptive, other species of *Yersinia* may grow along with some enteric microorganisms. Therefore, full identification tests must be carried out on suspect colonies.

Serratia liquefaciens, *Citrobacter freundii* and *Enterobacter agglomerans* may give a colonial morphology resembling *Y. enterocolitica*. These organisms can be differentiated from *Y. enterocolitica* by biochemical tests such as growth on nutrient and MacConkey agars, indole and urease production and for acid reactions from sucrose, cellobiose, amygdalin, melibiose, rhamnose and raffinose. Tests should be carried out at 30°C rather than 37°C.^{4,5}

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of *Yersinia enterocolitica* is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the device. The precision of Yersinia Selective Agar Base (Dehydrated) (CM0653) was demonstrated by an overall pass rate of 100% obtained for the product over 1.5 years of testing (March 2020 – November 2021; 10 batches). This shows that the performance is reproducible.


The devices are tested in-house as part of the QC process. For target organisms, when using 10-100 cfu inoculum of *Yersinia enterocolitica* (ATCC® 27729 and NCTC 10460) and 10³-10⁴ cfu inoculum of *Yersinia enterocolitica* (ATCC® 9610) and incubating the device at 30°C ± 2°C for 18-24 hours, the user can recover organisms with colony size and morphology as listed in this document. For non-target organisms, when using 10-100 cfu inoculum of *Citrobacter freundii* (ATCC® 8090), 10⁴-10⁵ cfu inoculum of *Proteus mirabilis* (ATCC® 29906) and 10⁴-10⁶ cfu of *Escherichia coli* (ATCC® 25922), *Pseudomonas aeruginosa* (ATCC® 27853) and *Enterococcus faecalis* (ATCC® 29212) and incubating the device at 30°C ± 2°C for 18-24 hours, the user can recover organisms with colony size and morphology as listed in this document.

Bibliography

1. Public Health England. 2015. 'Identification of *Yersinia* species'. UK Standards for Microbiology Investigations. ID 21 Issue 3. Accessed 11 February 2022. <https://www.gov.uk/government/publications/smi-id-21-identification-of-yersinia-species-from-faeces>
2. Public Health England. 2020. 'Gastroenteritis'. UK Standards for Microbiology Investigations. S 7 Issue 2. Accessed 11 February 2022. <https://www.gov.uk/government/publications/smi-s-7-gastroenteritis-and-diarrhoea>
3. European Centre for Disease Prevention and Control. 2021. 'Yersiniosis-Annual epidemiological report for 2019'. ECDC. Accessed 11 February 2022. <https://www.ecdc.europa.eu/en/publications-data/yersiniosis-annual-epidemiological-report-2019>
4. Mair, N. S. and Fox E. 1986. 'Yersiniosis: Laboratory Diagnosis, Clinical Features and Epidemiology'. Pub. Hlth Lab. Ser. London.
5. Pai, C. H., Sorger S., Lafleur L., Lackman L. and Marks M. I. 1979. 'Efficacy of cold enrichment techniques for recovery of *Yersinia enterocolitica* from human stools'. Journal of Clinical Microbiology 9 (6): 712-715. <https://doi.org/10.1128/jcm.9.6.712-715.1979>

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date

	Keep away from sunlight
	Consult instructions for use
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
	Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union
	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

For technical assistance please contact your local distributor.

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Literatura

- Public Health England. 2015. 'Identification of *Yersinia* species'. UK Standards for Microbiology Investigations. ID 21 Issue 3. Accessed 11 February 2022. <https://www.gov.uk/government/publications/smi-id-21-identification-of-yersinia-species-from-faeces>
- Public Health England. 2020. 'Gastroenteritis'. UK Standards for Microbiology Investigations. S 7 Issue 2. Accessed 11 February 2022. <https://www.gov.uk/government/publications/smi-s-7-gastroenteritis-and-diarrhoea>
- European Centre for Disease Prevention and Control. 2021. 'Yersiniosis-Annual epidemiological report for 2019'. ECDC. Accessed 11 February 2022. <https://www.ecdc.europa.eu/en/publications-data/yersiniosis-annual-epidemiological-report-2019>
- Mair, N. S. and Fox E. 1986. 'Yersiniosis: Laboratory Diagnosis, Clinical Features and Epidemiology'. Pub. Hlth Lab. Ser. London.
- Pai, C. H., Sorger S., Lafleur L., Lackman L. and Marks M. I. 1979. 'Efficacy of cold enrichment techniques for recovery of *Yersinia enterocolitica* from human stools'. Journal of Clinical Microbiology 9 (6): 712-715. <https://doi.org/10.1128/jcm.9.6.712-715.1979>

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- Public Health England. 2015. 'Identification of *Yersinia* species'. UK Standards for Microbiology Investigations. ID 21 Issue 3. Accessed 11 February 2022. <https://www.gov.uk/government/publications/smi-id-21-identification-of-yersinia-species-from-faeces>
- Public Health England. 2020. 'Gastroenteritis'. UK Standards for Microbiology Investigations. S 7 Issue 2. Accessed 11 February 2022. <https://www.gov.uk/government/publications/smi-s-7-gastroenteritis-and-diarrhoea>
- European Centre for Disease Prevention and Control. 2021. 'Yersiniosis-Annual epidemiological report for 2019'. ECDC. Accessed 11 February 2022. <https://www.ecdc.europa.eu/en/publications-data/yersiniosis-annual-epidemiological-report-2019>
- Mair, N. S. and Fox E. 1986. 'Yersiniosis: Laboratory Diagnosis, Clinical Features and Epidemiology'. Pub. Hlth Lab. Ser. London.
- Pai, C. H., Sorger S., Laffleur L., Lackman L. and Marks M. I. 1979. 'Efficacy of cold enrichment techniques for recovery of *Yersinia enterocolitica* from human stools'. Journal of Clinical Microbiology 9 (6): 712-715. <https://doi.org/10.1128/jcm.9.6.712-715.1979>

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
Bibliografie

1. Public Health England. 2015. 'Identification of *Yersinia* species'. UK Standards for Microbiology Investigations. ID 21 Issue 3. Accessed 11 February 2022. <https://www.gov.uk/government/publications/smi-id-21-identification-of-yersinia-species-from-faeces>
2. Public Health England. 2020. 'Gastroenteritis'. UK Standards for Microbiology Investigations. S 7 Issue 2. Accessed 11 February 2022. <https://www.gov.uk/government/publications/smi-s-7-gastroenteritis-and-diarrhoea>
3. European Centre for Disease Prevention and Control. 2021. 'Yersiniosis-Annual epidemiological report for 2019'. ECDC. Accessed 11 February 2022. <https://www.ecdc.europa.eu/en/publications-data/yersiniosis-annual-epidemiological-report-2019>
4. Mair, N. S. and Fox E. 1986. 'Yersiniosis: Laboratory Diagnosis, Clinical Features and Epidemiology'. Pub. Hlth Lab. Ser. London.
5. Pai, C. H., Sorger S., Laffleur L., Lackman L. and Marks M. I. 1979. 'Efficacy of cold enrichment techniques for recovery of *Yersinia enterocolitica* from human stools'. Journal of Clinical Microbiology 9 (6): 712-715. <https://doi.org/10.1128/jcm.9.6.712-715.1979>

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Regatul Unit

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Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.
Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Storbritannien