www.thermofisher.com

Selenite Cystine Broth Base

[REF] CM0699B

Intended Use

Selenite Cystine Broth Base (CM0699B) is a selective enrichment medium for the isolation of *Salmonella* species from faecal and urine samples. Can also be used for testing water and food samples.

Selenite Cystine Broth Base (CM0699B) is used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having salmonellosis.

The devices are for professional use only, are not automated, nor are they companion diagnostics.

Summary and Explanation

Salmonella species are inherently pathogenic to humans. These bacteria can be ingested from contaminated food, or water - via the faecal-oral route – and will invade the gastrointestinal tract to cause enteric infection, with diarrhoea being the most common symptom. *Salmonella* infections can become invasive in infants, older adults and can be serious in immunocompromised individuals.¹

Salmonella species are listed as one of four key global causes of diarrhoeal disease worldwide.^{2,3} In Europe, salmonella infection – salmonellosis - is the second most common foodborne disease,⁴ and an estimated 550 million people fall ill with diarrhoeal disease every year from unsafe food.³

Principle of Method

Selenite Cystine Broth Base (CM0699B) facilitates the isolation of *Salmonella* spp. found in clinical faecal specimens. Tryptone provides essential amino acids, vitamins, and minerals for bacterial growth. Lactose is added as a fermentable carbohydrate source, sodium phosphate is added as a buffering agent. L-Cystine is added as a reducing agent to enhance the recovery of *Salmonella* species. As *Salmonella* species are motile, sampling from the upper portion of the broth may help to reduce overgrowth of non-motile microorganisms which settle to the bottom.

Typical Formula

	<u>grams per litre</u>
Tryptone	5.0
Lactose	4.0
Di-sodium phosphate	10.0
L-Cystine	0.01

Materials Provided

CM0699B: 500g of Dehydrated Selenite Cystine Broth

Materials Required but Not Supplied

- Inoculating loops, swabs, collection containers
- Incubators
- Quality control organisms
- Petri dish
- Sodium Biselenite (LP0121)

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- IF SWALLOWED: Clean mouth with water and drink afterwards plenty of water. Get medical attention if symptoms occur.
- IF IN EYES: Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Get medical attention immediately if symptoms occur.
- IF ON SKIN: Wash off immediately with plenty of water for at least 15 minutes. Get medical attention if symptoms occur.
- IF INHALED: Remove to fresh air. Get medical attention immediately if symptoms occur
- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimen should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 24 and S 7.

Procedure

Suspend 19g in 1 litre of distilled water to which 4g of Sodium Biselenite (LP0121) have been added. Warm to dissolve. Mix well and distribute into final containers. Sterilize by heating at 100°C for 10 minutes. DO NOT AUTOCLAVE. DO NOT OVERHEAT.

Interpretation

Once the medium is reconstituted:

- Straw colonies indicate indicate *Salmonella* Nottingham or *Salmonella* Enteritidis.
- Red colonies with black centre indicates *Salmonella* Typhimurium.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 hours @ 37°C, aerobic

Positive Controls Inoculate 10ml quantities of medium to achieve 1-10 colony-forming units/ml (cfu/ml) of <i>Salmonella</i> spp., add 10 ² to 10 ⁴ cfu/ml of <i>Escherichia coli</i> ATCC®25922. Incubate broths at 37°C for 18 hours. Subculture onto XLD Medium (CM0469B) and incubate plates at 37°C for 24 hours.	
<i>Salmonella</i> Typhimurium ATCC® 14028™	1-2 mm red colonies with black centre
Recovery of each organism shall be comparable to the standard. <i>Escherichia coli</i> should be easily distinguishable from <i>Salmonella</i> spp. and not overgrow the target organism. A satisfactory result is represented by recovery of <i>Salmonella</i> Typhimurium ATCC®14028 equal to or greater than a 2 log (10) increase.	
Inoculation with pure cultures Inoculate 10ml quantities of medium to achieve 1-10 colony-forming units/ml (cfu/ml) of the test organism. Incubate broths at 37°C for 18 hours. Subculture onto Tryptone Soya Agar (CM0131B) and incubate plates at 37°C for 24 hours.	
<i>Salmonella</i> Enteritidis ATCC® 13076™	1-3mm straw colonies
<i>Salmonella</i> nottingham NCTC 7832	1-3mm straw colonies
A satisfactory result is represented by recovery of <i>Salmonella</i> strains equal to or greater than a 2 log (10) increase.	

Limitations

Non-target organisms that are not affected by the selective components may grow in the device once reconstituted.

Atypical organisms may give anomalous reactions.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of salmonella strains is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the device. The precision of Selenite Cystine Broth Base (CM0699B) devices was demonstrated by an overall pass rate of 100% obtained for the product over ten years of testing (10.01.2011 – 01.07.2021; 10 batches). This shows that the performance is reproducible.

Selenite Cystine Broth Base (CM0699B) device is tested in-house as part of the QC process since the products were launched in 1996. For target organisms, when using 10ml quantities of medium to achieve 1-10 colony-forming units/ml of *Salmonella* Typhimurium ATCC® 14028™,

Salmonella Enteritidis ATCC® 13076™ or *Salmonella* nottingham NCTC 7832 and incubating the device at 37°C for 18 hours, the user can recover organisms with colony size and morphology as listed in this document.

Bibliography

1. CDC. 2020. 'Salmonellosis (Nontyphoidal) - Chapter 4 - 2020 Yellow Book | Travelers' Health | CDC.' 2020. <https://wwwnc.cdc.gov/travel/yellowbook/2020/travel-related-infectious-diseases/salmonellosis-nontyphoidal>.
2. Kirk, Martyn D, Sara M Pires, Robert E Black, Marisa Caipo, John A Crump, Brecht Devleesschauwer, Dörte Döpfer, et al. 2015. 'World Health Organization Estimates of the Global and Regional Disease Burden of 22 Foodborne Bacterial, Protozoal, and Viral Diseases, 2010: A Data Synthesis.' <https://doi.org/10.1371/journal.pmed.1001921>.
3. WHO. 2018. 'Salmonella (Non-Typhoidal).' World Health Organisation. [https://www.who.int/news-room/fact-sheets/detail/salmonella-\(non-typhoidal\)](https://www.who.int/news-room/fact-sheets/detail/salmonella-(non-typhoidal)).
4. European Food Safety Authority. 2014. 'The European Union Summary Report on Trends and Sources of Zoonoses, Zoonotic Agents and Food-Borne Outbreaks in 2012.' *EFSA Journal* 12 (2). <https://doi.org/10.2903/J.EFSA.2014.3547>.

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union

	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier
Made in the United Kingdom	Made in the United Kingdom

©2022 Thermo Fisher Scientific Inc. All rights reserved.
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
NCTC and NCTC catalogue marks are a trademark of National Collection of Type Cultures.
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Hampshire,
Basingstoke, RG24 8PW, United Kingdom

For technical assistance please contact your local distributor.

Revision information

Version	Date of modifications introduced
2.0	2022-12-05

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
Bibliographie

1. CDC. 2020. 'Salmonellosis (Nontyphoidal) - Chapter 4 - 2020 Yellow Book | Travelers' Health | CDC.' 2020. <https://wwwnc.cdc.gov/travel/yellowbook/2020/travel-related-infectious-diseases/salmonellosis-nontyphoidal>.
2. Kirk, Martyn D, Sara M Pires, Robert E Black, Marisa Caipo, John A Crump, Brecht Devleesschauwer, Dörte Döpfer, et al. 2015. 'World Health Organization Estimates of the Global and Regional Disease Burden of 22 Foodborne Bacterial, Protozoal, and Viral Diseases, 2010: A Data Synthesis.' <https://doi.org/10.1371/journal.pmed.1001921>.
3. WHO. 2018. '*Salmonella* (Non-Typhoidal).' World Health Organisation. [https://www.who.int/news-room/fact-sheets/detail/salmonella-\(non-typhoidal\)](https://www.who.int/news-room/fact-sheets/detail/salmonella-(non-typhoidal)).
4. European Food Safety Authority. 2014. 'The European Union Summary Report on Trends and Sources of Zoonoses, Zoonotic Agents and Food-Borne Outbreaks in 2012.' *EFSA Journal* 12 (2). <https://doi.org/10.2903/J.EFSA.2014.3547>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Royaume-Uni

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Bibliographie



1. CDC. 2020. 'Salmonellosis (Nontyphoidal)' - Chapter 4 - 2020 Yellow Book | Travelers' Health | CDC.' 2020. <https://wwwnc.cdc.gov/travel/yellowbook/2020/travel-related-infectious-diseases/salmonellosis-nontyphoidal>.
2. Kirk, Martyn D, Sara M Pires, Robert E Black, Marisa Caipo, John A Crump, Brecht Devleesschauwer, Dörte Döpfer, et al. 2015. 'World Health Organization Estimates of the Global and Regional Disease Burden of 22 Foodborne Bacterial, Protozoal, and Viral Diseases, 2010: A Data Synthesis.' <https://doi.org/10.1371/journal.pmed.1001921>.
3. WHO. 2018. 'Salmonella (Non-Typhoidal).' World Health Organisation. [https://www.who.int/news-room/fact-sheets/detail/salmonella-\(non-typhoidal\)](https://www.who.int/news-room/fact-sheets/detail/salmonella-(non-typhoidal)).
4. European Food Safety Authority. 2014. 'The European Union Summary Report on Trends and Sources of Zoonoses, Zoonotic Agents and Food-Borne Outbreaks in 2012.' *EFSA Journal* 12 (2). <https://doi.org/10.2903/J.EFSA.2014.3547>.

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Oxoid Limited, Wade Road, Basingstoke,
 Hampshire, RG24 8PW, UK

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Bibliografia

1. CDC. 2020. 'Salmonellosis (Nontyphoidal)' - Chapter 4 - 2020 Yellow Book | Travelers' Health | CDC.' 2020. <https://wwwnc.cdc.gov/travel/yellowbook/2020/travel-related-infectious-diseases/salmonellosis-nontyphoidal>.
2. Kirk, Martyn D, Sara M Pires, Robert E Black, Marisa Caipo, John A Crump, Brecht Devleesschauwer, Dörte Döpfer, et al. 2015. 'World Health Organization Estimates of the Global and Regional Disease Burden of 22 Foodborne Bacterial, Protozoal, and Viral Diseases, 2010: A Data Synthesis.' <https://doi.org/10.1371/journal.pmed.1001921>.
3. WHO. 2018. 'Salmonella (Non-Typhoidal).' World Health Organisation. [https://www.who.int/news-room/fact-sheets/detail/salmonella-\(non-typhoidal\)](https://www.who.int/news-room/fact-sheets/detail/salmonella-(non-typhoidal)).
4. European Food Safety Authority. 2014. 'The European Union Summary Report on Trends and Sources of Zoonoses, Zoonotic Agents and Food-Borne Outbreaks in 2012.' EFSA Journal 12 (2). <https://doi.org/10.2903/J.EFSA.2014.3547>.

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Oxoid Limited, Wade Road, Basingstoke,
Hampshire, RG24 8PW, Reino Unido

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© 2022 Thermo Fisher Scientific Inc. Todos los derechos reservados.
Todas las demás marcas comerciales son propiedad de Thermo Fisher Scientific Inc. y sus filiales.